

Carbon Fiber AFOs

- Allard USA 2-4
- Euro International 4
- Fillauer 5
- Orthomerica® 6
- Ossur Americas® 6
- Spinal Solutions 7
- Trulife 7-8

Prefabricated Polypropylene AFOs

- AliMed® 9, 11
- Orthomerica 9
- Ossur Americas 10, 12
- ProCare (DJO) 10
- RCAI® 10
- Trulife 10-11

Static Positioning AFOs

- CoreLINE™ 12
- Anatomical Concepts 13-15
- DeRoyal® 15
- LEEDer Group 16
- OCSI/Neuroflex 16-17
- Orthomerica 17
- ProCare (DJO) 17
- RCAI® 18-20

Walking Boots

- Non-Pneumatic with Plastic Uprights 21-22
- Pneumatic with Plastic Uprights 22-25
- Non-Pneumatic with Metal Uprights 25-28
- Pneumatic with Metal Uprights 28-30
- Walking Boot Compensator 30
- Diabetic Walking Boots 31

Dorsal Plantar Fasciitis Splints

- Active Innovations 32
- AliMed 32
- Darco® 32
- Med Spec 33
- Ossur Americas 33

Plantar Fasciitis/Posterior Night Splint

- CoreLINE 34
- Apex 34
- AliMed 34-35
- Corflex 36
- Orthomerica 36
- Ossur Americas 36
- ProCare (DJO) 36-37
- RCAI 37

Hinged Ankle Splints

- Active Innovations 37
- Breg® 38
- Ossur Americas 38
- Swede-O 38

Air/Gel Ankle Stirrups

- Air Ankle Stirrups 39-41
- Gel/Foam Ankle Stirrups 41-42

Multiligamentous Ankle Supports

- Ankle Braces 42-49

Soft Ankle Wraps/Straps

- Aircast (DJO) 50
- Bauerfeind 50-51
- Bledsoe Brace™ 51
- Medi USA 52
- ProCare (DJO) 52

Miscellaneous

- AliMed 53
- DeRoyal 53
- DonJoy (DJO) 53
- DuroMed 53
- Ossur Americas 54
- WBC 54

AFO Socks

- Comfort Products 55
- Knit-Rite® 55-56

Functional Electrical Stimulators

- Innovative Neurotronics® 57-58
- Prizm Medical 59

Allard USA - ToeOFF® Standard

- The ToeOFF is the original dynamic response floor reaction ankle foot orthosis
- Suitable for many different indications and patient groups
- Keeps the foot up during swing phase and gives a soft heel strike, stability in stance and good toe-off
- The ToeOFF should always be worn with a soft interface between the leg and the anterior shell

Indications: Foot drop in combination with no spasticity to moderate spasticity; Limb proprioception deficit and mild proximal deficit; Partial foot amputations

Size	Foot plate Length	Height	ToeOFF Standard		ToeOFF w/D Ring	
			Left	Right	Left	Right
X-Small	8 ¹ / ₄ "	14"	283751010	283752010	283701010	283702010
Small	9"	15"	283751011	283752011	283701011	283702011
Medium	9 ⁵ / ₈ "	16"	283751012	283752012	283701012	283702012
Large	10 ⁵ / ₈ "	17"	283751013	283752013	283701013	283702013
X-Large	11 ¹ / ₄ "	17"	283751014	283752014	283701014	283702014



ToeOFF Short

- The ToeOFF Short is specially developed to accommodate individuals with shorter stature slightly broader-shaped tibia, and wider foot at the metatarsal area

Indications: Foot drop in combination with no spasticity to moderate spasticity; Limb proprioception deficit and mild proximal deficit

Size	Foot plate Length	Height	ToeOff Short	
			Left	Right
Small	8 ¹ / ₂ "	11 ⁵ / ₈ "	283801011	283802011
Medium	9 ¹ / ₄ "	12 ⁵ / ₈ "	283801012	283802012
Large	10"	13 ³ / ₈ "	283801013	283802013



TOEOFF ACCESSORIES

ComfortKIT

- Interface for ToeOFF and BlueRocker
- Designed for the patient that requires added cushioning and a more intimate fit due to the shape of the tibia, skin condition or activity level

Size	Product No.
X-Small	283980010
Small	283980011
Medium	283980012
Large/X-Large	283980013



SoftKIT™ Interfaces

- Every ToeOFF and BlueRocker should have a soft interface between the orthosis and the tibial crest
- SoftKIT is a padding system that includes two narrow neoprene pads to create a pressure-free channel to help relieve any potential discomfort from bony prominences
- The kit includes pre-cut pressure-sensitive MikroFIX tabs to easily affix the pad to ToeOFF, BlueRocker or ToeOFF Short
- Extra pad included



Size	Product No.	
	ToeOFF/BlueRocker	ToeOFF Short
X-Small	283900010	—
Small	283900011	283910011
Medium	283900012	283910012
Large/X-Large	283900013	283910013

ToeOff T-Strap

- Pulls the talocrual and subtalar joints laterally to position them back underneath the tibia
- Distal end attaches to bottom of foot plate – proximal “T” wraps around ankle and attaches to lateral strut
- Easily adjusted to customize to height of medial malleolus
- Plastazote inner core may be gently heated with dry heat to improve shape and fit

Indications: Severe pes valgus that cannot be controlled from ground up with foot orthotics or patient’s coming out of traditional plastic AFO’s



Size	Product No.	
	Left	Right
Universal	287771000	287772000

Foot Orthotic Shells (F.O.S.)

- Proper foot positioning is crucial to optimize biomechanical function and maximize ToeOFF product durability
- Fabricated with EFM (Extra Firm Material) and Microcell Puff Interface, both Ethylene Vinyl Acetate (EVA) foam
- EVA shock absorber is heat moldable and has shape retention and abrasion resistance

Indications: Used to optimize foot biomechanical function for mild to moderate pes planus/cavus



Size	Length	Ypsilon/ToeOFF/BlueRocker Foot plate Length	Product No.
X-Small	8 ³ / ₄ ”	8 ¹ / ₄ ”	680000010
Small	9 ¹ / ₂ ”	9”	680000011
Medium	10”	9 ⁵ / ₈ ”	680000012
Large	10 ³ / ₄ ”	10 ⁵ / ₈ ”	680000013
X-Large	11 ¹ / ₂ ”	11 ¹ / ₄ ”	680000014
2X-Large	12”	—	680000015

BlueRocker™

- The BlueRocker is identical in shape and design as the ToeOFF but offering more rigid orthotic control
- The extra stability will improve both balance and posture and give the wearer greater security, especially individuals with weak quadricep muscles
- The BlueRocker should always be worn with a soft interface between the leg and the anterior shell

Indications: Foot drop due to neurological disorders such as stroke, MS, CP, Muscular Dystrophy and Post Polio; weakness or impairment in multiple leg groups; Severe ankle instability; Proximal neuromuscular weakness/deficits; Can be used in conjunction with a socket and toe filler for management of partial amputations

Size	Foot plate Length	Height	Product No.	
			Left	Right
X-Small	8 ¹ / ₄ "	14"	284051010	284052010
Small	9"	15"	284051011	284052011
Medium	9 ⁵ / ₈ "	16"	284051012	284052012
Large	10 ⁵ / ₈ "	17"	284051013	284052013
X-Large	11 ¹ / ₄ "	17"	284051014	284052014

Ypsilon™

- Designed for the stable ankle to provide dynamic toe-off assistance while allowing natural ankle movements
- Provides opportunity for functional or potentially functional muscles, tendons, and ligaments to function and strengthen
- Individuals with limited energy or physiological endurance may require less energy to overcome ground reaction forces and obtain a smoother gait

Size	Foot plate Length	Height	Ypsilon Model	
			Left	Right
Small	9"	13"	284001011	284002011
Medium	9 ⁵ / ₈ "	13 ³ / ₈ "	284001012	284002012
Large	10 ⁵ / ₈ "	13 ³ / ₄ "	284001013	284002013

Euro International - PeroMax AFO

- The flexible heel and forefoot area allow the best possible physiologic and dynamic gait pattern for the patient
- Y-shaped calf collar allows patients with a bigger calf circumference to wear the AFO comfortably
- Weight limit up to 275 lbs
- Includes 2 calf-pads

Indications: Peroneal paralysis; Weakness of foot lever function; Ankle joint arthritis; Postoperative stabilization after Achilles tendon rupture

Size	Height	US Shoe Size		Product No.	
		Men's	Women's	Left	Right
1	15"	6-7	7 ¹ / ₂ -8 ¹ / ₂	197P4/L1	197P4/R1
2	15 ¹ / ₂ "	7-8 ¹ / ₂	8 ¹ / ₂ -10	197P4/L2	197P4/R2
3	15 ³ / ₄ "	8 ¹ / ₂ -11 ¹ / ₂	10-13	197P4/L3	197P4/R3
4	16"	11 ¹ / ₂ -14	13-15	197P4/L4	197P4/R4



BlueRocker



Ypsilon



PeroMax AFO

Fillauer - Dynamic Walk Dorsiflexion Assist AFO

- Utilizes medial and lateral PEEK rods for active motion
- Provides dorsi-assist with subtalar inversion and eversion
- Fits easily into regular shoes
- Flexible for multi-plane action
- Enables athletic activity such as moderate running

Size	Product No.				Padding Kit Product No.
	Regular Ankle		Wide Ankle		
	Left	Right	Left	Right	
Small	700111121	700111122	700111321	700111322	700100220
Medium	700111141	700111142	700111341	700111342	700100240
Large	700111161	700111162	700111361	700111362	700100260

Sizing Chart						
Size	US Shoe Sizes		Length	Height	Ankle Widths	
	Men's	Women's			Regular	Wide
Small	4-6 ¹ / ₂	5 ¹ / ₂ -8	9.1"	12.6"	2.7"	2.9"
Medium	7-9	8 ¹ / ₂ -10 ¹ / ₂	9.8"	13.4"	2.9"	3.2"
Large	10 ¹ / ₂ -14	12-15 ¹ / ₂	10.6"	14.6"	3.2"	3.5"



Dynamic Walk Single Side AFO

- Made of lightweight moldable carbon fiber composite material
- Easy donning and doffing
- Fits easily into regular or narrow shoes
- Choose between PEEK rods on medial or lateral side
- Provides dorsi-assist with subtalar inversion and eversion
- Flexible for multi-plane action
- Enables athletic activity such as moderate running

Size	Lateral Product No.		Medial Product No.		Padding Kit Product No.
	Left	Right	Left	Right	
Small	700114121	700114122	700113121	700113122	700100220
Medium	700114141	700114142	700113141	700113142	700100240
Large	700114161	700114162	700113161	700113162	700100260

Sizing Chart				
Size	US Shoe Sizes		Length	Height
	Men's	Women's		
Small	4-6 ¹ / ₂	5 ¹ / ₂ -8	9.1"	12.6"
Medium	7-9	8 ¹ / ₂ -10 ¹ / ₂	9.8"	13.4"
Large	10 ¹ / ₂ -14	12-15 ¹ / ₂	10.6"	14.6"



Orthomerica - Carbon Supra-Lite™ AFO

- Low-profile, dorsiflexion assist design
- Increased rigidity in vertical aspect
- Weight limit up to 220 lbs

Contraindicated: Severe ankle-foot deformities; Severe spasticity; Fluctuation edema; M-L instability; Extreme activity



Size	US Shoe Size		Brace Height	Product No.	
	Men's	Women's		Left	Right
Small	5-7	up to 6	10 ¹ / ₂ "	3520	3525
Medium	7 ¹ / ₂ -10	6 ¹ / ₂ +	11 ¹ / ₂ "	3521	3526
Large	10 ¹ / ₂ -13 ¹ / ₂	—	13 ¹ / ₂ "	3522	3527
X-Large	13	—	15 ¹ / ₂ "	3523	3528

Ossur Americas - AFO Dynamic

- Carbon fiber design provides excellent strength-to-weight ratio and dynamic motion
- The tapered carbon heel stores energy and manages forces at heel strike
- Full-length toe lever provides full support while foot is loaded, as well as reduced shock to the contralateral limb
- Weight limit up to 265 lbs

Indications: Drop foot secondary to CVA; MS, Charcot Marie Tooth disease; Other neurological conditions; Mild knee instability



Size	US Shoe Size		Foot Size	Foot plate	Calf Height*	Heel Height	Weight	Product No.	
	Men's	Women's						Left	Right
X-Small	3-5	4-6 ¹ / ₂	8"-9"	8 ³ / ₄ "	12"-13 ¹ / ₂ "	0.4"	2.7 oz	O-102021	O-102011
Small	5 ¹ / ₂ -7 ¹ / ₂	7-9	9"-10"	9 ¹ / ₂ "	12 ¹ / ₂ "-14 ¹ / ₄ "	0.6"	3.5 oz	O-102022	O-102012
Medium	8-10 ¹ / ₂	9 ¹ / ₂ -12	10"-10 ³ / ₄ "	10 ¹ / ₄ "	13 ¹ / ₄ "-15"	0.6"	4.2 oz	O-102024	O-102014
Large	11-14	—	10 ³ / ₄ "-11 ¹ / ₂ "	11"	14 ¹ / ₄ "-15 ³ / ₄ "	0.6"	5 oz	O-102026	O-102016

*the height at the thickest part of the calf

AFO Light

- Energy return and propulsion assistance of the limb at toe off
- Carbon fiber design provides excellent strength-to-weight ratio and dynamic motion
- Stiffness gradient at the heel and the toe correspond to ground reaction forces as they are applied to the foot
- Weight limit up to 220 lbs

Indications: Mild to moderate drop foot; Severe ankle-foot deformities; Severe spasticity; Fluctuation edema; M-L instability, or extreme activity



Size	US Shoe Size		Foot Size	Foot Plate	Calf Height*	Weight	Product No.	
	Men's	Women's					Left	Right
Small	up to 7 ¹ / ₂	7-9	8 ³ / ₄ "-10 ¹ / ₄ "	9 ¹ / ₂ "	12 ¹ / ₂ "-14 ¹ / ₄ "	2.8 oz	O-121022	O-121012
Medium	8-10 ¹ / ₂	9 ¹ / ₂ -12	9 ¹ / ₂ "-11"	10 ¹ / ₄ "	13 ¹ / ₂ "-15"	3.15 oz	O-121024	O-121014
Large	11-14	—	11"	11"	14 ¹ / ₄ "-15 ³ / ₄ "	3.5 oz	O-121026	O-121016
X-Large	12-15	—	11"-12 ¹ / ₂ "	11 ³ / ₄ "	15"-16 ¹ / ₂ "	3.85 oz	O-121028	O-121018

*the height at the thickest part of the calf

Spinal Solutions - EZ Stride

- True floor reaction; minimizes contact on the anterior tibia
- Contact at proximal tibia and foot plate only
- Flexible calf band conforms to patients
- Malleolus relief is anatomically positioned
- Lateral strut is fabricated to maintain 1/4" offset thereby eliminating potential contact with patient

Size	Brace Height	US Shoe Size		Product No.		Replacement Liners
		Men's	Women's	Left	Right	
X-Small	13"	3 ¹ / ₂ -7	5-8 ¹ / ₂	EZ-XSL-LAT	EZ-XSR-LAT	EZ-XS-FL
Small	14"	7 ¹ / ₂ -10	9-12 ¹ / ₂	EZ-SL-LAT	EZ-SR-LAT	EZ-S-FL
Medium	15"	10 ¹ / ₂ -11 ¹ / ₂	12-13	EZ-ML-LAT	EZ-MR-LAT	EZ-M-FL
Large	16"	11 ¹ / ₂ +	13+	EZ-LL-LAT	EZ-LR-LAT	EZ-L-FL



Trulife - Matrix/Matrix Max

- Energy storing, composite foot plate and strut
- Lateral strut is anterior to malleolus
- Open calcaneal design
- Custom-fit, anterior shell and foot plate
- Matrix Max is lightweight, carbon composite frame with enhanced stability

Indications: Matrix: Low activity patients with simple, unilateral drop foot; Maximum weight of 250 lb; **Matrix Max:** Low to high activity patients with unilateral or bilateral drop foot; Suitable for partial foot amputations

Size	Height	Foot Measurement	Product No.	
			Left	Right
Matrix				
Small	15 ¹ / ₂ "-16 ¹ / ₂ "	up to 9 ⁴ / ₅ "	FH200L02	FH200R02
Medium	16"-17"	9 ² / ₅ "-10 ³ / ₅ "	FH200L03	FH200R03
Large	17"-18"	9 ⁴ / ₅ "-11 ² / ₅ "	FH200L04	FH200R04
Matrix Max				
Small	15 ¹ / ₂ "-16 ¹ / ₂ "	up to 9 ⁴ / ₅ "	FH250L02	FH250R02
Medium	16"-17"	9 ² / ₅ "-10 ³ / ₅ "	FH250L03	FH250R03
Large	17"-18"	9 ⁴ / ₅ "-11 ² / ₅ "	FH250L04	FH250R04
X-Large	17"-18"	10 ¹ / ₅ "-12 ¹ / ₅ "	FH250L05	FH250R05
Accessories for Matrix and Matrix Max				
Description			Product No.	
Replacement Strap Set			FH250B	
Small Replacement Anterior Shell			FH250A02	
Medium Replacement Anterior Shell			FH250A03	
Large Replacement Anterior Shell			FH250A04	
Small Replacement Padding			FH250P02	
Medium Replacement Padding			FH250P03	
Large Replacement Padding			FH250P04	



Matrix



Matrix Max

Matrix Max Raw

- Dynamic foot plate and strut
- Lightweight carbon composite
- Ready for practitioner to create a custom AFO, KAFO or partial foot prosthesis
- 10 pack contains 5 left and 5 right

Indications: Drop foot, CMT, MS, ALS, Off-loading and PTB

Description	Height	Foot Length	Product No.	
			Left	Right
Matrix Max Raw, single	14 ¹ / ₂ "	11 ¹ / ₂ "	FH225-L	FH225-R
Matrix Max Raw, 10/pack	—	—	FH226	



Matrix/Matrix Max Split Toe (ST)

- Features a revolutionary split foot plate to improve patient balance and increase stability
- Foot plate acts like a natural spring
- Height-adjustable, anterior shell

Indications: Matrix Split Toe: Hallux amputations; Low activity patients with simple, unilateral drop foot; Matrix Max Split Toe: Hallux amputations; Low to high activity patients with unilateral or bilateral drop foot

Size	Height	Foot Length	Product No.	
			Left	Right
Matrix Split Toe				
Small	15 ¹ / ₂ "-16 ¹ / ₂ "	up to 9 ⁴ / ₅ "	FH210L02	FH210R02
Medium	16"-17"	9 ² / ₅ "-10 ³ / ₅ "	FH210L03	FH210R03
Large	17"-18"	9 ⁴ / ₅ "-11 ² / ₅ "	FH210L04	FH210R04
Matrix Max Split Toe				
Small	15 ¹ / ₂ "-16 ¹ / ₂ "	up to 9 ⁴ / ₅ "	FH215L02	FH215R02
Medium	16"-17"	9 ² / ₅ "-10 ³ / ₅ "	FH215L03	FH215R03
Large	17"-18"	9 ⁴ / ₅ "-11 ² / ₅ "	FH215L04	FH215R04
X-Large	17"-18"	10 ¹ / ₅ "-12 ¹ / ₅ "	FH215L05	FH215R05



Matrix Split Toe



Matrix Max Split Toe

Matrix Spiral

- Spiral spring resists plantar flexion and allows gradual dorsiflexion
- Lateral strut is anterior to malleolus for more comfortable fit
- Anterior band allows efficient loading of the spring and enhanced function
- Foot plate flexibility and toe spring allow for smooth transition to toe off

Indications: High activity patients with foot drop as a result of CVA, MS, Post-Polio, CMT or Peroneal Nerve Injury

Size	Height	Foot plate Length	Product No.	
			Left	Right
Small	14 ¹ / ₂ "	9 ⁴ / ₅ "	FH280L02	FH280R02
Medium	15 ⁴ / ₅ "	9 ¹ / ₂ "-10 ³ / ₅ "	FH280-03	FH280R03
Large	16 ¹ / ₂ "	9 ⁴ / ₅ "-11 ² / ₅ "	FH280L04	FH280R04



Prefabricated Polypropylene AFO's

- Provides dorsiflexion assist
- Heat moldable for a custom fit

Alimed	Freedom Swedish AFO						
	<ul style="list-style-type: none"> • Injection-molded polyethylene splint • Provides static dorsiflexion assistance and lateral stability for the entire foot-ankle area 						
			<u>White Product No.</u>		<u>Black Product No.</u>		
	Gender	US Shoe Size	Left	Right	Left	Right	
Men's	up to 12	64331	64330	64335	64334		
Women's	up to 9 ¹ / ₂	64329	64328	64333	64332		
Alimed	C-90 Superior Posterior Leaf Spring						
	<ul style="list-style-type: none"> • Injection-molded for durability • Provides full range of plantar and dorsiflexion 						
		<u>US Shoe Size</u>		Height	<u>Product No.</u>		
	Size	Men's	Women's		Left	Right	
	Small	4-6	6-8	11 ¹ / ₂ "	6603	6602	
Medium	7-8	8-10	12 ¹ / ₄ "	6605	6604		
Large	9-12	10-11	13 ³ / ₄ "	6607	6606		
Orthomerica	AFO						
	<ul style="list-style-type: none"> • Posterior ribs for increased support 						
		<u>US Shoe Size</u>		Height	<u>Product No.</u>		
	Size	Men's	Women's		Left	Right	
	Small	up to 6	—	12"	3540	3545	
	Medium	5-7	up to 6	13 ¹ / ₂ "	3541	3546	
	Large	7 ¹ / ₂ -10	6 ¹ / ₂ +	14"	3542	3547	
	X-Large	10 ¹ / ₂ -13	—	14"	3543	3548	
	Leaf Spring Orthosis						
	<ul style="list-style-type: none"> • Low-profile design fits in most normal lace-up shoes 						
		<u>US Shoe Size</u>		Height	<u>Product No.</u>		
	Size	Men's	Women's		Left	Right	
Small	5-7	up to 6	12"	3510	3515		
Medium	7 ¹ / ₂ -10	6 ¹ / ₂ +	13 ¹ / ₂ "	3511	3516		
Large	10 ¹ / ₂ -13	—	13 ¹ / ₂ "	3512	3517		
X-Large	13+	—	14"	3513	3518		
Orthomerica	Supra-Lite AFO						
	<ul style="list-style-type: none"> • Leaf spring design to correct mild to moderate foot drop • Injection molded • Special heel cutout accommodates a wide range of sizes 						
		<u>US Shoe Size</u>		Height	<u>Product No.</u>		
	Size	Men's	Women's		Left	Right	
	Small	5-7	up to 6	10 ¹ / ₂ "	3500	3505	
Medium	7 ¹ / ₂ -10	6 ¹ / ₂ +	11 ¹ / ₂ "	3501	3506		
Large	10 ¹ / ₂ -13	—	13 ¹ / ₂ "	3502	3507		
X-Large	13+	—	15 ¹ / ₂ "	3503	3508		

Prefabricated Polypropylene AFO's (continued)

Prefabricated Polypropylene AFO's (continued)						
Ossur Americas	AFO Leaf Spring					
	<ul style="list-style-type: none"> Absence of heel section for improved comfort 					
	Size		Height	Product No.		
				Left	Right	
Small/Medium		10 ³ / ₄ "	O-110227	O-110127		
Large		12"	O-110230	O-110130		
						
ProCare (DJO)	Posterior Leaf Spring - Mild to moderate foot drop					
	Size	US Shoe Size		Height	Product No.	
		Men's	Women's		Left	Right
	Small	6-8	6 or less	12"	79-97733	79-97743
	Medium	8 ¹ / ₂ -10	6 ¹ / ₂ -8	13 ¹ / ₂ "	79-97735	79-97745
	Large	10 ¹ / ₂ -12	8 ¹ / ₂ -10	13 ¹ / ₂ "	79-97737	79-97747
X-Large	12+	10	14"	79-97738	79-97748	
						
RCAI	Leaf Spring Orthosis					
	<ul style="list-style-type: none"> Provides functional alignment for flaccid extremities and resistance to mild extensor spasticity Set in 8° dorsiflexion to provide spring assist during ambulation 					
	Size	Foot Length		Height	Product No.	
		Patient	Product		Left	Right
	Small	8"-9 ³ / ₄ "	7 ¹ / ₂ "	14"	50LS-AFO-S-L	50LS-AFO-S-R
	Medium	9 ¹ / ₂ "-10"	8"	15"	50LS-AFO-M-L	50LS-AFO-M-R
	Large	10 ¹ / ₂ "-12"	8 ¹ / ₂ "	16"	50LS-AFO-L-L	50LS-AFO-L-R
	X-Large	12"-15"	9"	17"	50LS-AFO-XL-L	50LS-AFO-XL-R
	Rigid AFO					
	<ul style="list-style-type: none"> Provides functional alignment for flaccid extremities and resistance to mild extensor spasticity as well as mediolateral support for ankle instability 					
Size	Foot Length		Height	Product No.		
	Patient	Product		Left	Right	
Small	8"-9 ³ / ₄ "	7"	13"	60RAFO-S-L	60RAFO-S-R	
Medium	9 ¹ / ₂ "-10"	7 ¹ / ₂ "	14"	60RAFO-M-L	60RAFO-M-R	
Large	10 ¹ / ₂ "-12"	8"	15"	60RAFO-L-L	60RAFO-L-R	
X-Large	12"-15"	8 ¹ / ₂ "	15"	60RAFO-XL-L	60RAFO-XL-R	
						
Trulife	Semi-Solid AFO - For moderate to severe foot drop, MS, and CP					
	Size	Calf Circum.	Foot Plate Length	Height	Product No.	
					Left	Right
	Small	12 ⁴ / ₅ "	9 ¹ / ₂ "	12 ³ / ₅ "	115AFO2N316F1L	115AFO2N316F1R
	Medium	14"	10 ¹ / ₄ "	13 ³ / ₄ "	115AFO2N316F2L	115AFO2N316F2R
	Large	15 ¹ / ₄ "	10 ³ / ₄ "	15"	115AFO2N316F3L	115AFO2N316F3R
	Semi-Solid AFO - Metatarsal Head Foot Plate					
	Size	Calf Circum.	Foot Plate Length	Height	Product No.	
					Left	Right
	Small	12 ⁴ / ₅ "	6 ³ / ₅ "	12 ³ / ₅ "	115AFO2N316M1L	115AFO2N316M1R
Medium	14"	7 ¹ / ₁₀ "	13 ³ / ₄ "	115AFO2N316M2L	115AFO2N316M2R	
Large	15 ¹ / ₄ "	7 ³ / ₅ "	15"	115AFO2N316M3L	115AFO2N316M3R	
						

Prefabricated Polypropylene AFO's (continued)							
Trulife	Rigid AFO - For moderate to severe foot drop, MS, and CP						
	Size	Calf Circum.	Foot Plate Length	Height	Product No. Left Right		
	Small	12 ⁷ / ₈ "	9 ¹ / ₂ "	12 ⁵ / ₈ "	115AFO1N316F1L 115AFO1N316F1R		
	Medium	14"	10 ¹ / ₄ "	13 ⁴ / ₅ "	115AFO1N316F2L 115AFO1N316F2R		
	Large	15 ¹ / ₄ "	10 ⁴ / ₅ "	15"	115AFO1N316F3L 115AFO1N316F3R		
	Posterior Leaf Spring - Full Length						
	Size	Calf Circum.	Foot Plate Length	Height	Product No. Left Right		
	Small	12 ⁷ / ₈ "	9 ¹ / ₂ "	12 ⁵ / ₈ "	115AFO5N316F1L 115AFO5N316F1R		
	Medium	14"	10 ³ / ₁₀ "	13 ⁴ / ₅ "	115AFO5N316F2L 115AFO5N316F2R		
	Large	15 ¹ / ₃ "	10 ⁴ / ₅ "	15"	115AFO5N316F3L 115AFO5N316F3R		
	Posterior Leaf Spring with Metatarsal Head Foot Plate						
	Size	Calf Circum.	Foot Plate Length	Height	Product No. Left Right		
	Small	12 ⁷ / ₈ "	6 ³ / ₅ "	12 ⁵ / ₈ "	115AFO5N316M1L 115AFO5N316M1R		
Medium	14"	7 ¹ / ₁₀ "	13 ⁴ / ₅ "	115AFO5N316M2L 115AFO5N316M2R			
Large	15 ¹ / ₃ "	7 ³ / ₅ "	15"	115AFO5N316M3L 115AFO5N316M3R			

AliMed - FREEDOM Adjustable Foot Drop Brace

- Unique closure allows for user to customize level of support
- Quick-release attachment buckle permits easy transfer of the brace from shoe to shoe (one extra shoe strap included)
- Lightweight neoprene cuff and hook and loop closure make this brace easy to apply and adjust for a proper fit
- Brace must be worn with lace-up or hook and loop closure shoes



Description	Circumference*	Product No.
Foot Drop Brace	up to 10 ¹ / ₂ "	66220
Replacement Foot Strap 3/pack	—	66221

*Measure lower leg just above ankle joint at narrowest points

Ossur Americas - Foot-Up™

- Foot-up is a lightweight ankle-foot orthosis that offers dynamic support for drop foot or similar complaints
- Provides visible improvement in gait by supporting the foot the moment it is raised
- Cushioned ankle wrap made from breathable, 3-layer material that can be worn comfortably for long periods of time
- Plastic inlay fits discreetly between the tongue and laces of the shoe and is practically invisible when worn

Indications: Dynamic support for drop foot or similar complaints for which support of dorsiflexion is desirable



Size	Ankle Circumference	Product No.
Medium	7"-8 ¹ / ₄ "	07810-1
Large	8 ¹ / ₂ "-10 ¹ / ₄ "	07810-2
X-Large	10 ¹ / ₂ "-13"	07810-3
Accessories		
Plastic Inlay	Universal Size	17810
Small Shoeless Wrap	up to 7 ⁹ / ₁₀ "	27810-1
Medium Shoeless Wrap	7 ⁹ / ₁₀ "-9 ² / ₅ "	27810-2
Large Shoeless Wrap	9 ² / ₅ "-11"	27810-3
X-Large Shoeless Wrap	11"+	27810-4

STATIC POSITIONING AFOs

CoreLINE Static Positioning AFO

- Designed with wider toe post than other commercially available models to address a greater range of shapes and deformities.
- This design increases patient comfort and prevents edge pressures
- Kydex plastic shell allows re-forming and increases breakage resistance
- Universal sizing, without straps ensures easy fit and lower inventory

Indications: Prevents heel ulcers; Prevents contractures; Recumbent positioning



Fleece

Velfoam

Size	Calf Circum.	Product No.	
		Fleece	Velfoam
Regular	up to 16"	CL0330	CL0331
Large	16"+	CL0335	CL0336

**Anatomical Concepts - Bariatric APU®
Articulating Ankle/Foot Orthosis**

- Extends from 12¹/₂" to 13³/₄" in length and from 14" to 15" in height (can be trimmed shorter if needed)
- Heavy duty ³/₁₆" Aluminum APU® (Adjustable Posterior Upright) heel suspension bar is standard
- The adjustable calf strap extends from 21" to 32"
- The Kodel® liner will fit across a dorsum measurement up to 23"
- The width of the foot plate is 5"

Description	Foot Length	Calf Height	Calf Circum.	Product No.
Bariatric APU - Teal Kodel	9 ¹ / ₄ "-13 ¹ / ₄ "	12 ³ / ₄ "-13 ³ / ₄ "	26"	752SKT
Liner - Teal	—	—	—	L70SKT



Heel strap positioned behind the connecting bar (for display purposes only)

650APU® Adjustable Posterior Upright

- Adjustable length lessens the need for multiple sizes
- Features a distinctive and protective black powder coated metal upright and single jointed metal upright that provides a ROM of 30° plantar flexion and 30° of dorsiflexion
- Permanent alignment marks are in 5° intervals to facilitate exact prescribed adjustments

Size	Description	Product No.	
		White Liner	Black Liner
Adult	Adjustable Posterior Upright/full Kodel Liner	652SKG	652SKB
Adult	Replacement Kodel Liner	L62SKG	L62SKB

Sizing Chart for 650APU			
Size	Foot Length	Calf Height	Calf Circum.
Adult	6 ³ / ₄ "-12"	12"-13"	20"

652SKB



Dynamic Dorsi-Assist

- Dynamically stretch the gastroc muscles, offsetting equinus contracture
- Allows for virtually unrestricted range of plantar or dorsiflexion and inversion or eversion of the ankle/foot complex
- Two elastic straps attached at the distal foot section and proximal calf section can be adjusted to compensate and/or offset varus and valgus tendencies
- Maintains constant dorsi-assistive force and can be used for ambulation

Size	Foot Length	Foot Plate	Calf Height	Calf Circum.	Product No.	
					Full Kodel	Pad/Strap
Infant	3 ¹ / ₂ "-6 ¹ / ₂ "	2 ³ / ₄ "	6 ³ / ₈ "-7 ³ / ₈ "	13 ¹ / ₂ "	454RKG	—
Pediatric	4 ³ / ₄ "-8 ¹ / ₂ "	3 ¹ / ₂ "	9"-10"	15"	554SKG	—
Adult	6 ³ / ₄ "-12"	3 ³ / ₄ "	12"-13'	20"	654SKG	654PKG



PRAFO®

- Interchangeable custom-contoured heel connector bar helps control plantar and dorsiflexion by allowing infinite, measurable adjustments
- Calf and foot segments can be custom modified
- Adjustable anti-rotation bar helps unwanted extremity rotation
- Liner can easily be changed without tools

Description	Product No.		
	Adult	Infant	Pediatric
PRAFO with full Kodel liner	650SKT	450RKG	550SKG
PRAFO with Sheepskin liner	650SFY	—	—
PRAFO Pad and Strap style	650PKG	—	—
PRAFO with heavy duty bar	680SKT	—	—
PRAFO with stainless steel bar	690SKT	—	—
PRAFO with wrap over liner	650RKG	—	—
PRAFO Accessories			
Replacement Kodel liner	L60SKT	L40SKG	L50SKG
Replacement Sheepskin liner	L60SFY	—	—
Replacement Pad and Strap liner	L60PKG	—	—
Replacement Calf Strap	C60SOW	C40SOB	C50SOW

Sizing Chart for PRAFO			
Size	Foot Length	Calf Height	Calf Circum.
Infant	3 ¹ / ₂ "-6 ¹ / ₂ "	6 ³ / ₈ "-7 ³ / ₈ "	13 ¹ / ₂ "
Pediatric	4 ³ / ₄ "-8 ¹ / ₂ "	9"-10"	15"
Adult	6 ³ / ₄ "-12"	12"-13"	20"

PRAFO with full Kodel liner



PRAFO Pad and Strap style



650EV PRAFO® EV, with Full Kodel® Liner

- Offers positive heel suspension; accommodates and helps control varus and valgus conditions of the ankle/foot complex
- Calibrated settings on the bar allow for a more accurate documentation of patient's progress
- Adjustable anti-rotation bar helps to control unwanted extremity rotation

Size	Description	Product No.
Infant	With full Kodel Bootie	451RKG
Pediatric	With full Kodel Bootie	551SKG
Adult	With full Kodel Bootie	651SKT

Sizing Chart for 650EV PRAFO			
Size	Foot Length	Calf Height	Calf Circum.
Infant	3 ¹ / ₂ "-6 ¹ / ₂ "	6 ³ / ₈ "-7 ³ / ₈ "	13 ¹ / ₂ "
Pediatric	4 ³ / ₄ "-8 ¹ / ₂ "	9"-10"	15"
Adult	6 ³ / ₄ "-12"	12"-13"	20"



PKA - PRAFO® KAFO Attachment

- Reduces patient visits for adjustments to the ankle/foot segments of KAFOs, and improves patient compliance
- Bar and PKA system stops migration and rotation of your KAFO with practically no worry of skin impingement
- Foot lengths: Adult - 6³/₄"-12"; Pediatric - 4³/₄"-8¹/₂"; Infant - 3¹/₂"-6¹/₂"

Description	Product No.		
	Adult	Pediatric	Infant
PRAFO KAFO Attach	250SKT	150SKG	050SKG
PRAFO KAFO Attach/EV	251SKT	151SKG	—
PRAFO KAFO Attach/APU	252SKG	—	—
PRAFO KAFO Attach/APU/PS	252PKG	—	—
PRAFO KAFO Attach/DDA	254SKG	—	—



252PK6

DeRoyal - Ankle Contracture Boot

- Static AFO used for positioning and pressure reduction
- Structured to completely suspend the heel
- “No strap” closures make it easy for patients to apply and remove the boot
- Anti-rotation bar, padded toe piece and mesh laundry bag included
- Universal fit for left and right foot

Indications: Hip rotation; Foot drop; Not recommended for severe plantar flexion contractures, open wounds, severe spasticity or thrombosis

Description	Material	Product No.			
		Small	Medium	Large	X-Large
Ankle Contracture Boot	Velfoam	4302B	4302C	4302D	4302E
	Fleece	4307B	4307C	4307D	4307E
Ankle Contracture Boot with Boot Sole	Velfoam	4301B	4301C	4301D	4301E
	Fleece	4306B	4306C	4306D	4306E
Ankle Contracture Boot with Boot Sole and Extra Liner	Velfoam	4300B	4300C	4300D	4300E
	Fleece	4305B	4305C	4305D	4305E



Vel-Foam



Fleece

Ankle Contracture Boot Sizing Chart

Size	Calf Circumference (Measure 8" up from bottom of heel)	Foot Circumference (Measure around widest part of foot)
Small (B)	10"-14"	7"-9"
Medium (C)	14"-16"	9"-11"
Large (D)	16"-21"	11"-13"
X-Large (E)	16"-21"	13"+

LEEDer Group - LG-Kydex-Pro Ankle Foot Orthosis

- Professional MP Boot: Ambulation or Bed
- Easy application: wide closures for less pressure points (no straps)
- Adjustable anti-rotation bar prevents hip rotation and adjustable toe post

Indications: Non-fixed plantar flexion contracture; Plantar fasciitis; Equinovarus contracture limitation; Heel float; Foot drop; Achilles tendonitis; Heel decubitus

Size	Child Shoe Size	Calf Circum.	Product No.	
			Fleece	Velfoam
Pediatric	12-6	—	LBMP	LBMPNF
Regular	—	up to 16"	LBMR	LBMRNF
Large	—	16"-20"	LBML	LBMLNF



OCSI/Neuroflex - SoftPro™ Ambulating™ AFO

- Bacti-Ban™ keeps covers cleaner and more odor-resistant
- Rocker bottom design provides initial weight bearing and assisted gait training
- Non-removable ClickStep™ transfer pad provides audible cue with each heel strike
- Anti-rotation bar

Indications: Plantar Flexion

Size	Length of Foot	Max. Calf Circum.	Product No.	
			Smooth	Fleece
Small	6"-8"	15"	SP20-50	SP20-00
Medium	8"-10"	17"	SP20-51	SP20-01
Large	10"-11"	19"	SP20-52	SP20-02
X-Large	11"-12"	23"	SP20-53	SP20-03



SoftPro™ Gait Trainer™ AFO

- Removable ClickStep™ Transfer Pad provides audible cue with heel strike
- Rocker bottom design facilitates initial weight bearing and assisted gait training

Indications: Plantar flexion; Inversion; Eversion

Size	Length of Foot	Product No.	
		Smooth	Fleece
Small	6"-8"	SP30-50	SP30-00
Medium	8"-10"	SP30-51	SP30-01
Large	10"-11"	SP30-52	SP30-02
X-Large	11"-12"	SP30-53	SP30-03



SoftPro™ In-Bed AFO

- Effectively treats plantar flexion
- Non-skid bottom for assisted transfers

Size	Length of Foot	Max. Calf Circum.	Product No.	
			Smooth	Fleece
Small	6"-8"	15"	SP10-50	SP10-00
Medium	8"-10"	17"	SP10-51	SP10-01
Large	10"-11"	19"	SP10-52	SP10-02
X-Large	11"-12"	23"	SP10-53	SP10-03



Orthomerica - HealPRO™

- HealPRO's design reduces pressure on the toes from bed linens, increasing comfort and skin protection
- Reduced calf pressure due to proximal trim line of Coolfoam wrap
- Reduced heel pressure in supine position due to heel offset
- Internal or external rotation control with 90 degree stop

Size	Measured Foot Length	Maximum Calf Circum.	US Shoe Size		Coolfoam Lined	Fleece Lined
			Men's	Women's		
Small	8"-9"	16"	3-5	3-6	3530	3530F
Medium	9"-10"	16"	5-8	6-9	3531	3531F
Large	10"-11"	18"	8-11	9-12	3532	3532F
X-Large	11"-12"	20"	11-14	12-15	3533	3533F



ProCare (DJO) - Podous Boot

- Latex free; lightweight, quality plastic shell with poly/pile liner
- Adjustable toe extension helps protect toes from bedding
- Anti-rotation bar may be repositioned to assist with internal or external rotation
- Fits left or right

Indications: Helps treat and prevent lower extremity disorders associated with trauma or immobility including pressure necrosis, ankle contractures and post-operative hip abduction or rotation

Size	US Shoe Size		Product No.
	Men's	Women's	
Regular	up to 10	up to 11	79-90550
X-Large	10 ¹ / ₂ +	11 ¹ / ₂ +	79-90551
Liners 3/pack			
Regular	—	—	79-90560
X-Large	—	—	79-90561



RCAI - Corrxit Foam Foot Orthosis with Ambulatory Attachment

- The heel design of the splint eliminates pressure points and allows for visual examination and circulation of air over the heel
- The Toe Guard reduces bed linen skin irritation by separating the bedding from contact with the foot
- The Stabilizer Bar effectively prevents rolling of the patient’s leg and controls internal or external rotation
- The low profile Ambulatory Attachment (shown) is easily applied and removed without the use of tools



Size	Measure Heel to Great Toe	Circumference		Product No.
		Max. Foot	Max. Calf	
Average	8 ¹ / ₂ -9 ³ / ₄ "	13"	17"	03CX-A
Large	9 ³ / ₄ -11"	13"	17"	03CX-L
Bariatric	9 ³ / ₄ -11"	16"	20"	03CX-B

Dorsiflexion Multi Podus

- Provides adjustable static stretch of the plantar flexors
- Medial/lateral straps attached to the sides of the brace allow the foot and ankle to be controlled

Indications: Diminished foot and ankle range-of-motion; Pressure sores: Hip rotation



Size	Measure Foot Length	Product No.
Pediatric	4 ¹ / ₂ "-6"	11DFMP-P
Youth	6"-7"	11DFMP-Y
Small	7"-8"	11DFMP-S
Medium	8"-9"	11DFMP-M
Large	9"-10"	11DFMP-L
X-Large	10"-11"	11DFMP-XL

Phase II MultiPodus® System

- Adjustable positioning from 40° plantar flexion to 10° dorsiflexion in infinite increments, provides controlled static stretch correcting plantar flexion contractures
- Swivel rotator bar controls hip and leg rotation
- Inversion/eversion strap supplied simulates a T-strap on a short leg brace



Size	Measure Foot Length	Max. Calf Circum.	Product No.	
			Kodel	Breathable Foam
Pediatric	4 ¹ / ₂ "-6"	10"	20MP-A-P	21MPO-A-P
Youth	6"-7"	13"	20MP-A-Y	21MPO-A-Y
Small	7"-8"	16"	20MP-A-S	21MPO-A-S
Average	8"-9 ¹ / ₂ "	16"	20MP-A-A	21MPO-A-A
Large	9 ¹ / ₂ "-11"	20"	20MP-A-L	21MPO-A-L
X-Large	11"-12 ³ / ₄ "	20"	20MP-A-XL	21MPO-A-XL

**The Original Multi Podus[®] System
(with transfer attachment)**

- Systematic treatment of foot drop, hip rotation and pressure sores
- Dynamic flex action supplies continuous counter force pressure to the plantar surface of the foot
- Swivel rotator bar positioned to the side, controls hip and leg rotation
- Adjustable toe post relieves pressure on the toes, and can be positioned to the side in the treatment of lateral pressure sores



Size	Measure Foot Length	Max. Calf Circum.	Product No.
Toddler*	4"-5"	10"	10MP-T
Child*	5"-6"	10"	10MP-C
Youth*	6"-7"	13"	10MP-Y
Small	7"-8"	16"	10MP-S
Medium	8"-9"	18"	10MP-M
Large	9"-10"	20"	10MP-L
X-Large	10"-11"	22"	10MP-XL
2X-Large*	11"-12 ¹ / ₂ "	24"	10MP-XXL

*Transfer attachment not available

MPO 2000[®] (with transfer attachment)

- Features a breathable anti-microbial liner that wicks away moisture from the skin
- Includes a transfer attachment providing a nonskid standing surface
- Easily applied and removed with hook and loop straps
- Rotator bar controls hip and leg rotation



Size	Measure Foot Length	Max. Calf Circum.	Product No.
Infant*	3"-4"	8"	01MP-I
Toddler*	4"-5"	10"	01MP-T
Child*	5"-6"	10"	01MP-C
Youth*	6"-7"	13"	01MP-Y
Small	7"-8"	15"	01MP-S
Medium	8"-9"	17"	01MP-M
Large	9"-10"	17"	01MP-L
X-Large	10"-11"	17"	01MP-XL
2X-Large*	11"-12 ¹ / ₂ "	17"	01MP-XXL

*Transfer attachment not available

MPO Active 2000® (with ambulatory attachment)

- Includes a removable ambulatory attachment that converts from a bed positioning device to a limited ambulation device secured with a simple “twist and lock” motion
- Allows normal plantar and dorsiflexion and acts as a spring assist

Size	Measure Foot Length	Max. Calf Circum.	Product No.
Small	7"-8"	15"	02MP-S
Medium	8"-9"	17"	02MP-M
Large	9"-10"	17"	02MP-L
X-Large	10"-11"	19"	02MP-XL



MPO Active® (with ambulatory attachment)

- Features a removable ambulatory attachment that converts the brace from a bed positioning device to a limited ambulation device secured with a simple “twist and lock” motion
- Washable synthetic sheepskin liner

Size	Measure Foot Length	Max. Calf Circum.	Product No.
Small	7"-8"	16"	09MP-S
Medium	8"-9"	16"	09MP-M
Large	9"-10"	18"	09MP-L
X-Large	10"-11"	20"	09MP-XL



Non-Pneumatic Walking Boots with Plastic Uprights

- Provides comfort and stability for fractures, soft tissue injuries or severe sprains of the ankle or foot
- Rocker bottom sole for less impact at heel strike

CoreLINE	CoreLINE Walker					
	• Rocker bottom and shock absorbing sole provides stable gait and reduces impact on heel while walking					
		US Shoe Size		Product No.		
	Size	Men's	Women's	Short	Tall	
	X-Small	2-4	3 ¹ / ₂ -5 ¹ / ₂	CL0210	CL0200	
	Small	4 ¹ / ₂ -7	6-8	CL0212	CL0202	
Darco	FX Pro Walker					
	• Low profile design with ventilated padded liner for comfort					
	• Anatomically contoured "spatula" stirrups					
	• Pressure moldable and shock absorbing insole					
		US Shoe Size		Product No.		
	Size	Men's	Women's	Short	Standard	
X-Small	2-4 ¹ / ₂	3 ¹ / ₂ -6	FXS0	FX0		
Small	5-7 ¹ / ₂	6 ¹ / ₂ -9	FXS1	FX1		
Medium	8-10 ¹ / ₂	9 ¹ / ₂ -12	FXS2	FX2		
Large	11-13	12 ¹ / ₂ -15	FXS3	FX3		
X-Large	13 ¹ / ₂ +	15 ¹ / ₂ +	FXS4	FX4		
DeRoyal	Pacesetter II					
	• Low profile rocker sole and wide foot plate					
		US Shoe Size		Product No.		
	Size	Men's	Women's	Short	Standard	
	Small	3-6 ¹ / ₂	4-7 ¹ / ₂	CSG53520-A	CSG53510-A	
	Medium	7-8 ¹ / ₂	8-10	CSG53520-B	CSG53510-B	
Large	9-10 ¹ / ₂	10 ¹ / ₂ -12 ¹ / ₂	CSG53520-D	CSG53510-D		
X-Large	11-16	13+	CSG53520-C	CSG53510-C		
DonJoy (DJO)	Motion Control (MC) Walker					
	• Provides foot/ankle immobilization on 0°, 10° or 20° of plantar flexion					
	• Range of motion in 7.5° increments from 45° plantar flexion to 30° dorsiflexion					
		US Shoe Size		Product No.		
	Size	Men's	Women's			
Small	up to 6	up to 7	11-0459-2-06136			
Medium	6-10	7-11	11-0459-3-06136			
Large	10+	11+	11-0459-4-06136			
Ossur Americas	Equalizer Walker					
	• Contoured strut design allows walker frame to conform to patients anatomy					
	• Promotes a natural gait and increases stability					
		US Shoe Size		Product No.		
	Size	Men's	Women's	Low Top	Standard	
	X-Small	2-4	3 ¹ / ₂ -5 ¹ / ₂	W0300BLKB	W0200BLKB	
Small	4 ¹ / ₂ -7	6-8	W0500BLKB	W0400BLKB		
Medium	7 ¹ / ₂ -10 ¹ / ₂	8 ¹ / ₂ -11 ¹ / ₂	W0700BLKB	W0600BLKB		
Large	10 ¹ / ₂ -12 ¹ / ₂	11 ¹ / ₂ -13 ¹ / ₂	W0900BLKB	W0800BLKB		
X-Large	12 ¹ / ₂ +	13 ¹ / ₂ +	W1100BLKB	W1000BLKB		

Non-Pneumatic Walking Boots with Plastic Uprights (continued)					
ProCare (DJO)	Nextep™ Contour Walker <ul style="list-style-type: none"> Lightweight, durable, one-piece wrap-around shell Nylon foam liner with removable toe cover 				
		<u>US Shoe Size</u>			
	Size	Men's	Women's	Product No.	
	X-Small	3 ¹ / ₂ -7 ¹ / ₂	4 ¹ / ₂ -8 ¹ / ₂	79-95063	
	Small	7-10	8-11 ¹ / ₂	79-95065	
	Medium	10-13	11-14	79-95067	
	Nextep Contour 2 Walker <ul style="list-style-type: none"> Lightweight, durable, one-piece wrap-around posterior shell Brushed nylon foam liner, low profile, cushioned non-skid sole 				
		<u>US Shoe Size</u>			
	Size	Men's	Women's	Product No.	
	X-Small	3 ¹ / ₂ -7 ¹ / ₂	4 ¹ / ₂ -8 ¹ / ₂	79-95073	
	Small	7-10	8-11 ¹ / ₂	79-95075	
	Medium	10-13	11-14	79-95077	
Nextep Contour "Shortie" Walker <ul style="list-style-type: none"> Lightweight, durable, with shorter one-piece wrap-around shell Brushed nylon foam liner, cushioned non-skid sole 					
	<u>US Shoe Size</u>				
Size	Men's	Women's	Product No.		
Small	3 ¹ / ₂ -7 ¹ / ₂	4 ¹ / ₂ -8 ¹ / ₂	79-95083		
Medium	7-10	8-11 ¹ / ₂	79-95085		
Large	10-13	11-14	79-95087		

Pneumatic Walking Boots with Plastic Uprights						
<ul style="list-style-type: none"> Feature aircells for compression, stabilization and comfort Provide comfort and stability for fractures, soft tissue injuries or sprains of the ankle or foot Rocker bottom sole for less impact at heel strike 						
CoreLINE	CoreLINE Pneumatic Walker <ul style="list-style-type: none"> Rocker bottom and shock absorbing sole provides stable gait and reduces impact on heel while walking 					
		<u>US Shoe Size</u>		<u>Product No.</u>		
	Size	Men's	Women's	Short		Tall
	X-Small	2-4	3 ¹ / ₂ -5 ¹ / ₂	CL0190		CL0180
	Small	4 ¹ / ₂ -7	6-8	CL0192		CL0182
	Medium	7 ¹ / ₂ -10 ¹ / ₂	8 ¹ / ₂ -11 ¹ / ₂	CL0194		CL0184
Large	10 ¹ / ₂ -12 ¹ / ₂	11 ¹ / ₂ -13 ¹ / ₂	CL0196	CL0186		
X-Large	12 ¹ / ₂ +	13 ¹ / ₂ +	CL0198	CL0188		
Aircast (DJO)	XP Walker - Extra Pneumatic <ul style="list-style-type: none"> 4 custom inflated aircells for individual fit and support 					
		<u>US Shoe Size</u>				
	Size	Men's	Women's	Product No.		
	Pediatric	up to 4	up to 5	01P-P		
	Small	4-7	5-8	01P-S		
	Medium	7-10	8-11	01P-M		
Large	10-13	11-15	01P-L			
X-Large	13+	15+	01P-XL			

Pneumatic Walking Boots with Plastic Uprights (continued)

Aircast (DJO)	FP Walker - Foam Pneumatic				
	<ul style="list-style-type: none"> • 2 custom inflated distal aircells • Lightweight, low profile design: rigid foot, flexible around calf for strength, fit and comfort 				
		US Shoe Size			
	Size	Men's	Women's	Product No.	
	Pediatric	up to 4	up to 5	01F-P	
	Small	4-7	5-8	01F-S	
	Medium	7-10	8-11	01F-M	
	Large	10-13	11-15	01F-L	
	X-Large	13+	15+	01F-XL	
	SP Walker - Short Pneumatic				
	<ul style="list-style-type: none"> • Shorter height makes it ideal when limb length or calf circumference prohibits use of a standard height walking brace 				
		US Shoe Size			
	Size	Men's	Women's	Product No.	
Pediatric	up to 4	up to 5	01A-P		
Small	4-7	5-8	01A-S		
Medium	7-10	8-11	01A-M		
Large	10-13	11-15	01A-L		
X-Large	13+	15+	01A-XL		
XP Diabetic Walker					
<ul style="list-style-type: none"> • Dual density insole - Impax™ Grid technology - to reduce shear stress and eliminate pressure points • Semi-rigid shell combined with full aircell coverage to maximize plantar unloading, regulate shear stress, and provide protection and immobilization • Regulated aircell inflation with hand bulb pressure gauge 					
	US Shoe Size				
Size	Men's	Women's	Product No.		
Small	4-7	5-8	01PD-S		
Medium	7-10	8-11	01PD-M		
Large	10-13	11-15	01PD-L		
X-Large	13+	15+	01PD-XL		
Replacement Bulbs			01PG		
DeRoyal	Pacesetter Air Walker				
	<ul style="list-style-type: none"> • Non-welded air bladder reduces chance of over inflation • Wide foot plate 				
		US Shoe Size			
	Size	Men's	Women's	Product No.	
	Small	3-6½	4-7½	15780005	
	Medium	7-8½	8-10	15780006	
Large	9-10½	10½-12½	15780067		
X-Large	11-16	13+	15780008		

Pneumatic Walking Boots with Plastic Uprights (continued)

DeRoyal	Inline Pneumatic Walker - Open Heel					
	<ul style="list-style-type: none"> • Wide foot plate for patient comfort and stability and threaded sole for patient safety • Heavy duty flared uprights 					
		<u>US Shoe Size</u>		<u>Product No.</u>		
	Size	Men's	Women's	Low Top	Standard	
	Pediatric	8-12 (child)		15480002	15470002	
	Small	3 ¹ / ₂ -7	5-8 ¹ / ₂	15480005	15470005	
	Medium	7 ¹ / ₂ -8 ¹ / ₂	9-10	15480006	15470006	
	Large	9-10 ¹ / ₂	10 ¹ / ₂ -12	15480067	15470067	
	X-Large	11+	12 ¹ / ₂ +	15480008	15470008	
	Inline Pneumatic Walker - Closed Heel					
<ul style="list-style-type: none"> • Wide foot plate for patient comfort and stability and threaded sole for patient safety • Heavy duty flared uprights 						
	<u>US Shoe Size</u>		<u>Product No.</u>			
Size	Men's	Women's	Low Top	Standard		
Pediatric	8-12 (child)		15500002	15490002		
Small	3 ¹ / ₂ -7	5-8 ¹ / ₂	15500005	15490005		
Medium	7 ¹ / ₂ -8 ¹ / ₂	9-10	15500006	15490006		
Large	9-10 ¹ / ₂	10 ¹ / ₂ -12	15500067	15490067		
X-Large	11+	12 ¹ / ₂ +	15500008	15490008		
Ossur Americas	Equalizer Air Walker					
	<ul style="list-style-type: none"> • Patented pneumatic system easily adjusts to achieve desired compression 					
		<u>US Shoe Size</u>		<u>Product No.</u>		
	Size	Men's	Women's	Low Top	Standard	
	X-Small	2-4	3 ¹ / ₂ -5 ¹ / ₂	AW0300BLKB	AW0200BLKB	
	Small	4 ¹ / ₂ -7	6-8	AW0500BLKB	AW0400BLKB	
	Medium	7 ¹ / ₂ -10 ¹ / ₂	8 ¹ / ₂ -11 ¹ / ₂	AW0700BLKB	AW0600BLKB	
	Large	10 ¹ / ₂ -12 ¹ / ₂	11 ¹ / ₂ -13 ¹ / ₂	AW0900BLKB	AW0800BLKB	
	X-Large	12 ¹ / ₂ +	13 ¹ / ₂ +	AW1100BLKB	AW1000BLKB	
	Equalizer Pre-Inflated Air Walker					
	<ul style="list-style-type: none"> • Pre-inflated pneumatic cells in liner 					
		<u>US Shoe Size</u>		<u>Product No.</u>		
	Size	Men's	Women's	Low Top	Standard	
	X-Small	2-4	3 ¹ / ₂ -5 ¹ / ₂	PW0300B	PW0200B	
	Small	4 ¹ / ₂ -7	6-8	PW0500B	PW0400B	
Medium	7 ¹ / ₂ -10 ¹ / ₂	8 ¹ / ₂ -11 ¹ / ₂	PW0700B	PW0600B		
Large	10 ¹ / ₂ -12 ¹ / ₂	11 ¹ / ₂ -13 ¹ / ₂	PW0900B	PW0800B		
X-Large	12 ¹ / ₂ +	13 ¹ / ₂ +	PW1100B	PW1000B		
Rebound Air Walker						
<ul style="list-style-type: none"> • Flex Edge over-molding conforms to anatomy to reduce formation of pressure points 						
	<u>US Shoe Size</u>		<u>Product No.</u>			
Size	Men's	Women's	Low Top	Standard		
X-Small	2-4	3 ¹ / ₂ -5 ¹ / ₂	B-242900061	B-242900001		
Small	4 ¹ / ₂ -7	6-8	B-242900062	B-242900002		
Medium	7 ¹ / ₂ -10 ¹ / ₂	8 ¹ / ₂ -11 ¹ / ₂	B-242900063	B-242900003		
Large	10 ¹ / ₂ -12 ¹ / ₂	11 ¹ / ₂ -13 ¹ / ₂	B-242900064	B-242900004		
X-Large	12 ¹ / ₂ +	13 ¹ / ₂ +	B-242900065	B-242900005		

Pneumatic Walking Boots with Plastic Uprights (continued)					
ProCare (DJO)	Nextep Contour Air Walker • A secure, full wrap-around shell provides maximum support for ligament and tendon repair				
	US Shoe Size				
	Size	Men's	Women's	Product No.	
	Small	3 ¹ / ₂ -7 ¹ / ₂	4 ¹ / ₂ -8 ¹ / ₂	79-95273	
	Medium	7-10	8-11 ¹ / ₂	79-95275	
	Large	10-13	11-14	79-95277	
	Nextep Contour 2 Air Walker • Lightweight, durable, one-piece wrap-around posterior shell • Brushed nylon/foam liner, low profile, cushioned non-skid sole				
	US Shoe Size				
	Size	Men's	Women's	Product No.	
	Small	3 ¹ / ₂ -7 ¹ / ₂	4 ¹ / ₂ -8 ¹ / ₂	79-95173	
	Medium	7-10	8-11 ¹ / ₂	79-95175	
	Large	10-13	11-14	79-95177	
	Nextep "Shortie" Air Walker • Lightweight, durable, with shorter one-piece wrap-around shell • Brushed nylon/foam liner, cushioned non-skid sole				
US Shoe Size					
Size	Men's	Women's	Product No.		
Small	3 ¹ / ₂ -7 ¹ / ₂	4 ¹ / ₂ -8 ¹ / ₂	79-95443		
Medium	7-10	8-11 ¹ / ₂	79-95445		
Large	10-13	11-14	79-95447		

Non-Pneumatic Walking Boots with Metal Uprights						
<ul style="list-style-type: none"> Uprights are malleable to accommodate ankle edema or bulky dressing Provides comfortable stability for acute ankle sprains, soft tissue and stress fractures to the lower leg and ankle 						
CoreLINE	CoreLINE Walker with Metal Uprights • Provides maximum comfort and stability for fractures and soft tissue injuries					
	US Shoe Size					
	Size	Men's	Women's	Product No.		
				Short		Tall
	X-Small	2-4	3 ¹ / ₂ -5 ¹ / ₂	CL0010		CL0020
	Small	4 ¹ / ₂ -7	6-8	CL0012		CL0022
	Medium	7 ¹ / ₂ -10 ¹ / ₂	8 ¹ / ₂ -11 ¹ / ₂	CL0014		CL0024
Large	11-12 ¹ / ₂	11 ¹ / ₂ -13 ¹ / ₂	CL0016	CL0026		
X-Large	13-15	13 ¹ / ₂ +	CL0018	CL0028		

Non-Pneumatic Walking Boots with Metal Uprights (continued)

AdjustaFit

- Rocker design allows patients to simulate a near normal gait without abnormal loads on any joints

Size	US Shoe Size		Product No.	
	Men's	Women's	Short	Standard
X-Small	1-2 ¹ / ₂	1-3 ¹ / ₂	AL056001	AL053001
Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	AL056003	AL053003
Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	AL056005	AL053005
Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	AL056007	AL053007
X-Large	12 ¹ / ₂ -17	13+	AL056009	AL053009



Hi-Top Boot

- Low profile rocker bottom helps eliminate pelvic tilt and alleviate low back pain
- Contains 2 layers of pads that contour to the ankle shape and support the ankle for total control

Size	US Shoe Size		Product No.	
	Men's	Women's	With Cuff	With Bootie
X-Small	1-2 ¹ / ₂	1-3 ¹ / ₂	AL036001	AL036021
Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	AL036003	AL036023
Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	AL036005	AL036025
Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	AL036007	AL036027
X-Large	12 ¹ / ₂ -17	13+	AL036009	AL036029

Hi-Top Boot EZ Set (Hinged Version)

- Designed for range-of-motion from 20° of dorsiflexion to 40° of plantar flexion in 10° increments

Size	US Shoe Size		Product No.	
	Men's	Women's	With Cuff	With Bootie
X-Small	1-2 ¹ / ₂	1-3 ¹ / ₂	AL072001	AL072021
Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	AL072003	AL072023
Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	AL072005	AL072025
Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	AL072007	AL072027
X-Large	12 ¹ / ₂ -17	13+	AL072009	AL072029



Hi-Top Boot - Achilles

- Boot has four 10° removable wedges allowing the foot to be plantar flexed while maintaining good ground contact
- Contains three layers of pads that contour to the ankle shape and support the ankle for total control

Size	US Shoe Size		Product No.
	Men's	Women's	
Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	AL036013
Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	AL036015
Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	AL036017
X-Large	12 ¹ / ₂ -17	13+	AL036019

Mid-Calf Fixed Boot

- Designed for range-of-motion from 20° of dorsiflexion to 40° of plantarflexion in 10° increments
- Contains two layers of pads that contour to the ankle shape and support the ankle for total control

Size	US Shoe Size		Product No.
	Men's	Women's	
X-Small	1-2 ¹ / ₂	1-3 ¹ / ₂	AL055001
Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	AL055003
Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	AL055005
Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	AL055007
X-Large	12 ¹ / ₂ -17	13+	AL055009



Bledsoe Brace

Non-Pneumatic Walking Boots with Metal Uprights (continued)

Bledsoe Brace	ProGait					
	<ul style="list-style-type: none"> • Straps engage and lock to brace uprights, holding the leg steady • Computer-designed rocker tread allows a near normal walking gait and low heel provides lift 					
		<u>US Shoe Size</u>		<u>Product No.</u>		
	Size	Men's	Women's	Short	Standard	
	X-Small	1-2 ¹ / ₂	1-3 ¹ / ₂	AL052001	AL032001	
	Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	AL052003	AL032003	
	Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	AL052005	AL032005	
	Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	AL052007	AL032007	
	X-Large	12 ¹ / ₂ -17	13+	AL052009	AL032009	
	ProGait EZ Set (Hinged Version)					
<ul style="list-style-type: none"> • Designed for range-of-motion from 20° of dorsiflexion to 40° of plantarflexion in 10° increments 						
	<u>US Shoe Size</u>		<u>Product No.</u>			
Size	Men's	Women's	Short	Standard		
X-Small	1-2 ¹ / ₂	1-3 ¹ / ₂	AL075001	AL074001		
Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	AL075003	AL074003		
Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	AL075005	AL074005		
Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	AL075007	AL074007		
X-Large	12 ¹ / ₂ -17	13+	AL075009	AL074009		
J Walker - Full bootie foam liner						
	<u>US Shoe Size</u>		<u>Product No.</u>			
Size	Men's	Women's				
X-Small	1-2 ¹ / ₂	1-3 ¹ / ₂	BL010001			
Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	BL010003			
Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	BL010005			
Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	BL010007			
X-Large	12 ¹ / ₂ -17	13+	BL010009			
Ossur Americas	Range Of Motion Walker					
	<ul style="list-style-type: none"> • Joint can be adjusted in 7.5° increments from 0°-45° 					
		<u>US Shoe Size</u>		<u>Product No.</u>		
	Size	Men's	Women's			
Small	4 ¹ / ₂ -7	6-8	RW0400			
Medium	7 ¹ / ₂ -10 ¹ / ₂	8 ¹ / ₂ -11 ¹ / ₂	RW0600			
Large	10 ¹ / ₂ -12 ¹ / ₂	11 ¹ / ₂ -13 ¹ / ₂	RW0800			
ProCare (DJO)	MaxTrax Walker					
	<ul style="list-style-type: none"> • Soft nylon/foam liner is constructed to ensure no seams are in critical post-surgical areas • Designed with a wider foot bed for comfort and stability following lower leg, foot and ankle trauma and post-op procedures • Swiveling D-rings allow for criss-cross or parallel strapping 					
		<u>US Shoe Size</u>		<u>Product No.</u>		
	Size	Men's	Women's	Short		Standard
	X-Small	—	up to 4	79-95342		79-95322
	Small	up to 5	4 ¹ / ₂ -6	79-95343		79-95323
Medium	5 ¹ / ₂ -10	6 ¹ / ₂ -11	79-95345	79-95325		
Large	10 ¹ / ₂ -13 ¹ / ₂	11 ¹ / ₂ -14 ¹ / ₂	79-95347	79-95327		
X-Large	14-17	—	79-95348	79-95328		

WALKING BOOTS WITH METAL UPRIGHTS

Non-Pneumatic Walking Boots with Metal Uprights (continued)

ProCare (DJO)	MaxTrax ROM Walker					
	<ul style="list-style-type: none"> • ROM settings between 45° plantarflexion and 30° dorsiflexion in 7.5° increments • Uprights will lock in fixed positions of 0°, 7.5°, 15°, 22.5°, and 30° plantar and dorsiflexion • Ideal following Achilles tendon repair for controlled rehab 					
		US Shoe Size				
	Size	Men's	Women's	Product No.		
	Pediatric/X-Small	—	up to 4	79-95332		
	Small	up to 5	4 ¹ / ₂ -6	79-95333		
Medium	5 ¹ / ₂ -10	6 ¹ / ₂ -11	79-95335			
Large	10 ¹ / ₂ -13 ¹ / ₂	11 ¹ / ₂ -14 ¹ / ₂	79-95337			
X-Large	14-17	—	79-95338			
Trulife	CAM Walker II					
	<ul style="list-style-type: none"> • Controlled ankle motion adjustability ranging from neutral to 45° • Unique ankle joint allows easy access for adjustments 					
		US Shoe Size				
	Size	Men's	Women's	Product No.		
	Small	5-7	7-9	0991302		
	Medium	8-10	10-12	0991303		
	Large	11-13	13-15	0991304		
	Pin CAM Walker					
	<ul style="list-style-type: none"> • Adjusts from 45° plantarflexion to 30° dorsiflexion in 7.5° increments • Adjustable straps and soft liner accommodate volumetric soft tissue changes 					
		US Shoe Size				
	Size	Men's	Women's	Product No.		
	Small	up to 6	up to 7	0990902		
	Medium	6-10	7-11	0990903		
	Large	10+	11+	0990904		
	Value Walker					
<ul style="list-style-type: none"> • Comfort, function and low cost • Molded of high strength reinforced nylon 						
	US Shoe Size					
Size	Men's	Women's	Product No.			
Small	up to 6	up to 7	0991002			
Medium	6-10	7-11	0991003			
Large	10+	11+	0991004			

Pneumatic Walking Boots with Metal Uprights

- Provides comfortable stability for acute ankle sprains, soft tissue and stress fractures to the lower leg and ankle
- Adjustable air chambers to accommodate changes in edema during rehab and maximize patient compliance

CoreLINE	Pneumatic Walker with Metal Uprights					
	<ul style="list-style-type: none"> • Low natural heel height, with rocker bottom and slight toe flexion promotes natural gait 					
		US Shoe Size		Product No.		
	Size	Men's	Women's	Short		Standard
	X-Small	2-4	3 ¹ / ₂ -5 ¹ / ₂	CL0030		CL0040
	Small	4 ¹ / ₂ -7	6-8	CL0032		CL0042
Medium	7 ¹ / ₂ -10 ¹ / ₂	8 ¹ / ₂ -11 ¹ / ₂	CL0034	CL0044		
Large	11-12 ¹ / ₂	11 ¹ / ₂ -13 ¹ / ₂	CL0036	CL0046		
X-Large	13-15	13 ¹ / ₂ +	CL0038	CL0048		

Pneumatic Walking Boots with Metal Uprights (continued)

Bledsoe Brace	Hi Top Boot <ul style="list-style-type: none"> Four cast interlocking 10° urethane wedges provide fast change of foot angle Air ankle/heel pad allows for ability to adjust ankle pressure for increased comfort and provides pneumatic ankle support 					
			<u>US Shoe Size</u>		<u>Product No.</u>	
	Size	Men's	Women's			
	Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	AL236013		
	Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	AL236015		
	Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	AL236017		
	X-Large	12 ¹ / ₂ -17	13+	AL236019		
	J Walker Boot <ul style="list-style-type: none"> The low heel height mimics a normal shoe for a balanced gait and reduces pressure on the knee and back 					
			<u>US Shoe Size</u>		<u>Product No.</u>	
	Size	Men's	Women's	Standard	With MC	
X-Small	1-2 ¹ / ₂	1-3 ¹ / ₂	BL210001	BL211001		
Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	BL210003	BL211003		
Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	BL210005	BL211005		
Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	BL210007	BL211007		
X-Large	12 ¹ / ₂ -17	13+	BL210009	BL211009		
ProGAIT EZ Set with Air <ul style="list-style-type: none"> Air ankle/heel pad allows for ability to adjust ankle pressure for increased comfort and provides pneumatic ankle support Designed for range-of-motion from 20° of dorsiflexion to 40° of plantar-flexion in 10 degree increments 						
		<u>US Shoe Size</u>		<u>Product No.</u>		
Size	Men's	Women's	Short	Standard		
X-Small	1-2 ¹ / ₂	1-3 ¹ / ₂	AL275001	AL274001		
Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6	AL275003	AL274003		
Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	AL275005	AL274005		
Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	AL275007	AL274007		
X-Large	12 ¹ / ₂ -17	13+	AL275009	AL274009		
Breg	Vectra Premium Air Walker <ul style="list-style-type: none"> Features a unique fold-over design that houses the integrated air pump and release valve Dual-axis rocker bottom 					
			<u>US Shoe Size</u>		<u>Product No.</u>	
	Size	Men's	Women's	Short	Standard	
	X-Small	up to 2 ¹ / ₂	up to 3 ¹ / ₂	28451	28441	
	Small	2 ¹ / ₂ -5 ¹ / ₂	3 ¹ / ₂ -6 ¹ / ₂	28452	28442	
	Medium	5 ¹ / ₂ -9	6-9 ¹ / ₂	28453	28443	
	Large	9-12 ¹ / ₂	9 ¹ / ₂ -13	28454	28444	
X-Large	12 ¹ / ₂ -17	13+	28455	28445		

Pneumatic Walking Boots with Metal Uprights (continued)

ProCare (DJO)	MaxTrax Air Walker					
	<ul style="list-style-type: none"> • Integrated pump designed to make the liner simple to inflate and adjust • Adjustable medial/lateral aircell system 					
		US Shoe Size		Product No.		
	Size	Men's	Women's	Short		Standard
	X-Small	—	up to 4	79-95422		79-95412
	Small	up to 5	4 ¹ / ₂ -6	79-95423		79-95413
	Medium	5 ¹ / ₂ -10	6 ¹ / ₂ -11	79-95425	79-95415	
	Large	10 ¹ / ₂ -13 ¹ / ₂	11 ¹ / ₂ -14 ¹ / ₂	79-95427	79-95417	
	X-Large	14-17	—	79-95428	79-95418	
	ProCare (DJO)	MaxTrax Air ROM				
<ul style="list-style-type: none"> • Integrated pump designed to make the liner simple to inflate and adjust • ROM settings between 45° plantarflexion and 30° dorsiflexion in 7.5° increments • Uprights will lock in fixed positions of 0°, 7.5°, 15°, 22.5° and 30° plantar- and dorsiflexion 						
		US Shoe Size		Product No.		
Size		Men's	Women's			
X-Small		—	up to 4	79-95432		
Small		up to 5	4 ¹ / ₂ -6	79-95433		
Medium		5 ¹ / ₂ -10	6 ¹ / ₂ -11	79-95435		
Large	10 ¹ / ₂ -13 ¹ / ₂	11 ¹ / ₂ -14 ¹ / ₂	79-95437			
X-Large	14-17	—	79-95438			
Trulife	Air CAM					
	<ul style="list-style-type: none"> • Controlled, ankle motion adjustability ranging from neutral to 45° • Unique ankle joint allows for easy adjustments 					
		US Shoe Size		Product No.		
	Size	Men's	Women's			
	Small	5-7	7-9	0992302		
Medium	8-10	10-12	0992303			
Large	11-13	13-15	0992304			

WALKING BOOT COMPENSATOR

Procare (DJO) Evenup™

- Simple, fast, easy to apply, the Evenup uses an innovative design which helps to equalize a patient's limb length and reduce body strain while walking in a cast or walker
- The Evenup comes in three different sizes and is adjustable to two thicknesses

Size	US Shoe Size		Product No.
	Men's	Women's	
Small	7-8 ¹ / ₂	6 ¹ / ₂ -9 ¹ / ₂	79-95483
Medium	9-11 ¹ / ₂	9 ¹ / ₂ -12 ¹ / ₂	79-95485
Large	12+	12 ¹ / ₂ +	79-95487



Bledsoe Brace - Conformer

- Outperforms the Total Contact Cast in reducing pressure
- Proprietary Auto-mold foam holds shape to the buttress of the foot and prevents movement that can breakdown the skin
- Dual density mid-sole foam promotes perimeter loading

Indications: Ulcerative/pre-ulcerative conditions on the bottom of the foot for sensate and non-sensate patients



Size	US Shoe Size		Product No.	
	Men's	Women's	Low Profile	High Profile
04	4-5 ¹ / ₂	5-6 ¹ / ₂	AL041104	AL041204
06	6-7 ¹ / ₂	7-8 ¹ / ₂	AL041106	AL041206
08	8-9 ¹ / ₂	9-10 ¹ / ₂	AL041108	AL041208
10	10-11 ¹ / ₂	11-12 ¹ / ₂	AL041110	AL041210
12	12-13 ¹ / ₂	13-14 ¹ / ₂	AL041112	AL041212
14	14-15 ¹ / ₂	15-16 ¹ / ₂	AL041114	AL041214

Darco - Body Armor® Walker II

- Rigid shell provides uniform compression which aids to reduce edema and secure immobilization
- Ski-boot closure allows for a secure and comfortable fit
- Flexible upper conforms to the leg and provides more uniform compression of the foot and ankle
- Removable Plastazote insoles allow for targeted off-loading when required



Size	US Shoe Size		Product No.	
	Men's	Women's	Low Profile	High Profile
X-Small	0-3	0-4	BATS0	BAT0
Small	3 ¹ / ₂ -6	4 ¹ / ₂ -6 ¹ / ₂	BATS1	BAT1
Medium	6 ¹ / ₂ -9	7-10	BATS2	BAT2
Large	9 ¹ / ₂ -13	10 ¹ / ₂ -13	BATS3	BAT3
X-Large	13 ¹ / ₂ -17	—	BATS4	BAT4

Active Innovations - Dorsal Night Splint

- Unique small anterior dorsal clamshell designed holds the foot in neutral position to prevent contracture, supports the arch and allows for proper healing
- Neutral positioning helps to alleviate morning pain and discomfort associated with plantar fasciitis
- Anti-slip pad on bottom facilitates easy, safe night-time ambulation
- New triple stitched stretch resistant strap

Size	US Shoe Size		Product No.
	Men's	Women's	
Small	5-10	5-9 ^{1/2}	DNS220
Large	10 ^{1/2} -16	10-16	DNS240



AliMed - FREEDOM Dorsal PF Night Splint

- Less bulky, cooler and lower profile than plantar-based splints
- Wide neoprene calf and toe straps with gel padding on top of foot helps reduce pressure points
- Flexible shell and plastic-free bottom allow safer and easier nighttime ambulation

Size	US Shoe Size		Product No.
	Men's	Women's	
Medium	5-9	6-10	66113
Large	9 ^{1/2} -14	10 ^{1/2} -15	66114



Darco - Body Armor® Night Splint

- Toe strap dorsiflexes the first MPJ to provide a specific and sustained stretch to the plantar fascia and a sustained low load stretch to the flexor tendons, Achilles tendon and calf muscles
- Loading the windlass mechanism naturally brings the foot/ankle to a near 90 degree angle, thus eliminating the “lever effect” and the resulting pressure on the apex of the ankle present in other dorsal models
- Elastic connector ring allows for natural flexion of the foot and toes during sleep, returning to the dorsiflexed position as muscles relax
- Low profile design allows for maximum user comfort and compliance



Product No. **BADS**

Med Spec® - Phantom™ Dorsal Night Splint

- The fully padded internal dorsal stays may be adjusted to increase or decrease dorsiflexion to the desired level
- Non-skid material sewn to sole provides better traction when patients need to ambulate during the night
- Elastic material sewn under midfoot provides additional tension on plantar fascia
- Low-profile bilateral design fits left or right foot

Size	US Shoe Size		Product No.
	Men's	Women's	
Small	4 ¹ / ₂ -7	5 ¹ / ₂ -8	264122
Medium	7 ¹ / ₂ -10	8 ¹ / ₂ -11	264124
Large	10 ¹ / ₂ -12 ¹ / ₂	11 ¹ / ₂ -13 ¹ / ₂	264125
X-Large	13-15	14-16	264126
2X-Large	15+	16+	264127



Ossur Americas - Exoform® Dorsal Night Splint

- Integrated support shell with adjustable aluminum stay maintains the therapeutic angle of the foot and lower leg at 85° to 90°
- Flex Edge™ overmolding eliminates pressure points and discomfort while protecting sleeping partners
- Hinged shell at the calf area prevents brace rotation
- Semi-rigid plantar pad reduces discomfort along the sides of the foot while providing an improved angle of dorsiflexion

Indications: Plantar Fasciitis; Achilles Tendonitis; Drop Foot; Post-Static Pain

Size	US Shoe Size		Product No.
	Men's	Women's	
Medium	5-9	6-10	W-50085
Large	9 ¹ / ₂ -14	10 ¹ / ₂ -15	W-50087



Airform® Night Splint

- Its light weight, low profile design makes it more comfortable
- Slips on and adjusts easily
- Fits left or right foot


Indications: Plantar Fasciitis; Achilles Tendonitis; Drop Foot; Post-static pain

Size	US Shoe Size		Product No.	
	Men's	Women's	Single	6 Pack
Medium	5-9	6-10	50045C	500456
Large	9 ¹ / ₂ -14	10 ¹ / ₂ -15	50047C	500476



Plantar Fasciitis/Posterior Night Splints

- Designed to gently stretch the plantar fascia and Achilles tendon and prevent ankle contractures
- Lightweight, low profile designs with adjustable straps
- Universal left or right

CoreLINE	Posterior Night Splint					 <p>Rear View</p>	
	• Dorsiflexion straps adjust from neutral (90°) up to 10° dorsiflexion						
	US Shoe Size						
	Size	Men's	Women's	Product No.			
Small	up to 7	up to 7 ¹ / ₂	CL0082				
Medium	7 ¹ / ₂ -9	8-9 ¹ / ₂	CL0084				
Large	9 ¹ / ₂ -11 ¹ / ₂	10-12 ¹ / ₂	CL0086				
Apex	Foot Pain (PF) Night Splint						
	• Includes 1/4" heel lifts with adhesive						
	Size		Product No.				
Universal		Men's	Women's				
		3102	3101				
AliMed	Foot Drop Night Splint						
	• Prevents heel cord and Achilles tendon shortening during prolonged bed rest						
	• Aliplast foam liner						
	• Foot plate set at angle slightly less than 90°						
	Size		Calf		Product No.		
			Instep	Metatarsal	Left		Right
	Infant	For infants (up to 18 months)			6771	6770	
	Junior	Children (up to 7 years)			6773	6772	
	Medium	10" or less	9" or less	7" or less	6775	6774	
	Large	10"-13"	9"-12"	7"-8 ³ / ₄ "	6777	6776	
X-Large	13"+	12"+	8 ³ / ₄ "+	6779	6778		
Soft Plantar Fasciitis Night Splint							
• Lightweight, compact and foldable and provides safe ambulation							
• Straps may be set in multiple configurations to set foot at 90° or in dorsiflexion							
Size		US Shoe Size		Product No.			
		Men's	Women's				
Small	—		5-9	64027			
Medium	7-10		9-11	64028			
Large	10-13		—	64029			

Plantar Fasciitis/Posterior Night Splints (continued)

AllMed

FREEDOM® *comfort*™ Night Splint

- Relieves inflammation of the tissue supporting the arch by passively stretching the plantar fascia Achilles tendon

Size	Product No.
Universal	64255



Early Fit Night Splint Standard

- Fixed at 90° for patients with less range-of-motion


Size	US Shoe Size		Product No.
	Men's	Women's	
Small	—	4-7 ^{1/2}	6633/NA/SM
Medium	6-8 ^{1/2}	8-10	6633/NA/MD
Large	9+	10 ^{1/2} +	6633/NA/LG



Early Fit Adjustable Ventilated Night Splint

- Adjustment in 7° increments from 40° plantarflexion to 10° dorsiflexion


Size	US Shoe Size		Product No.
	Men's	Women's	
Small	—	4-7 ^{1/2}	6652/NA/SM
Medium	6-8 ^{1/2}	8-10	6652/NA/MD
Large	9+	10 ^{1/2} +	6652/NA/LG



PF Night Splint

- Foot plate set in 5° dorsiflexion
- Strap at inset; liner included

Size	US Shoe Size		Product No.
	Men's	Women's	
Medium	—	4-7 ^{1/2}	6019/NA/SM
Large	6-8 ^{1/2}	8-10	6019/NA/MD
X-Large	9+	10 ^{1/2} +	6019/NA/LG



Plantar Fasciitis Night Splint II

- Adjusts in 2.5° increments to 10° dorsiflexion or plantarflexion
- Enlarged malleolus areas prevent tightness
- Heel is protected in a ballooned area

Size	US Shoe Size		Product No.
	Men's	Women's	
Medium	6-8	4-10	60973/NA/MD
Large	9+	10+	60973/NA/LG
X-Large	9+	10+	60973/NA/XL



Turnbuckle Ankle Orthosis

- Apply low-load torque at the ankle joint to control muscle imbalance and reduce contracture
- The splint features a pair of bilateral turnbuckles, controlling both dorsiflexion/plantarflexion and eversion/inversion simultaneously

Size	US Shoe Size		Product No.
	Men's	Women's	
Small	5-6	5-7	66426/NA/SM
Medium	6-9	8-11	66426/NA/MD
Large	10+	12+	66426/NA/LG



Plantar Fasciitis/Posterior Night Splints (continued)

Corflex	Ultra Stretch PF Night Splint <ul style="list-style-type: none"> • Polyethylene shell with padded cover • Features bilateral dorsiflexion straps and removable wedge for added toe flexion 					
	US Shoe Size					
	Size	Men's	Women's	Product No.		
	Small	5 ¹ / ₂ -6 ¹ / ₂	5-6	37-4191		
	Medium	7-8	6 ¹ / ₂ -7 ¹ / ₂	37-4192		
	Large	8 ¹ / ₂ -9 ¹ / ₂	8-9	37-4193		
X-Large	10-11	9 ¹ / ₂ -10 ¹ / ₂	37-4194			
Orthomerica	UFO™ - Universal Plantar Fasciitis Orthosis <ul style="list-style-type: none"> • Comfortably supports the affected foot in a controlled amount of dorsiflexion, applying a Low Load Prolonged Stretch (LLPS) to the soft tissues of plantar fascia and Achilles tendon • CoolFoam replaceable liner 					
	US Shoe Size					
	Size	Men's	Women's	Product No.		
	X-Small	up to 5	up to 6	3619		
	Small	5-7	6-8	3620		
	Medium	7-9	8-10	3621		
Large	9-11	10+	3622			
X-Large	11+	—	3623			
Ossur Americas	Form Fit Night Splint <ul style="list-style-type: none"> • Center release buckles simplify application • Adjusts from 10°-90° dorsiflexion • 3-dimensional padding throughout the shell and adjustable flexion straps that control the amount of stretch provided for effective treatment 					
	US Shoe Size					
	Size	Men's	Women's	Product No.		
	Small	up to 7	up to 7 ¹ / ₂	50003C		50023C
	Medium	7 ¹ / ₂ -10	8-10 ¹ / ₂	50005C		50025C
	Large	10 ¹ / ₂ +	11+	50007C		50027C
	Replacement Liners - for Form Fit Splints					
	US Shoe Size					
	Size	Men's	Women's	Product No.		
	Small	—	—	11760		
Medium	—	—	11761			
Large	—	—	11762			
Optional Foam Wedge <ul style="list-style-type: none"> • for increased plantar/dorsiflexion for Form Fit Splints 				10700		
ProCare (DJO)	Dorsiwedge Night Splint <ul style="list-style-type: none"> • Padded, contoured posterior shell provides passive dorsiflexion • Includes toe wedge for additional stretch of the plantar fascia 					
	US Shoe Size					
	Size	Men's	Women's	Product No.		
	Small	up to 6	up to 6 ¹ / ₂	79-81403		
	Medium	6 ¹ / ₂ -9 ¹ / ₂	7-10	79-81405		
	Large	10-12	10 ¹ / ₂ -12 ¹ / ₂	79-81407		
X-Large	12 ¹ / ₂ +	13+	79-81408			

Plantar Fasciitis/Posterior Night Splints (continued)					
ProCare (DJO)	Plantar Fasciitis Night Splint				
	<ul style="list-style-type: none"> • Lightweight, sturdy, one-piece flexible shell • Shipped with toe wedge for additional stretch on plantar fascia 				
		US Shoe Size			
	Size	Men's	Women's	Product No.	
	Small	up to 6	up to 7	79-97753	
	Medium	6 ¹ / ₂ -8 ¹ / ₂	7 ¹ / ₂ -9 ¹ / ₂	79-97755	
Large	9-11 ¹ / ₂	10-12 ¹ / ₂	79-97757		
X-Large	12+	13+	79-97758		
RCAI	Plantar Fasciitis Night Splint				
	<ul style="list-style-type: none"> • Padded, contoured posterior shell • Wide foam inset strap maintains proper foot alignment • Includes toe wedge for additional stretch of the plantar fascia 				
	Size	Foot Length	Height	Product No.	
	Youth	up to 6 ³ / ₄ "	12"	22PFNS-Y	
	X-Small	up to 7 ³ / ₄ "	12 ¹ / ₂ "	22PFNS-ES	
	Small	up to 8 ³ / ₄ "	13"	22PFNS-S	
	Medium	up to 9 ³ / ₄ "	13 ¹ / ₂ "	22PFNS-M	
	Large	up to 10 ³ / ₄ "	14 ¹ / ₂ "	22PFNS-L	
X-Large	up to 11 ³ / ₄ "	15"	22PFNS-XL		

HINGED ANKLE SPLINTS

Active Innovations - Active Ankle T2

- Solid U-shaped design provides inversion and eversion protection
- Bilateral hinge allows for dorsiflexion and plantar flexion movement
- Single strap adjusts for either high or low top shoes
- EVA padding system ensures snug fit and comfort

Indications: Prophylactic use as well as for minor ankle sprains and injuries

Size	US Shoe Size		Product No.	
	Men's	Women's	White	Black
Small	5-8 ¹ / ₂	6 ¹ / ₂ -9 ¹ / ₂	TRT120	TRT220
Medium	9-12	10-13	TRT130	TRT230
Large	12 ¹ / ₂ -16	13 ¹ / ₂ -16	TRT140	TRT240



Active Ankle Multi-Phase™

- Dual-density padding provides compression for edema control in acute phase
- Easily modifies for comfort in prophylactic phase
- Solid U-shape design eases the weight bearing on the ankle

Size	US Shoe Size		Product No.
	Men's	Women's	
X-Small	—	5-6	MPP210
Small	5-8 ¹ / ₂	6 ¹ / ₂ -9 ¹ / ₂	MPP220
Medium	9-12	10-13	MPP230
Large	12 ¹ / ₂ -16	13 ¹ / ₂ -16	MPP240



Breg - Ultra Ankle XT

- LiteFit padding system with lock-down strap
- Anti-slip heel
- Circumferential cuff along with a dynamic hinge that allows for normal plantar or dorsiflexion
- Universal left or right

Indications: High ankle sprains; Sub-acute ankle sprains; Chronic ankle instabilities; Prophylactic use

Size	US Shoe Size		Product No.
	Men's	Women's	
Small	5-9	6-10	11022
Medium	9 ¹ / ₂ -12	10 ¹ / ₂ -13	11023
Large	12 ¹ / ₂ -15	13 ¹ / ₂ +	11024



Ossur Americas - Rebound® Ankle Brace

- Offers more stability than soft ankle supports, and is perfect for transitioning out of a walker boot
- Application is simple and reduces plantar flexion, making it ideal for patients with limited range-of-motion
- The low-profile design is comfortable and fits easily in most shoes

Indications: Recommended for ankle injury rehabilitation of grade I, II, III and high ankle sprains and strains; Transition from walker boot or cast to daily activities or sports after an injury; Support for chronic instability following rehabilitation; Prophylactic use

*Rebound Ankle Brace
shown with optional
Stability Strap*



Size	US Shoe Size		Rebound Ankle Brace		with Stability Strap	
	Men's	Women's	Left	Right	Left	Right
Small	6-7 ¹ / ₂	7 ¹ / ₂ -9	B-231600012	B-231500012	B-231602512	B-231502512
Medium	8-12	9 ¹ / ₂ -13 ¹ / ₂	B-231600013	B-231500013	B-231602513	B-231502513
Large	12 ¹ / ₂ -16	—	B-231600014	B-231500014	B-231602514	B-231502514

Swede-O - Arch Lok®

- Modified orthotic foot plate provides critical foot/arch support and resists inversion/eversion
- Neutral arch position promotes healing of injured ligaments
- Three strap adjustment system provides customized fit

Size	US Shoe Size		Black Product No.		White Product No.	
	Men's	Women's	Left	Right	Left	Right
X-Small	—	4-6	23111	23011	03111	03011
Small	4-7	7-9	23112	23012	03112	03012
Medium	8-10	10-12	23113	23013	03113	03013
Large	11-13	13-15	23114	23014	03114	03014
X-Large	14+	16+	23115	23015	03115	03015



Air Ankle Stirrups						
<ul style="list-style-type: none"> • Air bladders provide compression and comfort • Useful for ankle sprains, ankle instabilities 						
CoreLINE	Air/Gel Stirrup Ankle Brace <ul style="list-style-type: none"> • Ankle wrap ensures control of swelling to speed recovery while minimizing ankle rotation • Air/Gel bladder may be used for cold therapy 			 <p><i>Internal view of air bladder</i></p>		
	Size	Height	Product No.			
	X-Small	6"	CL0090			
	Small	8 ³ / ₄ "	CL0092			
	Medium	9 ¹ / ₂ "	CL0094			
Large	10"	CL0096				
Aircast (DJO)	Air-Stirrup® Universe <ul style="list-style-type: none"> • Patented Duplex™ aircell system to enhance circulation and reduce swelling • Adjustable heel pad 					
	Size	Height	Product No.			
	One Size	9"	Single			12-Pack
			02E	02E12		
	Air-Stirrup Ankle Brace <ul style="list-style-type: none"> • Semi-rigid, anatomically designed shells for protection, comfort, and the prevention of inversion/eversion • Patented Duplex™ aircell system to enhance circulation and reduce swelling 					
	Size	Height	Product No.			
			Left			Right
	Pediatric	6"	02JL			02JR
	Small	8 ³ / ₄ "	02CL			02CR
	Medium	9"	02BL	02BR		
	Large	10 ¹ / ₂ "	02AL	02AR		
	Air-Stirrup Plus <ul style="list-style-type: none"> • Combines the clinically proven Air-Stirrup Ankle Brace with an ankle wrap 					
	Size	Height	Product No.			
			Left			Right
	Small	8 ⁷ / ₁₀ "	02CLP			02CRP
Medium	9"	02BLP	02BRP			
Large	10 ¹ / ₂ "	02ALP	02ARP			
Sport-Stirrup® <ul style="list-style-type: none"> • Air-Stirrup ankle brace with a narrower shell for less bulk in a sneaker or shoe 						
Size	Height	Product No.				
		Left			Right	
Universal	9"	02DL	02DR			

Air Ankle Stirrups (continued)

Air Ankle Stirrups (continued)						
Corflex	Air Stirrup					
	<ul style="list-style-type: none"> • Prevents eversion/inversion of ankle • Constructed of semi-rigid shells connected with soft, adjustable foot-piece 					
	Size	Product No.				
	Universal fit, up to 14" ankle circumference		75-1015			
Ossur Americas	Airform Inflatable Stirrup Ankle Brace					
	<ul style="list-style-type: none"> • Inflatable air bladders provide adjustable pneumatic compression • Adjustable, pivoting heel strap increases comfort and allows for full range-of-motion 					
	Size	Product No.				
		Left	Right			
	Youth/Pony	70200	70201			
	Adult	80200	80201			
	Airform Universal Inflatable Stirrup Ankle Brace - White					
	<ul style="list-style-type: none"> • Soft Flex-Edges provide a universal fit • Inflatable air bladders provide adjustable pneumatic compression 					
	Size	Product No.				
		Each	6/pack			
	Youth	80251	802516			
	Adult	80250	802506			
Airform Universal Inflatable Stirrup Ankle Brace - Black (not pictured)						
Size	Product No.					
Adult - each	80350					
Airform Pre-Inflated Stirrup Ankle Brace						
<ul style="list-style-type: none"> • Pre-inflated, pulsating air bladder with regulated chambers reduces bottoming out and decreases edema • Adjustable, pivoting heel strap 						
Size	Product No.					
	Left	Right				
Youth/Pony	70100	70101				
Adult	80100C	80101C				
ProCare (DJO)	Surround® Air Ankle					
	<ul style="list-style-type: none"> • Two full-length pre-inflated air chambers lined with soft foam help reduce swelling and enhance circulation • Air valve at proximal end of bladder allows for simple adjustment of air • Injection molded shells help minimize ankle rotation and inversion/eversion without restricting plantar or dorsiflexion 					
	Size	Height	Product No.			
			Left			Right
Trainer	9"	79-81715	79-81705			
Regular	10"	79-81717	79-81707			

Air Ankle Stirrups (continued)				
ProCare (DJO)	Surround Gel Ankle			
	<ul style="list-style-type: none"> • Rigid thermoplastic shells with adjustable heel strap to help prevent inversion/eversion • Air and gel bladder may be used for cold therapy to treat tender or swollen ankles 			
	Size	Height	Product No.	
	Small	8 ¹ / ₂ "	79-97863	
Trainer	9"	79-97865		
Regular	10"	79-97867		

Gel/Foam Ankle Stirrups				
Ossur Americas	Form Fit® with Honeycomb Ankle Brace			
	<ul style="list-style-type: none"> • Patented honeycomb gel pads • Flex-edges and adjustable pivoting heel strap 			
	Size	Product No.		
		Left	Right	
	Youth/Pony	40100	40101	
	Adult	30100C	30101C	
	Form Fit Foam Ankle Brace			
	<ul style="list-style-type: none"> • Patented 3-dimensional molded pads • Flex-edges and adjustable pivoting heel strap 			
	Size	Product No.		
		Left	Right	
Youth/Pony	20100C	20101C		
Adult	50100C	50101C		
Gel Ankle Brace				
<ul style="list-style-type: none"> • Safe hot/cold therapy • Patented lace-lock feature • Available in 50% surface gel or 80% surface gel 				
Size	Product No.			
	Left	Right		
Youth/Pony - 50% Gel	E200PLC	E200PRC		
Adult - 50% Gel	E5050LC	E5050RC		
Adult - 80% Gel	E8020LC	E8020RC		
ProCare (DJO)	Excelerator Ankle Brace			
	<ul style="list-style-type: none"> • Constructed of a three-layer foam system for added comfort and bilateral rigid shells to help prevent inversion and eversion 			
	Size	Height	Product No.	
	Small	8"	79-81323	
Large	9 ¹ / ₂ "	79-81327		

Gel/Foam Ankle Stirrups (continued)					
ProCare (DJO)	Universal Ankle Brace				
	<ul style="list-style-type: none"> • Rigid thermoplastic one piece outer shell helps prevent inversion/eversion • Cushioned liner incorporates neoprene forefoot closure • Elastic proximal strap provides additional compression 				
	Size	Height	Product No.		
	Universal	10"	79-81330		
	Surround FLOAM™ Ankle				
<ul style="list-style-type: none"> • Provides stability to the ankle by minimizing ankle rotation • Comprised of two polypropylene outer shells connected with soft, adjustable foot piece • FLOAM™ bladder provides uniform compression and a customized fit 					
Size	Height	Product No.			
Pediatric	6"	79-81193			
Medium	9"	79-81195			
Large	10"	79-81197			
Professional Products	"Robust" Stirrup Ankle Brace				
	<ul style="list-style-type: none"> • Semi-rigid shell material • Impact-resistant plastic • Patient can adjust distal pressure while wearing device 				
	Size	Height	Product No.		
One Size	10 ³ / ₁₀ "	Left	Right		
		00840-LEFT	00840-RIGHT		

Ankle Braces							
<ul style="list-style-type: none"> • For ankle support and the treatment of ankle sprains 							
Aircast (DJO)	A60 Ankle Support						
	<ul style="list-style-type: none"> • 60 degree stabilizer guards against rollover and helps prevent sprains 						
	Size	US Shoe Size				Product No.	
		Men's	Women's			Left	Right
	Small	up to 7	up to 8 ¹ / ₂	02TSL	02TSR		
	Medium	7 ¹ / ₂ -11 ¹ / ₂	9-13	02TML	02TMR		
	Large	12+	13 ¹ / ₂ +	02TLL	02TLR		
	Airsport Ankle Brace						
	<ul style="list-style-type: none"> • Foam-filled aircell that cushions and supports the ankle • ATF cross-strap for ankle stabilization • Integral forefoot and shin wraps for support 						
	Size	US Shoe Size				Product No.	
	Men's	Women's	Left			Right	
X-Small	up to 5	up to 5	02MXSL			02MXSR	
Small	5 ¹ / ₂ -7	5-8 ¹ / ₂	02MSL			02MSR	
Medium	7 ¹ / ₂ -11	9-12 ¹ / ₂	02MML	02MMR			
Large	11 ¹ / ₂ -13	13-14 ¹ / ₂	02MLL	02MLR			
X-Large	13 ¹ / ₂ +	15+	02MXLL	02MXLR			

Ankle Braces (continued)						
Bauerfeind	MalleoLoc® Ankle Brace <ul style="list-style-type: none"> Stimulates the tibialis anterior muscle Titanium color 					
	Size	Heel Width	Product No.			
			Left	Right		
	1	up to 2 ³ / ₈ "	12013013080701	12013013080601		
2	2 ³ / ₈ "+	12013013080702	12013013080602			
Bledsoe Brace	Axiom Ankle Brace <ul style="list-style-type: none"> Separate upper shells allow the material to naturally form to any ankle Foot plate material is heat formable Center pivot on the malleolus eliminates pressure and supports natural anatomical movement 					
	Size	US Shoe Size		Product No.		
		Men's	Women's	Left		Right
	Small	up to 7 ¹ / ₂	up to 9	AS020103		AS020203
	Medium	8-11 ¹ / ₂	9 ¹ / ₂ -13	AS020105	AS020205	
	Large	12+	13 ¹ / ₂ +	AS020107	AS020207	
	Wrapr Ankle Stabilizer <ul style="list-style-type: none"> Figure 8 strapping mimics traditional taping with gapping High ankle strap provides additional support and secures medial and lateral stirrup straps If shoe size is close to border between ankle brace sizes, choose the smaller size if foot is narrow; choose larger size if foot is wide 					
	Size	US Shoe Size		Product No.		
		Men's	Women's	Regular Laces		Speed Laces
	X-Small	6-7	7-8	SA702501		SA702001
Small	8-9	9-10	SA702503	SA702003		
Medium	9-11	10-12	SA702505	SA702005		
Large	11-13	12-14	SA702507	SA702007		
X-Large	14+	15+	SA702509	SA702009		
2X-Large	14W+	15W+	SA702511	SA702011		
Breg	Lace-Up Ankle <ul style="list-style-type: none"> Bilateral figure 8 strapping Padding over Achilles tendon 					
	Size	US Shoe Size		Product No.		
		Men's	Women's			
	2X-Small	up to 6	up to 5	97060		
	X-Small	6-7	5-6	97061		
	Small	7 ¹ / ₂ -9	6 ¹ / ₂ -8	97062		
	Medium	9 ¹ / ₂ -11	8 ¹ / ₂ -10	97063		
	Large	11 ¹ / ₂ -13	10 ¹ / ₂ -12	97064		
X-Large	13 ¹ / ₂ -15	12 ¹ / ₂ -14	97065			
2X-Large	15 ¹ / ₂ +	—	97066			

Ankle Braces (continued)

		Web Ankle Brace			
		• Bungee closure provides compression and easy adjustments			
Darco	Size	US Shoe Size		Product No.	
		Men's	Women's		
	X-Small	3-5	4-6	WB0	
	Small	6-7	6 ¹ / ₂ -9	WB1	
	Medium	7 ¹ / ₂ -10	9 ¹ / ₂ -11	WB2	
	Large	10 ¹ / ₂ -12	11 ¹ / ₂ -13	WB3	
X-Large	12 ¹ / ₂ -14	13 ¹ / ₂ -15	WB4		
					
		Drytex RocketSoc			
		• Patented design features calcaneal heel lock and peroneal forefoot strap to control abnormal plantar flexion and inversion without restricting range-of-motion			
		US Shoe Size		Product No.	
		Size	Men's	Women's	Left
DonJoy (DJO)	X-Small	4-5 ¹ / ₂	5-7 ¹ / ₂	11-0667-1-06000	11-0665-1-06000
	Small	6-7 ¹ / ₂	8-9 ¹ / ₂	11-0667-2-06000	11-0665-2-06000
	Medium	8-9 ¹ / ₂	10-11 ¹ / ₂	11-0667-3-06000	11-0665-3-06000
	Large	10-11 ¹ / ₂	12-13 ¹ / ₂	11-0667-4-06000	11-0665-4-06000
	X-Large	12-13 ¹ / ₂	—	11-0667-5-06000	11-0665-5-06000
	2X-Large	14-16	—	11-0667-6-06000	11-0665-6-06000
					
		Lace-Up Drytex RocketSoc			
		• Patented design features calcaneal heel lock and peroneal forefoot strap to control abnormal plantar flexion and inversion without restricting range-of-motion			
		US Shoe Size		Product No.	
		Size	Men's	Women's	Regular
DonJoy (DJO)	X-Small	4-5 ¹ / ₂	5-7 ¹ / ₂	11-0339-1-06000	11-0338-1-06000
	Small	6-7 ¹ / ₂	8-9 ¹ / ₂	11-0339-2-06000	11-0338-2-06000
	Medium	8-9 ¹ / ₂	10-11 ¹ / ₂	11-0339-3-06000	11-0338-3-06000
	Large	10-11 ¹ / ₂	12-13 ¹ / ₂	11-0339-4-06000	11-0338-4-06000
	X-Large	12-13 ¹ / ₂	—	11-0339-5-06000	11-0338-5-06000
	2X-Large	14-16	—	11-0339-6-06000	11-0338-6-06000
					
		Stabilizing Ankle			
		• Features non-stretch nylon figure 8 straps which lock the calcaneus into place to help control abnormal eversion and inversion			
		Ankle Circumference		Product No.	
		Size			
DonJoy (DJO)	X-Small	10"-11"		11-0451-1-06000	
	Small	11"-12"		11-0451-2-06000	
	Medium	12"-13"		11-0451-3-06000	
	Large	13"-14"		11-0451-4-06000	
	X-Large	14"-15"		11-0451-5-06000	
	2X-Large	15"-16"		11-0451-6-06000	
	3X-Large	16"-17"		11-0451-7-06000	
					

Ankle Braces (continued)

Freeman	Canvas Ankle Brace				
	<ul style="list-style-type: none"> • Heavy duty canvas construction with flannel lining • A floating tongue has two removable pre-bent aluminum stays included for additional support 				
	Size	Ankle Circumference	Product No.		
	Small	7 ¹ / ₂ "-8 ¹ / ₂ "	912-S		
	Medium	8 ¹ / ₂ "-9 ¹ / ₂ "	912-M		
Hely & Weber	Webly Ankle Orthosis				
	<ul style="list-style-type: none"> • Constructed of heavy gauge ballistic nylons • Limits inversion/eversion while not restricting plantar and dorsiflexion • Neoprene tongue and Achilles relief padding for comfort 				
	Size	US Shoe Size			Product No.
		Men's	Women's		
	2X-Small	—	5-6		304XXS
	X-Small	6-7	6-8		304XS
	Small	7 ¹ / ₂ -9	8-10		304S
	Medium	9 ¹ / ₂ -11	10-12	304M	
	Large	11 ¹ / ₂ -13	12-14	304L	
	X-Large	13 ¹ / ₂ -15	—	304XL	
	2X-Large	14-16	—	304XXL	
	Optional Stays for Webly Ankle				
	Description		Product No.		
Single Stays		308			
Pair of Stays		308A			
Med Spec	ASO Ankle Stabilizer				
	<ul style="list-style-type: none"> • The patented figure-eight stabilizing straps are angled to better contour the malleoli for enhanced comfort and support 				
	Size	Ankle Circumference (around the heel)	Product No.		
			Black		White
	2X-Small	9"-10"	264019		264009
	X-Small	10"-11"	264011		264001
	Small	11"-12"	264012		264002
	Medium	12"-13"	264014		264004
	Large	13"-14"	264015		264005
	X-Large	14"-15"	264016		264006
2X-Large	15"-16"	264017	264007		
3X-Large	16"+	264018	264008		

Ankle Braces (continued)

Med Spec

ASO EVO Ankle Stabilizer

- Figure-eight straps for enhanced comfort and support
- Dynamic cuff protects against high ankle sprains

Size	Ankle Circumference (around the heel)	Product No.	
		Black	White
X-Small	10"-11"	264091	264081
Small	11"-12"	264092	264082
Medium	12"-13"	264094	264084
Large	13"-14"	264095	264085
X-Large	14"-15"	264096	264086
2X-Large	15"-16"	264097	264087



EVO Speed Lacer Ankle Stabilizer

- Figure-eight straps for enhanced comfort and support
- Dynamic cuff protects against high ankle sprains
- Speed lace closure for faster and easier application

Size	Ankle Circumference (around the heel)	Product No.
X-Small	10"-11"	264111
Small	11"-12"	264112
Medium	12"-13"	264114
Large	13"-14"	264115
X-Large	14"-15"	264116
2X-Large	15"-16"	264117



ASO Flex-Hinge

- The internal plastic shell with hinge is built into a ballistic nylon boot and flexes to conform to the anatomy of the ankle
- The upper portion of the internal plastic shell is incorporated into the eyelet lace closure

Size	Ankle Circumference (around the heel)	Product No.
X-Small	10"-11"	264071
Small	11"-12"	264072
Medium	12"-13"	264074
Large	13"-14"	264075
X-Large	14"-15"	264076
2X-Large	15"-16"	264077



ASO Speed Lacer Ankle Stabilizer

- Figure-eight straps for enhanced comfort and support
- Speed Laces for faster and easier application

Size	Ankle Circumference (around the heel)	Product No.	
		Black	White
X-Small	10"-11"	223611	223601
Small	11"-12"	223612	223602
Medium	12"-13"	223614	223604
Large	13"-14"	223615	223605
X-Large	14"-15"	223616	223606
2X-Large	15"-16"	223617	223607



Ankle Braces (continued)					
Med Spec	ASO MAX <ul style="list-style-type: none"> • Figure-eight straps for enhanced comfort and support • Features Skinloc™ material which allows the braces to be worn without a sock 				
	Size	Ankle Circumference (around the heel)	Product No.		
	X-Small	10"-11"	264041		
	Small	11"-12"	264042		
	Medium	12"-13"	264044		
	Large	13"-14"	264045		
	X-Large	14"-15"	264046		
2X-Large	15"-16"	264047			
Ossur Americas	Gameday® Ankle Brace <ul style="list-style-type: none"> • Durable ballistic nylon shell with inner non-slip coating adds support by reducing slippage • Figure-8 straps mimic athletic taping for inversion/eversion support • Measure ankle circumference diagonally over the base of the heel and across the ankle bone 				
	Size	Ankle Circumference	Product No.		
	2X-Small	9"-10"	W10601		
	X-Small	10"-11"	W10602		
	Small	11"-12"	W10603		
	Medium	12"-13"	W10605		
	Large	13"-14"	W10607		
	X-Large	14"-15"	W10608		
	2X-Large	15"-16"	W10609		
	Form Fit® Ankle Brace <ul style="list-style-type: none"> • Stirrup shaped exoskeleton provides additional support • Speed laced design provides the option of a quick lace up • Measure ankle circumference diagonally over the base of the heel and across the ankle bone 				
	Size	Ankle Circumference	Product No.		
			Form Fit Ankle		w/ Figure 8 Straps
	X-Small	11"-12"	B-212010001		B-212000001
Small	12"-13"	B-212010002	B-212000002		
Medium	13"-14"	B-212010003	B-212000003		
Large	14"-15"	B-212010004	B-212000004		
X-Large	15"-16"	B-212010005	B-212000005		
ProCare (DJO)	Lace-Up Ankle Brace <ul style="list-style-type: none"> • Canvas stirrup design with removable, plastic medial/lateral stays provide support for acute ankle ligament sprains and chronic instabilities • Fits easily into any athletic shoe 				
	Size	Ankle Circumference	Product No.		
	X-Small	5 ¹ / ₂ "-7"	79-81312		
	Small	7"-8 ¹ / ₂ "	79-81313		
	Medium	8 ¹ / ₂ "-10"	79-81315		
	Large	10"-11 ¹ / ₂ "	79-81317		
X-Large	11 ¹ / ₂ "-13 ¹ / ₂ "	79-81318			

Ankle Braces (continued)

ProCare (DJO)		Stabilized Ankle Support				
		<ul style="list-style-type: none"> • Non-stretch nylon • Lace-up design and trimmable figure 8 strapping • Removable medial and lateral stays 				
Size	Ankle Circumference	Product No.				
X-Small	10"-11"	79-81352				
Small	11"-12"	79-81353				
Medium	12"-13"	79-81355				
Large	13"-14"	79-81357				
X-Large	14"-15"	79-81358				
2X-Large	15"-16"	79-81359				
3X-Large	16"-17"	79-81359-10				
Swede-O		Ankle Lok				
		<ul style="list-style-type: none"> • Removable side stabilizer inserts provide additional medial and lateral support • Internal U-shaped spiral stays provide extra support • ANKLE LOK® offset panel traps the laces between the inner and outer flap to hold laces tighter and longer 				
Size	US Shoe Size		Product No.			
	Men's	Women's	Black	White		
2X-Small	—	2-3	—	01210R		
X-Small	3-5 ¹ / ₂	4-6 ¹ / ₂	21211R	01211R		
Small	6-7 ¹ / ₂	7-9 ¹ / ₂	21212R	01212R		
Medium	8-10 ¹ / ₂	10-11 ¹ / ₂	21213R	01213R		
Large	11-12 ¹ / ₂	12-13 ¹ / ₂	21214R	01214R		
X-Large	13-14 ¹ / ₂	14-15 ¹ / ₂	21215R	01215R		
2X-Large	15+	—	21216R	01216R		
Swede-O		Ankle Lok with Padded Tongue - White				
		<ul style="list-style-type: none"> • Padding tongue provides easier application and comfort for fresh sprains • Removable side stabilizer inserts provide additional medial and lateral support • Internal U-shaped spiral stays provide extra support 				
Size	US Shoe Size		Product No.			
	Men's	Women's				
X-Small	3-5 ¹ / ₂	4-6 ¹ / ₂	01011R			
Small	6-7 ¹ / ₂	7-9 ¹ / ₂	01012R			
Medium	8-10 ¹ / ₂	10-11 ¹ / ₂	01013R			
Large	11-12 ¹ / ₂	12-13 ¹ / ₂	01014R			
X-Large	13-14 ¹ / ₂	14-15 ¹ / ₂	01015R			
2X-Large	15+	—	01016R			

Ankle Braces (continued)

Inner Lok 8®

- Figure-eight straps are pre-configured inside the brace
- Quick release strap makes brace removal quick

Size	US Shoe Size		Product No.	
	Men's	Women's	Black	White
X-Small	3-5 ^{1/2}	4-6 ^{1/2}	22611	02611
Small	6-7 ^{1/2}	7-9 ^{1/2}	22612	02612
Medium	8-10 ^{1/2}	10-11 ^{1/2}	22613	02613
Large	11-12 ^{1/2}	12-13 ^{1/2}	22614	02614
X-Large	13-14 ^{1/2}	14-15 ^{1/2}	22615	02615



Tarsal Lok

- Combines the support of a rigid ankle brace with the comfort of a lace-up brace
- Patented high impact built-in stabilizer design helps control the midtarsal joint

Size	US Shoe Size		Product No.
	Men's	Women's	
X-Small	3-5 ^{1/2}	4-6 ^{1/2}	22811
Small	6-7 ^{1/2}	7-9 ^{1/2}	22812
Medium	8-10 ^{1/2}	10-11 ^{1/2}	22813
Large	11-12 ^{1/2}	12-13 ^{1/2}	22814
X-Large	13-14 ^{1/2}	14-15 ^{1/2}	22815



X8

- Twin dual-purpose straps serve as both the figure eight straps and the top locking straps
- Straps are pre-positioned halfway through the figure-eight configuration for easier application and to ensure proper strap placement

Size	US Shoe Size		Product No.
	Men's	Women's	
X-Small	3-5 ^{1/2}	4-6 ^{1/2}	22711
Small	6-7 ^{1/2}	7-9 ^{1/2}	22712
Medium	8-10 ^{1/2}	10-11 ^{1/2}	22713
Large	11-12 ^{1/2}	12-13 ^{1/2}	22714
X-Large	13-14 ^{1/2}	14-15 ^{1/2}	22715



Swede-O

Aircast (DJO) - Air Heel™

- Made of lightweight breathable fabric and is simple to apply and adjust with one strap
- Two interconnected aircells located under the foot arch and in back of the heel apply pulsating compression with every step; helps reduce swelling and discomfort, enhances circulation

Indications: Plantar fasciitis; Achilles tendonitis; Heel pain



Size	US Shoe Size		Product No.
	Men's	Women's	
Small	up to 7	up to 8 ¹ / ₂	09AS
Medium	7 ¹ / ₂ -11	8 ¹ / ₂ -12 ¹ / ₂	09AM
Large	11 ¹ / ₂ +	13+	09AL

Bauerfeind - AchilloTrain® Achilles Tendon Support

- Anatomical contouring and compression promote proprioception and muscular stabilization
- Viscoelastic insert runs alongside the Achilles tendon for uniform application of pressure and targeted intermittent compression
- Separate heel cushion provided for unaffected leg



Size	Ankle Circum.	Product No.		
		Black	Titanium	Nature
Left Side				
1	6 ³ / ₄ "-7 ¹ / ₂ "	11011013070701	11011013080701	11011013010701
2	7 ¹ / ₂ "-8 ¹ / ₄ "	11011013070702	11011013080702	11011013010702
3	8 ¹ / ₄ "-9"	11011013070703	11011013080703	11011013010703
4	9"-9 ⁷ / ₈ "	11011013070704	11011013080704	11011013010704
5	9 ⁷ / ₈ "-10 ⁵ / ₈ "	11011013070705	11011013080705	11011013010705
Right Side				
1	6 ³ / ₄ "-7 ¹ / ₂ "	11011013070601	11011013080601	11011013010601
2	7 ¹ / ₂ "-8 ¹ / ₄ "	11011013070602	11011013080602	11011013010602
3	8 ¹ / ₄ "-9"	11011013070603	11011013080603	11011013010603
4	9"-9 ⁷ / ₈ "	11011013070604	11011013080604	11011013010604
5	9 ⁷ / ₈ "-10 ⁵ / ₈ "	11011013070605	11011013080605	11011013010605

MalleoTrain® Ankle Support

- Lightweight, breathable, knitted support provides stabilization through graduated compression
- Two anatomical silicone inserts massage the soft tissue, increasing circulation and reducing pain and swelling
- Promotes rapid absorption of edema

Size	Ankle Circum.	Product No.		
		Black	Titanium	Nature
Left Side				
1	6 ³ / ₄ "-7 ¹ / ₂ "	11011102070701	11011102080701	11011102010701
2	7 ¹ / ₂ "-8 ¹ / ₄ "	11011102070702	11011102080702	11011102010702
3	8 ¹ / ₄ "-9"	11011102070703	11011102080703	11011102010703
4	9"-9 ⁷ / ₈ "	11011102070704	11011102080704	11011102010704
5	9 ⁷ / ₈ "-10 ⁵ / ₈ "	11011102070705	11011102080705	11011102010705
6	10 ⁵ / ₈ "-11 ³ / ₈ "	11011102070706	11011102080706	11011102010706
Right Side				
1	6 ³ / ₄ "-7 ¹ / ₂ "	11011102070601	11011102080601	11011102010601
2	7 ¹ / ₂ "-8 ¹ / ₄ "	11011102070602	11011102080602	11011102010602
3	8 ¹ / ₄ "-9"	11011102070603	11011102080603	11011102010603
4	9"-9 ⁷ / ₈ "	11011102070604	11011102080604	11011102010604
5	9 ⁷ / ₈ "-10 ⁵ / ₈ "	11011102070605	11011102080605	11011102010605
6	10 ⁵ / ₈ "-11 ³ / ₈ "	11011102070606	11011102080606	11011102010606

**Bledsoe Brace - PFS Plantar Fasciitis Strap**

- Reduces plantar fasciitis pain through continuous elastic tension and pressure day or night
- Easily fits into most shoes for daytime wear

Indications: Plantar Fasciitis

Size	US Shoe Size		Product No.	
	Men's	Women's	Left	Right
Small	up to 8	up to 8 ¹ / ₂	ES007103	ES007203
Medium	8 ¹ / ₂ -11	9-11 ¹ / ₂	ES007105	ES007205
Large	11 ¹ / ₂ +	12+	ES007107	ES007207



Medi USA - Ankle Support

- Provides relief of the irritated joint and improves blood flow to the soft tissue areas surrounding the joint by producing a constant compression and massage action during patient movement
- Constructed of a lightweight, breathable material which makes it comfortable to wear, thus improving patient compliance and therapeutic results
- Compression class II

Indications: Inflammation or strains of the ankle; Muscle strains; Bruise; Prophylactic and therapeutic use in sport activities

Size	Ankle Circumference	Product No.
X-Small	6 ³ / ₄ "-8"	50101
Small	8"-8 ³ / ₄ "	50102
Medium	8 ³ / ₄ "-9 ¹ / ₂ "	50103
Large	9 ¹ / ₂ "-10 ¹ / ₄ "	50104
X-Large	10 ¹ / ₄ "-11"	50105
2X-Large	11"-11 ³ / ₄ "	50106

Measure circumference one inch above ankle bone



Levamed® Ankle Support with Silicone Inserts

- Anatomically tailored support with medial/lateral integrated silicon pads
- Knitted comfort-zone over the instep prevents constriction and pressure areas

Indications: Slight instability of the ankle joint; Football player's ankle, sprains; Joint effusions and swellings due to arthritis and osteoarthritis; Chronic, postoperative and post-traumatic inflammation

Size	Ankle Circumference	Product No.	
		Sand	Black
I	6 ³ / ₄ "-8"	50771	50851
II	8"-8 ³ / ₄ "	50772	50852
III	8 ³ / ₄ "-9 ¹ / ₂ "	50773	50853
IV	9 ¹ / ₂ "-10 ¹ / ₄ "	50774	50854
V	10 ¹ / ₄ "-11"	50775	50855
VI	11"-11 ³ / ₄ "	50776	50856

Circumference of ankle at narrowest point above ankle bone



ProCare (DJO) - Elastic Ankle Support

- Constructed of durable, breathable cotton elastic
- Pull on design with open heel for mild compression and support

Indications: Mild ankle sprains and strains; Edema; Osteoarthritis

Size	Ankle Circumference	Product No.
Small	7 ³ / ₄ "-8 ¹ / ₂ "	79-81123
Medium	8 ¹ / ₂ "-9 ¹ / ₂ "	79-81125
Large	9 ¹ / ₂ "-10 ¹ / ₂ "	79-81127
X-Large	10 ¹ / ₂ "-11 ¹ / ₂ "	79-81128
2X-Large	11 ¹ / ₂ "-12 ¹ / ₂ "	79-81129



AliMed - Heelift® Suspension Boot

- Relieves pressure and helps control foot drop
- Patented elevation pad disperses pressure evenly across calf
- Smooth tricot fabric covering on posterior side encourages mobility
- Ventilated with 1³/₄" air holes and exposed anterior side
- Two adjustable hook and loop straps

Size	Calf Circum.	Product No.	
		Smooth Interior	Convoluted Interior
Petite	6"-10"	66016	66015
Standard	8"-14"	63261	62884
Bariatric	12"-23"	66017	—



Convoluted Interior Shown

DeRoyal - Universal Bucks Traction Boot

- Convoluted foam construction helps prevent migration
- Medial and lateral stays provide additional stability
- 1" ventilation holes increase airflow
- Spreader bar attached
- Trimmable foam body

Indications: Hip fractures; Hip dislocations; Femur fractures

Product No. **8121-00**



DonJoy (DJO) - Dura*Soft® Surgical Foot/Ankle

- Universal fit; Ultra-wick fabric with antimicrobial protection
- Ankle wrap with nylon wraps and elastic straps
- Includes two ice inserts which can be frozen and exchanged to provide extended cold therapy

Product No. **11-0911-0-02000**



DuroMed - Heelbo® Heel and Elbow Standard and Premium Protectors

- Treat and help prevent dermal ulcers
- Allows skin to breathe, doesn't restrict blood flow, stays in place
- Standard features double-layer of foam padding minimizing pressure
- Sold as pair
- Measure at the largest circumference of arm/ankle/knee

Size	Description	Product No.
Small	8" Yellow	12037
Regular	8 ¹ / ₂ " Blue	12038
Large	9 ¹ / ₂ " White	12039
X-Large	11 ¹ / ₂ " Green	12040
2X-Large	12 ¹ / ₂ " Beige	12041



Ossur Americas - Universal Bucks Traction Boot

- Constructed of a durable foam pocket to cradle foot and lower leg
- Adjustable to suit all sizes
- Convenient loop lock closures
- Universal size

Product No. **204230**



WBC - AFO Pads

- Felt 1/4" thick
- Sold 10/pack
- Available in White

Width	Product No.
1"	SPS90001
1 1/2"	SPS9002
2"	SPS9003
1" x 2 1/2"	SPS9001W
1 1/2" x 3 3/8"	SPS9002W
2" x 3 3/4"	SPS9003W



AFO Pads

- Felt 1/8" thick
- Sold 10/pack
- Available in White

Width	Product No.	
	Slender	Wide
1"	AWH10	AWH10W
1 1/2"	AWH15	AWH15W
2"	AWH20	AWH20W

AFO Pads

- Felt 1/8" thick
- Sold as each
- Available in Black and Beige

Width	Product No.		
	Beige Slender	Beige Regular	Black Regular
1"	AFO1B	WBC6124MB	WBC6124ML
1 1/2"	AFO15B	WBC6124LB	WBC6124LL
2"	AFO2B	WBC6124XLB	—

Comfort Products - AFO Liner Socks

- Available in the following styles:
 - Traditional: All cotton, fully terry lined
 - CoolMax™: Wicks moisture away from skin
 - SilverMax™: Combines X-Static silver and CoolMax
- Available three pair (6 socks) per package

Description/ Size	Black Cotton/Nylon	Natural Cotton	White CoolMax	Silver SilverMax
Child's 4-6 ¹ / ₂	AFOLSB04	AFOLS04	AFOLSCM04	AFOLSX04
Child's 6-8 ¹ / ₂	AFOLSB06	AFOLS06	AFOLSCM06	AFOLSX06
Child's 12 ¹ / ₂	AFOLSB09	AFOLS09	AFOLSCM09	AFOLSX09
Women's 4-9 ¹ / ₂ Men's 3-8 ¹ / ₂	AFOLSB07	AFOLS07	AFOLSCM07	AFOLSX07
Women's 10-13 Men's 9-12	AFOLSB10	AFOLS10	AFOLSCM10	AFOLSX10
Women's 13+ Men's 12+	AFOLSB13	AFOLS13	AFOLSCM13	AFOLSX13



Knit-Rite - SmartKnit® AFO Interface Socks

- Seamless knitting and form fitting core-spun yarns help eliminate wrinkling, bunching and seam pressure
- Provides exceptional softness, moisture wicking and antimicrobial to help prevent odor in the sock
- "Heel-less" design ensures a perfect fit every time
- Sold as pair

Description	Width/Length	Product No.		
		White	Pink	Purple
Infant	2" x 9"	1SP1AFNF	1S21AFNF	1S31AFNF
Child Small	2 ¹ / ₄ " x 11"	1SP1AFCS	1S21AFCS	1S31AFCS
Child Regular	2 ¹ / ₂ " x 14"	1SP1AFCR	1S21AFCR	1S31AFCR
Child X-Long	2 ³ / ₄ " x 20"	1SP1AFCX	—	—
Adult Small	3" x 18"	1SP1AFAS	—	—
Adult Regular	3 ¹ / ₂ " x 20"	1SP1AFAR	—	—



AFO Interface Socks with X-Static® - The Silver Fiber™

- Sold as Pair

Description	Width/Length	Navy	Product No.	
			Gray	Black
Infant	2" x 9"	1NX1AFNF	1SX1AFNF	1BX1AFNF
Child Small	2 ¹ / ₄ " x 11"	1NX1AFCS	1SX1AFCS	1BX1AFCS
Child Regular	2 ¹ / ₂ " x 14"	1NX1AFCR	1SX1AFCR	1BX1AFCR
Child X-Long	2 ³ / ₄ " x 20"	—	1SX1AFCX	—
Adult Small	3" x 18"	—	1SX1AFAS	1BX1AFAS
Adult Regular	3 ¹ / ₂ " x 20"	—	1SX1AFAR	1BX1AFAR

SmartKnit® AFO Liner™

- Innovative protective sleeve covers the AFO, not the leg
- Shields the wearer's skin from direct contact with the AFO reducing sheering and chafing
- Provides breathable cushion for the foot without being hot or bulky
- Sold each



Description	Size (inches)	Flap Length	CoolMax White	Polyester Tan	X-Static Black
Child Small	2 ³ / ₄ x 2 ³ / ₄ x 7	5"	1SP1ABCS	—	—
Child Regular	3 x 3 x 9	5"	1SP1ABCR	1TP1ABCR	1BX1ABCR
Narrow Short	2 ¹ / ₂ x 2 ¹ / ₂ x 11	5"	1SP1ABNS	1TP1ABNS	1BX1ABNS
Narrow Long	2 ¹ / ₂ x 2 ¹ / ₄ x 14	5"	1SP1ABNL	1TP1ABNL	1BX1ABNL
Medium Short	3 ¹ / ₂ x 2 ³ / ₄ x 13	5"	1SP1ABMS	1TP1ABMS	1BX1ABMS
Medium Regular	3 ¹ / ₂ x 2 ³ / ₄ x 14	7 ¹ / ₂ "	1SP1ABMR	1TP1ABMR	1BX1ABMR
Wide Medium	5 x 3 ¹ / ₂ x 14 ¹ / ₂	7"	1SP1ABWM	1TP1ABWM	1BX1ABWM
X-Wide Medium	6 x 4 x 15	7"	1SP1AXWM	1TP1AXWM	1BX1AXWM

SmartKnit AFO Liner Sizing Chart				
Description	Overall Length	Upper Width	Mid-Calf Width	Ankle Width
Child Small	11"-14"	4"-8"	3 ¹ / ₂ "-6"	3 ¹ / ₂ "-7 ¹ / ₂ "
Child Regular	13 ¹ / ₂ "-17 ¹ / ₂ "	4"-8"	3 ¹ / ₂ "-6"	3 ¹ / ₂ "-7 ¹ / ₂ "
Narrow Short	16"-20 ¹ / ₂ "	2 ¹ / ₂ "-7"	2"-6"	2"-5"
Narrow Long	19"-24"	2"-6"	1 ³ / ₄ "-5"	1 ³ / ₄ "-4"
Medium Short	17"-22"	4"-8"	3 ¹ / ₂ "-6 ¹ / ₂ "	3 ¹ / ₂ "-6 ¹ / ₂ "
Medium Regular	19"-25"	4"-9"	3 ¹ / ₂ "-7"	3 ¹ / ₂ "-7 ¹ / ₂ "
Wide Medium	19"-26"	6"-10"	4"-10"	5"-10"
X-Wide Medium	22"-27"	8"-13"	6"-12"	7"-10"

Innovative Neurotronics - WalkAide® System

- Utilizes a patented tilt accelerometer to adjust stimulation timing, promoting improved gait pattern
- No heel sensor required; patients can wear any footwear - no orthopedic shoes or special shoes required, users can walk barefoot
- Each unit is wirelessly programmed and customized by the clinician using the WalkAnalyst software
- Adaptations for vision and hearing-impaired users
- 3 selections of electrode option
- 4 cuff size choices
- Certification: Clinicians wishing to fit the WalkAide must attend a 2-day educational training session. Please contact customer service or an SPS Sales Manager for training locations and dates

Indications: Foot drop secondary to upper motor neuron injuries, such as stroke; multiple sclerosis; Traumatic brain injury; Cerebral palsy; Incomplete spinal cord injury

Description	Product No.
WalkAide Patient Kit includes:	20-1000
WalkAide Control Unit	—
WalkAide Electrode Lead cable	20-0110
WalkAide Standard Electrodes 1.25", 4/pk	20-0350
WalkAide Cuff, Small 10"-14"	20-0210
WalkAide Cuff, Medium 14"-18"	20-0220
WalkAide Cuff, Large 18"-22"	20-0230
WalkAide Patient Heel Sensor - Left	20-0400
WalkAide Patient Heel Sensor - Right	20-0401
WalkAide Standard Electrodes, 1.25", 4/pk	20-0500
WalkAide Standard Electrodes, 1.25", Case (40)	20-0510
WalkAide 1.25" Premium Round Electrodes, 4/pk	20-0521
WalkAide 1.25" Premium Round Electrodes, 40/case	20-0531
WalkAide 1.875" Premium Round Electrodes, 4/pk	20-0540
WalkAide 1.875" Premium Round Electrodes, 40/case	20-0550
WalkAide 2x2 Blue Gel Electrodes, 4 per pack	20-0520
WalkAide 2x2 Blue Gel Electrodes, Case (40)	20-0530
WalkLink Clinician Kit includes:	20-0333
WalkLink Cable	20-0311
Clinical Heel Sensor	20-0320
Bluetooth Adapter	20-0350
MiniStim	20-0340
WalkAnalyst Software	20-0603
Optional Parts	
Disposable Cuff Liner (10/pk)	20-0290
Clinician Tool Case	20-0313
Netbook, with WalkAnalyst & Bluetooth Installed	20-0314



WalkAide® Bi-Flex Cuff

- Dual durometer construction: rigid side helps secure the cuff with easy one-handed donning and doffing, while soft side conforms to the leg for total electrode contact
- Visual indicator helps accurately align the WalkAide System for consistent positioning
- Electrode locators help ensure precise electrode placement every time for optimized stimulation
- Universal fit can be used on either left or right leg

Size	Circumference	Product No.
Pediatric	7 ¹ / ₂ "-10"	20-0240B
Small	10"-15"	20-0250
Medium	11"-17"	20-0260
Large	12"-19"	20-0270

Optional Parts for Bi-Flex Cuff

Size	Product No.	
	Liner	Strap
Small	20-0251	20-0252
Medium	20-0261	20-0262
Large	20-0271	20-0272

Electrode Locators

Description	Product No.	
	Black	Red
Locator for 1 ¹ / ₄ " Electrode	20-0201	20-0202
Locator for 2" Electrode	20-0203	20-0204
Locator for 1 ⁷ / ₈ " Electrode	20-0205	20-0206

Washable Liner promotes patient comfort and hygiene

Vented Design allows for increased air circulation and increased breathability



Unique Dual Durometer Construction Rigid side helps secure the cuff with easy one-handed operation. Soft side conforms to the leg for total electrode contact

Customizable Strap Kid-friendly optional ribbons*

Pediatric Bi-Flex Cuff

Electrode Locators help ensure precise electrode placement every time for optimized stimulation

Visual Indicator helps accurately align the WalkAide System for consistent positioning



Universal Fit can be used on either the left or right leg and is available in three convenient sizes

Easy-to-Use Latch securely holds the cuff to the leg for consistent reproducibility

Adult Bi-Flex Cuff

Prizm Medical - Micro-Z™ Mini Stimulator

- A clinically proven, Medicare approved, non-pharmaceutical alternative for the treatment of idiopathic pain, post-op pain and chronic intractable pain
- Lightweight (2 oz), convenient and easy to use
- Programmable for day or night use
- Wide selection of garments to facilitate use on upper or lower limbs to achieve pain relief
- The Silver-Thera Stocking™ or Silver-Thera Prosthetic Sock™ facilitate nocturnal treatment (cycles 20 min ON and 40 min OFF over 8 hours) to achieve pain relief and increase microcirculation for those with impaired vascular function

Indications: Pain relief; Increase microcirculation; Maintain or increase range-of-motion; Prevention or retardation of disuse atrophy; Relaxation of muscles spasms



Silver-Thera Stocking



Silver-Thera Prosthetic Sock

Silver-Thera Stocking and Silver-Thera Prosthetic Sock (above) are knitted garment electrodes which conduct energy (electrotherapy) to the foot/leg or residual limb respectively via the Micro-Z Mini Stimulator.

Description	Product No.
Micro-Z Mini 100 volt	MZM100
Accessories	
Electrodes 3" x 4" 2/pack	EMSP003
Conduct Mist Electrolyte spray 250 ml	CMS250
Thera-Cream conductive lotion	TCT004

Silver-Thera Sock Electrodes		
Size	Shrinker Size	Product No.
Narrow Short	B, 0, 1	STP5011
Regular Short	1, 2	STP5012
Regular Medium	1, 2	STP5013
Wide Medium	3+	STP5014

Silver-Thera Stocking Electrodes			
Size	US Shoe Size		Product No.
	Men's	Women's	
Small	—	up to 8 ¹ / ₂	STS2312
Medium	up to 11	up to 9	STS2313
Large	12-14	—	STS2314
X-Large	14+	—	STS2315

Silver-Thera Glove Electrodes		
Size	Glove Size	Product No.
Small	6 ¹ / ₂	STG2412
Medium	7	STG2413
Large	8	STG2414

