

**Shoes**

- Answer 2® ..... 147-149
- Apex® ..... 150-158
- Hush Puppies® ..... 158-160
- Mt. Emey® ..... 160-164
- New Balance® ..... 164-166
- SureFit® ..... 167-169
  - Women's ..... 167
  - Men's ..... 168-169

**Relief and Healing Shoes**

- Bauerfeind® ..... 169-170
- Darco® ..... 170-173
- DeRoyal ..... 173-175
- Hely & Weber™ ..... 175-176
- Ossur Americas® ..... 176-178
- Procare (DJO) ..... 178-181
- RCAI® ..... 181

**Arch Supports and Inserts**

- Acor® ..... 181-183
- Alimed® ..... 184-185
- Answer 2 ..... 185
- Apex ..... 186-188
- Bauerfeind ..... 189
- DonJoy (DJO) ..... 189
- Euro International ..... 190
- Fillauer ..... 190-191
- Hapad, Inc. .... 191
- Langer ..... 191-192
- PediFix® ..... 192-193
- Pedors ..... 194
- ProCare (DJO) ..... 194
- Silipos® ..... 194
- SoleTech® ..... 195
- Spenco® ..... 195-196
- Stable Step, Inc. .... 196-197
- UCO International ..... 198
- Wrymark ..... 198

**Diabetic Socks**

- Apex ..... 199
- BSN Medical ..... 199
- Comfort Products™ ..... 200
- Knit-Rite® ..... 200-202
- Silipos ..... 203

**Foot and Digital Care**

- Gel Cups ..... 204
- Recover Heel Cups ..... 205
- Posted Heel Cups ..... 205
- Tuli's ..... 206
- Leg Length Discrepancy Pads ..... 206
- M/L Lifts/Sports Wedges ..... 206
- Viscoelastic Varus/Valgus Wedges ..... 206
- Apex ..... 207
- Silipos ..... 207

# Shoe Width Conversion

The following conversion charts are provided as a general “rule of thumb” to assist you when ordering shoes available through SPS.

To improve fitting results, it is recommended that you use a “Ritz” measuring stick (p/n 1141) or “Brannock” device (p/n 1131M or 1131W) to more accurately determine the correct shoe width and size for your patient. Please contact your SPS Sales Manager or Customer Service for assistance.

Men’s Shoe Width Cross Reference	
N	B
M	D
W	E
WW	2E
WWW	3E-4E (depending upon the manufacturer)

Women’s Shoe Width Cross Reference	
N	2A
M	B
W	D
WW or XW	2E
WWW or XXW	3E




Place patient’s heel at rear edge of the appropriate (right or left) heel cup.



HEEL to TOE: Take measurement of longest toe on the center scale. As pictured, measurement would be recorded as an 8, since longest toe length is between 7½ & 8.



HEEL to BALL (Arch Length): Position heel to ball device snug against the first metatarsal head. Record the arch length. As pictured, the measurement for the left foot is 8+.



WIDTH: Slide bar against the outside of the foot. Locate the heel to toe length on the sloped scale and read the width immediately below. Heel to toe length was between 7½ & 8. The corresponding width is between an “A” and a “B”. Order a “B” width or “N”.

**Answer 2 Men's Synchronized Walking Shoe**

- Can be modified to increase severity at toe-off by adding on another rocker sole
- Designed to meet stability needs while standing at mid-stance
- Provides increased propulsion to ease smoothly into roll-over and toe-off

**Indications:** Metatarsalgia; Hallux Limitus; Minor ankle pain; Toe amputation; Maintenance of healed ulcer or wearers of AFO/KAFO



Description	Sizes	Widths	Product No.
Black Lace	8 1/2-13	M, W, XW	<b>565-1</b>

**Ordering Example**  
 Model ..... **565-1**  
 Size ..... **11**  
 Width ..... **XW**  
**565-111XW**

**Answer 2 Men's Walking Shoe**

- Full grain genuine calf leather uppers & leather lining
- Added depth to accommodate AFOs, KAFOs, and foot orthotics
- Added width at the toe and extra girth at the ankle
- Extended medial counters provide support and stability

Description	Sizes	Widths	Product No.
<b>Men's Answer 2 Walking Shoe</b>			
Black Lace	6-13, 14, 15	M, W, XW	<b>557-1</b>
White/Gray Lace	6-13, 14, 15	M, W, XW	<b>557-3</b>
Black Hook and Loop	6-13, 14, 15	M, W, XW	<b>558-1</b>
White/Gray Hook and Loop	6-13, 14, 15	M, W, XW	<b>558-3</b>
<b>Men's Answer 2 Casual Shoe</b>			
Black Lace	6-13, 14, 15	M, W, XW	<b>555-1</b>
Brown Lace	6-13, 14, 15	M, W, XW	<b>555-2</b>
Black Hook and Loop	6-13, 14, 15	M, W, XW	<b>556-1</b>
Brown Hook and Loop	6-13, 14, 15	M, W, XW	<b>556-2</b>



557-1  
Lace



558-3  
Hook & Loop



555-1  
Lace



556-1  
Hook & Loop

**Ordering Example**  
 Model ..... **556-1**  
 Size ..... **12**  
 Width ..... **XW**  
**556-112XW**

Shoes

**Answer 2 Men's Athletic Boot**

- Leather boots with long counter for maximum stability
- Added depth to accommodate AFOs, KAFOs, and foot orthotics
- Symmetrical cut line for easy sole modification

Description	Sizes	Widths	Product No.
<b>Men's Answer 2 Athletic Boot</b>			
Black Lace	6-13, 14, 15	M, W, XW	<b>552-1</b>
Black Hook and Loop	6-13, 14, 15	M, W, XW	<b>551-1</b>
<b>Men's Answer 2 Athletic Shoe</b>			
Black Lace	6-13, 14, 15	M, W, XW	<b>554-1</b>
White/Navy Lace	6-13, 14, 15	M, W, XW	<b>554-3</b>
Black Hook and Loop	6-13, 14, 15	M, W, XW	<b>553-1</b>
White/Navy Hook and Loop	6-13, 14, 15	M, W, XW	<b>553-3</b>



552-1  
Lace



551-1  
Hook & Loop



554-3  
Lace



553-1  
Hook & Loop

Ordering Example	
Model . . . . .	<b>552-1</b>
Size . . . . .	<b>10</b>
Width . . . . .	<b>XW</b>
<b>552-110XW</b>	

**Answer 2 Women's Athletic Boot**

- Leather shoe with long counter provides maximum stability
- Accommodates AFO's and KAFO's
- Symmetrical cut straight last, easy sole modification

Description	Sizes	Widths	Product No.
Black Hook and Loop	5-11	M, W, XW	<b>441-1</b>
Black Lace	5-11	M, W, XW	<b>442-1</b>



441-1  
Hook & Loop



442-1  
Lace

Ordering Example	
Model . . . . .	<b>441-1</b>
Size . . . . .	<b>10</b>
Width . . . . .	<b>W</b>
<b>441-110W</b>	

**Answer 2 Women's Casual Shoe**

- Soft genuine leather uppers & leather lining
- Added depth to accommodate AFOs, KAFOs, and foot orthotics
- Added width at the toe and extra girth at the ankle
- Extended medial counters provide support and stability
- Removable three-layer cushioning inserts
- Symmetrical cut line for easy sole modification



445-1  
Lace



446-4  
Hook & Loop

Description	Sizes	Widths	Product No.		
			Casual	Walking	Athletic
Black Lace	5-11	M, W, XW	445-1	447-1	444-1
Black Hook & Loop	5-11	M, W, XW	446-1	448-1	443-1
White/Gray Lace	5-11	M, W, XW	447-3	—	444-3
White/Gray Hook & Loop	5-11	M, W, XW	448-3	—	443-3
Bone Lace	5-11	M, W, XW	445-4	—	—
Bone Hook and Loop	5-11	M, W, XW	446-4	—	—




447-3  
Lace



448-1  
Hook & Loop



444-1  
Lace



443-3  
Hook & Loop

Shoes

**Ordering Example**  
 Model . . . . . 444-1  
 Size . . . . . 9  
 Width . . . . . W  
**444-19W**

**Apex - Edge™ Runner**

- Cobra support pod prevents over-pronation and provides midfoot stability
- Aeromax™ mesh upper for breathability
- Q last for superior stability
- TPU footbridge controls motion and stabilizes midsole
- Innovative heel cradle midsole for rearfoot stability, unmatched shock attenuation and pressure relief to the calcaneus



Q560M

Description	Product No.
Men's Black	<b>Q560M</b>
Men's Silver/Blue	<b>Q562M</b>
Men's Gun Metal/Yellow	<b>Q563M</b>
Women's Black	<b>Q560W</b>
Women's White/Light Blue	<b>Q562W</b>
Women's Silver/Raspberry	<b>Q568W</b>



Q560W

Edge Runner Sizing Chart		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E)	7-13, 14
Women's	Medium (B/C), Wide (D/E)	5-11

**Ordering Example**  
 Model ..... **Q560M**  
 Width ..... **M**  
 Size ..... **7**  
**Q560MM07**

**Women's Xspress™ Fitness Runner**

- Breathable Airlex™ mesh upper promotes maximum airflow
- S last for superior cushioning and support
- Skylight™ TPU external heel counter cups the rearfoot for extra stability
- Innovative "Donut" heel design unloads stress from the calcaneus



Color	Widths	Sizes	Product No.
Black/Coral	Medium (B/C), Wide (D/E)	5-11	<b>S680W</b>
White/Navy	Medium (B/C), Wide (D/E)	5-11	<b>S681W</b>
Silver/Blue	Medium (B/C), Wide (D/E)	5-11	<b>S682W</b>

**Ordering Example**  
 Model ..... **S680W**  
 Width ..... **M**  
 Size ..... **11**  
**S680WM11**

### Sierra Trail Runner

- SmartMask™ toe guard & external heel counter
- Nylon Cordura™ uppers with nubuck leather
- Sculpted EVA midsole with EGX™ footbridge
- Scotchguard™ protection for water resistance
- V last for stability

Color	Product No.	
	Men's	Women's
Brown	<b>V751M</b>	<b>V751W</b>
Gray	<b>V753M</b>	<b>V753W</b>

Sierra Trail Runner Sizing Chart		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	6 1/2-13, 14, 15
Women's	Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4 1/2-11, 12, 13



**Ordering Example**  
 Model ..... **V751M**  
 Width ..... **M**  
 Size ..... **13**  
**V751MM13**

### Stealth™ Runner

- Highly breathable Airlex™ sandwich mesh
- Multi-density midsole with supporting EGX™ upper footbridge
- Medial Support Pod to minimize over pronation
- Molded, external heel counter for extra rearfoot control

Color	Product No.	
	Men's	Women's
Blue	<b>V551M</b>	<b>V551W</b>



Shoes

### Rhino & Reina Runner

- Open air sandwich mesh uppers for breathability
- Strategically placed flex grooves for traction during gait cycle
- Carboplast® footbridge for motion control & stability
- 2 removable layers provide fitting flexibility

Description	Product No.
Men's Blue	<b>X532M</b>
Women's Periwinkle	<b>X532W</b>

Stealth, Rhino and Reina Runner Sizing Chart		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	6 1/2-13, 14, 15
Women's	Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4 1/2-11, 12, 13



**Ordering Example**  
 Model ..... **X532M**  
 Width ..... **W**  
 Size ..... **11**  
**X532MW11**

**Boss™ Runner**

- X last for superior motion control
- 2 removable layers provide fitting flexibility
- Mozaic® pressure relief layer for customized comfort
- Carboplast® footbridge for motion control and stability
- Strategically placed flex groves for traction during gait cycles

Description	Product No.
Men's White/Blue	<b>X521M</b>
Men's White/Blue	<b>X522M</b>
Women's White/Periwinkle	<b>X521W</b>
Women's White/Pale Blue	<b>X522W</b>
Women's Silver/Sea Blue	<b>X527W</b>

Boss Runner Sizing Chart		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	6½-13, 14, 15
Women's	Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4½-11, 12, 13



X521M



X522M



X527W

**Ordering Example**

Model ..... **X521M**  
 Width..... **W**  
 Size..... **11**  
**X521MW11**

**Double Strap Walker**

- Full-grain leather upper
- X last for superior stability and motion control
- Mozaic pressure relief layer for customized comfort
- Two removable layers provide fitting flexibility

Description	Product No.
Men's Black	<b>X920M</b>
Men's White	<b>X926M</b>
Men's Tan	<b>X929M</b>
Women's White	<b>X926W</b>

Double Strap Walker Sizing Chart		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	6½-13, 14, 15
Women's	Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4½-11, 12, 13



X926W

**Ordering Example**

Model ..... **X920M**  
 Width..... **W**  
 Size..... **7**  
**X920MW07**



### Men's Triple Strap Walker

- X last for superior motion control
- 2 removable layers provide fitting flexibility
- Mozaic pressure relief layer for customized comfort
- Full grain leather upper for comfort and breathability
- Hook and loop closure
- Available in Black

Widths	Sizes	Product No.
Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	6 <sup>1</sup> / <sub>2</sub> -13, 14, 15	<b>X903M</b>



Ordering Example	
Model .....	<b>X903M</b>
Width .....	<b>W</b>
Size .....	<b>13</b>
<b>X903MW13</b>	

### Voyage Walker - Lace

- V last for superior stability
- 2 removable layers provide fitting flexibility
- Pressure relief layer for customized comfort
- Full grain leather upper for comfort and breathability

Description	Product No.
Men's White/Blue	<b>V854M</b>
Women's White/Periwinkle	<b>V854W</b>

Voyage Walker Sizing Chart		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	6 <sup>1</sup> / <sub>2</sub> -13, 14, 15
Women's	Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4 <sup>1</sup> / <sub>2</sub> -11, 12, 13



V854M

Shoes

### Athletic Walker - Lace

- X last for superior motion control
- 2 removable layers provide fitting flexibility
- Pressure relief layer for customized comfort
- Full grain leather upper for comfort and breathability

Color	Product No.	
	Men's	Women's
Black	<b>X801M</b>	<b>X801W</b>
White/Gray	<b>X821M</b>	<b>X821W</b>
White	<b>X826M</b>	<b>X826W</b>

Athletic Walker Sizing Chart		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	6 <sup>1</sup> / <sub>2</sub> -13, 14, 15
Women's	Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4 <sup>1</sup> / <sub>2</sub> -11, 12, 13



X821M

Ordering Example	
Model .....	<b>X801W</b>
Width .....	<b>W</b>
Size .....	<b>12</b>
<b>X801WW12</b>	

**Men's Lexington™ Dress Casual**

- Full-grain leather upper with Mozaic® insole for customized comfort and pressure relief
- 1/4" removable depth to accommodate many orthotics
- Firm heel counter & soft padded collar

Color	Product No.	
	Cap Toe Oxford	Moc Toe Oxford
Black	<b>LT600M</b>	<b>LT900M</b>
Brown	<b>LT610M</b>	<b>LT910M</b>

Lexington Oxford Sizing Chart		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	7-13, 14



**Ordering Example**  
 Model . . . . . **LT600M**  
 Width . . . . . **M**  
 Size . . . . . **7**  
**LT600MM07**

**WOMEN'S COMFORT CASUALS**  
**DOUBLE DEPTH SHOES**

**Julia**

- Full grain leather upper with textured accents
- Soft, padded collars help reduce heel slippage
- Roomy toe box design for forefoot freedom
- Advanced polyurethane outsoles for shock absorption and resilience
- Narrow available in black only

Color	Product No.
Black	<b>E360W</b>
Brown	<b>E361W</b>



**Eliza**

- Soft leather and suede upper
- Soft, padded collars help reduce heel slippage
- Roomy toe box design for forefoot freedom
- Advanced polyurethane outsoles for shock absorption and resilience
- Narrow available in black only

Color	Product No.
Black	<b>E730W</b>
Brown	<b>E731W</b>



**Ordering Example**  
 Model . . . . . **E731W**  
 Width . . . . . **W**  
 Size . . . . . **11**  
**E731WW11**

Eliza Sizing Chart		
Shoe	Widths	Sizes
Women's	Narrow (A/AA), Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4 1/2-11, 12, 13

**Frances**

- Full grain leather upper
- Single-strap closure
- Roomy toe box design for forefoot freedom
- Advanced polyurethane outsoles for shock absorption and resilience

Product No. **E820W**

Frances Sizing Chart		
Shoe	Widths	Sizes
Women's	Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4 <sup>1</sup> / <sub>2</sub> -11, 12, 13



**Ordering Example**  
 Model ..... **E820W**  
 Width..... **M**  
 Size..... **7**  
**E820WM07**

**Mary**

- Soft leather and suede upper with double strap closure
- Soft, padded collars help reduce heel slippage
- Roomy toe box design for forefoot freedom
- Advanced polyurethane outsoles for shock absorption and resilience
- Narrow available in black only

Color	Product No.
Black	<b>E830W</b>
Brown	<b>E831W</b>

Mary Sizing Chart		
Shoe	Widths	Sizes
Women's	Narrow (A/AA), Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4 <sup>1</sup> / <sub>2</sub> -11, 12, 13



**Ordering Example**  
 Model ..... **E830W**  
 Width..... **M**  
 Size..... **7**  
**E830WM07**

Shoes

**Ambulator® Biomechanical Footwear**

- Genuine leather uppers and SmartGrip™ polyurethane soles
- 1/2" removable hidden depth in 3 layers
- Anti-Shox® footbed
- Soft linings with leather for maximum comfort and breathability

Description	Closure	Product No.
<b>Men's Biomechanical Shoes</b>		
Black Leather Oxford	Lace	<b>BIO2000M</b>
Brown Leather Oxford	Lace	<b>BIO2100M</b>
Black Leather	Hook and Loop	<b>BIO3000M</b>
Brown Leather	Hook and Loop	<b>BIO3100M</b>
Black Leather Boot	Lace	<b>BIO4000M</b>
<b>Women's Biomechanical Shoes</b>		
Black Leather Oxford	Lace	<b>BIO2000W</b>
Black Leather	Hook and Loop	<b>BIO3000W</b>
Black Leather Mary Jane	Hook and Loop	<b>BIO6000W</b>
Brown Suede Mary Jane	Hook and Loop	<b>BIO6100W</b>
Taupe Leather Mary Jane	Hook and Loop	<b>BIO6400W</b>

<b>Biomechanical Ambulator Sizing Chart</b>		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	6 1/2-13, 14, 15, 16
Women's	Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4 1/2-11, 12, 13

**Ordering Example**

Model ..... **BIO2000M**  
 Width..... **M**  
 Size..... **7**  
**BIO2000MM07**



*BIO2000M/BIO2100M*



*BIO3000M*



*BIO4000M*



*BIO6000W*



*BIO2000W*

### Ambulator® Conform

- Genuine leather uppers with SmartGrip™ polyurethane soles
- Plastazote® foam lining molds to the foot for a customized fit
- 1/2" removable hidden depth in 2 layers
- Removable Conform Orthotic with Plastazote
- High toe boxes

Description	Product No.
<b>Men's Conform Shoes</b>	
Black Leather Double Strap	<b>1260M</b>
Brown Leather Double Strap	<b>1261M</b>
Black Leather Oxford	<b>1270M</b>
Brown Leather Oxford	<b>1271M</b>
<b>Women's Conform Shoes</b>	
Black Leather Double Strap	<b>1260W</b>
Taupe Leather Double Strap	<b>1264W</b>
Black Leather Oxford	<b>1270W</b>

Conform Ambulator Sizing Chart		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	6 1/2-13, 14, 15, 16
Women's	Medium (B/C), Wide (D/E), X-Wide (2E/3E)	4 1/2-11, 12, 13

### Men's Ambulator® Active Walkers

- Genuine leather uppers with moisture wicking foam backed fabric linings for cushioning & protection
- 1/2" removable hidden depth in 2 layers
- High toe boxes

Description	Product No.
Black Lace	<b>G7000M</b>
Black Double Strap	<b>G8000M</b>
White Lace	<b>G7200M</b>
White Double Strap	<b>G8200M</b>

Ambulator Active Walker Sizing Chart		
Shoe	Widths	Sizes
Men's	Medium (C/D), Wide (E/2E), X-Wide (3E/4E)	6 1/2-13, 14, 15, 16

**Ordering Example**

Model ..... **1260M**  
 Width..... **W**  
 Size..... **10**  
**1260MW10**



1260M/1261M



1270M/1271M



1260W/1264W



1270W



G8000M/G8200M



G7000M/G7200M

Shoes

**Ambulator® Stretchables**

- Seamless spandex uppers with Plastazote® foam lining molds to the foot for a customized fit
- Plastazote insole for additional protection
- 1/2" removable hidden depth in 2 layers
- High toe boxes
- Unisex shoe; specify gender, size & width when ordering

Description	Product No.
Black Double Strap	<b>1200</b>
Brown Single Strap	<b>1201</b>
Black Triple Strap Boot	<b>1209</b>

Widths	Sizes	
	Men's	Women's
<b>Stretchable Double Strap Sizing Chart</b>		
Medium	3 1/2-10 1/2, 11, 12	—
Wide	5-13, 14, 15, 16	4 1/2-11, 12, 13
X-Wide	6-13, 14, 15, 16	6-14
2X-Wide	—	7-13
<b>Stretchable Single Strap Sizing Chart</b>		
Medium	3 1/2-12	—
Wide	4 1/2-14	4 1/2-13
X-Wide	—	5 1/2-14
<b>Stretchable Triple Strap Sizing Chart (sold each)</b>		
Medium	2 1/2-12	—
Wide	4 1/2-14	4 1/2-13
X-Wide	5 1/2-16	6 1/2-14
2X-Wide	—	5 1/2-15



1200M/1200W - Double Strap



1201M/1201W - Single Strap



1209M /1209W- Triple Strap

**Hush Puppies - Women's Achieve**

- Casual walking shoe
- Flexible rubber sole
- Bounce technology for superior shock absorption and comfort
- Drilex® lining for odor and stain control

Color	Product No.
Black	<b>H502316</b>
Dark Brown	<b>H502317</b>

<b>Hush Puppies Achieve Sizing Chart</b>	
Sizes	Widths
7-10, 11	N
5-10, 11, 12	M
6-10, 11, 12	W
6-10, 11	EW



**Ordering Example**

Model ..... **H502317**  
 Size ..... **9**  
 Width ..... **N**  
**H50231709N**

**Women's Powerwalker and Powerwalker II**

- Walking shoe available in lace or hook and loop
- Flexible rubber sole
- Drilex® lining for odor and stain control

Color	Product No.	
	Powerwalker Lace	Powerwalker II Hoop & Loop
White	<b>H54025</b>	<b>H70293</b>
Black	<b>H54026</b>	<b>H70294</b>
Taupe	<b>H54027</b>	<b>H70295</b>

Powerwalker and Powerwalker II Sizing Chart	
Sizes	Widths
7-10, 11	N
5-10, 11, 12	M
6-10, 11, 12	W
6-10, 11	EW



*H54026  
Lace*



*H70295  
Hook & Loop*

**Ordering Example**  
 Model ..... **H54026**  
 Size ..... **10**  
 Width ..... **M**  
**H5402610M**

Shoes

**Men's Gil**

- Modified Oxford style leather casual shoe
- Hook and loop closure

Color	Product No.
Black	<b>H18800</b>
Brown	<b>H18801</b>
Beige	<b>H18802</b>

Hush Puppies Gil Sizing Chart	
Sizes	Widths
7-12, 13, 14, 15, 16	M
7-12, 13, 14, 15	W
7-12, 13, 14	EW



**Ordering Example**  
 Model ..... **H18802**  
 Size ..... **13**  
 Width ..... **W**  
**H1880213W**

**Men's Gus**

- Weightless design combined with dual density
- Elon® midsole and blended rubber outsole for durability

Color	Product No.
Black	<b>H18770</b>
Brown	<b>H18772</b>
Beige	<b>H18773</b>



H18770

**Men's Glen**

- ZeroG® system reduces weight without sacrificing comfort
- Elon midsole and blended rubber outsole for durability

Color	Product No.
Black	<b>H19072</b>
Brown	<b>H19078</b>



H19072

**Men's Burke**

- Stylish casual comfort shoe
- Strobel construction and a contoured bottom allow for flexibility and comfort
- Available in Black

Product No. **H101407**



H101407

**Hush Puppies Gus, Glen and Burke Sizing Chart**

Sizes	Widths
9-12, 13	N
7-12, 13, 14, 15, 16	M
7-12, 13, 14, 15	W
7-12, 13, 14	EW

**Ordering Example**

Model ..... **H18770**  
 Size ..... **12**  
 Width ..... **EW**  
**H1877012EW**

**Mt. Emey - Men's 9501 Dual-Depth Leather Dress Shoe (Lace)**

- Feather light shoes with 3 layers of removable inserts
- High Density Unit EVA Flex Outsole
- Available in Black

Product No. **9501BK**

**Men's 9501 Dual-Depth Shoe Sizing Chart**

Sizes	Widths
5-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17, 18	D, 2E
6-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17, 18, 19, 20	4E
7-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17, 18, 19, 20	6E
7-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17	9E



9501

**Ordering Example**

Model ..... **9501BK**  
 Size ..... **17**  
 Width ..... **D**  
**9501BK17D**



**Men's 801/802 Supra-Depth Shoe**

- Extra soft genuine Napa leather uppers
- Seamless lining
- Extra length counter support
- Direct injected Flex PU outside

Color	Product No.	
	801 (Lace)	802 (Hook & Loop)
Black	<b>801BK</b>	<b>802BK</b>
Brown	<b>801BR</b>	<b>802BR</b>

Men's 801/802 Supra-Depth Shoe Sizing Chart	
Sizes	Widths
5-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17	D
6-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17	4E, 6E, 9E



801 Lace



802  
Hook and Loop

**Ordering Example**  
 Model ..... **801BK**  
 Size ..... **5**  
 Width ..... **D**  
**801BK05D**

**Men's 708 Dual-Depth Leather (Lace) and 718 Dual-Depth Leather (Hook and Loop)**

- Genuine leather upper
- HD EVA sole
- Extended firm heel counter
- Stretchable vamp
- Bunion/bunionette flex area

Color	Product No.	
	708 (Lace)	718 (Hook & Loop)
Black	<b>708-BLK</b>	<b>718-BLK</b>
Beige	<b>708-BEIG</b>	<b>718-BEIG</b>

Men's 708 and 718 Dual-Depth Shoe Sizing Chart	
Sizes	Widths
5-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17	D
6-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17	4E
7-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17	6E



708 Lace



718  
Hook & Loop

**Ordering Example**  
 Model ..... **708-BLK**  
 Width ..... **D**  
 Size ..... **8**  
**708-BLK-D8**

Shoes

**Men's 728 Bunion Edema Shoes (Hook and Loop)**

- Seamless vamp lining
- Soft leather, fabric lining
- Re-heat-moldable toe-puff and long counter
- Padded collar; Available in Black

Product No. **728-E-BLK**

Men's 728 Dual-Depth Shoe Sizing Chart	
Sizes	Widths
5-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17	D
6-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17	4E
7-11 <sup>1</sup> / <sub>2</sub> , 12, 13, 14, 15, 16, 17	6E, 10E, 14E



**Ordering Example**

Model . . . . . **728-E-BLK**  
 Width . . . . . **6E**  
 Size . . . . . **12**  
**728-E-BLK-6E12**

**Women's 618 Dual-Depth Leather (Hook and Loop)**

- Lycro shoes with HD EVA/RBR sole
- Extended firm heel counter
- Stretchable vamp
- Bunion/bunionette flex area

Color	Product No.
Black	<b>618-BLK</b>
Beige	<b>618-BEIG</b>

Women's 618 Dual-Depth Shoe Sizing Chart	
Sizes	Widths
4 <sup>1</sup> / <sub>2</sub> -10 <sup>1</sup> / <sub>2</sub> , 11, 12, 13, 14, 15	B, D, 3E, 5E



**Ordering Example**

Model . . . . . **618-BLK**  
 Width . . . . . **3E**  
 Size . . . . . **12**  
**618-BLK-3E12**

**Women's 628 Dual-Depth Leather (Hook and Loop)**

- Leather counters, single strap
- Extended firm heel counter
- Stretchable vamp
- Bunion/bunionette flex area; Available in Black

Product No. **628-E-BLK**

Women's 628 Dual-Depth Shoe Sizing Chart	
Sizes	Widths
4 <sup>1</sup> / <sub>2</sub> -10 <sup>1</sup> / <sub>2</sub> , 11, 12, 13, 14, 15	B, D, 3E, 5E, 10E, 14E



**Ordering Example**

Model . . . . . **628-E-BLK**  
 Width . . . . . **3E**  
 Size . . . . . **11**  
**628-E-BLK-3E11**

**Women's 9205 Added-Depth Extrem-Light Mary Jane**

- Seamless linings
- Added width at the toe & extra support at the ankle
- Extended lateral counter provide support & stability
- High density injected EVA outsole with symmetrical cut line for easy sole modifications (with 2 layer of: 1/8" base layer + top 5/16" insert #308 inside shoes)
- B (M) width is defined by including 2 layer removable inserts; D (W) width is defined by taking out the bottom insert



Color	Product No.
Black	<b>9205-BLK</b>
Beige	<b>9205-BGE</b>
Brown	<b>9205-BRN</b>
Gray	<b>9205-GRY</b>
Jade	<b>9205-JADE</b>
Navy	<b>9205-NAVY</b>
Red	<b>9205-RED</b>
White	<b>9205-WHT</b>
Sky	<b>9205-SKY</b>
Wine	<b>9205-WINE</b>
Taupe	<b>9205-TPE</b>

Ordering Example	
Model . . . . .	<b>9205-BLK</b>
Width . . . . .	<b>3E</b>
Size . . . . .	<b>8</b>
<b>9205-BLK-3E8</b>	

Women's 9205 Shoe Sizing Chart	
Sizes	Widths
5-10 <sup>1</sup> / <sub>2</sub> , 11	B/D (M/W), 3E (WW), 5E (XW)

**Women's 9301 Dual-Depth Leather Comfort Shoes**

- Double depth ladies' feather light shoes with 3 layers of removable inserts
- Extra soft leather uppers
- Durable HD Pure EVA outsole, easy to work on, such as Build-up, Rocker etc.



Color	Product No.
Black	<b>9301-BLK</b>
Beige	<b>9301-BEIG</b>
White	<b>9301-WHT</b>

Ordering Example	
Model . . . . .	<b>9301-BLK</b>
Width . . . . .	<b>3E</b>
Size . . . . .	<b>10</b>
<b>9301-BLK-3E10</b>	

Women's 9301 Shoe Sizing Chart	
Sizes	Widths
4 <sup>1</sup> / <sub>2</sub> -10 <sup>1</sup> / <sub>2</sub> , 11, 12, 13, 14, 15	B, D, 3E
5-10 <sup>1</sup> / <sub>2</sub> , 11, 12, 13, 14, 15	5E, 7E

Shoes

**Women's 9106 and 9108 Supra Depth Comfort Shoe**

- Leather upper and lining
- Triple depth with broad oblique toe
- Deep toe box for severe deformities

Color	Product No.	
	9106 (Hook & Loop)	9108 (Lace)
Black	<b>9106BK</b>	<b>9108BK</b>
Beige	<b>9106BE</b>	<b>9108BE</b>

Women's 9106 and 9108 Comfort Shoe Sizing Chart	
Sizes	Widths
4 <sup>1</sup> / <sub>2</sub> -11, 12	2A
4 <sup>1</sup> / <sub>2</sub> -10 <sup>1</sup> / <sub>2</sub> , 11, 12, 13, 14, 15	B, D, 3E
5-10 <sup>1</sup> / <sub>2</sub> , 11, 12, 13, 14, 15	5E, 7E



9106  
Hook & Loop



9108  
Lace

**Ordering Example**  
 Model . . . . . **9106BE**  
 Size . . . . . **6**  
 Width . . . . . **2A**  
**9106BE062A**

**New Balance - 577 Women's Walking Shoe**

- Durable walking shoe with synthetic upper
- Absorb® heel for added comfort

Color	Lace	Hook & Loop
Black	<b>WW577BK</b>	<b>WW577VK</b>
White	<b>WW577WT</b>	<b>WW577VW</b>
Beige	—	<b>WW577VB</b>

Women's New Balance 577 Sizing Chart			
Color	Sizes	Widths	
		Lace	Hook & Loop
Black	6-11, 12	2A, 2E	—
Black	5-11, 12	B, D	2A, B, D, 2E
White	5-11, 12	2A, B, D, 2E	2A, B, D, 2E
Beige	5-11, 12	—	2A, B, D, 2E



**Ordering Example**  
 Model . . . . . **WW577VB**  
 Size . . . . . **11**  
 Width . . . . . **B**  
**WW577VB11B**

### 812 Women's Performance Walking Shoe

- Man-made synthetic material
- Light midsole and outsole cushion increases durability and comfort
- ROLLBAR® graphite technology for rearfoot stability
- Lace closure

Women's New Balance 812 Sizing Chart			
Color	Sizes	Widths	Product No.
Black	6-12, 13	2A, 2E, 4E	WW812BK
	5-12, 13	B, D	
White	6-12, 13	2A, 2E, 4E	WW812WT
	5-12, 13	B, D	
Beige	5-12, 13	B, D	WW812BE
	6-12, 13	2A, 2E	



**Ordering Example**  
 Model ..... **WW812WT**  
 Size ..... **13**  
 Width ..... **D**  
**WW812WT13D**

### 812 Men's Performance Walking Shoe

- Light midsole and outsole cushion for durability and comfort
- ROLLBAR® graphite technology for rearfoot stability

Color	Product No.	
	Lace	Hook & Loop
Black	MW812BK	MW812VK
White	MW812WT	MW812VW
Bone	MW812BE	MW812VB

Men's New Balance 812 Sizing Chart			
Color	Sizes	Widths	
		Lace	Hook & Loop
Black	9-12, 13, 14, 15, 16	B	B
	7-12, 13, 14, 15, 16	D, 2E, 4E, 6E	D, 2E, 4E
White	9-12, 13, 14, 15, 16	2A, B	D, 2E, 4E
	7-12, 13, 14, 15, 16	D, 2E, 4E, 6E	D, 2E, 4E
Beige	7-12, 13, 14, 15	D, 2E, 4E	D, 2E, 4E



MW812BK



MW812VW

**Ordering Example**  
 Model ..... **MW812BK**  
 Size ..... **10**  
 Width ..... **B**  
**MW812BK10B**

Shoes

**927 Women's**

- ROLLBAR® graphite technology provides stability and minimal flex
- Designed for over pronation, motion control with a straight last

Color	Product No.	
	Lace	Hook & Loop
Black	<b>WW927BK</b>	—
White	<b>WW927WT</b>	<b>WW927VW</b>

Women's New Balance 927 Sizing Chart	
Sizes	Widths
6-11, 12	2A
5-11, 12, 13	B, D, 2E, 4E



WW927WT



WW927VW  
not available  
in 2A and 4E

**Ordering Example**  
 Model . . . . . **WW927WT**  
 Size . . . . . **9**  
 Width . . . . . **B**  
**WW927WT09B**

**927 Men's**

- Full leather design
- Straight last for over-pronation and motion control
- ROLLBAR® graphite technology provides excellent stability
- Black hook and loop also available (sizes D, 2E, 4E and 6E only)

Color	Product No.	
	Lace	Hook & Loop
Black	<b>MW927BK</b>	<b>MW927VK*</b>
White	<b>MW927WT</b>	—

\*MW927VK not available in B width

Men's New Balance 927 Sizing Chart	
Sizes	Widths
9-12, 13, 14, 15, 16	B
7-13, 14, 15, 16	D, 2E, 4E, 6E



MW927WT

**Ordering Example**  
 Model . . . . . **MW927WT**  
 Size . . . . . **10**  
 Width . . . . . **2E**  
**MW927WT102E**

## SUREFIT - WOMEN'S SHOES

### Monaco

- Leather construction with stitch detailing
- Classic Mary Jane design
- Hook and loop closure

### Acapulco

- Light flexible shoe constructed of Nubuck upper with leather lining for comfort and durability
- Hook and loop closure

### Sydney

- Mary Jane style with hook and loop closure
- Accommodates claw, hammer toes and bunions for that special fit



Acapulco



Monaco



Sydney

Color	Sizes	Widths	Product No.		
			Monaco	Acapulco	Sydney
Black	5-12	M, W, WW	<b>S656-1</b>	<b>S606-1</b>	<b>S636-1</b>
Brown	5-12	M, W, WW	<b>S656-2</b>	—	—
Camel	5-12	M, W, WW	—	<b>S606-2</b>	—
Navy	5-12	M, W, WW	—	<b>S606-7</b>	—
Gray	5-12	M, W, WW	—	<b>S606-8</b>	—

Ordering Example	
Model .....	<b>S636-1</b>
Size .....	<b>8</b>
Width .....	<b>WW</b>
<b>S636-108WW</b>	

### Manila (Lace) and Athens (Hook & Loop)

- Leather lining with light perforations to keep cool
- Great shoe for everyday walking

Color	Sizes	Widths	Product No.	
			Manila	Athens
Black	5-12, 13	M, W, WW	<b>S125-1</b>	<b>S136-1</b>
White	5-12, 13	M, W, WW	<b>S125-3</b>	<b>S136-3</b>
Beige	5-12, 13	M, W, WW	<b>S125-4</b>	<b>S136-4</b>
Taupe	5-12, 13	M, W, WW	<b>S125-5</b>	<b>S136-5</b>



Manila



Athens

Ordering Example	
Model .....	<b>S125-1</b>
Size .....	<b>10</b>
Width .....	<b>W</b>
<b>S125-110W</b>	

Shoes

**SUREFIT - MEN'S SHOES**

**Copenhagen**

- Overlay of leather overlay for a traditional gentleman's look
- Design features low counters and ample toe room
- Hook and loop closure



**London**

- Full grain leather upper for durability and easy cleaning
- Soft leather liner for added comfort
- Hook and loop closure



Color	Sizes	Widths	Product No.	
			Copenhagen	London
Black	7-13, 14, 15	M, W, WW	<b>S206-1</b>	<b>S196-1</b>
Brown	7-13, 14, 15	M, W, WW	<b>S206-2</b>	—
Camel	7-13, 14, 15	M, W, WW	—	<b>S196-2</b>

Ordering Example	
Model .....	<b>S206-1</b>
Size .....	<b>12</b>
Width .....	<b>M</b>
<b>S206-112M</b>	

**Rome**

- Classic Oxford look made of Napa leather with a wide toe box
- Added depth can be achieved when insole layers are removed
- Hook and loop closure



Color	Sizes	Widths	Product No.
Black	7-13, 14, 15, 16	M, W, WW	<b>S116-1</b>
Brown	7-13, 14, 15, 16	M, W, WW	<b>S116-2</b>

Ordering Example	
Model .....	<b>S116-2</b>
Size .....	<b>9</b>
Width .....	<b>WW</b>
<b>S116-209WW</b>	



**Barbados Boat Shoe**

- Classic boat shoe style with added comfort
- Lace closure



Barbados



Milan

**Milan**

- Handcrafted leather with whip stitching and rounded seams
- Exceptional casual comfort shoe
- Lace closure

**Nassau**

- Supple soft upper and lightweight soles absorb and reduce shock
- Spandex panels capture and protect the foot



Nassau



Zurich

**Zurich Boot**

- Durable oiled leather boot
- Comfortable fit with padded collars
- Firm heel counter
- Lace closure

Color	Sizes	Widths	Product No.			
			Barbados	Milan	Nassau	Zurich
Black	7-13, 14, 15	M, W, WW	—	S215-1	S146-1	S155-1
Brown	7-13, 14, 15	M, W, WW	S175-2	S215-2	S146-2	S155-R
Cordovan	7-13, 14, 15	M, W, WW	—	S215-C	—	—

Relief and  
Healing Shoes

**RELIEF AND HEALING SHOES**

**Bauerfeind - GloboPed® Heel Relief Shoe**

- Off-loading diabetic shoe reduces convalescent time for early patient mobility
- 10° incline built into rear puts foot in plantar flexion to offload heel area
- Three adjustable closure straps with cutout in heel and forefoot for breathability and enhanced healing
- Black denim-like fabric



Size	Foot Length inches	US Shoe Size		Product No.
		Men's	Women's	
Small	8 <sup>5</sup> / <sub>8</sub> -9 <sup>3</sup> / <sub>8</sub>	—	6-8 <sup>1</sup> / <sub>2</sub>	36628000010000
Medium	9 <sup>3</sup> / <sub>8</sub> -10 <sup>3</sup> / <sub>8</sub>	6-9	9-11	36628000020000
Large	10 <sup>3</sup> / <sub>8</sub> -11 <sup>3</sup> / <sub>8</sub>	10-12	12-13	36628000030000
X-Large	11 <sup>3</sup> / <sub>8</sub> -12 <sup>3</sup> / <sub>8</sub>	13-15	—	36628000040000

**GloboPed Forefoot Relief Shoe**

- Reduced convalescent time for early patient mobility
- Three straps for optimal fit, soft inside lining with no seams
- Optional shield is made from heat moldable ABS and prevents undesired contact with foreign objects
- Black denim-like fabric

Size	Foot Length inches	US Shoe Size		Product No.
		Men's	Women's	
Small	8 <sup>5</sup> / <sub>8</sub> -9 <sup>3</sup> / <sub>8</sub>	—	6-8 <sup>1</sup> / <sub>2</sub>	<b>36627000010000</b>
Medium	9 <sup>7</sup> / <sub>8</sub> -10 <sup>3</sup> / <sub>8</sub>	6-9	9-11	<b>36627000020000</b>
Large	10 <sup>5</sup> / <sub>8</sub> -11 <sup>3</sup> / <sub>8</sub>	10-12	12-13	<b>36627000030000</b>
X-Large	11 <sup>3</sup> / <sub>8</sub> -12 <sup>3</sup> / <sub>8</sub>	13-15	—	<b>36627000040000</b>

Protective Shield	Product No.
Small/Medium	<b>36624900020000</b>
Large/X-Large	<b>36624900030000</b>



**Darco - Body Armor® Cast Shoe**

- Shields the cast from moisture, shock and everyday wear
- EVA unibody design provides shock absorption
- Bungee closure system is easy to use and adjust; secures the cast in the shoe
- Built-in toe guard protects toes and prevents debris from entering the toe area

Size	US Shoe Size		Product No.
	Men's	Women's	
X-Small	—	up to 5	<b>BACS0</b>
Small	—	5 <sup>1</sup> / <sub>2</sub> -8	<b>BACS1</b>
Medium	6-8	8 <sup>1</sup> / <sub>2</sub> -10	<b>BACS2</b>
Large	8 <sup>1</sup> / <sub>2</sub> -11	10 <sup>1</sup> / <sub>2</sub> -13	<b>BACS3</b>
X-Large	11 <sup>1</sup> / <sub>2</sub> -14	13 <sup>1</sup> / <sub>2</sub>	<b>BACS4</b>



**Softie™ Shoe**

- Semi-rigid metatarsal shank provides superior control of forefoot motion to ensure postoperative healing remains undisturbed
- Better protection after osteotomies.
- Tri-laminated upper made of light weight tricot-foam material is breathable and comfortably conforms to bony abnormalities
- The hook and loop straps are easy to fasten and unfasten with one hand

Size	US Shoe Size		Product No.	
	Men's	Women's	Men's	Women's
Small	6-8	4-6	<b>STM1B</b>	<b>STW1B</b>
Medium	8 <sup>1</sup> / <sub>2</sub> -10	6 <sup>1</sup> / <sub>2</sub> -8	<b>STM2B</b>	<b>STW2B</b>
Large	10 <sup>1</sup> / <sub>2</sub> -12	8 <sup>1</sup> / <sub>2</sub> -10	<b>STM3B</b>	<b>STW3B</b>
X-Large	12 <sup>1</sup> / <sub>2</sub> -14	—	<b>STM4B</b>	—



### HeelWedge™ and OrthoWedge Shoe

- Square-toe design acts as a bumper to protect the toes when k-wires are present and provides better left/right fit
- Ankle strap seats foot firmly in shoe to reduce heel slippage and eliminate friction
- Removable forefoot closure provides security while eliminating buckle pressure
- HeelWedge is designed to relieve pressure from the heel areas and OrthoWedge is designed to relieve pressure from the forefoot and toe areas



HeelWedge



OrthoWedge

Size	US Shoe Size		Product No.	
	Men's	Women's	HeelWedge	OrthoWedge
X-Small	—	4-7	<b>HQ0B</b>	<b>OQ0B</b>
Small	6-8	7 <sup>1</sup> / <sub>2</sub> -10	<b>HQ1B</b>	<b>OQ1B</b>
Medium	8 <sup>1</sup> / <sub>2</sub> -10	10 <sup>1</sup> / <sub>2</sub> -13	<b>HQ2B</b>	<b>OQ2B</b>
Large	10 <sup>1</sup> / <sub>2</sub> -12	13+	<b>HQ3B</b>	<b>OQ3B</b>
X-Large	12 <sup>1</sup> / <sub>2</sub> -14	—	<b>HQ4B</b>	<b>OQ4B</b>

### Wound Care Shoe System

- Four multi-density insoles
- Ultra-soft Plastazote® Lining provides a virtually friction-free contact surface that won't exacerbate existing wounds
- Genuine leather upper can be modified or cut away to provide even greater levels of pressure relief
- Sold in pairs

Size	US Shoe Size		Product No.
	Men's	Women's	
X-Small	—	4-5 <sup>1</sup> / <sub>2</sub>	<b>WCS0B</b>
Small	6-7 <sup>1</sup> / <sub>2</sub>	6-8 <sup>1</sup> / <sub>2</sub>	<b>WCS1B</b>
Small/Medium	8-8 <sup>1</sup> / <sub>2</sub>	9-9 <sup>1</sup> / <sub>2</sub>	<b>WCS7B</b>
Medium	9-9 <sup>1</sup> / <sub>2</sub>	10-10 <sup>1</sup> / <sub>2</sub>	<b>WCS2B</b>
Large	10-11 <sup>1</sup> / <sub>2</sub>	11-12 <sup>1</sup> / <sub>2</sub>	<b>WCS3B</b>
X-Large	12-13 <sup>1</sup> / <sub>2</sub>	13-14 <sup>1</sup> / <sub>2</sub>	<b>WCS4B</b>
2X-Large	14-16	15-17	<b>WCS5B</b>



Relief and  
Healing Shoes

### Original Med-Surg™ Shoe

- Patented metatarsal shank provides superior control of forefoot motion to ensure post operative healing remains undisturbed
- Better protection after osteotomies
- Reinforced padded heel grips the foot and reduces slippage
- Breathable mesh upper light weight material allows for proper ventilation and conforms to bony abnormalities



Size	US Shoe Size			Product No.		
	Pediatric	Men's	Women's	Pediatric	Men's	Women's
Pediatric	12-1	—	—	<b>MSPN</b>	—	—
Small	—	6-8	4-6	—	<b>MSM1N</b>	<b>MSW1N</b>
Medium	—	8 <sup>1</sup> / <sub>2</sub> -10	6 <sup>1</sup> / <sub>2</sub> -8	—	<b>MSM2N</b>	<b>MSW2N</b>
Large	—	10 <sup>1</sup> / <sub>2</sub> -12	8 <sup>1</sup> / <sub>2</sub> -10	—	<b>MSM3N</b>	<b>MSW3N</b>
X-Large	—	12 <sup>1</sup> / <sub>2</sub> -14	—	—	<b>MSM4N</b>	—

## Med-Surg™ Shoe

- Square-toe design acts as a bumper to protect the toes when k-wires are present and provides better left/right fit
- MetaShank™ protection provides more rigid control under the metatarsal heads
- Strapless closure eliminates buckle pressure on the forefoot and can expand to accommodate dressing up to 17 inches in circumference
- Fully compatible with the PegAssist™ off-loading insole



Size	US Shoe Size			Product No.		
	Pediatric	Men's	Women's	Pediatric	Men's	Women's
Pediatric	12-1	—	—	<b>MQPB</b>	—	—
Small	—	6-8	4-6	—	<b>MQM1B</b>	<b>MQW1B</b>
Medium	—	8 <sup>1</sup> / <sub>2</sub> -10	6 <sup>1</sup> / <sub>2</sub> -8	—	<b>MQM2B</b>	<b>MQW2B</b>
Large	—	10 <sup>1</sup> / <sub>2</sub> -12	8 <sup>1</sup> / <sub>2</sub> -10	—	<b>MQM3B</b>	<b>MQW3B</b>
X-Large	—	12 <sup>1</sup> / <sub>2</sub> -14	—	—	<b>MQM4B</b>	—

## APB™ All-Purpose Boot

- Square-toe design acts as a bumper to protect the toes when k-wires are present and provides better left/right fit
- Ankle strap seats foot firmly in shoe to reduce heel slippage and eliminate friction
- Removable forefoot closure provides security while eliminating buckle pressure
- Unique metatarsal shank under the MP joints increases support and rigidity under the forefoot



Size	US Shoe Size		Product No.
	Men's	Women's	
X-Small	—	4-7	<b>APQ0B</b>
Small	6-8	7 <sup>1</sup> / <sub>2</sub> -10	<b>APQ1B</b>
Medium	8 <sup>1</sup> / <sub>2</sub> -10	10 <sup>1</sup> / <sub>2</sub> -13	<b>APQ2B</b>
Large	10 <sup>1</sup> / <sub>2</sub> -12	13+	<b>APQ3B</b>
X-Large	12 <sup>1</sup> / <sub>2</sub> -14	—	<b>APQ4B</b>

## SlimLine™ Cast Boot

- Square-toe design acts as a protective bumper and gives the patient more room to accommodate larger casts and bandages
- Forefoot closure can be adjusted to accommodate casts or bandages and can be switched between left/right closure
- Higher ankle strap adds stability
- Elasticized strap allows for fixed ankle ambulation
- Rocker sole allows for a smooth roll throughout the gait pattern



Size	US Shoe Size		Product No.
	Men's	Women's	
X-Small	—	up to 5	<b>SLQ0B</b>
Small	—	5 <sup>1</sup> / <sub>2</sub> -8	<b>SLQ1B</b>
Medium	6-8	8 <sup>1</sup> / <sub>2</sub> -10	<b>SLQ2B</b>
Large	8 <sup>1</sup> / <sub>2</sub> -11	10 <sup>1</sup> / <sub>2</sub> -13	<b>SLQ3B</b>
X-Large	11 <sup>1</sup> / <sub>2</sub> -14	13 <sup>1</sup> / <sub>2</sub>	<b>SLQ4B</b>

**Body Armor™ Toe Guard**

- Hard ABS construction is durable and withstands blows and strikes to the toe region in everyday situations
- Moldable side wings can be bent/molded to conform to any foot size
- Strong adhesive backing ensures secure, lasting fit
- Round-toe or square-toe versions available



Description	Quantity	Product No.	
		Men's	Women's
Square	6 pair per box	<b>BATGQ2</b>	<b>BATGQ1</b>
Round	each	<b>BATGS2</b>	<b>BATGS1</b>

**PegAssist™**

- Removable pegs allow for localized off-loading of wounds and ulcerations of the foot
- Poron® cover eliminates the incidence of ring edema and edge abrasion
- Stabilizer board prevents adjacent pegs from collapsing



Size	US Shoe Size		With Toe Cover		Insole Only
	Men's	Women's	Men's	Women's	
Small	6-8	4-6	<b>PTQM1</b>	<b>PTQW1</b>	<b>PQ1</b>
Medium	8 <sup>1</sup> / <sub>2</sub> -10	6 <sup>1</sup> / <sub>2</sub> -8	<b>PTQM2</b>	<b>PTQW2</b>	<b>PQ2</b>
Large	10 <sup>1</sup> / <sub>2</sub> -12	8 <sup>1</sup> / <sub>2</sub> -10	<b>PTQM3</b>	<b>PTQW3</b>	<b>PQ3</b>
X-Large	12 <sup>1</sup> / <sub>2</sub> -14	—	<b>PTQM4</b>	—	<b>PQ4</b>

Relief and  
Healing Shoes

**DeRoyal - Comfo-EZE® Rocker Cast Shoe**

- Constructed of durable navy canvas material
- Non-slip rubber sole for patient safety; rocker bottom promotes rolling gait
- Heel strap helps prevent foot migration
- Open toe design



Size	Length	Width	Product No.
Pediatric	5 <sup>1</sup> / <sub>2</sub> "	2 <sup>1</sup> / <sub>2</sub> "	<b>2004-00</b>
X-Small	7"	3 <sup>1</sup> / <sub>2</sub> "	<b>2004-01</b>
Small	8 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	<b>2004-02</b>
Medium	10"	4 <sup>3</sup> / <sub>4</sub> "	<b>2004-03</b>
Large	11 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>4</sub> "	<b>2004-04</b>
X-Large	13"	5 <sup>1</sup> / <sub>2</sub> "	<b>2004-05</b>

## Deluxe Wooden Post-Op Shoe

- Rigid wood sole with rubber bottom for support
- Perforated suede upper for comfort and durability
- Universal left or right

Size	US Shoe Size		Product No.	
	Men's	Women's	Men's	Women's
<b>Hook and Loop Closure</b>				
Child	2-4 <sup>1</sup> / <sub>2</sub>	2-3 <sup>1</sup> / <sub>2</sub>	<b>2007-01</b>	<b>2007-06</b>
Small	5-6 <sup>1</sup> / <sub>2</sub>	4-6	<b>2007-02</b>	<b>2007-07</b>
Medium	7-8 <sup>1</sup> / <sub>2</sub>	6 <sup>1</sup> / <sub>2</sub> -8	<b>2007-03</b>	<b>2007-08</b>
Large	9-11	8 <sup>1</sup> / <sub>2</sub> -10	<b>2007-04</b>	<b>2007-09</b>
X-Large	11 <sup>1</sup> / <sub>2</sub> -13	10 <sup>1</sup> / <sub>2</sub> -12	<b>2007-05</b>	<b>2007-10</b>
<b>Lace Up Closure</b>				
Child	2-4 <sup>1</sup> / <sub>2</sub>	2-3 <sup>1</sup> / <sub>2</sub>	<b>2006-01</b>	<b>2006-06</b>
Small	5-6 <sup>1</sup> / <sub>2</sub>	4-6	<b>2006-02</b>	<b>2006-07</b>
Medium	7-8 <sup>1</sup> / <sub>2</sub>	6 <sup>1</sup> / <sub>2</sub> -8	<b>2006-03</b>	<b>2006-08</b>
Large	9-11	8 <sup>1</sup> / <sub>2</sub> -10	<b>2006-04</b>	<b>2006-09</b>
X-Large	11 <sup>1</sup> / <sub>2</sub> -13	10 <sup>1</sup> / <sub>2</sub> -12	<b>2006-05</b>	<b>2006-10</b>



Hook and Loop Closure



Lace-up Closure

## Canvas Post-Op Shoe

- Rigid rubber sole provides controlled forefoot motion
- Canvas upper for durability
- Hook and loop closure
- Universal left or right

Size	US Shoe Size			Product No.		
	Child	Men's	Women's	Child	Men's	Women's
Child	6-8	—	—	<b>2041-00</b>	—	—
Small	—	6-8	4-6	—	<b>2041-01</b>	<b>2041-07</b>
Medium	—	8 <sup>1</sup> / <sub>2</sub> -10	6 <sup>1</sup> / <sub>2</sub> -8	—	<b>2041-02</b>	<b>2041-08</b>
Large	—	10 <sup>1</sup> / <sub>2</sub> -12	8 <sup>1</sup> / <sub>2</sub> -10	—	<b>2041-03</b>	<b>2041-09</b>



## Canvas Cast Shoe

- Constructed of navy canvas material
- Padded canvas footbed; Non-slip rubber sole
- Heel strap helps prevent foot migration
- Hook and loop closure for easy application

Size	Length	Width	Product No.
Child	5 <sup>1</sup> / <sub>2</sub> "	2 <sup>1</sup> / <sub>2</sub> "	<b>2035-01</b>
Small	8 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	<b>2035-02</b>
Medium	10"	4 <sup>3</sup> / <sub>4</sub> "	<b>2035-03</b>
Large	11 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>4</sub> "	<b>2035-04</b>
X-Large	13"	5 <sup>1</sup> / <sub>2</sub> "	<b>2035-05</b>



**Vinyl Cast Boot**

- Constructed of durable vinyl material
- Rocker bottom promotes rolling gait
- Open toe with closed heel design prevents foot migration
- Hook and loop closure for easy application

Size	Length	Width	Product No.
X-Small	7"	3 <sup>1</sup> / <sub>2</sub> "	<b>2027-01</b>
Small	8 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	<b>2027-02</b>
Medium	10"	4 <sup>3</sup> / <sub>4</sub> "	<b>2027-03</b>
Large	11 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>4</sub> "	<b>2027-04</b>
X-Large	13"	5 <sup>1</sup> / <sub>2</sub> "	<b>2027-05</b>



**Aegis Post-Op Shoe**

- Semi-rigid rubber sole for support with reinforced padded heel counter
- Hook and loop closure for easy application
- Universal left or right

Size	US Shoe Size		Product No.	
	Men's	Women's	Men's	Women's
Small	6-8	4-6	<b>2044-01</b>	<b>2044-07</b>
Medium	8 <sup>1</sup> / <sub>2</sub> -10	6 <sup>1</sup> / <sub>2</sub> -8	<b>2044-02</b>	<b>2044-08</b>
Large	10 <sup>1</sup> / <sub>2</sub> -12	8 <sup>1</sup> / <sub>2</sub> -10	<b>2044-03</b>	<b>2044-09</b>
X-Large	12 <sup>1</sup> / <sub>2</sub> -14	—	<b>2044-04</b>	—




Relief and  
Healing Shoes

**Hely & Weber - Molded Sole Cast Shoe**

- Canvas upper for strength; open toe
- Two loop lock closures secure cast to shoe
- Molded sole does not breakdown when using fiberglass casting tape
- Adjustable hook and loop strap on back for a snug fit around heel

Size	A	B	Product No.
Pediatric	7 <sup>1</sup> / <sub>4</sub> "	3 <sup>11</sup> / <sub>16</sub> "	<b>364PED</b>
Small	9 <sup>3</sup> / <sub>8</sub> "	4 <sup>5</sup> / <sub>16</sub> "	<b>364S</b>
Medium	11"	4 <sup>7</sup> / <sub>8</sub> "	<b>364M</b>
Large	12 <sup>1</sup> / <sub>2</sub> "	5 <sup>3</sup> / <sub>8</sub> "	<b>364L</b>
X-Large	13 <sup>9</sup> / <sub>16</sub> "	5 <sup>3</sup> / <sub>8</sub> "	<b>364XL</b>



## Wooden Post-Op Shoe

- Nylon mesh upper with open toe design
- Padded heel and soft tricot liner for greater comfort
- Wooden bottom with slip-resistant sole
- Loop & lock closures with movable hook/loop tabs provide unlimited strap adjustability

Size	US Shoe Size		Product No.
	Men's	Women's	
Female Small	—	5-7 <sup>1</sup> / <sub>2</sub>	<b>393FS</b>
Female Medium	7-8	7 <sup>1</sup> / <sub>2</sub> -9	<b>393FM</b>
Female Large	8 <sup>1</sup> / <sub>2</sub> -9 <sup>1</sup> / <sub>2</sub>	9 <sup>1</sup> / <sub>2</sub> -10 <sup>1</sup> / <sub>2</sub>	<b>393FL</b>
Male Small	8 <sup>1</sup> / <sub>2</sub> -9 <sup>1</sup> / <sub>2</sub>	9 <sup>1</sup> / <sub>2</sub> -10 <sup>1</sup> / <sub>2</sub>	<b>393MS</b>
Male Medium	10-11	11-12	<b>393MM</b>
Male Large	11 <sup>1</sup> / <sub>2</sub> -13 <sup>1</sup> / <sub>2</sub>	12 <sup>1</sup> / <sub>2</sub> +	<b>393ML</b>



## Flexible Sole Post-Op Shoe

- Nylon mesh upper with open toe design
- Padded heel and soft tricot liner for greater comfort
- Urethane rubber bottom with slip-resistant sole
- Hook and loop closures with movable hook/loop tabs provide unlimited strap adjustability

Size	US Shoe Size		Product No.	
	Men's	Women's	Men's	Women's
X-Small	—	3-5	<b>394XS</b>	
Female Small	—	5-7 <sup>1</sup> / <sub>2</sub>	<b>394FS</b>	
Female Medium	7-8	7 <sup>1</sup> / <sub>2</sub> -9	<b>394FM</b>	
Female Large	8 <sup>1</sup> / <sub>2</sub> -9 <sup>1</sup> / <sub>2</sub>	9 <sup>1</sup> / <sub>2</sub> -10 <sup>1</sup> / <sub>2</sub>	<b>394FL</b>	
Male Small	8 <sup>1</sup> / <sub>2</sub> -9 <sup>1</sup> / <sub>2</sub>	9 <sup>1</sup> / <sub>2</sub> -10 <sup>1</sup> / <sub>2</sub>	<b>394MS</b>	
Male Medium	10-11	11-12	<b>394MM</b>	
Male Large	11 <sup>1</sup> / <sub>2</sub> -13 <sup>1</sup> / <sub>2</sub>	12 <sup>1</sup> / <sub>2</sub> +	<b>394ML</b>	
X-Large	14-15	—	<b>394XL</b>	



## Ossur Americas - Mesh Top Post-Op Shoe

- Breathable mesh upper with reinforced heel counter and padded heel collar, provide more stability overall to the rear of the foot while reducing slippage
- Semi-rigid sole provides superior protection
- Extended strap and dual side closure accommodates extra bulky dressings

Size	US Shoe Size		Product No.	
	Men's	Women's	Men's	Women's
Small	4-6	5-6	<b>17013</b>	<b>17003</b>
Medium	6 <sup>1</sup> / <sub>2</sub> -8 <sup>1</sup> / <sub>2</sub>	7-8	<b>17015</b>	<b>17005</b>
Large	9-11	9+	<b>17017</b>	<b>17007</b>
X-Large	11 <sup>1</sup> / <sub>2</sub> -14	—	<b>17018</b>	—





**DH Offloading Post-Op Style Shoe**

- Sold in single shoe; fits either left or right foot
- Designed to accommodate bulky dressings
- Super soft nylon upper with reinforced heel counter provides extra comfort and stability
- Patented removable toe piece aids in wound accessibility for daily wound site inspection



Size	US Shoe Size		Shoe	Replacement Insoles
	Men's	Women's		
Small	4-6	6 <sup>1</sup> / <sub>2</sub> -8 <sup>1</sup> / <sub>2</sub>	<b>10340C</b>	<b>12365</b>
Medium	6 <sup>1</sup> / <sub>2</sub> -8 <sup>1</sup> / <sub>2</sub>	9-11	<b>10341C</b>	<b>12366</b>
Large	9-11	12+	<b>10342C</b>	<b>12367</b>
X-Large	11 <sup>1</sup> / <sub>2</sub> -14	13+	<b>10343C</b>	<b>12368</b>

**Soft Top Post-Op Shoe**

- Breathable cotton-poly blend upper provides added comfort
- Designed specifically for recovery after soft tissue procedures, post trauma or to relieve the geriatric or insensitive foot
- Contact closure system designed for ease of use and provides even torsion and easy entry and removal



Relief and  
Healing Shoes

Size	US Shoe Size		Product No.	
	Men's	Women's	Men's	Women's
Small	4-6	5-6	<b>18013</b>	<b>18003</b>
Medium	6 <sup>1</sup> / <sub>2</sub> -8 <sup>1</sup> / <sub>2</sub>	7-8	<b>18015</b>	<b>18005</b>
Large	9-11	9+	<b>18017</b>	<b>18007</b>
X-Large	11 <sup>1</sup> / <sub>2</sub> -14	—	<b>18018</b>	—

**Blue Nylon Soft Top Post-Op Shoe**

- Available with Lace or hook and loop closure



Lace Closure



Hook and Loop Closure

Size	US Shoe Size		Product No.	
	Men's	Women's	Men's	Women's
<b>Lace Closure</b>				
Pediatric	8 <sup>1</sup> / <sub>2</sub> "	—	<b>116PEDB</b>	—
Small	10 <sup>1</sup> / <sub>2</sub> "	9 <sup>1</sup> / <sub>2</sub> "	<b>116MSB</b>	<b>116WSB</b>
Medium	11 <sup>1</sup> / <sub>2</sub> "	10"	<b>116MMB</b>	<b>116WMB</b>
Large	12 <sup>1</sup> / <sub>2</sub> "	10 <sup>1</sup> / <sub>2</sub> "	<b>116MLB</b>	<b>116WLB</b>
<b>Hook and Loop Closure</b>				
Pediatric	8 <sup>1</sup> / <sub>2</sub> "	—	<b>117PEDB</b>	—
Small	10 <sup>1</sup> / <sub>2</sub> "	9 <sup>1</sup> / <sub>2</sub> "	<b>117MSB</b>	<b>117WSB</b>
Medium	11 <sup>1</sup> / <sub>2</sub> "	10"	<b>117MMB</b>	<b>117WMB</b>
Large	12 <sup>1</sup> / <sub>2</sub> "	10 <sup>1</sup> / <sub>2</sub> "	<b>117MLB</b>	<b>117WLB</b>

# RELIEF AND HEALING SHOES

## Rocker Bottom Cast Shoes

- Canvas upper, contact closure, open toe
- Cushioned rocker bottom cast sandal

Size	Insole Length	Product No.	
		Blue	Black
X-Small	7"	<b>308XS</b>	<b>B308XS</b>
Small	7 <sup>1</sup> / <sub>2</sub> "-9"	<b>308SM</b>	<b>B308SM</b>
Small/Medium	9"-9 <sup>3</sup> / <sub>4</sub> "	<b>308SM/MD</b>	<b>B308SM/MD</b>
Medium	10"-10 <sup>1</sup> / <sub>2</sub> "	<b>308MD</b>	<b>B308MD</b>
Medium/Large	10 <sup>1</sup> / <sub>2</sub> "-11"	<b>308MD/LG</b>	<b>B308MD/LG</b>
Large	11"-12"	<b>308LG</b>	<b>B308LG</b>
X-Large	12"-13 <sup>3</sup> / <sub>4</sub> "	<b>308XLG</b>	<b>B308XLG</b>



## Vinyl Rocker Bottom

- Vinyl Upper
- Rocker bottom

Size	Insole Length	Product No.	
		Open Toe	Closed Toe
X-Small	7"	<b>309XS</b>	<b>310XS</b>
Small	7 <sup>1</sup> / <sub>2</sub> -9"	<b>309SM</b>	<b>310SM</b>
Small/Medium	9-9 <sup>3</sup> / <sub>4</sub> "	<b>309SM/MD</b>	<b>310SM/MD</b>
Medium	10-10 <sup>1</sup> / <sub>2</sub> "	<b>309MD</b>	<b>310MD</b>
Medium/Large	10 <sup>1</sup> / <sub>2</sub> -11"	<b>309MD/LG</b>	<b>310MD/LG</b>
Large	11-12"	<b>309LG</b>	<b>310LG</b>
X-Large	12-13 <sup>3</sup> / <sub>4</sub> "	<b>309XLG</b>	<b>310XLG</b>



## ProCare (DJO) - Off-Loading Diabetic Shoe

- Designed to relieve pressure on ulcerated areas of the foot
- Impax Grid tri-laminate insole helps redistribute load away from ulcers
- Squared toe box and wide foot bed leave ample room for dressing
- Removable reinforced toe cover protects from external aggression
- Fits left or right

Size	US Shoe Size		Product No.
	Men's	Women's	
Small	4-6	6 <sup>1</sup> / <sub>2</sub> -8 <sup>1</sup> / <sub>2</sub>	<b>79-81513</b>
Medium	6 <sup>1</sup> / <sub>2</sub> -8 <sup>1</sup> / <sub>2</sub>	9-11	<b>79-81515</b>
Large	9-11	11 <sup>1</sup> / <sub>2</sub> -12 <sup>1</sup> / <sub>2</sub>	<b>79-81517</b>
X-Large	11 <sup>1</sup> / <sub>2</sub> -14	13+	<b>79-81518</b>



### Squared Toe Post-Op Shoe

- Allows dressings to be kept clean and toes dry following post-operative procedures and forefoot trauma
- Square toe design acts as a bumper, providing additional comfort and room for the patient
- Rigid E.V.A. rocker sole and contoured heel
- High ankle strapping that maintains foot in proper location
- Universal size fits left or right, male or female

Size	US Shoe Size		Product No.
	Men's	Women's	
X-Small	3-5	4-6	<a href="#">79-81232</a>
Small	5 <sup>1</sup> / <sub>2</sub> -7	6 <sup>1</sup> / <sub>2</sub> -8	<a href="#">79-81233</a>
Medium	7 <sup>1</sup> / <sub>2</sub> -9	8 <sup>1</sup> / <sub>2</sub> -10	<a href="#">79-81235</a>
Large	9 <sup>1</sup> / <sub>2</sub> -12	10 <sup>1</sup> / <sub>2</sub> -13	<a href="#">79-81237</a>
X-Large	12 <sup>1</sup> / <sub>2</sub> +	13 <sup>1</sup> / <sub>2</sub> +	<a href="#">79-81238</a>



### Deluxe Cast Boot Weatherproof

- Closed toe and heel protects from harsh weather
- Vinyl upper with rubber rocker sole to help cushion cast during ambulation
- Forefoot and heel contact closure for use with lower leg casted fractures
- Latex free

Size	Length x Toe Width	Product No.
X-Small	7 <sup>1</sup> / <sub>4</sub> " x 3 <sup>1</sup> / <sub>2</sub> "	<a href="#">79-81022</a>
Small	9 <sup>1</sup> / <sub>4</sub> " x 4 <sup>1</sup> / <sub>4</sub> "	<a href="#">79-81023</a>
Medium	11" x 4 <sup>3</sup> / <sub>4</sub> "	<a href="#">79-81025</a>
Large	12 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>4</sub> "	<a href="#">79-81027</a>
X-Large	13 <sup>1</sup> / <sub>2</sub> " x 6"	<a href="#">79-81028</a>



### Beach Boy Cast Boot

- Sandal style boot with rubber rocker sole
- Sturdy elastic straps with contact closure
- Open air design ideal for warmer climates
- For use with lower leg casted fractures
- Latex free

Size	Length x Toe Width	Product No.
X-Small	7 <sup>1</sup> / <sub>2</sub> " x 3 <sup>3</sup> / <sub>4</sub> "	<a href="#">79-81072</a>
Small	9" x 4 <sup>1</sup> / <sub>4</sub> "	<a href="#">79-81073</a>
Medium	10 <sup>1</sup> / <sub>2</sub> " x 4 <sup>1</sup> / <sub>2</sub> "	<a href="#">79-81075</a>
Large	12" x 5"	<a href="#">79-81077</a>



**Rocker Cast Boot**

- Sandal style with durable canvas upper
- Non-slip rocker sole to allow for more natural gait during ambulation
- Forefoot and heel hook and loop closure for ease of application
- Latex free
- Available in Black, pediatric sizes available in pediatric print

Size	Length	Toe Width	Product No.
3X-Small	4 <sup>1</sup> / <sub>2</sub> "	2 <sup>1</sup> / <sub>2</sub> "	79-81110
2X-Small	5 <sup>1</sup> / <sub>2</sub> "	3"	79-81111
X-Small	7 <sup>1</sup> / <sub>4</sub> "	3 <sup>1</sup> / <sub>2</sub> "	79-81112
Small	9 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>4</sub> "	79-81113
Small/Medium	10 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	79-81114
Medium	11"	4 <sup>3</sup> / <sub>4</sub> "	79-81115
Medium/Large	12"	5"	79-81116
Large	12 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>4</sub> "	79-81117
X-Large	13 <sup>1</sup> / <sub>2</sub> "	6"	79-81118



**Rocker Cast Shoe with Flex Sole**

- Durable canvas upper with flexible rocker sole
- Contact closure forefoot and heel straps for ease of application
- Latex free

**Indications:** Lower leg casted fractures

Size	Length	Toe Width	Product No.
X-Small	7 <sup>1</sup> / <sub>4</sub> "	3 <sup>1</sup> / <sub>2</sub> "	79-81222
Small	9 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>4</sub> "	79-81223
Medium	11"	4 <sup>3</sup> / <sub>4</sub> "	79-81225
Large	12 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>4</sub> "	79-81227
X-Large	13 <sup>1</sup> / <sub>2</sub> "	6"	79-81228



**Med/Surg Shoe**

- Soft foam mesh upper with padded tricot lining
- Semi-flexible skid resistant sole and comfortable contoured heel
- Large overlapping toe area accommodates for post-operative dressing and swelling
- Loop lock contact closure straps

Size	US Shoe Size		Product No.	
	Men's	Women's	Men's	Women's
X-Small	5-7	2-4	79-81132	79-81142
Small	7-9	4-6	79-81133	79-81143
Medium	9-11	6-8	79-81135	79-81145
Large	11-13	8-10	79-81137	79-81147
X-Large	13+	—	79-81138	—



**Post-Op Shoe**

- Low-profile foam/nylon mesh upper with tricot lining
- Rigid E.V.A. rocker sole and contoured heel for added patient comfort
- Protective tongue provides maximum comfort
- Available in loop-lock design

**Indications:** Post-foot surgery; Metatarsal injuries; Soft tissue procedures

Size	US Shoe Size		Product No.	
	Men's	Women's	Men's	Women's
X-Small	5-7	2-4	79-90182	79-90192
Small	7-9	4-6	79-90183	79-90193
Medium	9-11	6-8	79-90185	79-90195
Large	11-13	8-10	79-90187	79-90197
X-Large	13+	—	79-90188	—



**RCAI - Ambu-Assist**

- Encourages early ambulation and gait training accelerating the rehabilitative process
- Standard Ambu-Assist will not work with Phase II system

Size	Product No.	
	Standard	Phase II
Toddler/Child	12AA-P	20AA-A-P
Small/Youth	12AA-S	20AA-A-Y
Medium	12AA-M	20AA-A-A
Large/X-Large	12AA-L	20AA-A-L
2X-Large	12AA-XXL	—



Arch Supports  
and Inserts

**ARCH SUPPORTS AND INSERTS**

**Acor - Sole Defense Medical Active Daily Living**

- X-Static®-lined Impression Puff™ + P-Cell® and Microcel Puff®
- Weight Range: 100-300+ lbs
- Activity Level: mall and active walking

Size	US Shoe Size		Product No.
	Men's	Women's	
A	4	5-6	40000
B	5-6	7-8	40001
C	7-8	9-10	40002
D	9-10	11-12	40003
E	11-12	13-14	40004
F	13-14	15-16	40005



## QUIKfit® Lite

- Materials: P-Cell® + Microcel Puff®
- Weight Range: 100-225 lbs
- Activity Level: Minimal and community walking
- Replace ? with the appropriate die cut letter

Product No.

QF-07?



## QUIKfit Pro 9P6

- Biomedical multi layer
- P-Cell® + PORON Medical® Urethane + Microcel Puff®
- Weight Range: 100-300 lbs
- Activity Level: Shopping, Mall Walking and Active Walking
- Replace ? with Euro size from chart
- Whole sizes only

Product No.

9P6-L-?



## QUIKfit Pro 3P6

- Biomedical multi layer
- Plastazote® + PORON Medical® Urethane + Microcel Puff®
- Weight Range: 100-300+ lbs
- Activity Level: Shopping, Mall Walking and Active Walking
- Replace ? with Euro size from chart
- Whole sizes only

Product No.

3P6L?



QuikFit® Lite Sizing Chart - Oblique Toe					
Euro Size	U.S. Men	U.S. Women	Med. Width	Wide Width	X-Wide Width
37.5	3½	5½	A	A	B
38	4	6	A	B	B
38.5	4½	6½	B	C	C
39	5	7	B	C	D
39.5	5½	7½	C	D	D
40	6	8	C	D	E
40.5	6½	8½	D	E	F
41	7	9	E	E	G
41.5	7½	9½	E	F	G
42	8	10	F	G	H
42.5	8½	10½	G	H	I
43	9	11	H	I	I
43.5	9½	11½	I	I	J
44	10	12	I	J	K
44.5	10½		J	K	K
45	11		J	K	L
45.5	11½		K	L	L
46	12		K	L	L
47	13		L	M	M
48	14		M	M	M
Shoe Size			Die Cut Letter		

Fits the Broadway™ and other oblique last footwear

QuikFit® Lite Sizing Chart - Standard Toe					
Euro Size	U.S. Men	U.S. Women	Med. Width	Wide Width	X-Wide Width
37	3	5	N	O	O
37.5	3½	5½	N	O	P
38	4	6	O	P	P
38.5	4½	6½	P	P	Q
39	5	7	P	Q	Q
39.5	5½	7½	Q	Q	Q
40	6	8	Q	R	R
40.5	6½	8½	R	R	S
41	7	9	R	S	S
41.5	7½	9½	S	S	S
42	8	10	T	T	T
42.5	8½	10½	T	T	T
43	9	11	U	U	U
43.5	9½	11½	U	U	V
44	10	12	V	V	V
44.5	10½		V	W	W
45	11		W	W	W
45.5	11½		W	X	X
46	12		X	X	X
47	13		Y	Y	Y
48	14		Y	Z	Z
Shoe Size			Die Cut Letter		

Fits all ACOR comfort standard last shoes and most popular standard last footwear

**QUIKformables®**

- Custom foot orthoses fabricating blanks that cut practitioner’s time in half
- Heat, vacuum form, and grind
- Available in any material combination
- Sold in pairs



Size	US Shoe Size		Product No.
	Men’s	Women’s	
<b>Pink Plastazote + Blue PORON® Urethane + 35 dur Puff with 55 dur MultiCork™ Post</b>			
Small	4-6	6-8	<b>QFAB-0301-S</b>
Medium	7-9	9-11	<b>QFAB-0301-M</b>
Large	10-12	12-14	<b>QFAB-0301-L</b>
X-Large	13-15	15-17	<b>QFAB-0301-XL</b>
<b>Pink Plastazote + Blue PORON® Urethane + 35 dur Puff with 35 dur Puff Post</b>			
Small	4-6	6-8	<b>QFAB-0302-S</b>
Medium	7-9	9-11	<b>QFAB-0302-M</b>
Large	10-12	12-14	<b>QFAB-0302-L</b>
X-Large	13-15	15-17	<b>QFAB-0302-XL</b>
<b>35 dur Puff + Blue PORON® Urethane + 35 dur Puff with 55 dur MultiCork™ Post</b>			
Small	4-6	6-8	<b>QFAB-0401-S</b>
Medium	7-9	9-11	<b>QFAB-0401-M</b>
Large	10-12	12-14	<b>QFAB-0401-L</b>
X-Large	13-15	15-17	<b>QFAB-0401-XL</b>
<b>20 dur P-Cell® + Blue PORON® Urethane + 35 dur Puff with 35 dur MultiCork™ Post</b>			
Small	4-6	6-8	<b>QFAB-1502-S</b>
Medium	7-9	9-11	<b>QFAB-1502-M</b>
Large	10-12	12-14	<b>QFAB-1502-L</b>
X-Large	13-15	15-17	<b>QFAB-1502-XL</b>

**AliMed - FREEDOM<sup>®</sup> Basic Foot Orthosis BFO Half Sole**

- Offers mild cushioning and moderate control
- Neutral orthotic base readily accepts glued modifications
- Aliplast™ upper layer and supportive Plastazote™ layer
- Flexible and durable orthosis that can be used off-the-shelf or modified with a wedge or posting

Size	US Shoe Size			Product No.
	Child	Men's	Women's	
0 (Child's)	2-5 <sup>1</sup> / <sub>2</sub>	—	—	<b>6820</b>
1	—	—	5-6	<b>6821</b>
2	—	7-8	6-8	<b>6822</b>
3	—	8-9	8-9 <sup>1</sup> / <sub>2</sub>	<b>6823</b>
4	—	9-10 <sup>1</sup> / <sub>2</sub>	10	<b>6824</b>
5	—	10 <sup>1</sup> / <sub>2</sub> -11	11	<b>6825</b>
6	—	11-12	—	<b>6826</b>
7	—	12-13	—	<b>6827</b>
8	—	13-14	—	<b>6828</b>
9	—	14-15	—	<b>6829</b>
10	—	15-16	—	<b>68290</b>



**Poron<sup>®</sup> 4000 Insoles**

- Insole features a semirigid BFO half sole base providing cushioning and control
- Full length 1/8" top layer of Poron 4000 offers excellent shock absorption for plantar surface
- Blue Ultralure<sup>®</sup> fabric top affords comfort and protection against shear forces
- Accepts glued corrections to base for rear and forefoot correction

US Shoe Size		Product No.
Men's	Women's	
7-8	6-8	<b>6842</b>
8-9	8-9 <sup>1</sup> / <sub>2</sub>	<b>6843</b>
9-10 <sup>1</sup> / <sub>2</sub>	10	<b>6844</b>
10 <sup>1</sup> / <sub>2</sub> -11	11	<b>6845</b>



**Glass Composite Orthotics**

- Ultra thin (1.5 mm) and lightweight fits comfortably into shoe
- Glass-fiber reinforced thermoplastic composite orthotics can be adjusted for fit and functional control
- Accepts posting for added stability
- Retains rigidity and shape after heat modification

Description	Sizes	Product No.
Men's	7-8	<b>62463</b>
Men's	9-10	<b>62464</b>
Men's	11-12	<b>62465</b>
Women's	5-6	<b>62459</b>
Women's	7-8	<b>62461</b>
Women's	9-10	<b>62462</b>





**Ultimate Orthotic, Posted**

- Provides superior rigid arch support to help reduce pain associated with plantar fasciitis
- Rear foot posted (neutral)
- Soft Microdry® cover is durable and cushioning and helps reduce friction, heat and moisture

Size	US Shoe Size		Product No.
	Men's	Women's	
A	—	5-6 <sup>1</sup> / <sub>2</sub>	<b>63344</b>
B	—	7-8 <sup>1</sup> / <sub>2</sub>	<b>63345</b>
E	7-8 <sup>1</sup> / <sub>2</sub>	9-10 <sup>1</sup> / <sub>2</sub>	<b>63346</b>
F	9-10 <sup>1</sup> / <sub>2</sub>	—	<b>63347</b>
G	11-12 <sup>1</sup> / <sub>2</sub>	—	<b>63348</b>
H	13-14	—	<b>63349</b>



**Answer 2 Men's Shoe Inserts**

- Thermozote® cover, ShoreA 45 durometer density base
- Straight lasted, mild arch support
- Available sizes: 5-13, 14, 15
- 3 widths (M, WW, XW)
- Replace 1st ? with size and 2nd ? with width
- Available in Pink/Sky Blue

Product No. **309-??**



**Answer 2 Women's Shoe Inserts**

- Thermozote® cover, ShoreA 45 durometer density base
- Straight lasted, mild arch support
- Available sizes: 5-11
- 3 widths (M, WW, XW)
- Replace 1st ? with size and 2nd ? with width
- Available in Pink/Sky Blue

Product No. **310-??**



**Multi-Density Plastazote Insert**

- 1/8 inch thermozote and 1/4 EVA base
- Close-cell materials for better moisture control for diabetics
- Multi-density inserts help prevent ulceration, provides offload and shear reduction
- Available sizes: 5-12, 13-20 (full sizes only)
- Available widths: **Men's:** D, 2E, 4E, 6E, 9E; **Women's:** B, D, 3E, 5E, 7E
- Replace 1st ? with width and 2nd ? with size

Product No. **308-INSRT-??**



## Apex - Conform Insoles with Gel

- Comfortable, lightweight insoles that help absorb shock
- Gel technology for extra protection and shock absorption
- The Plastazote® foam top cover conforms to accommodate pressure “hot” spots
- Ideal for diabetic and arthritic patients
- Replace ? with size

Description	Sizes	Product No.
Men's	6-16	<a href="#">40665M?</a>
Women's	4-13	<a href="#">40665W?</a>



## Anatomical ThermoThotics™ (Polyetherdiol Base)

- Conformable 1/8" Plastazote top cover
- Base made from ED-9 Polyetherdiol urethane, a state-of-the-art heat moldable base compound that absorbs shock forces and cushions the plantar aspect of the foot
- Supports the arch and unloads stress at the forefoot while sitting naturally in the shoe
- Replace ? with size

Description	Sizes	Product No.
Men's	6-16	<a href="#">T40664M?</a>
Women's	4-13	<a href="#">T40664W?</a>



## Anatomical ThermoThotics (ThermoSKY® EVA Base - 45 Dur. Shore A)

- Multi-density orthotics feature a conformable Plastazote 1/8" top cover
- ThermoSKY® EVA 3/16" base is light-weight, durable, and heat moldable
- Designed to support the longitudinal arch and unload stress at the ball-of-the-foot, while sitting naturally in the shoe
- Replace ? with size

Description	Sizes	Product No.
Men's	6-16	<a href="#">40668M?</a>
Women's	4-13	<a href="#">40668W?</a>



## Contour ThermoThotics (Polyetherdiol Base)

- Features a conformable 1/8" Plastazote top cover
- Base made from ED-9 Polyetherdiol urethane that absorbs shock forces and cushions the plantar aspect of the foot
- Designed to provide support, stability and protection
- Replace ? with size

Description	Sizes	Product No.
Men's	6-16	<a href="#">T40777M?</a>
Women's	4-13	<a href="#">T40777W?</a>



**Contour ThermoThotics (ThermoSKY EVA Base)**

- Multi-density orthotics feature a conformable Plastazote 1/8" top cover
- ThermoSKY EVA 3/16" base is lightweight, durable and heat moldable
- Base designed to provide support, stability and protection
- Replace ? with size



Description	Sizes	Product No.
Men's	6-16	<b>1584M?</b>
Women's	4-13	<b>1584W?</b>

**Walker/Runner Orthotics**

- Walker/Runner with gel to absorb shock forces
- Patented heel cup
- Longitudinal arch support
- Soft metatarsal cushion to comfort the ball-of-foot




US Shoe Sizes		Product No.	
Men's	Women's	Men's	Women's
8/9	5/6	<b>40663M089</b>	<b>40663W056</b>
10/11	7/8	<b>40663M1011</b>	<b>40663W078</b>
12/13	9/10	<b>40663M1213</b>	<b>40663W910</b>
14/15	11/12	<b>40663M1415</b>	<b>40663W1112</b>

Arch Supports  
and Inserts

**Lynco Conform 200 Series**

- Recommended to protect and comfort the diabetic foot and the hyper sensitive arthritic foot
- Constructed with a Plastazote® top cover that customizes to the foot during the gait cycle and accommodates pressure "hot spots" for maximum protection
- Features a molded Pedic™ core and a Plastazote X-Firm bottom for resilience and shock absorption
- Replace ? with size




Description	Sizes	Product No.
Men's	6-16	<b>L200M?</b>
Women's	4-13	<b>L200W?</b>


## Lynco Sports 400 Series

- Designed specifically for comfort, support and weight redistribution in athletic, walking and comfort footwear.
- These tri-density orthotics feature a Pro-Shox® top cover for superior cushioning and shock absorption
- Pro-Shox's high-tech polyester fibers wick moisture away from the foot and provide a cool, friction free surface to prevent blisters and other foot problems
- Replace ? with size

Description	Sizes	Product No.
Men's	6-16	<b>L400M?</b>
Men's	6-16	<b>L420M?</b>
Men's	6-16	<b>L425M?</b>
Women's	4-13	<b>L400W?</b>
Women's	4-13	<b>L420W?</b>
Women's	4-13	<b>L425W?</b>



L400



L425

## Lynco Sports 400 Series with Metatarsal Pads

- Full length athletic orthotic featuring a Pedic™ rubber core and plastazote bottom
- Designed for medium and high arch and relief from related problems such as Metatarsalgia and Morton's Neuroma
- Absorbs shock and shear forces
- Replace ? with size

Description	Sizes	Product No.
Men's	6-16	<b>L405M?</b>
Women's	4-13	<b>L405W?</b>



## Lynco Dress 500 Series

- Tri-density orthotics feature a full grain leather top cover, a molded Pedic™ core for stability and resilience and a suede bottom to prevent sliding inside your shoes
- Cushions the foot, absorbs shock and fits comfortably in most footwear styles

Description	Sizes	Product No.
Men's	6-16	<b>L505M</b>
Women's	4-13	<b>L505W</b>



**Bauerfeind - ViscoPed® S Full Length Insoles**

- Full length insoles provide superior shock load reduction
- Silicone with softer blue areas for cushioning, latex free
- Ventilation holes help to dissipate heat, easily trimmed for precise fit, comfortable and long-lasting

**Indications:** Foot, joint and back pain; Forefoot and toe deformities; Diabetic foot; Pressure redistribution; Arthralgia; Metatarsalgia; Foot pad atrophy

Size	US Shoe Size		Product No.
	Men's	Women's	
0	3-4 <sup>1</sup> / <sub>2</sub>	4-5	13817320000000
1	5-6	5 <sup>1</sup> / <sub>2</sub> -7	13817320000001
2	6 <sup>1</sup> / <sub>2</sub> -7 <sup>1</sup> / <sub>2</sub>	7 <sup>1</sup> / <sub>2</sub> -9	13817320000002
3	8-9	9 <sup>1</sup> / <sub>2</sub> -10	13817320000003
4	9 <sup>1</sup> / <sub>2</sub> -10 <sup>1</sup> / <sub>2</sub>	10 <sup>1</sup> / <sub>2</sub> -11 <sup>1</sup> / <sub>2</sub>	13817320000004
5	11-12 <sup>1</sup> / <sub>2</sub>	12-13 <sup>1</sup> / <sub>2</sub>	13817320000005
6	13-14	14-15	13817320000006



**ViscoPed® Full Length Insoles**

- Full length, dual density, insoles
- Made of 100% medical grade silicone with softer blue areas for maximum cushioning, latex free
- Raised metatarsal pad aids in realignment of metatarsal heads, can be easily trimmed for precise fit

**Indications:** Metatarsalgia; Morton's Neuroma; Shock load reduction

Size	US Shoe Size		Product No.
	Men's	Women's	
1	5-6	6 <sup>1</sup> / <sub>2</sub> -7 <sup>1</sup> / <sub>2</sub>	13017300000001
2	6 <sup>1</sup> / <sub>2</sub> -7 <sup>1</sup> / <sub>2</sub>	8 <sup>1</sup> / <sub>2</sub> -9	13017300000002
3	8-8 <sup>1</sup> / <sub>2</sub>	9 <sup>1</sup> / <sub>2</sub> -10	13017300000003
4	9-10 <sup>1</sup> / <sub>2</sub>	10 <sup>1</sup> / <sub>2</sub> -12	13017300000004
5	11 <sup>1</sup> / <sub>2</sub> -12	13-14	13017300000005



Arch Supports  
and Inserts

**DonJoy (DJO) - Arch Rival™ Cavus Foot Orthotic**

- Treats Subtle Cavus Foot by promoting natural foot function
- Recessed area under 1st metatarsal head accommodates the plantar-flexed 1st ray
- X-Static anti-microbial material helps to prevent blisters and odor
- Sold by the pair

Size	US Shoe Size		Product No.
	Men's	Women's	
A	—	4 <sup>1</sup> / <sub>2</sub> -6	11-2015-1
B	5 <sup>1</sup> / <sub>2</sub> -7	6 <sup>1</sup> / <sub>2</sub> -8	11-2015-2
C	7 <sup>1</sup> / <sub>2</sub> -9	8 <sup>1</sup> / <sub>2</sub> -10	11-2015-3
D	9 <sup>1</sup> / <sub>2</sub> -11	10 <sup>1</sup> / <sub>2</sub> -12	11-2015-4
E	11 <sup>1</sup> / <sub>2</sub> -13	12 <sup>1</sup> / <sub>2</sub> -14	11-2015-5
F	13 <sup>1</sup> / <sub>2</sub> -15	14 <sup>1</sup> / <sub>2</sub> +	11-2015-6
G	15 <sup>1</sup> / <sub>2</sub> -17	—	11-2015-7



# ARCH SUPPORTS AND INSERTS

## Euro International - ReflexPlus Diabetic Insole

- Maximizes the foot protection and comfort of diabetic patients
- Made of tri-laminated materials
- The insole is thermoformable and heat-moldable

Size	US Shoe Size		Product No.
	Men's	Women's	
Small	6-7	8-9	4-EU-270/S
Medium	8-9	10-11	4-EU-270/M
Large	10-11	12-13	4-EU-270/L



## Fillauer - Nyloplex Arch Support Blanks

- Not pre-formed

Size	Thickness	Length		Product No.
		in	cm	
1	2.0 mm	5 <sup>1</sup> / <sub>2</sub>	14.0	073007
2	2.0 mm	5 <sup>3</sup> / <sub>4</sub>	14.7	073015
3	2.5 mm	6 <sup>1</sup> / <sub>16</sub>	15.4	073023
4	2.5 mm	6 <sup>3</sup> / <sub>8</sub>	16.2	073031
5	3.0 mm	6 <sup>9</sup> / <sub>16</sub>	16.9	073049
6	3.0 mm	6 <sup>15</sup> / <sub>16</sub>	17.6	073056
7	3.5 mm	7 <sup>1</sup> / <sub>4</sub>	18.4	073064
8	3.5 mm	7 <sup>1</sup> / <sub>2</sub>	19.1	073072
9	3.5 mm	7 <sup>13</sup> / <sub>16</sub>	19.8	073080
10	4.0 mm	8	20.5	073098
11	4.0 mm	8 <sup>3</sup> / <sub>8</sub>	21.3	073106
12	4.0 mm	8 <sup>5</sup> / <sub>8</sub>	22.0	073114



## Comfil® Thermo Formable Composite (TFC)

- Combines the strength of carbon composite with the molding properties of thermoplastics
- Flexible insole is 1.5 mm thick, best for shaping double curved surfaces (ex: heel)
- Rigid insole is 2.5 mm thick, best if used as a rolling insole for diabetic patients; Made for shaping single curves

Insoles (sold as pair)	Unidirectional Rigid	Bidirectional Flexible
35-37 mm	700030141	700030121
38-40 mm	700030142	700030122
41-43 mm	700030143	700030123
44-46 mm	700030144	700030124

Arch Support Blanks (sold each)	Bidirectional Rigid	Bidirectional Flexible
0-23 cm	700500	700540
23-26 mm	700520	700550
26-32 mm	700530	700560



**Pre-Fab Kits**

- TFC Solid AFO/Rigid Toe Kits

Size	Product No.
Small	700000
Medium	700005
Large	700010



**Hapad, Inc. - Comf-Orthotic 3/4 Length Insole**

- Supports the arch to limit strain on the plantar fascia while cushioning the resulting painful heel
- Recommended for many other common painful foot problems associated with flat feet, arthritis and the normal aging process
- Recommended for dress and casual style shoes without a removable insole

Shoe Size	Product No.
<b>Women's Comf-Orthotic Insoles</b>	
4 <sup>1</sup> / <sub>2</sub> -5 <sup>1</sup> / <sub>2</sub>	COWXS
6-7	COWS
7 <sup>1</sup> / <sub>2</sub> -8 <sup>1</sup> / <sub>2</sub>	COWM
9-10	COWL
10 <sup>1</sup> / <sub>2</sub> -11 <sup>1</sup> / <sub>2</sub>	COWXL
<b>Men's Comf-Orthotic Insoles</b>	
5-6 <sup>1</sup> / <sub>2</sub>	COWXSEW
7-8 <sup>1</sup> / <sub>2</sub>	COMS
9-9 <sup>1</sup> / <sub>2</sub>	COMM
10-11	COMMB
11 <sup>1</sup> / <sub>2</sub> -12 <sup>1</sup> / <sub>2</sub>	COML
13-14 <sup>1</sup> / <sub>2</sub>	COMXL
15-16 <sup>1</sup> / <sub>2</sub>	COMXXL
<b>X-Wide Comf-Orthotic Insoles</b>	
7-8 <sup>1</sup> / <sub>2</sub> EEE	COSEW
9-9 <sup>1</sup> / <sub>2</sub> EEE	COMEW
10-11 EEE	COMEWB
11 <sup>1</sup> / <sub>2</sub> -12 <sup>1</sup> / <sub>2</sub> EEE	COLEW
13-14 <sup>1</sup> / <sub>2</sub> EEE	COLEW



Arch Supports  
and Inserts

**Langer - PPT Longitudinal Arch Pads**

- 6 pair per bag

Size	Length x Width x Thickness	Product No.
Small	3 <sup>1</sup> / <sub>2</sub> " x 2" x 1 <sup>1</sup> / <sub>2</sub> "	22400
Medium	3 <sup>3</sup> / <sub>4</sub> " x 2 <sup>1</sup> / <sub>2</sub> " x 1 <sup>1</sup> / <sub>2</sub> "	22410
Large	4 <sup>1</sup> / <sub>4</sub> " x 2 <sup>1</sup> / <sub>2</sub> " x 1 <sup>1</sup> / <sub>2</sub> "	22420



## PPT Metatarsal Bars

- 6 pair per bag

Size	Product No.
Small	22500
Medium	22510
Large	22520



## PPT Metatarsal Pads

- 6 pair per bag

Size	Product No.
Small	22600
Medium	22610
Large	22620



## PediFix® - GelStep® Thin Insole with Soft Reliefs

- Provides extra shock absorption under the heel and forefoot, with soft metatarsal pad to minimize pressure on metatarsal heads
- Anti-microbial cover absorbs odor, prevents bacterial growth and allows foot to slide easily into the shoe

Size	US Shoe Size		Product No.
	Men's	Women's	
Small	4-5 <sup>1</sup> / <sub>2</sub>	5-7	7250SC
Medium	6-7 <sup>1</sup> / <sub>2</sub>	7 <sup>1</sup> / <sub>2</sub> -9 <sup>1</sup> / <sub>2</sub>	7251SC
Large	8-9	10-12	7252SC
X-Large	9 <sup>1</sup> / <sub>2</sub> -11	12+	7253SC



## TL-2100® Preform Orthotics

- Super thin orthotics fit comfortably in most footwear
- Deep heel cups, longitudinal arch, no metatarsal raise
- Genuine TL-2100® carbon graphite composite for rigid control
- Use as is or customize with top covers, rearfoot posts, etc.

Size	Product No.	
	Men's	Women's
4	—	4860-L-4
5	—	4860-L-5
6	4860-M-6	4860-L-6
7	4860-M-7	4860-L-7
8	4860-M-8	4860-L-8
9	4860-M-9	4860-L-9
10	4860-M-10	4860-L-10
11	4860-M-11	4860-L-11
12	4860-M-12	—
13	4860-M-13	—






**GelStep 3/4 Length Dual Durometer Deep Insole**

- High medial and lateral walls cradle soft rear-foot tissues
- “Blue zone” area of softer silicone provides spot relief under heel area for maximum comfort
- Package contains 1 pair

Size	US Shoe Size		Product No.	
	Men’s	Women’s	Covered	Uncovered
Small	4-5 <sup>1</sup> / <sub>2</sub>	5-7	<b>7000SC</b>	<b>7000S</b>
Medium	6-7 <sup>1</sup> / <sub>2</sub>	7 <sup>1</sup> / <sub>2</sub> -9 <sup>1</sup> / <sub>2</sub>	<b>7001SC</b>	<b>7001S</b>
Large	8-9	10-12	<b>7002SC</b>	<b>7002S</b>
X-Large	9 <sup>1</sup> / <sub>2</sub> -11	12+	<b>7003SC</b>	<b>7003S</b>
2X-Large	11 <sup>1</sup> / <sub>2</sub> -12 <sup>1</sup> / <sub>2</sub>	—	<b>7004SC</b>	<b>7004S</b>



**GelStep Replacement Insole with Soft Heel and Met Zones**

- Relieves heel pain and improves stability
- Optional anti-microbial cover absorbs odor, prevents bacterial growth and allows foot to slide easily into the shoe
- Package contains 1 pair

Size	US Shoe Size		Product No.	
	Men’s	Women’s	Covered	Uncovered
Small	4-5 <sup>1</sup> / <sub>2</sub>	5-7	<b>7050SC</b>	<b>7050S</b>
Medium	6-7 <sup>1</sup> / <sub>2</sub>	7 <sup>1</sup> / <sub>2</sub> -9 <sup>1</sup> / <sub>2</sub>	<b>7051SC</b>	<b>7051S</b>
Large	8-9	10 - 12	<b>7052SC</b>	<b>7052S</b>
X-Large	9 <sup>1</sup> / <sub>2</sub> -11	12+	<b>7053SC</b>	<b>7053S</b>
2X-Large	11 <sup>1</sup> / <sub>2</sub> -12 <sup>1</sup> / <sub>2</sub>	—	<b>7054SC</b>	<b>7054S</b>
3X-Large	13+	—	<b>7055SC</b>	<b>7055S</b>



Arch Supports  
and Inserts

**GelStep Thin Dress Insole with Soft Wide Metatarsal Pad**

- For absorbing shock, enhanced comfort and forefoot pressure relief while wearing high heeled and tight-fitting shoes
- Anti-microbial top cover

Size	US Shoe Size		Product No.
	Men’s	Women’s	
Small	4-5 <sup>1</sup> / <sub>2</sub>	5-7	<b>7029SC</b>
Medium	6-7 <sup>1</sup> / <sub>2</sub>	7 <sup>1</sup> / <sub>2</sub> -9 <sup>1</sup> / <sub>2</sub>	<b>7030SC</b>
Large	8-9	10-12	<b>7031SC</b>



## Pedors - 3P Inserts (Unisex)

- Easily heat moldable to ensure contact with the plantar surface of the foot and is extremely durable
- The bottom layer is  $\frac{3}{16}$ " and consists of a 40 durameter EVA
- The middle layer is  $\frac{1}{16}$ " Poron
- The top layer is  $\frac{1}{8}$ " Polyethylene foam

US Shoe Size		Product No.
Men's	Women's	
3	5	<b>3P-S05000</b>
4	6	<b>3P-S06000</b>
5	7	<b>3P-S07000</b>
6	8	<b>3P-S08000</b>
7	9	<b>3P-S09000</b>
8	19	<b>3P-S10000</b>
9	11	<b>3P-S11000</b>
10	12	<b>3P-S12000</b>
11	13	<b>3P-S13000</b>
12	14	<b>3P-S14000</b>
13	15	<b>3P-S15000</b>
14	16	<b>3P-S16000</b>



## ProCare (DJO) - Achilles Wedge

- 4" neoprene padded foam wedge
- Fits into most closed heel walkers

Product No. **79-95130**



## Silipos - WonderSport™

- Maximum support and shock absorption using medical grade silicone
- Cushions both heel arch and metatarsal area
- For metatarsalgia, fat pad atrophy, diabetic feet, elongated metatarsals, Morton's neuroma
- Clear with green dots
- Sold as a pair

US Shoe Size		Product No.
Men's	Women's	
6-7	—	<b>4508</b>
8-9	—	<b>4509</b>
10-11	—	<b>4510</b>
12-13	—	<b>4511</b>
14+	—	<b>4512</b>
—	4-5	<b>4513</b>
—	6-7	<b>4514</b>
—	8-9	<b>4515</b>
—	10+	<b>4516</b>



**SoleTech - Oversize OrthoSole**

- Rubber replacement sole

Size	US Shoe Size		Product No.
	Men's	Women's	
Small	6-8	Up to 10	<b>98081008</b>
Medium	9-11	10+	<b>98081010</b>
Large	11-12	—	<b>98081012</b>
X-Large	13-14	—	<b>98081014</b>



**Dual Density Diabetic Insoles**

- Designed with a top layer of 1/8" lightweight, heat moldable pink plastazote type foam and a firm base of 3/16" Shore "A" EVA to provide total contact with the patient's foot
- Heat moldable with an external heat source

Size	US Shoe Size		Product No.
	Men's	Women's	
X-Small	—	5-7	<b>E510101</b>
Small	6-7	7 1/2-9	<b>E510102</b>
Medium	7 1/2-8 1/2	9 1/2-12	<b>E510103</b>
Large	9-11	—	<b>E510104</b>
X-Large	11 1/2-13	—	<b>E510105</b>
2X-Large	14-15	—	<b>E510106</b>



Arch Supports  
and Inserts

**Spenco - PolySorb® Cross Trainer Insoles**

- Provides extra cushioning and support
- Contoured heel cup and arch support position feet to reduce the impact stress from hard surfaces
- Designed to return energy and help you push off for the next movement

Size	US Shoe Size		Product No.
	Men's	Women's	
1	—	5-6	<b>3803401</b>
2	6-7	7-8	<b>3803402</b>
3	8-9	9-10	<b>3803403</b>
4	10-11	11-12	<b>3803404</b>
5	12-13	—	<b>3803405</b>
6	14-15	—	<b>3803406</b>



## Slip-In Comfort Insoles

- Provides long lasting heel-to-toe comfort

Size	US Shoe Size		Product No.
	Men's	Women's	
0	—	3-4	<a href="#">40-212-00</a>
1	—	5-6	<a href="#">40-212-01</a>
2	6-7	7-8	<a href="#">40-212-02</a>
3	8-9	9-10	<a href="#">40-212-03</a>
4	10-11	11-12	<a href="#">40-212-04</a>
5	12-13	—	<a href="#">40-212-05</a>
6	14-15	—	<a href="#">40-212-06</a>



## Arch Supports

- Provides long lasting heel to toe comfort with an added soft arch cushion for extra cushion
- Available in full length or in a <sup>3</sup>/<sub>4</sub> length cushion extending from the heel to the base of the toes

Size	US Shoe Size		Product No.	
	Men's	Women's	<sup>3</sup> / <sub>4</sub> Length	Full Length
0	—	3-4	<a href="#">43-158-00</a>	<a href="#">43-042-00</a>
1	—	5-6	<a href="#">43-158-01</a>	<a href="#">43-042-01</a>
2	6-7	7-8	<a href="#">43-158-02</a>	<a href="#">43-042-02</a>
3	8-9	9-10	<a href="#">43-158-03</a>	<a href="#">43-042-03</a>
4	10-11	11-12	<a href="#">43-158-04</a>	<a href="#">43-042-04</a>
5	12-13	—	<a href="#">43-158-05</a>	<a href="#">43-042-05</a>
6	14-15	—	<a href="#">43-158-06</a>	<a href="#">43-042-06</a>



## Stable Step, Inc - Powerstep® ProTech Full-Length

- Polypropylene shell encased in a unique double-layered Poron®/EVA cushioning for shock
- Stability is further enhanced by “neutral” forefoot/rearfoot posts and a well-contoured heel cradle

Size	US Shoe Size		Product No.
	Men's	Women's	
AS	3-3 <sup>1</sup> / <sub>2</sub>	5-5 <sup>1</sup> / <sub>2</sub>	<a href="#">PA1S</a>
A	4-4 <sup>1</sup> / <sub>2</sub>	6-6 <sup>1</sup> / <sub>2</sub>	<a href="#">PA1</a>
B	5-5 <sup>1</sup> / <sub>2</sub>	7-7 <sup>1</sup> / <sub>2</sub>	<a href="#">PB1</a>
C	6-6 <sup>1</sup> / <sub>2</sub>	8-8 <sup>1</sup> / <sub>2</sub>	<a href="#">PC1</a>
D	7-7 <sup>1</sup> / <sub>2</sub>	9-9 <sup>1</sup> / <sub>2</sub>	<a href="#">PD1</a>
E	8-8 <sup>1</sup> / <sub>2</sub>	10-10 <sup>1</sup> / <sub>2</sub>	<a href="#">PE1</a>
F	9-9 <sup>1</sup> / <sub>2</sub>	11-11 <sup>1</sup> / <sub>2</sub>	<a href="#">PF1</a>
G	10-10 <sup>1</sup> / <sub>2</sub>	12-12 <sup>1</sup> / <sub>2</sub>	<a href="#">PG1</a>
H	11-11 <sup>1</sup> / <sub>2</sub>	13-13 <sup>1</sup> / <sub>2</sub>	<a href="#">PH1</a>
IJ	12-13 <sup>1</sup> / <sub>2</sub>	14-15 <sup>1</sup> / <sub>2</sub>	<a href="#">PIJ1</a>
K	14-14 <sup>1</sup> / <sub>2</sub>	16-16 <sup>1</sup> / <sub>2</sub>	<a href="#">PK1</a>



**Powerstep ProTech Classic Plus**

- Full-length double-layered orthotics made of resilient inner foam and a soft EVA outer foam for greater comfort
- Medical-grade polypropylene supports with heel cradle for flexibility and comfort

Size	US Shoe Size		Product No.
	Men's	Women's	
AS	3-3 <sup>1</sup> / <sub>2</sub>	5-5 <sup>1</sup> / <sub>2</sub>	<b>PCASM</b>
A	4-4 <sup>1</sup> / <sub>2</sub>	6-6 <sup>1</sup> / <sub>2</sub>	<b>PCA</b>
B	5-5 <sup>1</sup> / <sub>2</sub>	7-7 <sup>1</sup> / <sub>2</sub>	<b>PCB</b>
C	6-6 <sup>1</sup> / <sub>2</sub>	8-8 <sup>1</sup> / <sub>2</sub>	<b>PCC</b>
D	7-7 <sup>1</sup> / <sub>2</sub>	9-9 <sup>1</sup> / <sub>2</sub>	<b>PCD</b>
E	8-8 <sup>1</sup> / <sub>2</sub>	10-10 <sup>1</sup> / <sub>2</sub>	<b>PCE</b>
F	9-9 <sup>1</sup> / <sub>2</sub>	11-11 <sup>1</sup> / <sub>2</sub>	<b>PCF</b>
G	10-10 <sup>1</sup> / <sub>2</sub>	12-12 <sup>1</sup> / <sub>2</sub>	<b>PCG</b>
H	11-11 <sup>1</sup> / <sub>2</sub>	13-13 <sup>1</sup> / <sub>2</sub>	<b>PCH</b>
IJ	12-13 <sup>1</sup> / <sub>2</sub>	14-15 <sup>1</sup> / <sub>2</sub>	<b>PCIJ</b>
K	14-14 <sup>1</sup> / <sub>2</sub>	16-16 <sup>1</sup> / <sub>2</sub>	<b>PCK</b>



**Powerstep ProTech 3/4**

- Polypropylene shell encased in a unique double-layered Poron®/EVA cushioning for shock absorption and patient comfort
- Stabilizing Heel Cradle with opening filled with Poron cushion

Size	US Shoe Size		Product No.
	Men's	Women's	
AS	4-4 <sup>1</sup> / <sub>2</sub>	6-6 <sup>1</sup> / <sub>2</sub>	<b>P2ASM</b>
A	5-6 <sup>1</sup> / <sub>2</sub>	7-8 <sup>1</sup> / <sub>2</sub>	<b>P2A</b>
B	7-8 <sup>1</sup> / <sub>2</sub>	9-10 <sup>1</sup> / <sub>2</sub>	<b>P2B</b>
C	9-9 <sup>1</sup> / <sub>2</sub>	11-11 <sup>1</sup> / <sub>2</sub>	<b>P2C</b>
D	10-10 <sup>1</sup> / <sub>2</sub>	12-12 <sup>1</sup> / <sub>2</sub>	<b>P2D</b>
E	11-11 <sup>1</sup> / <sub>2</sub>	13-13 <sup>1</sup> / <sub>2</sub>	<b>P2E</b>
F	12-13	14-15	<b>P2F</b>



Arch Supports  
and Inserts

**Powerstep Slim-Tech® 3/4**

- Prefabricated orthotic provides versatility of 3/4 length support in slimmer, thinner profile
- Fits in most casual, slip-on and dress shoes

Size	US Shoe Size		Product No.
	Men's	Women's	
AS	4-4 <sup>1</sup> / <sub>2</sub>	6-6 <sup>1</sup> / <sub>2</sub>	<b>AS2STB</b>
A	5-6 <sup>1</sup> / <sub>2</sub>	7-8 <sup>1</sup> / <sub>2</sub>	<b>A2STB</b>
B	7-8 <sup>1</sup> / <sub>2</sub>	9-10 <sup>1</sup> / <sub>2</sub>	<b>B2STB</b>
C	9-9 <sup>1</sup> / <sub>2</sub>	11-11 <sup>1</sup> / <sub>2</sub>	<b>C2STB</b>
D	10-10 <sup>1</sup> / <sub>2</sub>	12-12 <sup>1</sup> / <sub>2</sub>	<b>D2STB</b>
E	11-11 <sup>1</sup> / <sub>2</sub>	13-13 <sup>1</sup> / <sub>2</sub>	<b>E2STB</b>
F	12-13	14-15	<b>F2STB</b>



## UCO International - QuadTech™ Soft Plus

- Four-layer, full-length insole with total contact top layer and 45° Shore A base layer
- Contoured inner layer and THK-1 reinforcement strengthens arch and resists compression
- Heat adjustable at 250° F
- Replace ? with metric sizes: 34-48

Product No.

21/0898-?



## Wrymark - Dynaflex Composite Plates

- Designed to control pathologies involved with the range-of- motion of the metatarsal phalangeal joints of the foot
- Resistive energy is stored within the material and released at the terminal stages of the gait, facilitating a more normal toe-off (“energy storage”)
- Replace xx with plate length and ? with L (left) or R (right)



Available in Contoured (A) or Flat (B) configurations in three alternative flexions

Description	Flex	Firm	Rigid
A. Contoured Foot plate	DFC Xxx?	DFCMxx?	DFCRxx?
B. Flat Foot plate	DFFXxx?	DFFMxx?	DFFRxx?

Dynaflex Sizing Chart		
US Men's	US Women's	Plate Length
—	5-6	21 <sup>1</sup> / <sub>2</sub> cm
6-7	7-8	23 <sup>1</sup> / <sub>2</sub> cm
8-9	9-10	25 <sup>1</sup> / <sub>2</sub> cm
10-11	11-12	27 <sup>1</sup> / <sub>2</sub> cm
12-13	—	29 <sup>1</sup> / <sub>2</sub> cm
14-15	—	31 <sup>1</sup> / <sub>2</sub> cm

**Ordering Example**

Model . . . . . **DFCX**

Plate Length . . . . . **21<sup>1</sup>/<sub>2</sub>**

Side . . . . . **R**

**DFCX215R**

**Apex - Seamfree Diabetic Socks**

- Seamless socks designed to assure a smooth fit without pressure
- Thin, strong yarn to assure proper room in your shoe
- Spandex for a soft stay-up design

Size	US Shoe Size		Product No.			
	Men's	Women's	Black	Brown	White	Navy
Small	3-7	4-8	S700S	S710S	S720S	S750S
Medium	8-11	9-12	S700M	S710M	S720M	S750M
Large	12+	13+	S700L	S710L	S720L	S750L



**BSN Medical - SensiFoot Diabetic Sock**

- Flat, soft, low-profile toe seam reduces pressure and irritation on toes
- Extra padding in the foot, heel and toe reduces friction and provides extra comfort and protection
- Acrylic multi-fiber yarns wick away moisture to keep feet comfortable and dry

Size	Sock Style	Product No.			
		White	Black	Navy	Brown
X-Small	Knee	110830	110865	110860	110855
Small	Knee	110831	110866	110861	110856
Medium	Knee	110832	110867	110862	110857
Large	Knee	110833	110868	110863	110858
X-Large	Knee	110834	110869	110864	110859
X-Small	Crew	110835	110850	110845	110840
Small	Crew	110836	110851	110846	110841
Medium	Crew	110837	110852	110847	110842
Large	Crew	110838	110853	110848	110843
X-Large	Crew	110839	110854	110849	110844
X-Small	Mini-Crew	110875	—	—	—
Small	Mini-Crew	110876	—	—	—
Medium	Mini-Crew	110877	—	—	—
Large	Mini-Crew	110878	—	—	—
X-Large	Mini-Crew	110879	—	—	—

Sizing Chart					
Gender	X-Small	Small	Medium	Large	X-Large
Men's	3 1/2"-5 1/2"	6"-8"	8 1/2"-10"	10 1/2"-12"	12 1/2"-14"
Women's	4 1/2"-6 1/2"	7"-9"	9 1/2"-11"	11 1/2"-13"	—



Knee Hi Sock



Crew Sock



Mini Crew Sock

Diabetic Socks

# DIABETIC SOCKS

## Comfort Products - Diabetic Socks

- Designed for those with diabetes or circulation problems
- Available in unbleached cotton as prescribed by many doctors
- Non-restrictive top extends to mid calf
- Fully cushioned throughout for easy walking
- Low profile seam in toe area reduces risk of irritation
- Reversible – can be worn with either the cushion or smooth side against foot
- Sold in units of 6 pair



Size/Seam Color	US Shoe Size		Product No.					
	Men's	Women's	Natural	Black	Brown	Gray	Navy	Beige
Small/Red Seam	3-8 <sup>1</sup> / <sub>2</sub>	4-9 <sup>1</sup> / <sub>2</sub>	DSW	DSWBL	DSWBR	DSWGR	DSWN	DSWBE
Medium/White Seam	9-12	10-13	DSM	DSMBL	DSMBR	DSMGR	DSMN	DSMBE
Large/Green Seam	12+	13+	DSG	DSGBL	DSGBR	DSGGR	DSGN	DSGBE

## Knit-Rite - SmartKnit® Seamless Diabetic Socks

- Maximum diabetic protection and everyday comfort
- Patented seamless knitting and form fitting core-spun yarns help eliminate wrinkling, bunching and seam pressure
- Made with moisture wicking fibers, including X-Static®, for a drier sock environment; antimicrobial to inhibit odor in the sock
- Sold as a pair



Size	Style	US Shoe Size		Coolmax White	X-Static Black	X-Static Gray	X-Static Navy
		Men's	Women's				
Small	Mini-Crew	5-7	5-8	72122	—	72622	—
Medium	Mini-Crew	8-10	9-11	72132	—	72632	—
Large	Mini-Crew	11-13	12-14	72142	—	72642	—
X-Large	Mini-Crew	14+	—	72152	—	72652	—
Small	Crew	5-7	5-8	71122	71722	71622	71422
Medium	Crew	8-10	9-11	71132	71732	71632	71432
Large	Crew	11-13	12-14	71142	71742	71642	71442
X-Large	Crew	14+	—	71152	71752	71652	71452
Small	Over-the-calf	5-7	5-8	73122	73722	—	73422
Medium	Over-the-calf	8-10	9-11	73132	73732	—	73432
Large	Over-the-calf	11-13	12-14	73142	73742	—	73442
X-Large	Over-the-calf	14+	—	73152	73752	—	73452
Small	Wide crew	5-7	5-8	71162	71762	—	—
Medium	Wide crew	8-10	9-11	71172	71772	—	—
Large	Wide crew	11-13	12-14	71182	71782	—	—
X-Large	Wide crew	14+	—	71192	71792	—	—



**Therasock® Care Sox Plus**

- Designed for those with serious edema or larger extremities, stretches up to 23 inches in circumference
- Hi-bulk acrylic moisture-wicking fibers help keep feet dry
- Also work well as a covering over bandages, dressings or casts
- 91% acrylic, 8% nylon; 1% spandex
- Sold as a pair

Size	US Shoe Size		White Crew	Product No.		
	Men's	Women's		Black Crew	Brown Crew	Navy Crew
Small	4-6	5-7	11122	11222	11322	11422
Medium	6-9	7-10	11132	11232	11332	11432
Large	9-12	10-13	11142	11242	11342	11442
X-Large	12-15	13+	11152	11252	11352	11452



**Therasock® Wide Sock System**

- Double layer construction reduces friction against the feet
- White inner liner protects skin from dyes and additives
- Ideal for people with wide feet
- Sold as a pair
- 72% acrylic, 28% nylon (outer layer)
- 70% Driwright, 30% nylon (inner layer)

Size	US Shoe Size (Wide)		White Crew	Black Crew
	Men's	Women's		
Medium	6-9	7-10	31132	31232
Large	9-12	10-13	31142	31242
X-Large	12-15	13+	31152	31252



Diabetic Socks

**Therasock® Comfort System Lite Diabetic Socks**

- Comfort and protection in a soft, lightly padded sock
- Exterior seam design protects sensitive toes from irritation; allows for completely seam-free feel on the inside
- Lightweight pads cushion the bottoms of the feet
- Latex free material; 76% acrylic, 19% nylon, 5% spandex

Size	US Shoe Size		White Mini Crew	White Crew	Black Crew	Brown Crew	Navy Crew
	Men's	Women's					
Small	4-6	5-7	42122	41122	41222	—	—
Medium	6-9	7-10	42132	41132	41232	41332	41432
Large	9-12	10-13	42142	41142	41242	41342	41442
X-Large	12-15	13+	42152	41152	41252	41352	41452



**Therасock® Comfort System Plus  
(Maximum Padding)**

- Extra padding and a seamless interior
- Moisture-wicking fibers help keep feet dry and reduces odor
- Seamless interior gives the sock a seamless feel and reduces the friction across the toes
- 86% acrylic, 10% nylon, 4% spandex (outer layer)
- 70% Driwright, 30% nylon (inner layer)
- Sold as a pair



Size	US Shoe Size		Product No.				Over the Calf	
	Men's	Women's	White Mini Crew	White Crew	Black Crew	White	Black	
Small	4-6	5-7	52122	51122	51222	—	—	
Medium	6-9	7-10	52132	51132	51232	53132	53232	
Large	9-12	10-13	52142	51142	51242	53142	53242	
X-Large	12-15	13+	52152	51152	51252	53152	53252	

**Therасock® Double Sock System**

- Made with a special “sock-within-a-sock” design that helps protect sensitive feet from colored dyes, additives, friction and shear
- Feet stay dry and healthy with moisture-wicking fibers
- Outer layer: 72% acrylic, 24% nylon, 4% lycra
- Inner Layer: 70% Driwright, 30% nylon
- Sold as a pair



Size	US Shoe Size		Product No.			
	Men's	Women's	White Mini Crew	White Crew	Black Crew	Navy Crew
Small	4-6	5-7	22122	21122	21222	—
Medium	6-9	7-10	22132	21132	21232	21432
Large	9-12	10-13	22142	21142	21242	21442
X-Large	12-15	13+	22152	21152	21252	21452

**Sleep Socks**

- Non-skid sleep sock
- Ultra-soft interior provides extra warmth and comfort for cold feet
- Moisture-wicking fibers help keep feet dry
- One size fits all
- Available in Light Blue



Product No. **91602**

**Silipos - Arthritic/Diabetic Gel Sock**

- Single layer sock with 1/16" of gel throughout the sole
- Can be used with or without diabetic shoes
- Provides optimum protection for the neuropathic foot
- Helps reduce friction, abrasion, and shear forces while preventing callus formation, helps prevent foot ulcerations

Sock Size	US Shoe Size		Product No.	
	Men's	Women's	White	Black
7-9	—	up to 6	1701	1706
9-11	up to 9	6+	1702	1707
10-13	9+	—	1703	1708



**Soft Sock**

- Double layer sock with 1/4" of gel throughout the sole
- Fluid like properties of the gel help massage and comfort the foot while walking
- Requires extra depth footwear
- Available in White; Packaged in a polybag

**Indications:** Diabetic; Rheumatoid; Arthritic and sensitive feet


Sock Size	US Shoe Size		Product No.
	Men's	Women's	
S (7-9)	—	up to 7	10605
M (9-11)	up to 9	7 1/2 - 9 1/2	10615
L (11-13)	9 1/2 - 12 1/2	—	10625
XL (13-16)	13+	—	10635



Diabetic Socks

## Gel Heel Cups

- Cupped design transfers pressure away from the rear foot
- Absorbs shock relieving pain from heel spurs and plantar fasciitis
- Made of medical grade silicone, sold per pair

Apex	<b>Gel Heel Cups - "Cupped"</b>				
	• Medical grade silicone gel heel cup				
	<b>Size</b>	<b>Men's Sizing</b>	<b>Women's Sizing</b>	<b>Product No.</b>	
	Small	—	4-8 <sup>1</sup> / <sub>2</sub>	<b>50770S</b>	
	Medium	7-10 <sup>1</sup> / <sub>2</sub>	9-11 <sup>1</sup> / <sub>2</sub>	<b>50770M</b>	
	Large	11+	12+	<b>50770L</b>	
	<b>Gel Heel Cradles - "Cupped"</b>				
• Anti-Shox Gel Heel Cradles					
<b>Size</b>	<b>Men's Sizing</b>	<b>Women's Sizing</b>	<b>Product No.</b>		
Small	—	4-8 <sup>1</sup> / <sub>2</sub>	<b>50750S</b>		
Medium	7-10 <sup>1</sup> / <sub>2</sub>	9-11 <sup>1</sup> / <sub>2</sub>	<b>50750M</b>		
Large	11+	12+	<b>50750L</b>		
ProCare (DJO)	<b>Silicone Heel Cups</b>				
	• Gel heel cup, soft center				
	<b>Size</b>	<b>Men's Sizing</b>	<b>Women's Sizing</b>	<b>Product No.</b>	
	Child/Youth	—	—	<b>79-81102</b>	
Small/Medium	5-9	5 <sup>1</sup> / <sub>2</sub> -9 <sup>1</sup> / <sub>2</sub>	<b>79-81103</b>		
Large/X-Large	9 <sup>1</sup> / <sub>2</sub> +	10+	<b>79-81107</b>		
Silipos	<b>WonderSpur</b>				
	• Clear with oval blue dot				
	<b>Size</b>	<b>Men's Sizing</b>	<b>Women's Sizing</b>	<b>Product No.</b>	
	Small	2-4	4-5	<b>4500</b>	
	Medium	5-7	6-9	<b>4501</b>	
	Large	8-10	10-12	<b>4502</b>	
	X-Large	10+	12+	<b>4503</b>	
	<b>WonderCup</b>				
	• Clear with round red dot				
	<b>Size</b>	<b>Men's Sizing</b>	<b>Women's Sizing</b>	<b>Product No.</b>	
Small	2-4	4-5	<b>4504</b>		
Medium	5-7	6-9	<b>4505</b>		
Large	8-10	10-12	<b>4506</b>		
X-Large	10+	12+	<b>4507</b>		

Recover Heel Cups					
<ul style="list-style-type: none"> <li>Dual durometer for maximum comfort while standing, walking or running</li> <li>Sold per pair</li> </ul>					
Pedifix	<b>GelStep Medium - Recovery Heel Cup, Soft Center Spot</b>				
	Size	Men's Sizing	Women's Sizing	Product No.	
	Small	—	4-7	<b>5025S</b>	
	Medium	6-8 <sup>1</sup> / <sub>2</sub>	7 <sup>1</sup> / <sub>2</sub> -10	<b>5026S</b>	
	Large	9-11 <sup>1</sup> / <sub>2</sub>	10 <sup>1</sup> / <sub>2</sub> -12	<b>5027S</b>	
	<b>GelStep Medium - Recovery Heel Cup, Soft Spur Spot</b>				
	Size	Men's Sizing	Women's Sizing	Product No.	
Small	—	4-7	<b>5050S</b>		
Medium	6-8 <sup>1</sup> / <sub>2</sub>	7 <sup>1</sup> / <sub>2</sub> -10	<b>5051S</b>		
Large	9-11 <sup>1</sup> / <sub>2</sub>	10 <sup>1</sup> / <sub>2</sub> -12	<b>5052S</b>		

Posted Heel Cups					
<ul style="list-style-type: none"> <li>Slight medial posting redirect pressure away from painful area</li> <li>Made of medical grade silicone, sold per pair</li> </ul>					
Apex	<b>Gel Heel Cradles - "Posted"</b>				
	Size	Men's Sizing	Women's Sizing	Product No.	
	Small	—	4-8 <sup>1</sup> / <sub>2</sub>	<b>50760S</b>	
	Medium	7-10 <sup>1</sup> / <sub>2</sub>	9-11 <sup>1</sup> / <sub>2</sub>	<b>50760M</b>	
Large	11+	12+	<b>50760L</b>		
Bauerfeind	<b>ViscoSpot Heel Cushions</b>				
	• Latex free				
	Size	Men's Sizing	Women's Sizing	Product No.	
1	3 <sup>1</sup> / <sub>2</sub> -7 <sup>1</sup> / <sub>2</sub>	5-9	<b>13017201000501</b>		
2	7 <sup>1</sup> / <sub>2</sub> -12 <sup>1</sup> / <sub>2</sub>	9-14	<b>13017201000502</b>		
Pedifix	<b>GelStep Posted Heel Pad - Soft Spur Spot</b>				
	Size	Men's Sizing	Women's Sizing	Product No.	
	Small	—	4-7	<b>5075S</b>	
	Medium	6-8 <sup>1</sup> / <sub>2</sub>	7 <sup>1</sup> / <sub>2</sub> -10	<b>5076S</b>	
Large	9-11 <sup>1</sup> / <sub>2</sub>	10 <sup>1</sup> / <sub>2</sub> -12	<b>5077S</b>		

## Tuli's

- Waffle design reduces impact and vibration of walking and running
- Sold per pair

<b>Medi-Dyne</b>	<b>Heavy Duty Gel Heel Cups</b>			
	<b>Size</b>	<b>Weight</b>	<b>Product No.</b>	
	Regular	under 175 lbs	<b>10217</b>	
	Large	over 175 lbs	<b>10218</b>	

## Leg Length Discrepancy Pads

- Fabric top cover, non-slip
- Sold each

<b>AliMed</b>	<b>Leg Length Discrepancy Pads</b>					
	• Choose 6 mm ( <sup>1</sup> / <sub>4</sub> " ) or 9 mm ( <sup>3</sup> / <sub>8</sub> " )					
	<b>Size</b>	<b>Men's Sizing</b>	<b>Women's Sizing</b>	<b>Product No.</b>		
				6 mm	9 mm	
	Small	5-6	6 <sup>1</sup> / <sub>2</sub> -8 <sup>1</sup> / <sub>2</sub>	<b>6990</b>	<b>6993</b>	
	Medium	6 <sup>1</sup> / <sub>2</sub> -8 <sup>1</sup> / <sub>2</sub>	9-12	<b>6991</b>	<b>6994</b>	
	Large	9+	—	<b>6992</b>	<b>6995</b>	
<b>Fit Feet</b>	<b>Adjust-A-Lift</b>					
	• <sup>3</sup> / <sub>8</sub> " total thickness (3 removable <sup>1</sup> / <sub>8</sub> " thick layers)					
	<b>Size</b>	<b>Men's Sizing</b>	<b>Women's Sizing</b>	<b>Product No.</b>		
	Small	—	up to 8	<b>354501</b>		
	Medium	up to 11	8 <sup>1</sup> / <sub>2</sub> +	<b>354502</b>		
	Large	11 <sup>1</sup> / <sub>2</sub> +	—	<b>354503</b>		

## M/L Lifts/Sports Wedge

- Firm, lightweight made of nickelplast
- Available in 2 degree or 4 degree strips
- 2 per pack

<b>AliMed</b>	<b>M/L Lifts/Sports Wedge</b>		
	• Measurements: 3 <sup>1</sup> / <sub>2</sub> " W x 12" L		
	<b>Size</b>	<b>Product No.</b>	
	2 Degree	<b>6996</b>	
	4 Degree	<b>6997</b>	


## Viscoelastic Varus/Valgus Wedges

- Firm viscoelastic polymer allows for shock absorption along with pronation/supination control

<b>AliMed</b>	<b>Viscoelastic Varus/Valgus Wedges</b>				
	• 6 mm				
	<b>Size</b>	<b>Men's Sizing</b>	<b>Women's Sizing</b>	<b>Product No.</b>	
	Medium	6 <sup>1</sup> / <sub>2</sub> -8 <sup>1</sup> / <sub>2</sub>	9-12	<b>66437/NA/MD</b>	
	Large	9-14	—	<b>66437/NA/LG</b>	

**Apex - Bunion Night Splint**

- Keeps the big toe aligned and comforts pain caused by bunions
- The soft bumper pad cushions the big toe, while the hook and loop strap acts as a lever to gradually align the big toe, relieving discomfort
- Should be worn at night
- Sold individually



Size	US Shoe Size		Product No.	
	Men's	Women's	Left	Right
Small	—	4-7 <sup>1</sup> / <sub>2</sub>	3001SL	3001SR
Medium	6-9 <sup>1</sup> / <sub>2</sub>	8 -10 <sup>1</sup> / <sub>2</sub>	3001ML	3001MR
Large	10+	11+	3001LL	3001LR

**Silipos - Achilles Heel Pad™**

- A soft conforming 1/8" gel pad protects the posterior of the heel from pressure and friction
- Provides gentle compression and dissipation of pressure around an inflamed area
- Reduces pain caused by footwear pressure particularly on any irregular heel bumps



Size	Stretchable Ankle Circum.	Product No.
Small/Medium	7"-10"	10385
Large/X-Large	9"-12"	10395

Foot and  
Digital Care

**Malleolar Sleeve™**

- Two large discs of polymer gel molded on both medial and lateral malleolar reducing pressure and friction to the area
- Reduces abrasion and shear forces
- Assists with prevention of ulceration to malleolar
- One size

Product No. **15005**



