

Wrist and Hand - Thumb Spica

- Thumb Spicas/Supports 370-377

Wrist and Hand - Wrist Supports

- Cock-Up Splints 378-385
- Wrist Supports 386-388

Wrist and Hand - Resting/Positioning

- Alimed® 389-391
- DeRoyal® 391-392
- LEEDer Group 393
- Lenjoy 393-396
- OCSI/Neuroflex 396-397
- Orthomerica 398
- ProCare (DJO) 398
- RCAI® 398-401
- Truform® 401

Wrist and Hand - Range-of-Motion

- Alimed 402
- Becker Orthopedic 402
- DeRoyal 403
- Lenjoy 403

Wrist and Hand - Finger/Thumb Orthoses

- Alimed 404-405
- DeRoyal 405-409
- Hely & Weber™ 409
- Medline® 409

Arm and Elbow - Range-of-Motion

- Aircast (DJO) 410
- Alimed 410
- Bledsoe Brace™ 410-411
- Breg® 411-412
- DeRoyal 412-413
- DonJoy (DJO) 413
- Lenjoy 414
- OCSI/Neuroflex 414
- Ossur Americas® 415
- ProCare (DJO) 415
- Professional Products 415
- RCAI 416-417

Arm and Elbow - Tennis Elbow Supports

- Tennis Elbow Splints 418-419

Arm and Elbow - Slings

- AliMed 420
- Corflex 421
- DeRoyal 421-423
- Hely & Weber 423
- ProCare (DJO) 424

Support/Splint/Padding/Immobilizer

- Alimed 425
- Bauerfeind 426

- Corflex 426
- DeRoyal 427
- DonJoy (DJO) 428
- Hely & Weber 428-429
- Lenjoy 429-430
- Medi USA 430
- New Options Sports 431
- ProCare (DJO) 431
- RCAI 432
- Silipos® 432

Shoulder - Abduction Type

- CoreLINE 433
- Becker Orthopedic 433
- Bledsoe Brace 433-434
- Breg® 434
- Corflex 435-436
- DeRoyal 436
- DonJoy (DJO) 437-438
- Fillauer 439
- Hely & Weber 439
- Ossur Americas 440
- Orthomerica 440-441
- ProCare (DJO) 441
- Professional Products 441
- Truform 442
- Trulife 442

Shoulder - Shoulder Supports

- Alimed 443
- Bauerfeind 443
- Becker Orthopedic 444
- Corflex 444
- DonJoy (DJO) 444
- New Options Sports 445
- Orthomerica 445

Shoulder - Clavicle Straps

- Clavicle Straps 446-448






Shoulder - Shoulder Immobilizers

- Alimed 448
- DeRoyal 448-450
- Hely & Weber 451
- Ossur Americas 452-453
- ProCare (DJO) 452-453
- Professional Products 453-454
- Trulife 454






Fracture Bracing

- Alimed 455
- Corflex 455
- DeRoyal 456
- DonJoy (DJO) 457
- Orthomerica 457-460
- ProCare (DJO) 461-462
- RCAI 463
- Trulife 463-464

Thumb Spicas/Supports

CoreLINE	CoreLINE Thumb Spica 8"				
	<ul style="list-style-type: none"> • Durable neoprene construction has dual aluminum stays to support the palmar surface of wrist and extensor surface of the thumb • Contact closure straps and Lycra® lining provide added patient comfort • Provides wrist immobilization while allowing full hand functions • Offers immobilization of the wrist/thumb MP/CMC joints 				
	Size		Product No.		
Universal		CL0242			
AliMed	Comfort Wrist Immobilizer				
	<ul style="list-style-type: none"> • Soft padded 8" immobilizer with volar metal stay • Thumb is held in abduction and splinted with a short metal stay • Self-padded with hook and loop closure 				
	Size	Wrist Circumference	Product No.		
			Left	Right	
	Small	5 ³ / ₄ "-6 ¹ / ₂ "	51-215/NA/LS	51-215/NA/RS	
	Medium	6 ⁵ / ₈ "-7 ¹ / ₂ "	51-215/NA/LM	51-215/NA/RM	
	Large	7 ⁵ / ₈ "-8 ¹ / ₂ "	51-215/NA/LL	51-215/NA/RL	
	FREEDOM® Arthritis Support				
	<ul style="list-style-type: none"> • Double-thick, soft, foam-padded flannel protects inflamed and painful MP joints • First tab aligns index and middle fingers in neutral, separate tab aligns fourth finger • Flexible stiffener along ulnar side maintains natural wrist support 				
	Size	Wrist Circumference	Product No.		
			Left	Right	
	Small	5 ¹ / ₂ "-6 ¹ / ₂ "	5933	5932	
	Medium	6 ¹ / ₂ "-7 ¹ / ₄ "	5935	5934	
	Large	7 ¹ / ₄ "-8 ¹ / ₄ "	5937	5936	
	FREEDOM CMC ThumbFit™				
<ul style="list-style-type: none"> • Elastic strap lifts and supports CMC joint, reduces slipping and grinding • Breathable AliDry™ wicks moisture away • Ideal for any activity or sport requiring grip or pinch strength • Measure circumference of hand over MP joints at widest point 					
Size	Circumference of MP Joints	Product No.			
		Left	Right		
Small	6"-6 ⁷ / ₈ "	51251/NA/LS	51251/NA/RS		
Medium	7"-7 ⁷ / ₈ "	51251/NA/LM	51251/NA/RM		
Large	8"-8 ⁷ / ₈ "	51251/NA/LL	51251/NA/RL		
X-Large	9"-9 ⁷ / ₈ "	51251/NA/LXL	51251/NA/RXL		
FREEDOM comfort™ Thumb Wrap					
<ul style="list-style-type: none"> • Offers compression and supports with premium elastic 2-way stretch material • Extra long, malleable, removable metal stay extends from IP joint to wrist • Adjusts for custom thumb positioning and hook and loop closure at wrist and thumb 					
Size		Product No.			
Universal		77559			

Thumb Spicas/Supports (continued)

AllMed	FREEDOM ThumbKeeper™ with D-Ring <ul style="list-style-type: none"> Rigid 1/16" multiform thin shell lined with Aliplast can be customized Clenching D-ring closure IP joint can flex and the MP joint is maintained in preferred resting position Measure width of hand at MP joints from radical aspect of index finger to ulnar aspect of little finger 																																	
			<table border="1"> <thead> <tr> <th rowspan="2">Size</th> <th rowspan="2">Wrist Circumference</th> <th colspan="2">Product No.</th> </tr> <tr> <th>Left</th> <th>Right</th> </tr> </thead> <tbody> <tr> <td>Small</td> <td>2 1/4"-2 3/4"</td> <td>5921</td> <td>5920</td> </tr> <tr> <td>Medium</td> <td>2 3/4"-3 1/4"</td> <td>5923</td> <td>5922</td> </tr> <tr> <td>Large</td> <td>3 1/4"-3 3/4"</td> <td>5925</td> <td>5924</td> </tr> </tbody> </table>			Size	Wrist Circumference	Product No.		Left	Right	Small	2 1/4"-2 3/4"	5921	5920	Medium	2 3/4"-3 1/4"	5923	5922	Large	3 1/4"-3 3/4"	5925	5924											
	Size	Wrist Circumference	Product No.																															
			Left	Right																														
	Small	2 1/4"-2 3/4"	5921	5920																														
	Medium	2 3/4"-3 1/4"	5923	5922																														
	Large	3 1/4"-3 3/4"	5925	5924																														
	FREEDOM Low Profile ThumbKeeper <ul style="list-style-type: none"> Thumb can be used as a stable post in grasp Overlap palmar strap closure is easy to manipulate Easily customized; trimmable Measure width of hand at MP joints from radical aspect of index finger to ulnar aspect of little finger 																																	
		<table border="1"> <thead> <tr> <th rowspan="2">Size</th> <th rowspan="2">Wrist Circumference</th> <th colspan="2">Product No.</th> </tr> <tr> <th>Left</th> <th>Right</th> </tr> </thead> <tbody> <tr> <td>Small</td> <td>2 1/4"-2 3/4"</td> <td>5103</td> <td>5102</td> </tr> <tr> <td>Medium</td> <td>2 3/4"-3 1/4"</td> <td>5105</td> <td>5104</td> </tr> <tr> <td>Large</td> <td>3 1/4"-3 3/4"</td> <td>5107</td> <td>5106</td> </tr> </tbody> </table>		Size		Wrist Circumference	Product No.		Left	Right	Small	2 1/4"-2 3/4"	5103	5102	Medium	2 3/4"-3 1/4"	5105	5104	Large	3 1/4"-3 3/4"	5107	5106												
Size	Wrist Circumference	Product No.																																
		Left	Right																															
Small	2 1/4"-2 3/4"	5103	5102																															
Medium	2 3/4"-3 1/4"	5105	5104																															
Large	3 1/4"-3 3/4"	5107	5106																															
FREEDOM Thumb Spica <ul style="list-style-type: none"> Adjust shell with scissors to ensure that the distal edge of the thumb post doesn't interfere with the IP joint motion Non-rigid shells strap placement effectively restricts motion and offers supports 				 																														
		<table border="1"> <thead> <tr> <th rowspan="2">Size</th> <th rowspan="2">Circumference of MP Joints</th> <th colspan="2">Product No.</th> </tr> <tr> <th>Left</th> <th>Right</th> </tr> </thead> <tbody> <tr> <td colspan="4" style="text-align: center;">Beige</td> </tr> <tr> <td>Small/Medium</td> <td>up to 7 3/4"</td> <td>5783</td> <td>5782</td> </tr> <tr> <td>Medium/Large</td> <td>7 3/4"+</td> <td>5785</td> <td>5784</td> </tr> <tr> <td colspan="4" style="text-align: center;">Black</td> </tr> <tr> <td>Small/Medium</td> <td>up to 7 3/4"</td> <td>5483</td> <td>5482</td> </tr> <tr> <td>Medium/Large</td> <td>7 3/4"+</td> <td>5485</td> <td>5484</td> </tr> </tbody> </table>			Size	Circumference of MP Joints	Product No.		Left	Right	Beige				Small/Medium	up to 7 3/4"	5783	5782	Medium/Large	7 3/4"+	5785	5784	Black				Small/Medium	up to 7 3/4"	5483	5482	Medium/Large	7 3/4"+	5485	5484
Size	Circumference of MP Joints	Product No.																																
		Left	Right																															
Beige																																		
Small/Medium	up to 7 3/4"	5783	5782																															
Medium/Large	7 3/4"+	5785	5784																															
Black																																		
Small/Medium	up to 7 3/4"	5483	5482																															
Medium/Large	7 3/4"+	5485	5484																															
Corflex	Ultra Fit Wrist Splint with Adducted Thumb <ul style="list-style-type: none"> Features easy one handed application Provides comfortable thumb and wrist immobilization while preserving mobility of the other fingers Manufactured from durable lightweight material with removable malleable aluminum palmar and thumb stays, flexible dorsal stay and contact closure 																																	
			<table border="1"> <thead> <tr> <th rowspan="2">Size</th> <th rowspan="2">Wrist Circumference</th> <th colspan="2">Product No.</th> </tr> <tr> <th>Left</th> <th>Right</th> </tr> </thead> <tbody> <tr> <td>X-Small</td> <td>4"-5"</td> <td>73-5312</td> <td>73-5302</td> </tr> <tr> <td>Small</td> <td>5"-6"</td> <td>73-5313</td> <td>73-5303</td> </tr> <tr> <td>Medium</td> <td>6"-7"</td> <td>73-5314</td> <td>73-5304</td> </tr> <tr> <td>Large</td> <td>7"-8"</td> <td>73-5315</td> <td>73-5305</td> </tr> <tr> <td>X-Large</td> <td>8"-9"</td> <td>73-5316</td> <td>73-5306</td> </tr> </tbody> </table>			Size	Wrist Circumference	Product No.		Left	Right	X-Small	4"-5"	73-5312	73-5302	Small	5"-6"	73-5313	73-5303	Medium	6"-7"	73-5314	73-5304	Large	7"-8"	73-5315	73-5305	X-Large	8"-9"	73-5316	73-5306			
	Size	Wrist Circumference	Product No.																															
			Left	Right																														
	X-Small	4"-5"	73-5312	73-5302																														
	Small	5"-6"	73-5313	73-5303																														
Medium	6"-7"	73-5314	73-5304																															
Large	7"-8"	73-5315	73-5305																															
X-Large	8"-9"	73-5316	73-5306																															

Thumb Spicas/Supports (continued)

Comfort Cool Thumb Wrap

- Thin neoprene that conforms to hand contours
- Features tiny perforations for cool, comfortable wear
- Soft terry cloth lining

Size	Circumference of the MP Joints	Product No.	
		Left	Right
Small	6"-7"	348SL	348SR
Medium	7"-8"	348ML	348MR
Large	8"-9"	348LL	348LR



Premium Wrist and Thumb Support - 8"

- Breathable material with four-way stretch and D-ring closures
- Does not extend past the distal palmar crease allowing full finger range-of-motion

Size	Product No.	
	Left	Right
X-Small	350XSL	350XSR
Small	350SL	350SR
Medium	350ML	350MR
Large	350LL	350LR
X-Large	350XLL	350XLR
2X-Large	350XLL	350XLR



Thermo-Form Thumb Splint

- Heavy neoprene-type material and features a thermoplastic piece that can be easily molded directly on the patient using warm water (160°-170° F) for a custom fit
- Inside pocket holds molded thermoplastic piece that is removable when cool

Size	Circumference of the MP Joints	Product No.	
		Left	Right
Thermo-Form Thumb Splint - Short Length			
X-Small	up to 5 ¹ / ₂ "	359XSL	359XSR
Small	5 ¹ / ₂ "-6 ¹ / ₂ "	359SL	359SR
Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	359ML	359MR
Large	7 ¹ / ₂ "-8 ¹ / ₂ "	359LL	359LR
X-Large	8 ¹ / ₂ "-9 ¹ / ₂ "	359XLL	359XLR
Thermo-Form Thumb Splint - Mid Length			
X-Small	up to 5 ¹ / ₂ "	360XSL	360XSR
Small	5 ¹ / ₂ "-6 ¹ / ₂ "	360SL	360SR
Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	360ML	360MR
Large	7 ¹ / ₂ "-8 ¹ / ₂ "	360LL	360LR
X-Large	8 ¹ / ₂ "-9 ¹ / ₂ "	360XLL	360XLR
Thermo-Form Thumb Splint - Long Length			
X-Small	up to 5 ¹ / ₂ "	361XSL	361XSR
Small	5 ¹ / ₂ "-6 ¹ / ₂ "	361SL	361SR
Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	361ML	361MR
Large	7 ¹ / ₂ "-8 ¹ / ₂ "	361LL	361LR
X-Large	8 ¹ / ₂ "-9 ¹ / ₂ "	361XLL	361XLR








DeRoyal




Thumb Spicas/Supports (continued)

DeRoyal	Thumb Spica				
	<ul style="list-style-type: none"> • Lightweight splint constructed of perforated foam and a soft inner liner • May be trimmed for a custom fit 				
	Size	Circumference of the MP Joints	Product No.		
			Left	Right	
Small/Medium	2 ¹ / ₂ "-3 ¹ / ₂ "	340S/ML	340S/MR		
Medium/Large	3 ¹ / ₂ "-4 ¹ / ₄ "	340M/LL	340M/LR		
Hely & Weber	Modabber Thumb Orthosis				
	<ul style="list-style-type: none"> • Kuhl™ perforated neoprene for cool, comfortable compression • Aluminum stay and steel spiral stay are interchangeable and removable 				
	Size	Thumb Hole Size	Wrist Circum.	Product No.	
	X-Small	2 ¹ / ₈ "-3 ³ / ₈ "	4 ¹ / ₂ "-5 ¹ / ₂ "	5803XS	
	Standard	2 ¹ / ₂ "-3"	5"-9 ¹ / ₂ "	5803HEL	
	X-Large	3 ¹ / ₈ "-3 ⁵ / ₈ "	5"-9 ¹ / ₂ "	5803XL	
	Neoprene Thumb Orthosis				
	<ul style="list-style-type: none"> • Elasticized straps for comfortable compression • Supports the MP and CMC joints while allowing full finger motion 				
	Size		Product No.		
	Universal		3813		
	Phomfit™ Wrist Hand Thumb Orthosis				
	<ul style="list-style-type: none"> • Aluminum removable palmar stay for customer fit and adjustment • Stockinette tongue and dorsal hook and loop closure for simple application • Radial thumb stay is malleable 				
	Size	Wrist Circumference	Product No.		
			Left	Right	
	X-Small	4 ¹ / ₂ "-5 ¹ / ₂ "	447LXS	447RXS	
Small	5"-6 ¹ / ₂ "	447LS	447RS		
Medium	6"-7 ¹ / ₂ "	447LM	447RM		
Large	7"-8 ¹ / ₂ "	447LL	447RL		
X-Large	8"-9 ¹ / ₂ "	447LXL	447RXL		
Santa Barbara Thumb Splint					
<ul style="list-style-type: none"> • Kuhl perforated neoprene for cool, comfortable compression • Radial, palmar and dorsal packets contain malleable aluminum and spiral stays 					
Size	Thumb Hole Size	Wrist Circumference	Product No.		
Universal	2 ¹ / ₂ "-3 ³ / ₈ "	4 ¹ / ₂ "-5 ¹ / ₂ "	3842		
TITAN™ Thumb Orthosis					
<ul style="list-style-type: none"> • Soft orthopedic felt laminated to brush nylon exterior • Oblique thumb stay is malleable to allow for desired thumb position and rigid enough for excellent thumb stability • Provides abduction control, angle adjustability • Stockinette tongue and dorsal loop & lock closures 					
Size		Product No.			
		Left	Right		
Regular		455L	455R		
Long		459L	459R		






Thumb Spicas/Supports (continued)

Hely & Weber	Trimnable Thumb Orthosis				
	<ul style="list-style-type: none"> • 1/8" nylon covered neoprene provides heat and compression • Melco taping in thumb section for trimming • Wrap around design with built in assist strap for easy application 				
	Size	Wrist Circumference	Product No.		
	Universal	4 1/2"-5 1/2"	3810		
Med Spec	UNO WHT® Wrist Hand Thumb™				
	<ul style="list-style-type: none"> • Polyester outer shell with tricot covered foam lining • Stockinette tongue and dorsal loop & lock closures • Malleable oblique thumb stays allow desired thumb position and stability 				
	Size	Product No.			
	Universal	Left	Right		
		429L	429R		
Med Spec	Suede Thumb Support				
	<ul style="list-style-type: none"> • Malleable aluminum thumb splint may be adjusted for a custom fit 				
	Size	Length	Product No.		
	Short	7 1/2"	223180		
	Long	9"	223190		
	Tee Pee™ Thumb Support				
	<ul style="list-style-type: none"> • Three aluminum stays form the "Tee Pee" that provides a custom fit and safe and stable haven for the thumb • Each stay can be shaped and contoured with ease • Does not restrict wrist motion • Measure the linear distance from the dorsal crease of the thumb IP joint to the peak of the radial styloid process 				
	Size	Measurements	Product No.		
	X-Small	up to 3 1/2"	223081		
	Small	3 1/2"-4 1/8"	223082		
Medium	4 1/8"-4 5/8"	223084			
Large	4 5/8"-5 1/8"	223085			
X-Large	5 1/8"	223086			
Universal Ryno Lacer					
<ul style="list-style-type: none"> • Dorsal, palmar and radial aluminum stays for immobilization of the thumb while restricting flexion and extension on the wrist • Dorsal stay casing may be repositioned to fit various size patients 					
Size	Wrist Circumference	Product No.			
		Left	Right		
Universal Long	5 1/2"-8 1/2"	223520	223521		
Universal Short	5 1/2"-8 1/2"	223510	223511		





Thumb Spicas/Supports (continued)

North Coast Medical	Comfort Cool Thumb CMC Restriction Splint				
	<ul style="list-style-type: none"> • Wraparound supports are made of thin 1/16" perforated neoprene with a terry cloth liner • Reinforced stitching allows splint to be trimmed below the thumb MP joint • Contoured neoprene strap wraps around the thumb CMC joint and first metacarpal for direct CMC joint support 				
	Size	Circumference of the MP Joint	Product No.		
			Left	Right	
	Toddler	4 1/2"-5 1/2"	NC79546	NC79547	
	Youth	5 1/2"-6 1/2"	NC79548	NC79549	
	Small	6"-7"	NC79562	NC79563	
	Medium	7"-8"	NC79564	NC79565	
	Large	8"-9"	NC79566	NC79567	
	X-Large	9"-10"	NC79550	NC79551	
Small Plus	6 1/2"-7 1/4"	NC79556	NC79557		
Medium Plus	7 7/8"-8 1/4"	NC79558	NC79559		
Large Plus	8 1/4"-9 1/2"	NC79560	NC79561		
Orthomerica	European Thumb Spica				
	<ul style="list-style-type: none"> • Effective anatomical, flesh-colored splint for maintaining thumb in full extension while allowing finger and wrist functions following sprains or surgery • Includes Cotton/Lycra liner • Adjustable hook and loop wrist strap 				
	Size	Wrist Circum.	Brace Length	Product No.	
				Left Right	
	Small	under 7"	4"	1200.01 1205.01	
	Medium	7"-8"	5"	1201.01 1206.01	
Large	8"-9"	5"	1202.01 1207.01		
X-Large	9"+	6"	1203.01 1208.01		
Ossur Americas	FORMFIT® Thumb Spica				 <i>8" Thumb Spica</i>
	<ul style="list-style-type: none"> • 3-dimensional molding technology allows padding thickness to vary strategic locations for superior comfort and fit • Effective long-term immobilization and protection • Adjustable radial and palmar stays 				
	Size	Wrist Circumference	Product No.		
			Left	Right	
	FORMFIT Thumb Spica - 6"				
	Small	6 1/4"-7 1/4"	30304	31304	
	Medium	7 1/4"-8"	30504	31504	
	Large	8"-9"	30704	31704	
	X-Large	9"-10"	30804	31804	
	FORMFIT Thumb Spica - 8"				
	X-Small	5 1/4"-6 1/4"	3020C	3120C	
	Small	6 1/4"-7 1/4"	3030C	3130C	
	Medium	7 1/4"-8"	3050C	3150C	
	Large	8"-9"	3070C	3170C	
X-Large	9"-10"	3080C	3180C		

Thumb Spicas/Supports (continued)

Ossur Americas	Universal Thumb Spica				
	Size	Product No.			
		Left	Right		
	Universal	W-501000	W-501001		
ProCare (DJO)	Abducted ThumbSPICA				
	Size	Product No.			
			Left	Right	
	Universal	79-87101	79-87100		
	Quick-Fit® Wrist Thumb Orthosis (W.T.O.)				
	Size	Wrist Circumference	Product No.		
			Left	Right	
	Universal	up to 13"	79-87490	79-87480	
	ComfortFORM™ Wrist with Abducted Thumb				
	Size	Wrist Circumference	Product No.		
		Left	Right		
X-Small	up to 5 ¹ / ₂ "	79-87312	79-87302		
Small	5 ¹ / ₂ "-6 ¹ / ₂ "	79-87313	79-87303		
Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	79-87315	79-87305		
Large	7 ¹ / ₂ "-8 ¹ / ₂ "	79-87317	79-87307		
X-Large	8 ¹ / ₂ "+	79-87318	79-87308		
ThumbGUARD™					
Size	Wrist Circumference	Product No.			
		Left	Right		
X-Small	1 ³ / ₄ "-2 ¹ / ₄ "	79-87272	79-87262		
Small	2 ¹ / ₄ "-2 ³ / ₄ "	79-87273	79-87263		
Medium	2 ³ / ₄ "-3 ¹ / ₄ "	79-87275	79-87265		
Large	3 ¹ / ₄ "-3 ³ / ₄ "	79-87277	79-87267		
X-Large	3 ³ / ₄ "-4 ¹ / ₄ "	79-87278	79-87268		

Thumb Spicas/Supports (continued)

ProCare (DJO)	<p>ThumbSPICA™</p> <ul style="list-style-type: none"> • Lightweight perforated outer foam shell with a cool, cotton-terry liner • Elastic, circumferential, contact closure straps for proper positioning of splint • Small/Medium fits most women and smaller men; Large/X-Large fits most men and larger women 				
	Size	Length	Product No.		
			Left	Right	
	Small/Medium	9"	79-87114	79-87113	
Large/X-Large	9"	79-87118	79-87117		
Professional Products	<p>Universal Thumb-O-Prene</p> <ul style="list-style-type: none"> • Wraparound neoprene design provides a low-profile custom fit • Combines two stays, rigid and spiral plus a moveable elastic strap to address various pathologies of the thumb • May be wrapped between forefinger and thumb to help stabilize the CMC joint without restricting the use of the thumb or wrist 				
	Size		Product No.		
	Universal		79-82700		
Professional Products	<p>Wahoo™ II Wrist Support and Thumb Spica</p> <ul style="list-style-type: none"> • Abducted thumb and wrist support has three thumb stays that are designed to abduct the thumb and add virtually 360° support to the thumb joint • Two of the stays are actually removable so you can mold them to your thumb or remove them altogether for less rigidity in support • Fastens with a hook and loop closure 				
	Size		Product No.		
	Universal		11359-00		
RCAI	<p>Universal Thumb/Wrist Thumb Support</p> <ul style="list-style-type: none"> • Gently places the thumb in a position of function and limits extreme motion of the joints of the thumb • A bendable metal stay limits radial and palmar abduction of the CMC joint and flexion/extension of the MCP joint • Breathable airprene material 				 <p style="text-align: center;"><i>Wrist</i></p> <p style="text-align: center;"><i>Thumb Wrist</i></p>
	Size	Measure IP Crease to Wrist	Product No.		
			Thumb	Wrist Thumb	
	Youth	2 ¹ / ₂ "-3"	510-UTS-Y	520-UWTS-Y	
	X-Small	3"-3 ¹ / ₂ "	510-UTS-ES	520-UWTS-ES	
	Small	3 ¹ / ₂ "-4 ¹ / ₄ "	510-UTS-S	520-UWTS-S	
	Medium	4"-4 ¹ / ₂ "	510-UTS-M	520-UWTS-M	
Large	4 ¹ / ₂ "+	510-UTS-L	520-UWTS-L		

Cock-Up Splints

- Malleable palmar and dorsal stays provide anatomical fit and proper support
- Indications: Mild sprains and strains of the wrist; Tendonitis; Tendosynovitis

CoreLINE

Cock-Up Wrist Splint

- Provides full finger mobility with contoured aluminum palmar stay to maintain wrist in cock-up position
- Cool, soft nylon construction with single one-pull design fits wide range of wrist sizes
- Universal sizing for left or right hand

Size	Product No.
Universal	CL0145



Cock-Up Wrist Splint

- Non-removable aluminum dorsal stay
- Palmar stay can be removed and adjusted/formed for improved fit
- 8" long; Fits wrists up to 9.5" max circumference

Size	Wrist Circumference	Product No.	
		Left	Right
X-Small	up to 5"	CL0270	CL0280
Small	5"-6 ¹ / ₂ "	CL0272	CL0282
Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	CL0274	CL0284
Large	7 ¹ / ₂ "-8 ¹ / ₂ "	CL0276	CL0286
X-Large	8 ¹ / ₂ "+	CL0278	CL0288



Corflex

Ultra Fit Wrist Splint 8"

- Features easy one-hand application that provides comfortable support of the wrist while preserving finger mobility

Size	Wrist Circumference	Product No.	
		Left	Right
Universal	5"-9"	73-1001	73-1000
X-Small	4"-5"	73-1012	73-1002
Small	5"-6"	73-1013	73-1003
Medium	6"-7"	73-1014	73-1004
Large	7"-8"	73-1015	73-1005
X-Large	8"-9"	73-1016	73-1006



Ultra Fit Cool Wrist Splint 8"

- Features easy one-hand application that provides comfortable support of the wrist while preserving finger mobility
- Manufactured from durable, lightweight material, flexible dorsal stay and contact closure

Size	Wrist Circumference	Product No.	
		Left	Right
X-Small	4"-5"	73-1032	73-1022
Small	5"-6"	73-1033	73-1023
Medium	6"-7"	73-1034	73-1024
Large	7"-8"	73-1035	73-1025
X-Large	8"-9"	73-1036	73-1026



Cock-Up Splints (continued)

DeRoyal

Universal Cutaway Wrist and Wrist/Forearm Splint

- Perforated foam construction allows skin to breathe
- Loop lock closure at wrist provides compression and support
- Foam can be trimmed for desired fit

Size	Length	Product No.	
		Left	Right
Universal	8"	5040-04	5040-03
Universal	11"	5040-02	5040-01



Universal Foam Wrist and Wrist/Forearm Splint

- Perforated foam construction allows skin to breathe
- Moveable dorsal stay
- D-ring closures for easy application and removal

Size	Length	Product No.	
		Left	Right
Universal	8"	5066-83	5066-82
Universal	10"	BF5066-81	BF5066-80



Universal Wrist Splint

- Black foam material
- Moveable dorsal stay
- D-ring closures for easy application and removal

Size	Length	Product No.	
		Left	Right
Universal	7"	5067-02	5067-01



Black Wrist Splint

- Does not extend past distal palmar crease, allowing full finger range-of-motion

Size	Length	Product No.	
		Left	Right
Pediatric	8"	5073-12	5073-11
X-Small	8"	5073-06	5073-01
Small	8"	5073-07	5073-02
Medium	8"	5073-08	5073-03
Large	8"	5073-09	5074-04
X-Large	8"	5073-10	5074-05



Cock-Up Splints (continued)

Leatherette Wrist and Wrist/Forearm Splint

- Durable leatherette material with foam tricot lining
- Extra-wide loop lock strap at wrist for compression and support
- Does not extend past the distal palmar crease, allowing full finger range-of-motion

Size	Product No.			
	6" Left		6" Right	
Pediatric	8824-06		8824-01	
Size	8" Wrist		11" Wrist/Forearm	
	Left	Right	Left	Right
Small	8745-07	8745-02	8824-07	8824-02
Medium	8745-08	8745-03	8824-08	8824-03
Large	8745-09	8745-04	8824-09	8824-04
X-Large	8745-10	8745-05	8824-10	8824-05



Elastic Wrist Splint

- Cotton elastic construction provides compression
- Hook and loop closure for easy application and removal
- Does not extend past the distal palmar crease, allowing full finger range-of-motion

Size	6" Wrist Splint		8" Wrist Splint	
	Left	Right	Left	Right
X-Small	5015-06	5015-01	—	—
Small	5015-07	5015-02	5016-07	5016-02
Medium	5015-08	5015-03	5016-08	5016-03
Large	5015-09	5015-04	5016-09	5016-04
X-Large	5015-10	5015-05	5016-10	5016-05
2X-Large	—	—	5016-20	5015-15



Premium Universal Wrist Splint 8"

- Dynamic closure allows easy adjustment for an intimate fit

Size	Product No.
Universal (8")	5068-00



Universal Two-Panel Wrist and Wrist/Forearm Splint






- Blue canvas construction with foam tricot lining

Size	Product No.	
	8" Wrist	11" Wrist
Universal	1000275	1000285



DeRoyal

Cock-Up Splints (continued)

DeRoyal	Canvas Wrist and Wrist/Forearm Splint						
	<ul style="list-style-type: none"> Sturdy blue canvas material with foam tricot lining Extra-wide loop lock strap at wrist for compression and support Does not extend past the distal palmar crease, allowing full finger range-of-motion 						
	Size	4" Left	4" Right	6" Left	6" Right		
	Pediatric	1060115	1010115	1060025	1010025		
	Size	8" Left	8" Right	11" Left	11" Right		
	Small	1070115	1020115	1070025	1020025		
	Medium	1080115	1030115	1080025	1030025		
	Large	1090115	1040115	1090025	1040025		
	X-Large	1100115	1050115	1100025	1050025		
	DeRoyal	8" Suede Leatherette Wrist					
<ul style="list-style-type: none"> Extra-wide loop lock strap at wrist for compression and support Does not extend past the distal palmar crease, allowing full finger range-of-motion 							
Size		Product No.					
		Left			Right		
Small		12982105			12982205		
Medium		12982106			12982206		
Large		12982107			12982207		
X-Large		12982108			12982208		
DeRoyal		Wrist and Forearm Splint					
		<ul style="list-style-type: none"> Pre-molded aluminum design for anatomical fit and proper support Perforated design 					
	Size	Product No.					
		Left			Right		
	Youth	9101-05			9101-02		
	Adult	9101-06			9101-03		
	Freeman	7" Cock-Up Wrist Splint					
		<ul style="list-style-type: none"> Durable white canvas glove with soft flannel lining, elastic fly and removable 6" aluminum spoon Self adjusting two-cluster lacing attached to hook-and-pile straps ensures a proper fit along the entire length One pair of lightweight flexible stays added for extra stability 					
		Size	Wrist Circum.	Palm Circum.	Product No.		
					Left	Right	
2X-Small		4"-5 ¹ / ₄ "	5 ¹ / ₄ "-6 ¹ / ₄ "	8622-XXS	8621-XXS		
X-Small		4 ³ / ₄ "-6"	6"-7"	8622-XS	8621-XS		
Small		5 ¹ / ₄ "-6 ¹ / ₂ "	6 ³ / ₄ "-7 ¹ / ₂ "	8622-S	8621-S		
Medium		6"-7 ¹ / ₄ "	7 ¹ / ₄ "-8 ¹ / ₄ "	8622-M	8621-M		
Large		7"-8 ¹ / ₂ "	8 ¹ / ₄ "-9 ¹ / ₂ "	8622-L	8621-L		
X-Large		8 ¹ / ₄ "-9 ³ / ₄ "	9 ¹ / ₂ "-10 ³ / ₄ "	8622-XL	8621-XL		
Hely & Weber	UNO WHO Wrist Hand Orthosis						
	<ul style="list-style-type: none"> Stockinette tongue and dorsal loop and lock closures offer simple application Polyester outer shell with tricot covered foam lining for comfort and durability 						
	Size	Product No.					
	Left			Right			
Universal	428L			428R			

Cock-Up Splints (continued)

Suede Lacing Wrist Orthosis

- Perforated black suede laminated to moisture moving poly-felt
- Lacing design provides for easy secure closures
- Stockinette tongue eases application

Size	Wrist Circum.	Wrist		Wrist/Forearm	
		Left	Right	Left	Right
X-Small	4 ¹ / ₂ "-5 ¹ / ₂ "	434LXS	434RXS	—	—
Small	5"-6 ¹ / ₂ "	434LS	434RS	444LS	444RS
Medium	6"-7 ¹ / ₂ "	434LM	434RM	444LM	444RM
Large	7"-8 ¹ / ₂ "	434LL	434RL	444LL	444RL
X-Large	8"-9 ¹ / ₂ "	434LXL	434RXL	444LXL	444RXL



Lacing Wrist Orthosis



Lacing Wrist & Forearm

Phomfit™ Wrist Orthosis

- Perforated foam construction allows skin to breathe
- Moveable dorsal stay
- D-ring closures for easy application and removal

Size	Wrist Circum.	Short Product No.		Long Product No.	
		Left	Right	Left	Right
X-Small	4 ¹ / ₂ "-5 ¹ / ₂ "	436LXS	436RXS	437LXS	437RXS
Small	5"-6 ¹ / ₂ "	436LS	436RS	437LS	437RS
Medium	6"-7 ¹ / ₂ "	436LM	436RM	437LM	437RM
Large	7"-8 ¹ / ₂ "	436LL	436RL	437LL	437RL
X-Large	8"-9 ¹ / ₂ "	436LXL	436RXL	437LXL	437RXL



Short



Long

Universal Wrist Orthosis

- Comfortable orthopedic felt laminated to brushed nylon exterior
- Elastic carpal strap provides support and comfort

Size	Wrist Circum.	Product No.	
		Left	Right
Universal	4 ¹ / ₂ "-9 ¹ / ₂ "	438L	438R



TITAN™ Universal Wrist-Lacing Orthosis

- Soft, comfortable orthopedic felt laminated to brushed nylon exterior
- Adjustable lacing system handles a wide range of wrist sizes
- Stockinette tongue eases application

Size	Wrist		Wrist/Forearm	
	Left	Right	Left	Right
Universal	450L	450R	452L	452LR



Modabber™ Universal Wrist-Lacing Orthosis





- Patented patient assist single strap for easy application
- Kuhl perforated neoprene construction provides comfortable compression

Size	Short Product No.		Long Product No.	
	Left	Right	Left	Right
Universal	5819L	5819R	5818L	5818R








Hely & Weber




Cock-Up Splints (continued)

Med Spec	8" Wrist Lacer								
	<ul style="list-style-type: none"> Made with perforated suede and lined with polypropylene felt Uniform support through the entire wrist Blue palmar pad, anatomically contoured designed, and adjustable thumb strap provide a comfortable and supportive fit 								
	Size	Wrist Circum.	Black		Wrist/Forearm				
			Left	Right	Left	Right			
	X-Small	4 ¹ / ₂ "-5 ¹ / ₂ "	223921	223931	223221	223231			
	Small	5"-6 ¹ / ₂ "	223922	223932	223222	223232			
Medium	6"-7 ¹ / ₂ "	223924	223934	223224	223234				
Large	7"-8 ¹ / ₂ "	223925	223935	223225	223235				
X-Large	8"-9 ¹ / ₂ "	223926	223936	223226	223236				
Ossur Americas	Universal Wrist Lacer								
	<ul style="list-style-type: none"> Repositionable dorsal stay pocket allows the Universal Wrist Lacer to fit a wide range of patients Lined with polypropylene felt for comfort and has a single-pull lace closure for easy application 								
	Size	Wrist Circum.	Product No.						
			Left		Right				
	Universal (8")	5 ¹ / ₂ "-8 ¹ / ₂ "	223490		223491				
Universal (11")	5 ¹ / ₂ "-8 ¹ / ₂ "	223500		223501					
Ossur Americas	Form Fit Wrist Brace								
	<ul style="list-style-type: none"> Patented Form Fit 3-D molding technology allows us to vary padding thickness in strategic locations for superior comfort and fit Breathable Lycra-lined interior Lightweight durable fabric construction Contoured along distal palmar crease to allow for 90° M.P. flexion 								
	Size	Wrist Circum.	6" Wrist Brace		8" Wrist Brace			10" Wrist Brace	
			Left	Right	Left	Right		Left	Right
	X-Small	5 ¹ / ₄ "-6 ¹ / ₄ "	—	—	317082	317072		617082	617072
	Small	6 ¹ / ₄ "-7 ¹ / ₄ "	417083	417073	317083	317073		617083	617073
	Medium	7 ¹ / ₄ "-8"	417085	417075	317085	317075		617085	617075
	Large	8"-9"	417087	417077	317087	317077		617087	617077
	X-Large	9"-10"	417088	417078	317088	317078		617088	617078
	Exoform Wrist								
	<ul style="list-style-type: none"> Ergonomically contoured padding on the inside of the brace provides the best fit and most comfort Rigid plastic supports contour to the hand and wrist at the top and the bottom leading to greater support and immobilization Adjustable aluminum palmar stay allows for fine-tuning 								
Size	Wrist Circum.	Product No.							
		Left		Right					
X-Small	5 ¹ / ₄ "-6 ¹ / ₄ "	507082		507072					
Small	6 ¹ / ₄ "-7 ¹ / ₄ "	507083		507073					
Medium	7 ¹ / ₄ "-8"	507085		507075					
Large	8"-9"	507087		507077					
X-Large	9"-10"	507088		507078					

Cock-Up Splints (continued)

ProCare (DJO)	10" Deluxe Wrist and Forearm Supports						
	<ul style="list-style-type: none"> • Contoured design and durable perforated material provides improved patient fit and comfort • Loop-lock closure 						
	Size	Wrist Circum.	Product No.				
			Left		Right		
	X-Small	up to 5 ¹ / ₂ "	79-87052		79-87062		
	Small	5 ¹ / ₂ "-6 ¹ / ₂ "	79-87053		79-87063		
	Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	79-87055		79-87065		
	Large	7 ¹ / ₂ "-8 ¹ / ₂ "	79-87057		79-87067		
	X-Large	8 ¹ / ₂ "+	79-87058		79-87068		
	Quick-Fit® Wrist						
<ul style="list-style-type: none"> • Dorsal stay within adjustable stay pod can be moved proximal or distal for desired controlled range-of-motion • Designed to fit below palmar crease for full finger dexterity 							
Size	Wrist Circum.	Product No.					
		Left		Right			
Universal	up to 10"	79-87470		79-87460			
Universal/XL	10"+	79-87471		79-87461			
10" Quick-Fit Wrist and Forearm							
<ul style="list-style-type: none"> • Cool, soft felt/brush nylon construction with double-pull design for convenient application • Secondary dorsal stay may be removed to accommodate rehabilitation needs 							
Size	Wrist Circum.	Product No.					
		Left		Right			
Universal	up to 13"	79-87510		79-87500			
Universal/XL	13"+	79-87511		79-87501			
10" Universal Wrist and Forearm Supports							
<ul style="list-style-type: none"> • Breathable nylon/fiber with flannelette liner • Loop-lock closure 							
Size	Wrist Circum.	Product No.					
		Left		Right			
Universal	up to 9"	79-87060		79-87050			
Loop-Lock Cock-Up Splint							
<ul style="list-style-type: none"> • Sturdy canvas construction with loop-lock closure • Splint allows full finger function above palmar crease 							
Size	Wrist Circum.	Product No.					
		6" Left	6" Right	8" Left	8" Right		
X-Small/Ped.	up to 5 ¹ / ₂ "	79-87372	79-87362	79-87412	79-87402		
Small	5 ¹ / ₂ "-6 ¹ / ₂ "	79-87373	79-87363	79-87413	79-87403		
Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	79-87375	79-87365	79-87415	79-87405		
Large	7 ¹ / ₂ "-8 ¹ / ₂ "	79-87377	79-87367	79-87417	79-87407		
X-Large	8 ¹ / ₂ "+	79-87378	79-87368	79-87418	79-87408		

Cock-Up Splints (continued)

ProCare (DJO)	10" Foam Wrist Splint						
	• Foam tricot construction for breathable comfort						
			<u>Product No.</u>				
	Size	Wrist Circum.	Left		Right		
	X-Small	4 ¹ / ₂ "-5 ¹ / ₂ "	79-87432		79-87422		
	Small	5 ¹ / ₂ "-6 ¹ / ₂ "	79-87433		79-87423		
	Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	79-87435		79-87425		
	Large	7 ¹ / ₂ "-8 ¹ / ₂ "	79-87437		79-87427		
	X-Large	8 ¹ / ₂ "+	79-87438		79-87428		
	Lace-Up Wrist Support						
	• Perforated suede/flannel construction						
• Single pull laced closure for easy application							
		<u>Product No.</u>		<u>Product No.</u>			
Size	Wrist Circum.	7" Left	7" Right	10" Left	10" Right		
X-Small	up to 5 ¹ / ₂ "	79-87212	79-87202	79-87232	79-87222		
Small	5 ¹ / ₂ "-6 ¹ / ₂ "	79-87213	79-87203	79-87233	79-87223		
Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	79-87215	79-87205	79-87235	79-87225		
Large	7 ¹ / ₂ "-8 ¹ / ₂ "	79-87217	79-87207	79-87237	79-87227		
X-Large	8 ¹ / ₂ "+	79-87218	79-87208	79-87238	79-87228		
Truform	2364 Canvas Cock-Up Wrist Splint						
	• Durable 7" long canvas glove, flannel lined for comfort						
	• Soft elastic gusset						
	• Self-adjusting dual cluster lacing with criss-cross hook-and-pile fastening tabs						
	• Removable preshaped metal stay						
			<u>White</u>		<u>Gray</u>		
	Size	Width of Palm	Left	Right	Left	Right	
	2X-Small	3 ¹ / ₂ "-4 ¹ / ₂ "	2364/L-XXS	2364R-XXS	2364E/L-XXS	2364E/R-XXS	
X-Small	4 ¹ / ₂ "-5 ¹ / ₂ "	2364/L-XS	2364/R-XS	2364E/L-XS	2364E/R-XS		
Small	5 ¹ / ₂ "-6 ¹ / ₂ "	2364/L-S	2364/R-S	2364E/L-S	2364E/R-S		
Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	2364/L-M	2364/R-M	2364E/L-M	2364E/R-M		
Large	7 ¹ / ₂ "-8 ¹ / ₂ "	2364/L-L	2364/R-L	2364E/L-L	2364E/R-L		
X-Large	8 ¹ / ₂ "-9 ¹ / ₂ "	2364/L-XL	2364/R-XL	2364E/L-XL	2364E/R-XL		

Wrist Supports

- Indications: Frequently used for Carpal Tunnel Syndrome; Tenosynovitis; Post cast removal arthritis

WrisTimer PM™ Carpal Tunnel Support

- Foam padding; terry cloth covering
- Adjustable aluminum dorsal stay
- Hook and loop straps



Size	Product No.
Universal	510351

FREEDOM comfort™ Wrist Splint with MP Block


- Provides immobilization and protection of third, fourth and fifth MP and PIP joints
- Durable light weight fabric/foam laminate protects hand
- Ulnar stay positions MP and PIP joints in extension
- Hook and loop straps allow easy adjustment as hand volume changes



Size	Wrist Circumference	Product No.	
		Left	Right
Small/Medium	up to 8"	5910	5911
Large/X-Large	9"+	5912	5913

FREEDOM Short (5") Elastic Wrist Support

- Allows freedom to move and use fingers
- Has an adjustable, removable metal stay
- Elastic offers two varying degrees of compression: higher in the wrist, less on the hand



Size	Wrist Circumference	Product No.	
		Left	Right
X-Small	4 1/2"-5 1/2"	5760	5761
Small	5 1/2"-6 1/2"	5762	5763
Medium	6 1/2"-7 1/4"	5764	5765
Large	7 1/4"-8"	5766	5767

FREEDOM Long (8 1/2") Elastic Wrist Support

- Long forearm length optimally distributes pressure and hand weight without restricting hand function
- Elasticized cotton blend fabric creates compression in the wrist without feeling tight
- Adjustable removable aluminum stay
- Measure wrist circumference at base of hand



Size	Wrist Circumference	Product No.	
		Left	Right
2X-Small	0"-4 1/2"	5729	5728
X-Small	4 1/2"-5 1/2"	5731	5730
Small	5 1/2"-6 1/2"	5733	5732
Medium	6 1/2"-7 1/4"	5735	5734
Large	7 1/4"-8"	5737	5736
X-Large	8"-8 3/4"	5739	5738

VelPlush Universal Wrist Wrap Support

- Circumferential strap to help limit motion and removal splint for washing
- Hook and loop closure; Available in black








Size	Product No.	
	Left	Right
6"	W 60-L	W 60-R
8"	W68 L	W68 R

AIIMed

New Options
Sports

Wrist Supports (continued)

Ossur Americas	Exolite™ Carpal Tunnel Wrist				
	<ul style="list-style-type: none"> • 20% less circumferential bulk than other brands • Contoured for an anatomically correct fit • Preshaped into a neutral position for comfortable limited motion/flexibility 				
	Size	Wrist Circum.	Palm Circum.	Product No.	
				Left Right	
	Small	5 ¹ / ₄ "-6 ¹ / ₂ "	6 ³ / ₄ "-8"	517083 517073	
Medium	6 ¹ / ₂ "-7 ³ / ₄ "	8"-9 ¹ / ₄ "	517085 517075		
Large	7 ³ / ₄ "-10"	9 ¹ / ₄ "-10 ¹ / ₂ "	517087 517077		
ProCare (DJO)	CTS Wrist Supports				
	<ul style="list-style-type: none"> • Contoured cotton/elastic construction provides comfortable compression for proper support following mild sprains, strains and cast removal • Circumferential strap for improved wrist stabilization • Preformed aluminum stay provides anatomically correct fit 				
	Size	Wrist Circumference	Product No.		
			Left	Right	
	X-Small	up to 5 ¹ / ₂ "	79-87162	79-87152	
	Small	5 ¹ / ₂ "-6 ¹ / ₂ "	79-87163	79-87153	
Medium	6 ¹ / ₂ "-7 ¹ / ₂ "	79-87165	79-87155		
Large	7 ¹ / ₂ "-8 ¹ / ₂ "	79-87167	79-87157		
X-Large	8 ¹ / ₂ "+	79-87168	79-87158		
Professional Products	Jonah Neoprene Wrap Around Wrist Support				
	<ul style="list-style-type: none"> • Hook and loop compatible orthopedic neoprene • Removable palmer spoon • Universal sizing 				
	Size	Length	Product No.		
			Left	Right	
	Adult	7"	10017-00	10018-00	
	Ultimate™ V Wrist Supports				
	<ul style="list-style-type: none"> • Breathable high density hyper foam with black polyester cover or breathable moisture wicking orthopedic felt • EZY-On feature allows wrist support to be applied with one hand • EZY-Custom strap which allows equal tension regardless of forearm size 				
	Size	Product No.			
		Left	Right		
	8"	11479-15	11480-15		
8" with Thumb Spica	11481-15	11482-15			
10 ¹ / ₂ "	10692-15	10693-15			
Universal Ultimate Supreme V™ Wrist					
<ul style="list-style-type: none"> • Breathable high density hyper foam with black polyester cover or breathable moisture wicking orthopedic felt • Contoured design optimizes hand and wrist fit with MP function • Removable palmar stay for functional management of the wrist injury 					
Size	Product No.				
	Left	Right			
8"	11486-00	11487-00			
8" with Thumb Spica	11509-00	11510-00			
10 ¹ / ₂ "	11507-00	11508-00			

Wrist Supports (continued)

RCAI	Dorsal Carpal Tunnel Splint			
	<ul style="list-style-type: none"> • Lightweight, yet durable Kydex plastic, the splint can be heat molded for patient customization and positioning • Can be used as a wrist extension splint to support a weak or flaccid wrist • Available in black 			
	Size	Width of Palm	Product No.	
			Left	Right
	X-Small	2"-2 ¹ / ₂ "	44DCT-ES-L	44DCT-ES-R
	Small	2 ¹ / ₂ "-3"	44DCT-S-L	44DCT-S-R
Medium	3"-3 ¹ / ₂ "	44DCT-M-L	44DCT-M-R	
Large	3 ¹ / ₂ "-4"	44DCT-L-L	44DCT-L-R	
X-Large	4"-4 ¹ / ₂ "	44DCT-XL-L	44DCT-XL-R	
Truform	2360 Soft Fit Suede Finish Wrist Brace			
	<ul style="list-style-type: none"> • Designed to help relieve the discomfort associated with repetitive motion • Maintains the wrist in a neutral position to promote healing, while allowing normal activity • Opens flat for easy application; Hook and loop fastening economical and effective 			
	Size	Fits Wrist	Product No.	
			Left	Right
	X-Small	up to 4 ¹ / ₂ "	2360/L-XS	2360/R-XS
	Small	4 ¹ / ₂ "-5 ¹ / ₂ "	2360/L-S	2360/R-S
Medium	5 ¹ / ₂ "-6 ¹ / ₂ "	2360/L-M	2360/R-M	
Large	6 ¹ / ₂ "-7 ¹ / ₂ "	2360/L-L	2360/R-L	
X-Large	7 ¹ / ₂ "+	2360/L-XL	2360/R-XL	



AliMed - 3pp[®] Comforter[™] Splint

- Thick Kydex[™] shell positions the wrist in slight extension, holding the MP joints in neutral and allowing IP function
- Adjustable, padded finger and wrist stays for perfect alignment
- CMC joint is supported, allowing thumb to function freely
- Neoprene straps with finger holes facilitate splint application and strap adjustment



Size	Wrist Circum.	Product No.	
		Left	Right
Small	2 ⁵ / ₈ "-3 ¹ / ₈ "	51386/NA/LS	51386/NA/RS
Medium	3 ¹ / ₈ "-3 ¹ / ₂ "	51386/NA/LM	51386/NA/RM
Large	3 ¹ / ₂ "-3 ⁵ / ₈ "	51386/NA/LL	51386/NA/RL

Measure MP joints across dorsal surface from radial side of 2nd metacarpal to ulnar side of 5th metacarpal

FREEDOM[®] Ultimate Grip Splint

- Malleable aluminum core bends to fit required wrist and forearm angle
- For incremental finger extension, splint comes with two foam Ethafoam cylinders: one 1¹/₂" and one 1³/₄" diameter
- Three-part soft goods set is included to accommodate use with cylinder, or without for severe contractures
- One size with hook and loop closures and terry/foam laminate cover



Side	Product No.
Left	510841/NA/LF
Right	510841/NA/RT

ULTRApadded[™] Grip Orthosis

- Custom fit orthosis by bending foam-covered wire structural platform to accommodate any wrist angle
- Removable finger spacer and laundry bag included
- Max width at MP joints: 4¹/₄"; Max length from longest fingertip to mid palm: 5¹/₄"
- Available in Aqua



Product No. **510327**

Dorsal Resting Splint

- Allows for maximum range-of-motion without triggering flexor spasticity because contact is on dorsal surface of hand and wrist
- Secure three straps and splint is in place
- Integrated thumb piece allows thumb to rest naturally in slight flexion and opposition without risk of skin breakdown
- Malleable removable core



Side	Product No.
Left	513419
Right	513418

Grip Splint II

- Supports wrist and fingers in neutral position and thumb in functional position
- Heat-moldable Kydex® core with removable 1/4" Plastazote® foam padding
- One-piece cover in absorbent terry cloth or fleece
- One size fits left or right hand



Description	Product No.
Grip Splint II with Terry Cover	510328
Grip Splint II with Fleece Cover	510604

Adjustable Pil-O-Splint®

- Same as the standard Pil-O-Splint but adjustable
- Cushioning foam and breathable poly-cotton lining
- Two metal dorsal stays
- One size fits left or right

Product No. **51331**



Deluxe Functional Position Splint

- Classic hand positioning for post-stroke patients - thumb in full opposition, wrist in 30° extension, and MP joints at approx. 80°
- Made of adjustable thermoplastic with straps attached

Size	Measurements	Product No.	
		Left	Right
X-Small	fits 5 1/4"-6"	5180	5181
Small	fits 6 1/2"-7 1/4"	5182	5183
Medium	fits 7 1/4"-8 1/4"	5184	5185
Large	fits 8 1/4"-9 1/4"	5186	5187
Measure the circumference of the hand at the MP joints			



Ball Splint

- Uses a reflex inhibiting position to control spasticity and tone
- Adjustable thermoplastic
- Straps included

Size	Measurements	Product No.	
		Left	Right
Small	fits 6 1/2"-7 1/4"	5375	5378
Medium	fits 7 1/2"-8 1/4"	5376	5379
Large	fits 8 1/2"-9 1/4"	5377	5388
Replacement Strap Kit		51-101	
Measure the circumference of the hand at the MP joints			



Shown with strap set 51-101

FREEDOM® Omni Progressive WHT Orthosis

- Finger portion of splint offers a natural resting trough for each finger, helping to keep fingers appropriately aligned
- Wrist angle may be modified from 75° of flexion to 45° of extension; thumb piece can be adjusted to accommodate any degree of abduction or flexion
- Strong malleable aluminum core
- Sizing: Standard fits most women; Long fits larger women and most men



Size	Product No.	
	Left	Right
Standard	510741/NA/LF	510741/NA/RT
Long	513414/NA/LF	513414/NA/RT

DeRoyal - Hand Control Mitten with Loop Lock Closure

- Breathable coolknit or cotton construction with double-padded design (M2026V)
- Finger separators help prevent tactile movement
- Available with or without strap to limit range-of-motion if needed
- Sold per pair

Indications: Patients at risk of unplanned extubations and patients who disrupt medical treatment

Size	Style	Strap Attachment	Product No.
Universal	Coolknit	Optional - Tie	M2023V-A
Universal	Coolknit	None	M2040
Universal	Cotton	Optional - Tie	M2026V




Full Hand Splint

- Pre-molded aluminum design for anatomical fit and proper support
- Perforated design for air circulation
- X-ray lucent

Indications: Positioning of hand, fingers and thumb

Size	Product No.	
	Left	Right
Child	9151-04	9151-01
Youth	9151-05	9151-02
Adult	9151-06	9151-03




LMB Air-Soft™ Resting Hand Splint

- Supports fingers, thumb and wrist in functional/resting position with simple adjustment of the wrist and thumb
- Malleable aluminum frame allows for static progressive splinting
- Self wicking liner helps to control perspiration

Indications: Post-surgery/injury positioning; Wrist drop and finger flexor tightness; Not recommended when moderate to severe spasticity is required

Size	Measurements	Product No.	
		Left	Right
Small	2 ¹ / ₄ "-2 ³ / ₄ "	325BL	325BR
Medium	2 ⁷ / ₈ "-3 ¹ / ₄ "	325CL	325CR
Large	3 ³ / ₈ "-3 ³ / ₄ "	325DL	325DR



Pucci® Air Inflatable Hand Splints

- Wrist, hand and finger therapeutic orthoses apply low stretch therapy to contracted fingers, hand and wrist
- Features inflatable air bladders
- Fleece type liner helps reduce pressure points and provide comfortable fit
- Air-T version's air bladder has contoured area for thumb to rest in

Indications: Stroke; Spinal cord injuries; Head injuries; Geriatric patients

Size	Air-T		Air		Air Short	
	Left	Right	Left	Right	Left	Right
Pediatric	405PL	405PR	413PL	413PR	403PL	403PR
Adult	405L	405R	413L	413R	403L	403R



Air-T



Air



Air Short

LEEDer Group - 1-2-3 PRO-GRIP Orthosis

- Flex-E-Core Neuro-Dynamic Therapy designed for a “fisted” hand (clenched fingers) and moderate to severe wrist contracture
- Provides neuro-rehabilitation while the patient wears the device
- Provides gentle ulnar drift correction and is supplied with extra cover & finger separators

Indications: Moderate to severe flexion; Spasticity of wrist, hand and fingers



Size	Wrist to Elbow Crease	Product No.	
		Left	Right
Small	up to 8”	LGSL	LGSR
Medium	8”-9 ¹ / ₂ ”	LGML	LGMR
Large	9 ¹ / ₂ ”+	LGLL	LGLR
Air Grip	Universal Fit	LGAA	

1-2-3 PRO-REST Hand Splint

- Flex-E-Core Neuro-Dynamic Therapy designed for mild to moderate contractures of wrist, hand and fingers
- REST is adjustable without heat or tools
- Includes extra cover and finger separators
- Sized by distance between wrist and elbow crease



Size	Width of MP Joint	Wrist to Fingertips	Wrist to Elbow Crease	Fleece		Velfoam	
				Left	Right	Left	Right
Pediatric	up to 2 ¹ / ₂ ”	up to 6”	up to 6”	LRPL	LRPR	LRPLNF	LRPRNF
Small	2 ¹ / ₂ ”-3”	6”-7”	6”-8”	LRSL	LRSR	LRSLNF	LRSRNF
Medium	3”-3 ¹ / ₂ ”	7”-8”	8”-9 ¹ / ₂ ”	LRML	LRMR	LRMLNF	LRMRNF
Large	3 ¹ / ₂ ”-4”	8”-9 ¹ / ₂ ”	9 ¹ / ₂ ”	LRL	LRLR	LRLNF	LRLRNF

Lenjoy - Comfyrene Separate Finger Hand Orthosis

- Includes removable finger separators
- Adjust individual fingers and thumb
- Low-profile straps adhere anywhere
- Available in Camouflage and Purple

Indications: Dupuytren's; Swan neck; Boutonniere contractures; Sprains; Fractures



optional finger separator

Size	Wrist to Top of Middle Finger	Camouflage Product No.		Purple Product No.	
		Left	Right	Left	Right
Pediatric	4 ¹ / ₂ ”-6”	HSF-101-CP-PL-C-L	HSF-101-CP-PL-C-R	HSF-101-CP-PL-P-L	HSF-101-CP-PL-P-R
Adult	6”-8”	HSF-101-CP-A-C-L	HSF-101-CP-A-C-R	HSF-101-CP-A-P-L	HSF-101-CP-A-P-R

Comfyrene Hand Thumb Orthosis

- Supports the hand, wrist, fingers and thumb without stretching the thenar eminence
- Comes with disposable liners
- Available in Light Blue, Dark Blue and Purple

Indications: Adducted thumbs



Size	Wrist to Top of Middle Finger	Product No.		
		Light Blue	Dark Blue	Purple
Pediatric Small	2 ¹ / ₂ "-3 ¹ / ₂ "	HT-101-CP-PS-LB	—	—
Pediatric Medium	3"-4 ¹ / ₂ "	HT-101-CP-PM-LB	—	—
Pediatric Large	4 ¹ / ₂ "-6"	HT-101-CP-PL-LB	HT-101-CP-PL-DB	—
Adult	7"-8"	HT-101-CP-A-LB	—	—
Adult Small	6"-7"	—	—	—
Adult Large	8" +	—	—	HT-101-CP-AL-P

Comfy Adjustable Cone Hand Orthosis

- Developed to treat the patient with strong flexion synergy of the fingers and wrist
- Splint is worn on a lateral side of the arm and has an adjustable hinge at the wrist to help increase wrist extension
- Fingers positioned around a cone
- Finger separators can be placed on a cone to separate fingers and prevent pressure areas



Size	Wrist to Top of Middle Finger	Total Length	Product No.
Adult	7"-8"	12 ¹ / ₂ "	ACH-101
Adult with 2 Covers	7"-8"	12 ¹ / ₂ "	ACH-102

Comfy Hand Wrist Finger Orthosis

- Gradual extension of non-fixed wrist, hand or finger flexion
- Easy adjustability at all joints allows them to be used for multiple indications to optimize functionality
- Wings are bent upwards to prevent ulnar/radial deviation
- Available in Navy Blue



**Shown with optional finger separator*

Size	Wrist to Top of Middle Finger	Product No.
Pediatric Small	2 ¹ / ₂ "-3 ¹ / ₂ "	PH-101-S
Pediatric Medium	3"-4 ¹ / ₂ "	PH-101-M
Pediatric Large	4 ¹ / ₂ "-6"	PH-101-L
Adult	7"-8"	H-101
Accessories		
Cover for Comfy Hand		HC-1

Optional Comfy Finger Separator

- Separator has soft padded dividers that position fingers and prevent pressure areas
- Elastic strap allows for positioning and holds device in place
- FS-1 is not attached to the splint, but is an extra device that can be positioned to protect and prevent pressure areas between the or DIP joints in any angle
- Available in Navy Blue



Size	Measurements	Product No.
Pediatric Small	2 ¹ / ₂ "-3 ¹ / ₂ "	PFS1-S
Pediatric Medium	3"-4 ¹ / ₂ "	PFS1-M
Pediatric Large	4 ¹ / ₂ "-6"	PFS1-L
Adult	7"-8"	FS-1

Comfy Hand Thumb Orthosis

- Works the same as the Standard Comfy Hand Orthosis with extra tab for adducted thumb
- Splint is ambidextrous
- Optional removable finger separator (FS-1)
- Available in Navy Blue



Size	Wrist to Top of Middle Finger	Product No.
Pediatric Small	2 ¹ / ₂ "-3 ¹ / ₂ "	PHT-101-S
Pediatric Medium	3"-4 ¹ / ₂ "	PHT-101-M
Pediatric Large	4 ¹ / ₂ "-6"	PHT-101-L
Adult Small	6"-7"	HT-101-AS
Adult	7"-8"	HT-101

Comfy Grip Hand

- Adjustable orthosis used to support the wrist and gradually increases finger extension
- Frame is designed to slide on graduated finger rolls to increase finger extension
- Durable terry cloth absorbs moisture and prevents skin maceration
- Strap closes with hook and loop closure



Size	Wrist to Palmar Creases	Wrist Down Forearm	Total Length	Product No.	
				Left	Right
Small	2"	5"	8"	CGRIP-S-LT	CGRIP-S-RT
Medium	2 ³ / ₄ "	6 ¹ / ₄ "	9 ¹ / ₂ "	CGRIP-M-LT	CGRIP-M-RT
Large	3 ¹ / ₄ "	7"	12 ¹ / ₂ "	CGRIP-L-LT	CGRIP-L-RT

Comfy Air Hand Orthosis

- Adjustable air bladder that allows for gradual extension of the MP, DIP and POP joints, as well as an adjustable wrist for customized fit
- Use pump to inflate the bladder and increase finger extension

Product No.

HA-101



Hand Pump

Comfy Finger Extender Hand Orthosis

- Ideal for the severe flexor pattern of the wrist and fingers when used without the hand rolls
- Gradually build up finger extension by using graduated finger rolls and extending the spine of the splint for wrist extension
- Available in Navy Blue

Description	Wrist to Tip of Middle Finger	Total Length	Product No.
Adult with 1 Cover	7"-8"	12 ¹ / ₂ "	F-101
Adult with 2 Covers	7"-8"	12 ¹ / ₂ "	F-102



Covers

Comfy Opposition Thumb Hand Orthosis

- Adjustable "C" bar support to hold the thumb in opposition
- Provides gradual extension of the wrist, hand, finger and thumb
- Comes with patented Comfy frame and absorbent terry cloth cover
- Available in Navy Blue

Size	Wrist to Top of Middle Finger	Product No.	
		Left	Right
Pediatric Small	2 ¹ / ₂ "-3 ¹ / ₂ "	POPH-101-S-LT	POPH-101-S-RT
Pediatric Medium	3"-4 ¹ / ₂ "	POPH-101-M-LT	POPH-101-M-RT
Pediatric Large	4 ¹ / ₂ "-6"	POPH-101-L-LT	POPH-101-L-RT
Adult	7"-8"	OPH-101-L	OPH-101-R



OCSI/Neuroflex - AirPro™ Air Graduate Wrist Hand Finger Orthosis

- Conforms to the shape of the hand accommodating different degrees of finger contracture
- Removable cone shaped palmar bladder can be used independently of wrist support

Indications: Severely contracted fingers

Product No.

AP31AG-R/L



SoftPro™ Grip

- Palmar bar comfortable fits into the palm of a fist hand
- Two incrementally larger palmar rolls included to gradually increase extension of the MCP and PIP joints

Indications: Moderate to severe contractures of the wrist, hand and fingers

Size	Width	Length	Product No.	
			Left	Right
Small	2 ¹ / ₂ "-3"	5"-6"	SP31FG-SL	SP31FG-SR
Medium	3"-3 ¹ / ₂ "	6 ¹ / ₂ "-8"	SP31FG-ML	SP31FG-MR
Large	3 ¹ / ₂ "-4"	8"-9 ¹ / ₂ "	SP31FG-LL	SP31FG-LR



SoftPro™ Champ Resting Hand

- Treats mild to moderate contractures of the wrist/hand/fingers
- Maintains/supports the wrist/hand/fingers in functional alignment
- Unique strapping system provides optimal support
- Width is measurement of MP Joints; Length is measurement of wrist to fingertip

Size	Width	Length	Product No.	
			Left	Right
X-Small	up to 2 ¹ / ₂ "	up to 5"	SP31CHA-XSL	SP31CHA-XSR
Small	2 ¹ / ₂ -3"	5"-6"	SP31CHA-SL	SP31CHA-SR
Medium	3"-3 ¹ / ₂ "	6 ¹ / ₂ "-8"	SP31CHA-ML	SP31CHA-MR
Large	3 ¹ / ₂ "-4"	8"-9 ¹ / ₂ "	SP31CHA-LL	SP31CHA-LR



SoftPro™ Functional Resting Hand

- Treats moderate contractures of the wrist/hand/fingers
- Bend-to-fit; Easily molds without tools yet sturdy enough to maintain shape during wear
- Extra long wrist strap holds the wrist in the proper position and applies a gentle stretch
- Width is measurement of MP Joints; Length is measurement of wrist to fingertip


Size	Width	Length	Product No.	
			Left	Right
Small	2 ¹ / ₂ "-3"	5"-6"	SP31FFR-SL	SP31FFR-SR
Medium	3"-3 ¹ / ₂ "	6 ¹ / ₂ "-8"	SP31FFR-ML	SP31FFR-MR
Large	3 ¹ / ₂ "-4"	8"-9 ¹ / ₂ "	SP31FFR-LL	SP31FFR-LR



Orthomerica - Wrist-Hand Cock-Up Splint

- Anatomical thumb design for comfort and maximum control
- Maintains wrist and hand in slight extension
- Platform extends to the end of the fingers
- Kydex® material can be heat modified and trimmed

Indications: Arthritic hand positioning; Reconstructive surgery; Spasticity



Size	MP Width	Brace Length	Product No.	
			Left	Right
X-Small	up to 2 ¹ / ₂ "	11"	1019	1024
Small	3"	13"	1020	1025
Medium	3 ¹ / ₂ "	14"	1021	1026
Large	4 ¹ / ₂ "	14 ¹ / ₂ "	1022	1027

ProCare (DJO) - IMAK Pil-O-Splint

- Exceptional comfort and maximum pain relief during rest and sleep
- Soft cotton lining extends to fingertips to keep fingers from moving while sleeping
- Two rigid plastic dorsal stays immobilize the hand, wrist, and fingers in a neutral position
- Breathable foam covered with soft poly flannel promotes warmth and good circulation

Indications: Carpal Tunnel Syndrome; Forearm tendonitis; Arthritis, Other wrist and hand pain



Description	Wrist Circum.	Length	Product No.
Adjustable	12"+	12"	79-87550

RCAI - Adaptable Contour Hand Orthosis

- Made of a durable, foam covered aluminum, the Adaptable Contour Hand can be bent without heat or tools, offering progressive extension of the wrist
- Finger separators, sewn to the liner, provide finger abduction and discourage radial and ulnar deviation of the hand
- Breathable liner is made of an anti-microbial "wicking" material that keeps moisture away from the skin, maintaining good skin integrity

Size	Width of Palm	Product Length	Product No.	
			Left	Right
Small	2 ¹ / ₂ "-3"	13"	300-ACH-S-L	300-ACH-S-R
Medium	3"-3 ¹ / ₂ "	14 ¹ / ₂ "	300-ACH-M-L	300-ACH-M-R
Large	3 ¹ / ₂ "-3 ³ / ₄ "	16"	300-ACH-L-L	300-ACH-L-R



Adaptable Resting Hand

- Made of a durable, foam covered aluminum that can be bent without heat or tools, offering progressive extension of the wrist
- Finger separators, sewn to the liner, provide finger abduction and discourage radial and ulnar deviation of the hand
- Liner is made of an anti-microbial “wicking” material that keeps moisture away from the skin, maintaining good skin integrity

Size	Width of Palm	Product Length	Product No.	
			Left	Right
Small	2 ¹ / ₂ "-3"	11"	390-ARH-S-L	390-ARH-S-R
Medium	3"-3 ¹ / ₂ "	12"	390-ARH-M-L	390-ARH-M-R
Large	3 ¹ / ₂ "-3 ³ / ₄ "	14"	390-ARH-L-L	390-ARH-L-R



Adaptable Geriatric Hand

- Foam covered aluminum can be bent without heat or tools, offering progressive extension of the wrist and finger positioning
- Adjustable palmar web space abduction with the use of a foam roll for the severely adducted thumb
- Finger separators, sewn to the liner, provide finger abduction and discourage radial and ulnar deviation of the hand
- Liner is made of an anti-microbial “wicking” material that keeps moisture away from the skin, maintaining good skin integrity

Size	Width of Palm	Product Length	Product No.	
			Left	Right
X-Small	2"-2 ¹ / ₂ "	12"	320-AGH-ES-L	320-AGH-ES-R
Small	2 ¹ / ₂ "-3"	12"	320-AGH-S-L	320-AGH-S-R
Medium	3"-3 ¹ / ₂ "	13"	320-AGH-M-L	320-AGH-M-R
Large	3 ¹ / ₂ "-3 ³ / ₄ "	13"	320-AGH-L-L	320-AGH-L-R



Geriatric Hand with Finger Separators

- Lightweight, durable Kydex® plastic construction provides stable wrist and finger positioning
- Features adjustable palmar web space abduction with the use of a foam roll, for the severely adducted thumb
- A breathable anti-microbia wicking foam liner “wicks” away excess moisture and maintains skin integrity

Indications: Stroke; Arthritis; Hand injury; Treatment of finger, wrist and thumb contractures

Size	Measure Width of Palm	Length	Product No.	
			Left	Right
X-Small	2" - 2 ¹ / ₂ "	10"	32GHK-ES-L	32GHK-ES-R
Small	2 ¹ / ₂ " - 3"	13"	32GHK-S-L	32GHK-S-R
Medium	3" - 3 ¹ / ₂ "	14"	32GHK-M-L	32GHK-M-R
Large	3 ¹ / ₂ " - 3 ³ / ₄ "	15"	32GHK-L-L	32GHK-L-R
X-Large	3 ³ / ₄ " - 4 ¹ / ₄ "	16"	32GHK-XL-L	32GHK-XL-R



Contour Hand Orthosis Model 30CHK with Kydex®

- Can be positioned to meet progressive patient needs
- Anti-microbial liner included

Indications: Moderate to severely contracted hand and wrist



Size	Width of Palm	Product Length	Product No.	
			Left	Right
X-Small	2"-2 ¹ / ₂ "	7"	30CHK-ES-L	30CHK-ES-R
Small	2 ¹ / ₂ "-3"	8"	30CHK-S-L	30CHK-S-R
Medium	3"-3 ¹ / ₂ "	9"	30CHK-M-L	30CHK-M-R
Large	3 ¹ / ₂ "-3 ³ / ₄ "	9 ¹ / ₂ "	30CHK-L-L	30CHK-L-L

Resting Hand Orthosis

- Can be heat molded at low temperature
- Liner eliminates potential pressure points, while moisture is "wicked" away keeping the skin dry and comfortable
- Raised lateral palmar ridges of the orthosis discourage radial or ulnar deviation, and reinforce the orthosis to accommodate muscle tone and spasticity
- Finger separators provide abduction and skin integrity and promotes proper finger positioning

Indications: Stroke; Trauma; Injury; Rheumatoid arthritis, Arthroplasties; Burns and early wrist and finger contractures



Size	Measurements	Length	Without Separators		With Separators	
			Left	Right	Left	Right
Infant	3/4"-1 ¹ / ₂ "	5"	29RHO-I-L	29RHO-I-R	—	—
Toddler	1 ¹ / ₄ "-1 ³ / ₄ "	7"	29RHO-T-L	29RHO-T-R	—	—
Child	1 ³ / ₄ "-2"	8"	29RHO-C-L	29RHO-C-R	—	—
Youth	2"-2 ¹ / ₂ "	9"	29RHO-Y-L	29RHO-Y-R	—	—
X-Small	2"-2 ¹ / ₂ "	10"	29RHO-ES-L	29RHO-ES-R	—	—
Small	2 ¹ / ₂ "-3"	11"	29RHO-S-L	29RHO-S-R	39SRH-S-L	39SRH-S-R
Medium	3"-3 ¹ / ₂ "	11 ¹ / ₂ "	29RHO-M-L	29RHO-M-R	39SRH-M-L	39SRH-M-R
Large	3 ¹ / ₂ "-3 ³ / ₄ "	13"	29RHO-L-L	29RHO-L-R	39SRH-L-L	39SRH-L-R
X-Large	3 ³ / ₄ "-4 ¹ / ₄ "	14 ¹ / ₂ "	29RHO-XL-L	29RHO-XL-R	39SRH-XL-L	39SRH-XL-R

Replacement Liners for Resting Hand Orthosis

- Replace ? with size

Description	Product No.	
	Left	Right
Replacement Liner	29AK-?-L	29AK-?-R



Dorsal Resting Hand (Without Finger Separators)

- Dorsal arm platform frees the palmar area, allows for greater tone inhibition and enhances sensory feedback
- Made of lightweight, durable Kydex® plastic
- Anti-microbial foam liner eliminates potential pressure points and “wicks” moisture away from the skin

Indications: Treatment of spasticity and contracture due to CVA, closed head injuries, spinal cord injuries, and cerebral palsy



Size	Width of Palm	Product Length	Product No.	
			Left	Right
X-Small	2"-2 ¹ / ₂ "	12"	15DHO-ES-L	15DHO-ES-R
Small	2 ¹ / ₂ "-3"	13"	15DHO-S-L	15DHO-S-R
Medium	3"-3 ¹ / ₂ "	14 ¹ / ₂ "	15DHO-M-L	15DHO-M-R
Large	3 ¹ / ₂ "-3 ³ / ₄ "	16"	15DHO-L-L	15DHO-L-R
X-Large	3 ¹ / ₂ "-4 ¹ / ₄ "	16"	15DHO-XL-L	15DHO-XL-R

Truform - 1728 Wrist-Hand Orthosis

- Maintains wrist in a neutral position with fingers in extension and thumb in opposition
- Vacuum formed plastic shell with hook and loop encircling straps

Indications: Post stroke to prevent muscle contracture



Size	Palmar Width	Product No.	
		Left	Right
X-Small	up to 2 ¹ / ₂ "	1728/L-XS	1728/R-XS
Small	up to 2 ¹ / ₄ "	1728/L-S	1728/R-S
Medium	up to 3 ¹ / ₄ "	1728/L-M	1728/R-M
Large	up to 3 ³ / ₄ "	1728/L-L	1728/R-L
X-Large	up to 4"	1728/L-XL	1728/R-XL

AliMed - Turnbuckle Functional Position Splint


- Steel rod and turnbuckle mechanism provide a low-load stretch to the joint to gradually reduce a wrist flexion contracture while maintaining a functional hand position
- Rigid Kydex® splint shell is hinged at the wrist, with the turnbuckle connecting the hand and forearm portions of the splint
- Padded two-part Beta Pile II™ cover is comfortable against skin and is removable
- Sizing - approximate: B fits smaller women; C fits larger women and smaller men; D fits most men
- Adjusts from 80 degrees wrist flexion to 60 degrees wrist extension



Product No. **52395**

Becker Orthopedic - Model UM2.2 Dynamic WHO with Finger Extension


- Easily adjustable in several areas to ensure precision in fit and application of force
- Provides dynamic extension of DIP and PIP joints II-V and fixation of the wrist
- MCP extension block provides fixation of MCP joint II-V
- Replacement Bands (P/N UM-Bands)



Size	Measurements	Product No.	
		Left	Right
Small	2 ¹ / ₈ "-2 ³ / ₄ "	UM-2.2LS	UM-2.2RS
Medium	2 ³ / ₄ "-3 ³ / ₈ "	UM-2.2LM	UM-2.2RM
Large	3 ³ / ₈ "-4 ³ / ₈ "	UM-2.2LL	UM-2.2RL

Model UM3.3 Dynamic WHO with Wrist and Finger Extension

- Easily adjustable in several areas to ensure precision in fit and application of force
- Dynamic extension force to the wrist, DIP and PIP joints II-V
- Removable MCP extension block provides fixation of MCP joints II-V
- Replacement Bands (P/N UM-Bands)




Size	Measurements	Product No.	
		Left	Right
Small	2 ¹ / ₈ "-2 ³ / ₄ "	UM-3.3LS	UM-3.3RS
Medium	2 ³ / ₄ "-3 ³ / ₈ "	UM-3.3LM	UM-3.3RM
Large	3 ³ / ₈ "-4 ³ / ₈ "	UM-3.3LL	UM-3.3RL

DeRoyal - DeROM® Wrist

- Low load prolonged stretch to soft tissue to help restore joint's ROM
- Easy to set and release tension with integrated adjustment tool and tension release lever
- Tension Limiting Control (TLC) knob prohibits excessive force
- Range-of-motion stops for static progressive splinting
- Flexion and extension assist in single unit

Indications: Joint stiffness; ROM limitations; Soft tissue contractures



Size	Measurement	Product No.
Universal	Fits 2 ¹ / ₄ " - 3 ¹ / ₄ "	4219D
Replacement Softgood	—	4219DK

Lenjoy - Spring Loaded Goniometer Hand "Set ROM with Extension"

- Dynamic hand orthosis that gives a gentle push to extension of the wrist or can be set in a fixed static position
- Ability to set the range of flexion and extension desired at the wrist joint by setting patented Goniometer dial
- Re-adjustable custom fit of wrist, MP, PIP & DIP joints
- Goniometer dial and metal are portected with padded soft cover



Size	Wrist Measurement	Length	Product No.
Pediatric Small	2 ¹ / ₂ "-3 ¹ / ₂ "	7"	PSGH-101-S
Pediatric Medium	3"-4 ¹ / ₂ "	8"	PSGH-101-M
Pediatric Large	4 ¹ / ₂ "-6"	9 ¹ / ₂ "	PSGH-101-L
Adult Small	6"-7"	12 ¹ / ₂ "	SGH-101-AS
Adult	7"-8"	11"	SGH-101

AliMed - Inflatable Carrot

- Inflatable cone-shaped tube covered with velvety bacteriostatic cloth
- With the carrot deflated and the cover in place, the orthosis is pulled into the severely contracted hand using a placement wand
- Once inserted into the hand, the Carrot can be gradually inflated to reduce even the most severe contracture
- Max diameter: Length - 7¹/₄", Width - 5³/₄"



Description	Product No.
Carrot only	51594
Accessory Kit*	51595
Inflatable Carrot and Accessory Kit	51598

*Kit contains two cloth covers, inflator bulb and two wands

Joint Jack®

- Applies pressure on joint for progressive reduction of a flexion contracture
- Metal finger bed; cotton felt
- Padded screw tightens to increase tension on the finger joint

Size	Measurements	Product No.
Small	1 ¹ / ₂ "-2 ¹ / ₂ "	5147
Medium	2"-3"	51384
Large	3"-4"	51383
X-Large	3 ¹ / ₂ "-4 ¹ / ₂ "	5148
2X-Large	4"	5142

Measure length of finger from distal palmar crease



G-Force Boxer's Fracture Orthosis with MP Flexion

- Positions metacarpal joint in 75° of flexion during immobilization
- Allows index finger and thumb to perform minor tasks throughout the healing process
- Rigid shell lined with cushioning foam for immobilization and protection of fracture and surrounding joints

Indications: Uncomplicated fourth and fifth metacarpal fractures, including Boxer's fracture




Size	Standard Product No.		Long Product No.	
	Left	Right	Left	Right
Small	—	—	52202	52201
Medium	52208	52207	52204	52203
Large	52210	52209	52206	52205

Universal Digit Splint

- Aluminum finger strip
- Bendable aluminum stay
- Fits any finger
- Available in Light gray

Product No. **5398**



FREEDOM® comfort Boxer Fracture Prefab Orthosis with MP Extension

- T-Foam™ liner ensures comfort, support and conformity to the fracture volume
- Rigid Kydex® shell may be heated and reshaped to accommodate swelling
- Use for uncomplicated fourth and fifth metacarpal fractures
- Positions MP's in extension

Indications: Uncomplicated fourth and fifth metacarpal fractures, including Boxer's fracture



Size	Measurements	Product No.	
		Left	Right
Medium	7 -9"	5112	5111
Large	9"-10 ¹ / ₂ "	5114	5113


Measure circumference of hand at MP's

DeRoyal - Baseball Finger Splints

- Lightweight aluminum can be adjusted for proper fit; X-ray lucent
- Easily moldable to maintain desired position
- Packaged 12/case

Indications: Fractures of the distal phalanx

Size	Measurements	Product No.
Small	4 ¹ / ₄ "	9109-01
Medium	4 ³ / ₄ "	9109-02
Large	5 ³ / ₄ "	9109-03



Toad Finger Splints

- Lightweight aluminum can be adjusted for proper fit; X-ray lucent
- Malleable arms fold to hold finger in desired alignment
- Packaged 12/case

Indications: Distal and proximal interphalangeal joint positioning

Size	Measurements	Product No.
Small	2 ¹ / ₄ "	9110-01
Medium	2 ³ / ₄ "	9110-02
Large	2 ⁷ / ₈ "	9110-03



Fold Over Finger Splints

- Lightweight aluminum can be adjusted for proper fit
- X-ray lucent
- Packaged 12/case

Indications: Protects injured nails and fingertips; Baseball or transverse fractures of the finger

Size	Measurements	Product No.
Small	1"	9111-01
Medium	2"	9111-02
Large	3"	9111-03
X-Large	3 ¹ / ₂ "	9111-07



Four Prong Finger Splints

- Lightweight aluminum can be adjusted for proper fit
- X-ray lucent
- Available with or without foam padding
- Packaged 12/case

Indications: Provides protection and immobilization to injured fingers

Size	Measurements	Product No.	
		With Foam	Without Foam
Small	2"	9112-01	9112-04
Medium	3"	9112-02	9112-05
Large	3 ¹ / ₂ "	9112-03	9112-06



Foam Padded Finger Strips

- Aluminum with foam padding
- Easy to form and can be cut to desired length with scissors
- X-ray lucent
- Packaged 12/case

Indications: Phalangeal fractures

Size	Products No.
1/2" x 9"	9115-01
1/2" x 18"	9115-02
3/4" x 18"	9115-03
1" x 18"	9115-04
1 1/2" x 18"	9115-05
2" x 18"	9115-06



Plastalume Splints

- Blue plastic attached to white foam padding
- 12 in a box

Size	Product No.
1 1/4" Length	9119-01
2 1/4" Length	9119-02
3 1/4" Length	9119-03
4 1/4" Length	9119-04
3 1/4" Length, circular base	9119-05
4 1/4" Length, circular base	9119-06
5 1/4" Length, circular base	9119-07
6 1/4" Length, circular base	9119-08
7 1/4" Length, circular base	9119-09



STAX Finger Splint

- Includes an assortment of 30 splints packaged in convenient storage container
- Protects distal phalanx and extends the DIP joint
- Open tip permits examination of fingertip and nailbed

Indications: Fingertip injuries; Nailbed injuries; Mallet finger

Size	Quantity	Product No.
Complete Kit	1 kit	9122-00
Size 1	5 each	9121-01
Size 2	5 each	9121-02
Size 3	5 each	9121-03
Size 4	5 each	9121-04
Size 5	5 each	9121-05
Size 5 1/2	2 each	9121-06
Size 6	2 each	9121-07
Size 7	1 each	9121-08




Burnham Thumb Splint

- Pre-molded aluminum provides proper support for thumb
- Foam padding for patient comfort
- X-ray lucent
- Universal size fits most patients

Indications: Immobilizes wrist to alleviate muscle-tendon tensions that tend to distort fracture segments

Size	Product No.	
	Left	Right
Universal	9131-00	9130-00



LMB Spring Finger Extension Assist

- High quality Wire-Foam™ construction for a customized fit
- Versatile model extends PIP joint with a slight extension effect on MP joint
- Designed for patient comfort with an emphasis on distal pad curvature
- Force: 15° flexion ≈ approx. 8 oz.; 90° flexion ≈ approx. 4 lbs.

Indications: Finger flexion tightness; Boutonniere deformity (allows IP flexion)

Size	Measurement	Product No.
X-Small	2 ¹ / ₄ "	501AA
Small	2 ⁵ / ₈ "	501A
Medium	3"	501B
Large	3 ¹ / ₂ "	501C
X-Large	4"	501D



Finger Splint

- Lightweight aluminum can be easily adjusted for proper fit
- Features slight angle at tip to address a variety of phalangeal injuries
- X-ray lucent
- Packaged 12/case

Indications: Phalangeal and metacarpal fractures and finger dislocations

Measurements	Product No.
1 ¹ / ₂ " x 7 ⁷ / ₈ "	11201
3" x 7 ⁷ / ₈ "	11202
5 ¹ / ₂ " x 7 ⁷ / ₈ "	11203
7" x 7 ⁷ / ₈ "	11204
4" x 7 ⁷ / ₈ "	11205



Short Metacarpal Splint

- Pre-molded aluminum design for anatomical fit and proper support
- Perforated design for air circulation
- X-ray lucent

Indications: Phalangeal; Metacarpal fractures



Size	Product No.	
	Left	Right
Child	9136-02	9136-01
Youth	9134-02	9134-01
Adult	9135-02	9135-01

Hely & Weber - TKO The Knuckle Orthosis

- Provides pistoning for fractures or injuries to metacarpals, phalanges, MCP and IP joints
- Distal buddy strapping eliminates rotation of the finger
- Padded stays are malleable
- Reversible - one splint acts as a Boxer's Fracture Splint (ulnar gutter), Radial Gutter splint, and Middle & Ringer Finger splint

Description	Product No.	
	Left	Right
The Knuckle Orthosis	3848L	3848R



Radial Gutter



Splinted Middle/Ring Fingers



Boxer's Fracture

Medline - Fold Over Finger Cot

- Malleable aluminum lined with blue polyester
- Sold 12 per case

Size	Measurements	Product No.
Medium	3 1/4"	ORT32400M
Large	3 3/4"	ORT32400L



Aircast (DJO) - Mayo Clinic Elbow Brace

- Unique patented design provides static adjustable stretch of the elbow in flexion and extension
- Replaces the need for two braces for elbow fractures
- Eliminates the need for serial casting for elbow reconstructions and acute fracture dislocations
- Arc™ Forearm Rotational Brace may be used for additional pronation and supination control

Indications: Fracture dislocation; Bicep/tricep tendon rupture; Ulnar nerve transposition; Total elbow arthroplasty; Ligament repair; Radial fixation; Osteoarthritis; Epicondylitis



Description	Product No.	
	Left	Right
Mayo Clinic Elbow Brace	05EL	05ER

AliMed - Turnbuckle Elbow Orthosis

- Helps maintain the gains in ROM achieved through stretching and brace wear
- Turnbuckle can be quickly and easily adjusted with the splint in place for optimal positioning
- Two-piece T-Foam padded liner conforms for maximum comfort and is removable for hand washing
- Splint adjusts from 110° flexion to full extension



Size	Length	Forearm	Biceps	Product No.
A	11 ¹ / ₄ "	8"-10"	9"-11"	52396/NA/A
B	14 ¹ / ₂ "	8"-10"	9"-11"	52396/NA/B
BB	14 ¹ / ₂ "	10"-12 ¹ / ₂ "	11"-13 ¹ / ₂ "	52396/NA/BB
C	16 ¹ / ₄ "	10"-12 ¹ / ₂ "	11"-13 ¹ / ₂ "	52396/NA/C
D	16 ³ / ₄ "	12 ¹ / ₂ "-15"	13 ¹ / ₂ "-16"	52396/NA/D

Measure circumference at belly of forearm and biceps

Bledsoe Brace - Telescoping Elbow

- Push-button telescoping and malleable arm plates for perfect fit
- Locks the elbow in any position from -10° to 110° in 10° increments
- Easy to apply, no tools required

Size	Product No.	
	Left	Right
Universal	AE023100	AE023200



Extender Arm Brace

- Designed to handle post-op elbow arthroscopy, severe elbow sprains and muscle and ligament trauma
- Hinge is quick setting and limits extension from 0° to 90°, flexion from 0° to 135°, and locks 0° to 90° in 15° increments
- Detachable, formable wrist-arm splint for increased stability of the forearm-wrist-hand
- One-piece removal and single-hand patient reapplication



Description	Product No.	
	Left	Right
Extender Arm Brace	AE021100	AE021200

Simple Arm Brace

- Formable aluminum plates adjust to fit virtually any arm circumference
- Lateral hinge and arm design prevents rubbing or abrasions against the body
- Upper arm measurement from medial humeral epicondyle to axilla; Sizing: Small = 5¹/₂"-8", Large = 7"-11¹/₂"
- Lower arm measurement from ulnar styloid process to medial humeral epicondyle; Small = 7"-9¹/₂", Large = 8¹/₂"-12¹/₂"



Size	Product No.	
	Left	Right
Small	AE020101B	AE020101B
Large	AE020103B	AE020103B
Hand Extensions		
Small	AE013101	AE013201
Large	AE013103	AE013203

Breg - TScope Elbow®

- Lightweight, low-profile product provides telescoping forearm and humeral bars, range-of-motion control, and lock out ability at any angle
- Detachable shoulder sling included
- ROM control in 10° increments: extension from 10° to 110°, flexion control from 10° to 120°
- Weight: 13 oz

Indications: Stable fractures of the elbow; Distal humerus proximal radius or ulna; Tendon and ligament injuries or repairs (Tommy John surgery, distal biceps tendon repair); Chronic elbow injuries; Collateral ligament reconstructions; Elbow hyperextension; Range-of-motion control post injury; Tennis elbow release and triceps tendon repair



Size	Product No.	
	Left	Right
Regular	07251	07250
Long	07253	07252

HEX Elbow Brace

- “X” anterior straps with strap-lock buckle
- Adjustable hinges to control desired ROM (0°-40° extension/45°-90° flexion)
- Half Airmesh® for breathability and half neoprene for compression
- Additional elbow foam padding to protect olecranon process

Indications: Chronic elbow injuries; Elbow hyperextension; Elbow tendonitis; Post-elbow dislocations



Size	Arm Measurements	Product No.
Small	9"-10 ¹ / ₂ "	14482
Medium	10 ¹ / ₂ "-12"	14483
Large	12"-13 ¹ / ₂ "	14484
X-Large	13 ¹ / ₂ "-15"	14485
2X-Large	15"-16 ¹ / ₂ "	14486
Circumference taken at mid-bicep		

DeRoyal - Pucci® Inflatable Elbow Orthosis

- Gives continuous low-stretch therapy through inflation of upright air stay bladder
- Includes fabric covers, laundry bag, and hand bulb inflator
- Fits left or right

Indications: CVA; TBI; SCI treatment of non-fixed elbow contractures




Size	Circumference		Product No.
	Biceps	Forearm	
Small	6"-8 ¹ / ₂ "	5"-7"	4220B
Large	8 ¹ / ₂ "-13"	6 ¹ / ₂ "-11 ¹ / ₂ "	4220C

Transition Elbow Brace

- Breathable foam softgoods are trimmable for a customized fit
- Hinge can be easily adjusted in 10° increments without tools (Flexion: 120° to 0° and Extension: 90° to 0°)
- Loop lock strapping is trimmable to allow for superior fit
- Length adjustment from 19" to 16"
- Universal left or right

Indications: Post-operative rehabilitation




Size	Measurement	Product No.
Small	6"-8"	EB5000-01
Medium	8"-10"	EB5000-02
Large	10"+	EB5000-03
Measure from axis of elbow to axis of wrist		

DeROM® Elbow

- Dynamic splint provides low load prolonged stretch to soft tissue to restore joint’s ROM
- Flexion and extension assist in a single unit
- Easy to set and release tension with integrated adjustment tool and tension release lever
- Range-of-motion stops for static progressive splinting

Indications: Joint stiffness, ROM limitations; Soft tissue contractures

Size	Biceps Circumference	Product No.	
		Left	Right
C	9"-12"	4229CL	4229CR
D	12"+	4229DL	4229DR
Replacement Softgoods			
C		4229CLK	4229CRK
D		4229DLK	4229DRK



DonJoy (DJO) - ROM Elbow Deluxe

- 30° incremental settings
- Bilateral hinges, universal application
- Removable hand piece
- Prevents supination/pronation control

Indications: Post-op; Protected range-of-motion; When neutral position is desired

Size	Measurements	Product No.
Small	9"-11 ¹ / ₂ "	11-0520-2-06060
Medium	11 ¹ / ₂ "-14 ¹ / ₂ "	11-0520-3-06060
Large	14 ¹ / ₂ "-19"	11-0520-4-06060

Measure arm circumference at mid biceps



IROM Elbow

- 10° incremental settings (0°-120°)
- Allows for controlled flexion/extension
- Integrated straps for full circumferential compression
- Allows for supination and pronation of forearm and wrist

Indications: Stable supracondylar fractures; Stable radial fractures; Lateral/medial epicondylitis

Size	Measurement	Product No.	
		Left	Right
Small	up to 8"	11-0181-2-13066	11-0180-2-13066
Large	8+"	11-0181-4-13066	11-0180-4-13066

Measure wrist to elbow crease



Lenjoy - Comfy Goniometer Elbow “Set ROM”

- Treats non-fixed contractures of the elbow
- Works in the same way as serial casting, allowing progressive extension by gradually adjusting and re-adjusting in 10° increments
- Allen wrench conveniently stored under the cover in the frame of the orthotic
- One size fits most adults; Available in Navy Blue



Size	Product No.
Pediatric Small	PGE-101-S
Pediatric Medium	PGE-101-M
Pediatric Large	PGE-101-L
Adult	GE-101
Adult Small	GE-101-AS

Comfy Spring Loaded Goniometer Elbow “Set ROM with Extension”

- Dynamic splint works like the Goniometer elbow with added feature of having a spring that gently pushes
- Spring strength is fixed
- Allen wrench conveniently stored in frame under the cover
- One size fits most adults; Available in Navy Blue



Size	Product No.
Pediatric Small	PSGE-101-S
Pediatric Medium	PSGE-101-M
Pediatric Large	PSGE-101-L
Adult	SGE-101

OCSI/Neuroflex - DynaPro™ Flex Elbow

- Elbow cuff extenders
- Bacti-Ban™ keeps covers cleaner and more odor resistant
- Abnormal reflex arc therapy provided by the device assists in neuro-rehabilitation
- Axial hinge with legible goniometer facilitates brace fitting



Size	Circumference		Product No.
	Upper Arm	Lower Arm	
Small	8"-10"	5"-8"	DPFE20724-S
Medium	9"-11 1/2"	7 1/2"-10"	DPFE20724-M
Large	11"-14"	9"-12"	DPFE20724-L

Ossur Americas - Innovator X™ Post-Op Elbow

- OneTouch button enables easy release for slide-to-size length adjustments
- Easy to use Innovator hinge makes flexion and extension adjustment simple
- Sensil® padding on the strapping system reduces slipping and
- Malleable aluminum cuffs enable the fitter to customize the Innovator X to varying anatomy

Indications: Post elbow surgery; Elbow reconstruction; Ligament and tendon repairs; Surgical repair of a fracture; Biceps repair; Mild to severe soft tissue injuries



Description	Product No.
Left	W-823101
Right	W-823100
Elbow Arm Bar Kit	W-A82310

ProCare (DJO) - ElbowRANGER® Motion Control Splint

- Aluminum R.O.M. bilateral hinges allow adjustments in 30° increments with formed Kydex® shells
- May be used following early cast removal to stabilize fractures, and for post-operative immobilization
- Removable hand grip assists in the prevention of pronation/supination

Indications: Collateral ligament strains and ruptures; Ulnar Nerve Transposition



Size	Bicep Circumference	Product No.
Small	8"-11½"	79-94200
Medium	11½"-14½"	79-94210
Large	14½"-19"	79-94220

Professional Products - Advanced Rehab Design (ARD) Elbow Brace

- Black polyethylene bands with lightweight aircraft aluminum uprights and hypoallergenic liners
- Tool-less adjustable hinge
- 10° incremental adjustments from 10° to 90° in extension
- 10° incremental adjustments from 30° to 120° in flexion



Description	Product No.
ARD Post-Op Lateral Hinge, Left	1292-L
ARD Post-Op Lateral Hinge, Right	1292-R
ARD Post-Op Bilateral Hinge, Universal	1292-DH
Size	Bicep Circum.
X-Small	7-9"
Small-Medium	9"-11"
Large-X-Large	11"-14"

RCAI - Flex Cuff Respond® ROM Elbow Orthosis

- Lightweight and versatile, the Flex Cuff Elbow Orthosis features the Respond ROM range-of-motion hinge with flexion/extension stop settings, and static lock pem for immobilization or static progressive stretch
- Hinge arms positioned in media/lateral panels to prevent migration on the arm
- Length: 15"; Hinge Length: 12 1/2"

Size	Circumference	Product No.
Small	8"-10"	18ECO-S
Medium	10"-12"	18ECO-M
Large	12"-16"	18ECO-L
X-Large	16"-21"	18ECO-XL
Accessories		
Elbow Cap		18EC



Premier Elbow Corrective Orthosis

- Offers progressive static lock settings to permit gradual increases in elbow joint extension and positioning following surgical procedures
- "T" screws allow for degree of angle adjustment without the use of tools
- Thermoplastic humeral and forearm cuffs with lightweight nylon hinges can accommodate patient with 0° to moderate tone and mild spasticity
- Liner is constructed of one-piece breathable foam that wicks away moisture from the skin

Size	Mid Humerus Circum.	Product Length	Product No.
Small	9"-11"	12 1/2"	28ECO-S
Medium	11"-13"	13"	28ECO-M
Large	13"-15"	13"	28ECO-L




Ratchet POP Arm Brace

- Offers flexion/extension ROM control from 0-120° in 5° increments
- Static lock PEM® allows the brace to be used static progressively
- While the brace is on the arm, the angle of the limb can be increased in 10° increments from 0-120° without the use of settings or tools

Indications: Post surgery; Strains; Sprains; Instability; Ligament repairs; Soft tissue contractures, due to cerebral palsy, stroke, head trauma or spinal cord injury

Size	Mid Humerus Circum.	Product Length	Product No.	
			Left	Right
Small/Medium	11"-14"	15"	84RAPOP-AV-S/M-L	84RAPOP-AV-S/M-R
Large/X-Large	12"-16"	15"	84RAPOP-AV-LG/XL-L	84RAPOP-AV-LG/XL-R



Respond ROM® Elbow

- Provides 10° increment settings of flexion/extension elbow control
- The Respond ROM® hinge features a static lock pem that allows the brace to be used static progressively for positioning
- Elbow cap provides 3-point pressure leverage and elbow control
- Available in Royal Blue

Indications: Post surgical ROM limitation; Joint stiffness; Contracture from 0-120° in 5° increments



Size	Mid Humerus Circum.	Product Length	Product No.	
			Orthowick	Kodel
Small	9"-11"	14"	38ECO-S	34ECO-S
Medium	11"-13"	15"	38ECO-M	34ECO-M
Large	13"-15"	15"	38ECO-L	34ECO-L

Universal Arm Brace

- Bilateral hinges provide increased flexion and extension elbow control 0°-120°
- Precision length adjustment hinge arm assures a personalized fit
- Forearm design provides protection of the olecranon process






Indications: Elbow surgery; Tendon or muscle repair; Strains; Sprains; Stable humeral or forearm fractures





Size	Measure Arm Length	Max. Circum. Bicep	Max. Circum. Forearm	Product No.
Short	11"-13"	17 ¹ / ₂ "	14 ¹ / ₂ "	88UAB-SX
Average	14"-16"	19 ¹ / ₂ "	16 ¹ / ₂ "	88UAB-AV
Long	17"-20"	21 ¹ / ₂ "	18 ¹ / ₂ "	88UAB-LX

Tennis Elbow Splints

- Indications: Tennis elbow; Golfer's elbow (medial or lateral epicondylitis)

Aircast (DJO)	Pneumatic Arm Band			
	<ul style="list-style-type: none"> • A single pre-inflated air cell concentrates compression on the extensor muscle, not around the arm, providing more support, less constriction 			
	Size	Forearm Circumference	Product No. Beige Black	
Universal	8"-14"	05A 05A-B		
Bauerfeind	EpiPoint Elbow Strap			
	<ul style="list-style-type: none"> • Viscoelastic pad which can be rotated to be worn on the right or left forearm, exerting targeted compression • Elastic band with red warning section indicates whether compression is too tight • Titanium, latex free 			
	Size	Forearm Circumference	Product No.	
Universal	8 ³ / ₄ "-13"	12063501080000		
CHO-PAT	Premium Tennis Elbow Strap			
	<ul style="list-style-type: none"> • Reinforced, padded strap that molds to forearm shape • Hook and loop closure • Measure largest aspect of forearm circumference 			
	Size	Forearm Circumference	Product No.	
	Small	9"-10 ¹ / ₂ "	CP23	
	Medium	10 ¹ / ₂ "-12"	CP25	
Large	12"-13 ¹ / ₂ "	CP27		
X-Large	13 ¹ / ₂ "-15"	CP29		
DeRoyal	Tennis Elbow			
	<ul style="list-style-type: none"> • Support is low-profile increasing patient compliance • Inner pad can be strategically placed for maximum pain relief • D-ring closure for easy application 			
	Size	Length	Product No.	
Universal	2" x 16"	12870000		
DonJoy (DJO)	Universal Wrist/Forearm Splint			
	<ul style="list-style-type: none"> • Breathable nylon/fiber construction • Malleable dorsal and volar stays 			
	Size	Product No. Left Right		
	Universal	11-0221-9	11-0220-9	
	Universal Forearm Support/Forearm Strap			
<ul style="list-style-type: none"> • Provides compression, warmth and support • 1/8" neoprene 				
Size	Product No.			
Universal	11-0467-9			

Freeman	Neoprene Tennis Elbow Sleeve			
	<ul style="list-style-type: none"> • Pull-on design with adjustable loop-and-lock compressive strap • Support is perforated on medial side for ventilation • Sleeve should be work during athletic activity only 			
	Size	Forearm Circumference		Product No.
	X-Small	8"-9"		613-XS
	Small	9"-10"		613-S
	Medium	10"-11"		613-M
Hely & Weber	Tennis Elbow with Pressure Pad			
	<ul style="list-style-type: none"> • Low-profile design • 1/2" (3 mm) neoprene with 2-sided nylon outer covering • Pressure pad and elastic strap provide compression 			
	Size	Forearm Circumference		Product No.
	X-Small	8"-9"		3838XS
	Small	9"-10"		3838S
	Medium	10"-11"		3838M
	Large	11"-12"		3838L
X-Large	12"-14"	3838XL		
2X-Large	14"-16"	3838XXL		

AliMed - CVA Sling

- Has two large cuffs for hand and forearm support
- Strap is padded at the neck for maximum comfort
- Once adjusted, it is easy to apply or remove without assistance
- Adjusts to fit left or right arm

Size	Product No.
Small/Medium	5971
Medium/Large	5970



Harris Hemi-Arm Sling

- Weight of arm is broadly distributed across both shoulders
- Padded back panel prevents the sling from slipping
- Easy one-handed application
- Sizing: adjusts to any size user, one size fits all

Indications: Patients with hemiplegia

Description	Product No.	
	Black	White
Left	5877	5879
Right	5876	5878



GivMohr Sling

- Promotes functional positioning of upper extremity during standing and ambulation and naturally reduces subluxation
- Helps protect involved upper extremity from injury
- Available in Black

Size	Weight Range (lbs)	Height Range	Measurements	Product No.
X-Small	80-100	up to 5'0"	—	GMS05
Small	100-140	5'0"-5'4"	13"-14"	GMS10
Medium	135-185	5'3"-5'10"	15"-16"	GMS20
Large	180-235	5'10"-6'4"	17"-18"	GMS30
X-Large	235+	6'4"+	18"+	GMS40

Measure distance between right and left acromion of scapula



Corflex - Collar Cuff & Sling

- Tubular neck pad and flannel lined cuff for optimal patient comfort
- 1' durable strap with contact closure and wrist cuff
- Fits 6"-10" wrist circumference

Indications: Injuries to the arm, shoulder, wrist or hand

Product No. **21-2052**



Tricot Arm Sling

- Thumb loop to prevent wrist drop and hook and loop closure at top of sling envelope to provide a secure fit
- Soft yet durable tricot material

Indications: Injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position

Size	Measurements	Product No.
X-Small	5" x 12"	21-5060
Small	6" x 15"	21-5061
Medium	7" x 17"	21-5062
Large	8" x 19"	21-5063
X-Large	9" x 21"	21-5064



DeRoyal - Light Blue Arm Sling

- Light blue cotton material with deep pocket design
- Foam shoulder pad provides additional comfort
- Thumb loop helps prevent wrist drop
- Hook and loop closure

Indications: Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Strap	Strap Length	Product No.
Infant	5" x 9"	1" x 21"	10101
Child	6" x 11"	1" x 27"	10102
Small	7" x 13"	2" x 34"	10103
Medium	8" x 17"	2" x 37"	10104
Large	9" x 20"	2" x 37"	10105
X-Large	9" x 23"	2" x 42"	10106



Specialty Arm Sling

- Durable navy canvas construction with deep pocket design
- Thumb loop helps prevent wrist drop
- Adjustable slide buckle on strap for best fit
- Hook and loop closure with O-rings

Indications: Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Strap	Strap Length	Product No.	
			Foam Pad	Vinyl Pad
Infant	5" x 9"	1" x 21"	8004-01	8004-11
Child	6" x 11"	1" x 27"	8004-02	8004-12
Small	7" x 13"	2" x 34"	8004-03	8004-13
Medium	8" x 17"	2" x 37"	8004-04	8004-14
Large	9" x 20"	2" x 37"	8004-05	8004-15
X-Large	9" x 23"	2" x 42"	8004-06	8004-16



Arm Sling with Foam Strap

- Canvas or lightweight tietex material
- Deep pocket, closed pouch design
- Foam shoulder strap provides maximum comfort
- Hook closure fastens anywhere on foam strap for superior adjustment

Indications: Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains or anterior glenohumeral

Size	Pouch Size	Strap Length	Canvas	Tietex
Universal	8" x 19"	2" x 37"	—	TX8020-00
Small	7" x 13"	2" x 37"	8023-01	TX8023-21
Medium	8" x 17"	2" x 37"	8023-02	TX8023-22
Large	9" x 20"	2" x 37"	8023-03	TX8023-23
X-Large	9" x 22"	2" x 40"	8023-04	TX8023-24
2X-Large	10" x 23"	2" x 40"	8023-05	—



Arm Elevator/Sling with Ice Pockets

- Heavy duty canvas construction
- Includes two ice bags for cold therapy
- Converts easily to sling

Indications: Post upper extremity surgery and upper extremity edema

Size	Pouch Strap	Product No.
Small	10" x 14"	8011-01
Medium	10" x 16"	8011-02
Large	12" x 17"	8011-03



Universal Arm Sling

- Heavy-duty canvas material
- Easy fold-back universal design
- Foam shoulder pad provides additional comfort
- Hook and loop closure with O-rings

Indications: Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Size	Strap Length	Product No.
Child	5" x 9"	1" x 31"	8020-01
Adult	8" x 19"	2" x 37"	8020-00



Convuluted Foam Arm Elevator/Sling

- Convuluted foam liner with hook receivable outer lining
- Converts easily to sling with shoulder pad included
- Slide buckle for adjustment of shoulder strap
- Hook and loop closures for easy application
- Pouch size: 10" x 16"

Indications: Post upper extremity surgery; Upper extremity edema

Universal Fit

8788-03



CVA Sling

- Durable canvas material
- Tubular foam neck pad for patient comfort
- Strap features adjustable slide buckle closure

Indications: Stroke or paralysis affecting arm or shoulder

Size	Pouch Size	Strap Length	Product No.
Small/Medium	5 ¹ / ₂ " x 7 ¹ / ₂ "	2" x 30"	6157D02
Medium/Large	5 ¹ / ₂ " x 7 ¹ / ₂ "	2" x 38 ¹ / ₂ "	6157D03



Hely & Weber - Arm Sling

- Padded strap offers customized fitting and adjustment
- Deep pocket design accommodates the largest of casts and dressings
- Open end with thumb loop adds comfort and eliminates slippage

Size	Height x Length	Product No.
X-Small	7 ¹ / ₂ " x 11"	500XS
Small	7 ¹ / ₂ " x 13 ¹ / ₂ "	500S
Medium	7 ¹ / ₂ " x 15 ¹ / ₂ "	500M
Large	7 ¹ / ₂ " x 17 ¹ / ₂ "	500L
X-Large	7 ¹ / ₂ " x 20"	500XL



ProCare (DJO) - Deluxe Arm Sling with Pad

- Generous size envelope to support the arm
- Cotton/poly with fully padded web strap and contact closure
- Thumb loop helps prevent sling from migrating

Indications: Cast support; Injuries to the arm, wrist or hand

Size	Depth x Length	Product No.
Pediatric/2X-Small	6 ¹ / ₂ " x 9"	79-84001
Child/X-Small	7" x 11 ¹ / ₄ "	79-84002
Small	7 ¹ / ₂ " x 13 ¹ / ₂ "	79-84003
Medium	8" x 16"	79-84005
Large	8 ¹ / ₂ " x 18"	79-84007
X-Large	9" x 20 ¹ / ₄ "	79-84008



Clinic Arm Sling

- Cotton/poly envelope with web shoulder strap and slide buckle closure
- 6 per box

Indications: Cast support; Injuries to the arm, wrist or hand

Size	Depth x Length	Product No.
2X-Small	6" x 9"	79-84021
X-Small	6" x 11"	79-84022
Small	6" x 13"	79-84023
Medium	6" x 15 ¹ / ₂ "	79-84025
Large	6" x 18"	79-84027
X-Large	6" x 20"	79-84028



Universal Arm Sling

- Contact closure adjustment on envelope accommodates most patients
- Durable cotton/poly material, double D-ring and contact closure web shoulder strap for proper adjustment
- Length x Depth: 19" x 7"

Indications: Immobilization and support of the shoulder and elbow joints

Description	Product No.
With Foam Shoulder Pad	79-92070
Without Foam Shoulder Pad	79-84300



AliMed - Cubital Tunnel Splint

- Rigid, high-density, beige polyethylene shell with plush 3/8" thick terry cloth lining
- Immobilizes and positions the elbow at 122° so that the ulnar nerve is less prominent
- Restricts full extension and flexion to minimize stretching of ulnar nerve, and helps reduce pain and irritation due to contracture


Indications: Ulnar nerve transposition; Post-cast removal



Size	Fits	Product No.
Small/Medium	most women and smaller men	510503
Medium/Large	most men and larger women	510504

Cubital Tunnel Syndrome Support

- Helps prevent full extension and flexion at elbow to limit pain associated with cubital tunnel syndrome
- Elbow cutout helps position splint
- Soft beige tricot laminated to polyethylene foam
- Firm bilateral foam inserts



Size	Forearm	Biceps	Product No.
Small	9" and under	11" and under	510269/NA/SM
Medium	9"-10 ¹ / ₂ "	11 ¹ / ₂ "-13"	510269/NA/MD
Large	11"-12 ¹ / ₂ "	13 ¹ / ₂ "-15"	510269/NA/LG
X-Large	13"-15 ¹ / ₂ "	15 ¹ / ₂ "-17 ¹ / ₂ "	510269/NA/XL

Measure around proximal forearm and biceps

Pil-O-Splint® Elbow Support

- Limits excessive elbow flexion
- Removable interior plastic stay slides into hidden pocket to maintain full extension
- Padded cotton splint with two adjustable compression straps



Size	Forearm Circum.	Product No.
Standard	9 ¹ / ₂ "-11"	51329
X-Large	11"-15"	51996

Measure forearm circumference at midpoint

Bauerfeind - EpiTrain® Elbow Support

- Comfortable, lightweight, knitted support incorporates two silicone inserts providing intermittent compression to reduce pain, and swelling
- Indications:** Epicondylitis; Elbow sprains; Pre- and Post-op conditions; Tennis Elbow



Size	Circum.	Nature	Product No.	
			Black	Titanium
0	6 ¹ / ₂ "-7 ¹ / ₂ "	11061603010000	11061603070000	11061603080000
1	7 ¹ / ₂ "-8 ¹ / ₄ "	11061603010001	11061603070001	11061603080001
2	8 ¹ / ₄ "-9"	11061603010002	11061603070002	11061603080002
3	9"-9 ³ / ₄ "	11061603010003	11061603070003	11061603080003
4	9 ³ / ₄ "-10 ¹ / ₂ "	11061603010004	11061603070004	11061603080004
5	10 ¹ / ₂ "-11 ¹ / ₂ "	11061603010005	11061603070005	11061603080005
6	11 ¹ / ₂ "-12 ³ / ₄ "	11061603010006	11061603070006	11061603080006

Measurement taken 4³/₄" below epicondyle

Corflex - XR Cubital Tunnel Splint

- Innovative support that immobilizes and positions the elbow at 122° to help reduce pain and numbness associated with Cubital Tunnel Syndrome
- Extra rigid A.B.S. material with medial contouring on shell and plush cloth liner

Indications: Cubital Tunnel Syndrome; Post-op application; Post-cast removal



Size	Fits	Product No.
Small/Medium	most female and smaller framed males	37-2191
Medium/Large	most males and larger frame females	37-2192
Large/X-Large	large framed males & XL framed females	37-2193

Ventral Cubital Tunnel Splint

- Unique anterior design immobilizes and positions elbow at 122° while allowing easy access to dressing
- High density polyethylene shell with plush cloth liner

Indications: Cubital Tunnel Syndrome

Size	Fits	Product No.
Small/Medium	most females & smaller framed males	37-2291
Medium/Large	most males & larger framed females	37-2292
Large/X-Large	large framed males & XL framed females	37-2293



DeRoyal - Wrist and Elbow Splint

- Durable leatherette material
- Provides immobilization of elbow at 90°
- Aluminum frame provides protection to the elbow

Indications: Olecranon fractures; Radial head fractures; Hairline fractures of the radius, humerus and ulna; Ligamentous reconstruction procedures of the elbow

Size	Circumference		Product No.
	Biceps	Forearm	
Small	11 ¹ / ₂ "	9 ³ / ₄ "	6001-01
Medium	12 ¹ / ₂ "	10"	6001-02
Large	13 ¹ / ₂ "	11"	6001-03



Elbow Immobilizers

- Soft, perforated foam construction for patient comfort and breathability
- Metal T-bar stays for optimal immobilization and strength
- Prevents flexion of arm to prevent pulling NG tubes, tracheostomy tubes or I.V. lines
- Sold in pairs

Indications: Patients at risk of unplanned extubations; Patients who disrupt medical treatment

Size	Length x Width	Product No.
Infant	7 ¹ / ₂ " x 12 ¹ / ₂ "	M7035-S
Pediatric	7 ¹ / ₂ " x 17 ¹ / ₂ "	M7035-L
Adult	14" x 17"	M2035



Elbow Splint

- Pre-molded aluminum design for an anatomical fit and proper support
- Perforated design for air circulation; foam padded for patient comfort
- Universal fit for left and right arm
- X-ray lucent
- Aluminum bar can be bent to adjust angle

Indications: Elbow immobilization

Size	Product No.
Child	9118-01
Youth	9118-02
Adult	9118-03



DonJoy (DJO) - Elbow Guard

- Wraparound biceps/triceps portion with elastic strap for ease of application
- Bilateral polycentric hinges for flexion/extension and media/lateral control
- "X" strap for hyperextension control
- Elbow pad for protection in contact sports

Indications: Medial/lateral ligament instabilities

Size	Measurements	Product No.
Small	10"-11 ¹ / ₂ "	11-1003-2-06000
Medium	11 ³ / ₅ "-13"	11-1003-3-06000
Large	13 ¹ / ₁₀ "-14 ¹ / ₂ "	11-1003-4-06000
X-Large	14 ³ / ₅ "-16"	11-1003-5-06000
2X-Large	16 ¹ / ₁₀ "-17 ¹ / ₂ "	11-1003-6-06000
3X-Large	17 ³ / ₅ "-19"	11-1003-7-06000

Measurement taken is circumference of arm at mid-biceps



Spider® Pad Elbow

- Honeycomb design contours for comfortable movement while providing cushioning protection
- High impact-resistant foam; Low-profile and streamlined for protection during contact sports and job-related activities
- Available in Black

Size	Measurements	Product No.
X-Small	7"-8"	11-0392-1-06060
Small	8"-9"	11-0392-2-06060
Medium	9"-10"	11-0392-3-06060
Large	10"-11"	11-0392-4-06060
X-Large	11"-12"	11-0392-5-06060
2X-Large	12"-13"	11-0392-6-06060
3X-Large	13"-14"	11-0392-7-06060

Measurement taken at circumference of elbow joint



Hely & Weber - Hayes™ Padded Elbow Orthosis Wrap-Around

- Universal neoprene construction offers unlimited adjustments
- Memory Foam padding contours and cushions elbow while providing warmth and protection
- Loop & Lock Straps ease brace application while limiting inventory needs in a one-size-fits-all style

Indications: Tendonitis; Bursitis; Cubital Tunnel Syndrome; Arthritis

Product No. **3832**



Cubital Comfort Brace

- Two aluminum malleable stays positioned inside the elbow prevent excessive elbow flexion
- Aluminum stays are pre-bent at 20 degrees of full flexion
- Two stockinette arm bands aid in arm positioning

Indications: Treatment of Cubital Tunnel Syndrome



Size	Length*	Product No.
Short	up to 11 ¹ / ₂ "	3847S
Long	11 ¹ / ₂ "+	3847L

*Measurement taken from mid-forearm to mid-humerus

Lenjoy - Comfyrene Elbow Orthosis "Bend to Fit"

- Easy to use lower elbow splint
- Straps are held on using baby hook and loop that adhere at any location
- One size fits most adults
- Available in Dark Blue

Product No. **E-101-CP-ADULT-DB**



Comfy Combination Elbow Hand Orthosis

- Extends the elbow and hand simultaneously
- Includes hand roll attachment or full hand attachment, which can be rotated to accommodate pronation and supination
- Comes with one cover
- Available in Navy Blue; One size fits most adults

Indications: Spastic hand and elbow taking into account the flexion synergy of the upper extremity



Size	Elbow Crease to Wrist	Elbow Crease to Forearm	Elbow Crease to Shoulder	Product No.	
				Hand Roll	Full Hand
Pediatric Small	5"-7"	3 ¹ / ₄ "	3 ¹ / ₄ "	EH-101-HR-PS	EH-101-FH-PS
Pediatric Medium	7"-9"	3 ³ / ₄ "	3 ³ / ₄ "	EH-101-HR-PM	EH-101-FH-PM
Pediatric Large	9"-11"	4 ¹ / ₂ "	4 ¹ / ₂ "	EH-101-HR-PL	EH-101-FH-PL
Adult Small	11"-13"	5"	5"	EH-101-HR-AS	EH-101-FH-AS
Adult	13"-15"	5 ¹ / ₄ "	5 ¹ / ₄ "	EH-101-HR	EH-101-FH

Comfy Elbow “Bend to Fit”

- Cuffs are bendable to allow for optimal fit
- Use as a resting splint or for the gradual extension of non-fixed elbow contractures
- One size fits most adults
- Available in Navy Blue



Size	Elbow Crease To Wrist	Elbow Crease To Shoulder	Total Length	Circum. of Biceps	Product No.
Pediatric Small	3 ¹ / ₄ "	3 ¹ / ₄ "	6 ¹ / ₂ "	4"-6"	PE-101-S
Pediatric Medium	3 ³ / ₄ "	3 ³ / ₄ "	6 ¹ / ₂ "	6"-8"	PE-101-M
Pediatric Large	4 ¹ / ₂ "	4 ¹ / ₂ "	9"	8"-10"	PE-101-L
Adult Small	5"	5"	10"	9"-12"	E-101-AS
Adult	5 ¹ / ₄ "	5 ¹ / ₄ "	10 ¹ / ₂ "	11"-15"	E-101

Medi USA - Epicomed Elbow Support

- 3-dimensional knitting technique ensures a perfect fit
- Anatomically knit support with accu-pressure silicone inserts
- Comfort-zone knit at the elbow bend prevents constriction
- Clima-comfort fabric is highly elastic, breathable and moisture wicking
- Available in Sand



Size	Arm Circumference	Product No.
Small	8"-9"	82272
Medium	9"-10 ¹ / ₂ "	82273
Large	10 ¹ / ₂ "-11 ³ / ₄ "	82274
X-Large	11 ³ / ₄ "-12 ¹ / ₂ "	82275
2X-Large	12 ¹ / ₂ "-13 ³ / ₄ "	82276
Measure circumference at 2" below center of elbow joint		

New Options Sports - Hybrid™ Hyperextension Polycentric Hinged Elbow with “E-Z-E” Bicep Hook and Loop Closures

- 1/8" nylon 2-sides neoprene
- Covered, stainless steel hinges w/ hyperextension stop
- 2 loop and lock criss-cross straps to help prevent hyperextension
- Bicep "E-Z-E" Bicep hook and loop closures



Arm and Elbow

Size	Forearm Circumference	Product No.
X-Small	8"-9"	E15 XS
Small	9"-10"	E15 S
Medium	10"-11"	E15 M
Large	11"-12"	E15 L
X-Large	12"-14"	E15 XL
2X-Large	14"-16"	E15 XXL
3X-Large	16"-18"	E15 XXXL

ProCare (DJO) - Heel and Elbow Protectors

- Breathable interfacing with fiberfill lining
- Three layer comfort cushion helps to eliminate pressure points
- Hook and loop closure for ease of application and removal
- Poly/pile fiber with contact closure absorbs moisture and allows proper circulation of air



Description	Product No.
Bulk 18 pairs	79-81001
Individual 18 pairs	79-81002
Small	79-81003
Medium	79-81005
Large	79-81007

Heel/Elbow Protectors

- Choice of polyester pile or vented foam with contact closure
- One size fits all
- Sold in pairs



Size	Style	Product No.
Universal	Poly/Pile	79-81040
Universal	Foam	79-81050
Universal	Foam, 3 Pair	79-99740

RCAI - Elbow Stabilizer

- Removable plastic stay can be heat molded and repositioned to prevent full extension and flexion, as needed, to minimize stretching of the ulnar nerve
- Made of a breathable and anti-microbial “wicking” foam liner, which keeps the skin cool and dry, maintaining good skin integrity
- Adjustable compression straps assure proper individual support and elbow control

Indications: Epicondylitis; Cubital Tunnel Syndrome; Bursitis; Tendonitis



Size	Mid Humeral Circumference	Product Length	Product No.
Youth	6 ¹ / ₂ "-12"	9"	48ES-Y
Adult	9 ¹ / ₂ "-18"	12"	48ES-A
Adult Long	9 ¹ / ₂ "-18"	15"	48ES-A-LX

Elbow Sleeve

- Provides compression and therapeutic warmth to the elbow to help reduce inflammation and relieve stiffness and pain
- Adjustment tabs at the proximal and distal ends allow for a personalized fit
- Cutout at the elbow discourages migration on the limb
- Loop and lock compression strap assured added support

Indications: Tennis and Golfer's Elbow



Size	Forearm Circum.	Product No.
Small	8"-9"	630-NES-S
Medium	9"-10"	630-NES-M
Large	10"-12"	630-NES-L

Silipos - Elbow Slip Over

- Aids healing dry, cracked elbows
- Protective padding for Ulnar nerve or as an orthotic liner
- Mineral oil gel cushion contours to protect against pressure, friction and shear forces
- 5 mm thick gel pad

Size	Product No.
Small/Medium	15225
Large	15235
X-Large	15245



CoreLINE - Shoulder Sling with Pillow

- Breathable sling constructed of mesh material with quick-wicking liner for added patient comfort and support
- Body and shoulder straps provide easy, quick adjustment and application
- Quick release buckles allow easy donning/doffing

Indications: Post-operative immobilization; Capsular Shifts; Soft tissue strains/repairs

Size	Forearm Length	15° Abduction Pillow	Neutral Pillow
Small	11 ¹ / ₂ "-13"	CL0342	CL0352
Medium	13 ¹ / ₂ "-14 ¹ / ₂ "	CL0344	CL0354
Large	15"-16"	CL0346	CL0356
X-Large	16 ¹ / ₂ "-17 ¹ / ₂ "	CL0348	CL0358



Becker Orthopedic- Lerman Shoulder Holster Orthosis (SHO)

- Fully adjustable in all planes and is designed to enhance functional recovery, while effectively controlling shoulder movement
- May be used postoperatively or in conditions where the immobilization of the shoulder is the desired goal
- One size fits all - Adult

Color	Product No.	
	Left	Right
Black	185-L-B	185-R-B
White	185-L-W	185-R-W



Bledsoe Brace - Original ARC with Pillow

- One brace fits all size patients and is universal, right and left
- Traditional pillow design
- Perforated cool version available

Indications: Non-surgical or post-operative support and immobilization of the injured shoulder where glenohumeral positioning is key

Description	Product No.
Standard	AE050300
Cool	AE055300



ARC 2.0 Shoulder Brace

- Universal sling design fits all size patients
- Pistol grips adjusts with quick-pull tabs and keeps the arm from migrating forward out of the sling
- Optional underarm strap relieves pressure on the neck for larger patients and patients positioned in external rotation
- Aluminum support can be positioned in a variety of ways to accommodate numerous post-op requirements

Indications: Non-surgical or post-operative support and immobilization of the injured shoulder where glenohumeral positioning is key

Product No. **AE050400**



Breg - Slingshot® 3

- Comprised entirely of Airmesh for enhanced breathability
- Includes exercise ball
- 15° abduction, 90° neutral and 45° neutral
- Universal left or right

Indications: Rotator cuff repairs; Bankart lesions; SLAP lesions; Glenohumeral dislocations/subluxation; Posterior/Anterior capsule repairs; Total shoulder reconstructions; Global shoulder instability; Soft tissue repairs/strains

Size	Forearm Length	Product No.
Small	11 ¹ / ₂ "-13"	00042
Medium	13 ¹ / ₂ "-14 ¹ / ₂ "	00043
Large	15"-16"	00044
X-Large	16 ¹ / ₂ "-17 ¹ / ₂ "	00045



SlingShot® 2

- Features a comfortable, breathable Airmesh® sling and a 15° abduction pillow
- Quick release shoulder and waist strap buckles make this product easy to apply
- Includes exercise ball to stimulate circulation and a thumb rest to minimize migration

Indications: Rotator cuff repairs; Anterior repairs; Posterior dislocations; Capsular shifts; Global shoulder instability

Size	Forearm Length	Product No.
Small	11 ¹ / ₂ "-13"	08502
Medium	13 ¹ / ₂ "-14 ¹ / ₂ "	08503
Large	15"-16"	08504
X-Large	16 ¹ / ₄ "-17 ¹ / ₂ "	08505



Corflex - Shoulder Abduction Pillow with Sling

- Easy to apply design features pillow with malleable stay for proper positioning and padded neck strap
- Firm density foam covered in tricot material with tricot sling and exercise ball to stimulate circulation
- Pillow allows for 15° of abduction
- Firm Fit version (pictured) features a moisture wicking terry cloth

Indications: Treatment after Bankart repairs; Rotator cuff repairs; Shoulder arthroscopy; Separations and dislocations

Size	Hip		Product No.	
	Circumference	Forearm	Standard	Firm Fit
Small	24"-30"	9"-12"	23-1911	23-1931
Medium	30"-36"	12"-15"	23-1912	23-1932
Large	36"-42"	15"-18"	23-1913	23-1933
X-Large	42"-48"	18"-20"	23-1914	23-1934

Measurement of length of forearm taken wrist to elbow
NOTE: Sling may be purchased separately



Firm Fit

Neutral Shoulder Pillow with Sling

- Features contoured pillow that maintains shoulder in a natural state of external rotation
- Firm density foam covered in tricot material with 3D spacer fabric sling and exercise ball to stimulate circulation

Indications: Treatment after Bankart repairs; Rotator cuff repairs; Shoulder arthroscopy; Separations and dislocations

Size	Hip Circumference	Forearm	Product No.
Small	24"-30"	9"-12"	23-1921
Medium	30"-36"	12"-15"	23-1922
Large	36"-42"	15"-18"	23-1923
X-Large	42"-48"	18"-20"	23-1924

Measurement of length of forearm taken wrist to elbow



Ultra Shoulder Abduction Pillow with Sling

- Features extended pillow design with malleable stay for proper positioning, V-lock strapping system to prevent roll/shifting and padded neck strap for optimal patient comfort
- Firm density foam covered in tricot material with tricot sling and exercise ball to stimulate circulation

Indications: Treatment after Bankart repairs; Rotator cuff repairs; Shoulder arthroscopy; Separations and dislocations

Size	Hip Circumference	Forearm	Product No.
Small	24"-30"	9"-12"	23-1901
Medium	30"-36"	12"-15"	23-1902
Large	36"-42"	15"-18"	23-1903
X-Large	42"-48"	18"-20"	23-1904

Measurement of length of forearm taken wrist to elbow



Ultra ER Shoulder Abduction Pillow with Sling

- Removable wedge that provides 30° of external rotation, V-lock strapping system to prevent roll/shifting and padded neck strap for optimal patient comfort
- Firm density foam covered in tricot material with tricot sling and exercise ball to stimulate circulation

Indications: Protocols involving non-surgical immobilization in external rotation to help reduce the incidence of recurrent dislocation

Size	Hip	Forearm	Product No.
Small	24"-30"	9"-12"	23-1941
Medium	30"-36"	12"-15"	23-1942
Large	36"-42"	15"-18"	23-1943
X-Large	42"-48"	18"-20"	23-1944
Measurement taken from hip circumference and length of forearm from wrist to elbow			



DeRoyal - Shoulder P.A.D II

- Foam block provides 35° of abduction
- Quick release buckle on waist strap for ease of application and removal
- Exercise ball packaged with product
- Universal left or right

Indications: Surgical rotator cuff repairs; Anterior glenohumeral capsular reconstruction; Dislocation/subluxation of G-H joint; S/P glenohumeral acromioplasty; Bankart repairs; Soft tissue repairs; Shoulder sprains and strains

Size	Pouch Size	Product No.
Small	7" x 16"	11630005
Medium	7" x 17"	11630006
Large	7" x 18"	11630007
X-Large	7" x 19"	11630008



Shoulder Abduction Positioner

- Foam block provides 55° of abduction
- Adjustable wrist and humeral cuffs for proper positioning
- Foam shoulder pad for additional patient comfort
- Universal left or right

Indications: Acute shoulder instability; Surgical rotator cuff repairs; Anterior glenohumeral capsular reconstruction; Acute rotator cuff tendinopathy; Dislocation/subluxation of G-H joint; Acromioclavicular joint dislocations; S/P glenohumeral acromioplasty; Bankart repairs; Soft tissue repairs; Shoulder sprains and strains

Size	Dimensions	Product No.
Small	25 ¹ / ₄ " x 23 ³ / ₄ "	11650005
Medium	29 ¹ / ₂ " x 26 ³ / ₄ "	11650006
Large	30 ¹ / ₂ " x 36"	11650007



DonJoy (DJO) - UltraSling III

- Moisture wicking mesh fabric for greater comfort, particularly in long term use (15° of abduction with internal rotation)
- Promotes auxiliary air exchange to reduce risk of secondary infections and prevents post-operative internal rotation contractures
- AB Version: Large comfortable pillow gives greater abduction where required (45° or 60°)
- The De-Rotation Strap hooks to the UltraSling® III and prevents rotation by securely holding the patient's arm in neutral position
- Fits left or right arm

Indications: Rotator cuff repairs; Capsular shifts; Bankhart; Glenohumeral dislocations/subluxation; Soft tissue strains/repairs



Size	Measurements	Product No.	
		UltraSling III	UltraSling III AB
Small	up to 11"	11-0449-2	11-0450-2
Medium	11"-13"	11-0449-3	11-0450-3
Large	13"-15"	11-0449-4	11-0450-4
X-Large	15"+	11-0449-5	11-0450-5

Measure elbow crease to base of index finger

Ultrasling III ER

- Maintains the shoulder and arm in an externally rotated position
- Moisture wicking mesh fabric for greater comfort, particularly in longer term use
- 15° and 30° versions to suit your protocol
- Quick release buckles make it easy to take on and off



Size	Measurements	15° Product No.		30° Product No.	
		Left	Right	Left	Right
Small	up to 11"	11-1344-2-06000	11-1343-2-06000	11-1346-2-06000	11-1345-2-06000
Medium	11"-13"	11-1344-3-06000	11-1343-3-06000	11-1346-3-06000	11-1345-3-06000
Large	13"-15"	11-1344-4-06000	11-1343-4-06000	11-1346-4-06000	11-1345-4-06000
X-Large	15"+	11-1344-5-06000	11-1343-5-06000	11-1346-5-06000	11-1345-5-06000

Measure elbow crease to base of index finger

UltraSling II

- New moisture wicking mesh fabric for greater comfort
- Promotes auxiliary air exchange to reduce risk of secondary infections
- AB Version: Large comfortable pillow gives greater abduction where required (45° or 60°)
- Allows arm to be positioned in a variety of positions dependent on post-operative requirements (either more posterior or anterior)



Size	Measurements	Product No.	
		Ultra Sling II	AB Version
Small	up to 11"	11-0449-2-06000	11-0450-2-06000
Medium	11"-13"	11-0449-3-06000	11-0450-3-06000
Large	13"+	11-0449-4-06000	11-0450-4-06000

Measure elbow crease to base of index finger

UltraSling

- Cool breathable fabric
- Easy-open front panel encourages elbow exercises
- Encourages effective healing by allowing the shoulder and arm to remain in a neutral position
- 10° abduction pillow allows auxiliary air exchange to reduce the risk of secondary infections

Size	Measurements	Product No.
Small	up to 11"	11-0138-2-13130
Medium	11"-13"	11-0138-3-13130
Large	13"+	11-0138-4-13130

Measure elbow crease to base of index finger



S.C.O.I. Shoulder Brace

- Allows for 30°-150° of abduction/adduction range-of-motion hinge
- Universal, fits right or left
- Elbow range-of-motion hinge

Indications: Rotator cuff repairs; Capsular shifts; Prosthetic shoulder replacements; Bankhart procedures

Product No. **11-0100-9-13136**



Quadrant Shoulder Brace

- No auxiliary hardware
- Internal/external rotation
- Range-of-motion hinge for abduction/adduction, flexion/extension

Size	Product No.	
	Left	Right
Universal	11-0111-9-13136	11-0110-9-13136



Fillauer - E-Z Arm Abduction Orthosis

- “Worm gear” adjusts arm abduction angle from 20°-110° without removing any parts
- Waistband, forearm and humeral cuffs quickly position on slotted bars and lock by tightening set screws
- Malleable thoracic band adjustable to patient’s body contour
- Second snap-on shoulder strap optional for additional support

Product No. **109300**



Hely & Weber - Shoulder Houlder II

- NoSwet® mesh material combined with a soft moisture quick-wicking liner to provide your patient with comfort and support
- Body and shoulder straps provide quick adjustment, easy application, and padded comfort
- A posterior gusset comfortably holds the elbow in the desired position without providing unnecessary pressure to the area

Size	Height x Length*	w/ Abduction Pillow	w/ Neutral Pillow
Small	11 ¹ / ₂ " x 13"	514S	524S
Medium	13 ¹ / ₂ " x 14 ¹ / ₂ "	514M	524M
Large	15" x 16"	514L	524L
X-Large	16 ¹ / ₂ " x 17 ¹ / ₂ "	514XL	524XL

*Measure from olecranon to knuckles



Abduction Pillow



Neutral Pillow

Ossur Americas - SmartSling

- 3-in-1 design can provide immobilization (pictured), abduction and/or external rotation
- Includes abduction wedge and axilla pillow
- External rotation kit sold separately

Indications: Immobilization for rotator cuff repairs; Bankart repairs; Capsular shifts; Glenohumeral dislocations/subluxation; Soft tissue repairs/strains; Shoulder instabilities/injuries

Size	Length*	Product No.
Medium	13 ¹ / ₄ "-15 ³ / ₄ "	SSR
Large	16"-17 ¹ / ₂ "	SSL
External Rotation Kit	Universal	SSE

*Measure from olecranon to end of closed fist



Immobilization



Abduction




External Rotation

Orthomerica - Ultra Lite™ Airplane Splint

- Extremely lightweight plastic/aluminum, functional shoulder abduction splint

Indications: Rotator cuff repairs, involving humeral fractures; Post-op shoulder/arm surgery, where control of shoulder position and elbow flexion is desired; Tough shoulder immobilization problems

Size	Axilla to Medial Epicondyle	Axilla to Crest	Product No.
Small	6"-8 ¹ / ₂ "	7 ¹ / ₂ "-11"	1150
Medium/Large	8 ¹ / ₂ "-13"	11"-16"	1153




Newport Shoulder System

- Prefabricated functional shoulder system allows for full pelvic weight bearing, eliminating the need for shoulder straps
- Shoulder joint controls adduction/abduction, flexion/extension, and rotation of the shoulder and arm
- Scapular stabilization prevents anterior movement of the scapula so that connective tissue between the shoulder and the scapula is not stretched

Indications: Capsulitis

Size	Waist Circum.	Axilla to Crest	Axilla to Mid Palm	Product No.	
				Left	Right
Small	28"-34"	12"-15"	17"-26"	1180	1185
Medium	34"-40"	12"-15"	17"-26"	1181	1186
Large	34"-40"	15"-19"	19"-29"	1182	1187
X-Large	40"-46"	15"-19"	19"-29"	1183	1188



ProCare (DJO) - Shoulder Abduction Kit

- Large Comfortable pillow provides abduction where required (45° or 60°) to provide post-injury or post-surgical support
- Allows the shoulder and arm to remain in a neutral position
- Foam/nylon construction, universal design to fit most patients
- Universal fit for left/right

Product No. **79-84500**



Professional Products - The Six Shooter Abduction System

- Contoured non-slip foam provides constant abduction
- Deep arm pocket to improve immobilization
- 35° foam wedge provides constant abduction and comfort
- Universal left or right

Size	Forearm Measurement	Product No.
Small	9"-12"	2362-JA-C3-20
Medium	12"-15"	2362-JA-C3-30
Large	15"-18"	2362-JA-C3-40
X-Large	18"-20"	2362-JA-C3-50

Measure from lateral epicondyle to fifth metacarpo-phalangeal joint



Truform - Shoulder Abduction Support

- Carved foam block, 13" x 13" x 12", covered with stockinette and velfoam adjustment panels
- Reversible to maintain either 45° or 90° of abduction, contoured one side for positioning against rib cage
- Wrist and forearm cuffs, chest encircling strap with hook-and-pile fastening and adjustment
- One size fits up to 42" chest circumference; reversible and universally adjustable

Indications: Shoulder injury, separation or surgery; Post-stroke

Product No.

T84 0905



Trulife - Gunslinger II

- Variable abduction, adduction and forearm rotation based on height and position of the forearm component
- Extended forearm attachment provides increased variable abduction
- Forearm Abduction/Adduction: 75°; Shoulder Abduction/Adduction: 45°

Indications: Reconstruction surgery/shoulder replacement; Traumatic injuries/fractures to the proximal humerus and shoulder; Reduction of shoulder luxation; Post operative repair of rotator cuff; Bankart lesions

Size	Product No.	
	Left	Right
Universal	GS101L	GS101R



Gunslinger I

- Low density polyethylene shell with washable pads and liners
- Numbered, quick snap release buckles
- Forearm Abduction/Adduction: 75°

Indications: Reconstruction surgery/shoulder replacement; Traumatic injuries/fractures to the proximal humerus and shoulder; Reduction of shoulder luxation; Post operative repair of rotator cuff; Bankart lesions

Size	Product No.	
	Left	Right
Universal	GS100L	GS100R



AliMed - Hemi-Shoulder Sling

- Reduction joint-subluxation forces
- More lightweight and universal than conventional arm braces
- Shoulder saddle and humeral cuff are connected by three adjustable distraction straps that allow for full joint mobility while maintaining shoulder's integrity
- Elastic panel for muscle volume changes



Size	Measurement	Product No.	
		Left	Right
X-Small	9"	51-051	51-050
Small	11"	51-053	51-052
Medium	13"	51-055	51-054
Large	15"	51-057	51-056

Measure arm circumference at widest point of bicep

MFC™ Unilateral Shoulder Orthosis

- Helps prevent stretching of the tendons and ligaments that hold the humerus in position
- Can be worn under clothing
- Lightweight, breathable fabric



Size	Measurements	Product No.	
		Left	Right
Small	f9"-13"	5901/NA/LF	5901/NA/RT
Medium	12"-16"	5902/NA/LF	5902/NA/RT
Large	15"-20"	5903/NA/LF	5903/NA/RT

Measure the upper arm circumference at the largest point

Bauerfeind - OmoTrain® Shoulder Support

- Lightweight, knitted support guides and centers the glenohumeral joint promoting mobility while inhibiting rotation
- Adjustable, anatomically contoured silicone insert lies on shoulder providing intermittent compression to soft tissue

Indications: Shoulder pain, rotator cuff rehab, shoulder subluxation, post-surgery, osteoarthritis, and hemiparesis



Size	Loose Circumference at Thickest Part of Upper Arm	Product No.
0	7 ³ / ₄ "-8 ¹ / ₂ "	11071701080000
1	8 ¹ / ₂ "-9 ¹ / ₂ "	11071701080001
2	9 ¹ / ₂ "-10 ¹ / ₄ "	11071701080002
3	10 ¹ / ₄ "-11 ¹ / ₂ "	11071701080003
4	11 ¹ / ₂ "-12 ¹ / ₂ "	11071701080004
5	12 ¹ / ₂ "-14 ¹ / ₄ "	11071701080005
6	14 ¹ / ₄ "-15 ³ / ₄ "	11071701080006

Measure circumference of biceps while the arm is hanging loose

Becker Orthopedic - Wilmer® Carrying Orthosis

- Effectively neutralizes shoulder subluxation by using the weight of the forearm as a counterbalance to force the humerus into the shoulder complex
- Effectively suppresses edema formation in the forearm and hand

Size	Forearm Length	Forearm Circum.	Product No.	
			Left	Right
Normal	9"-11 ² / ₅ "	9"-11 ² / ₅ "	WCO-SNL	WCO-SNR
Large	9"-11 ² / ₅ "	11"-13 ² / ₅ "	WCO-SLL	WCO-SLR



Corflex - Cryo Pneumatic Shoulder

- Anatomical shoulder wrap designed to combine the benefits of pneumatic compression with the option of cold therapy; features pump bulb that detaches once desired compression is achieved
- Manufactured from plush foam laminate that is hook and loop compatible; optional gel pack with tricot lining; non-toxic gel
- Universal sizing

Indications: Post-surgical application; Acute injury as part of a physical therapy or athletic training program; Helps diminish swelling/edema; Reduce muscle spasms and decrease pain

Description	Product No.
0 Gel	72-7180
1 Gel	72-7181
2 Gels	72-7182
Additional Gels Sold Separately	
Single Shoulder Gel	72-7189



DonJoy (DJO) - Shoulder Stabilizer

- Provides abduction control, external rotation control and support for acromioclavicular separations
- The range-of-motion control can be customized with easy strap adjustments
- Latex-free black polyester Lycra material
- Available in bilateral version upon request

Indications: Glenohumeral dislocations; Rotator cuff tears; Acromioclavicular separation/subluxations

Size	Chest Circum.	Product No.	
		Left	Right
Small	30"-34"	11-1102-2-06000	11-1101-2-06000
Medium	34"-38"	11-1102-3-06000	11-1101-3-06000
Large	38"-41"	11-1102-4-06000	11-1101-4-06000
X-Large	41"-44"	11-1102-5-06000	11-1101-5-06000
2X-Large	44"-48"	11-1102-6-06000	11-1101-6-06000
3X-Large	48"-52"	11-1102-7-06000	11-1101-7-06000

Measurement is taken at chest circumference under arms and across nipples



**New Options Sports - S5 ARM-ADILLO™
Shoulder Stabilizer**

- 1/8" VelPlush™ perforated neoprene
- Two 1" x 3" hook straps originating from the bicep are attached to body of brace to resist abduction
- The two separate control straps are designed with alligator attachments on one end which allow them to be trimmed and re-attached to achieve desired length

Indications: Abduction, Grade I and II ACL separation; Pre- and post-surgery for Grade III separation; Arm and shoulder stabilization; Elbow subluxation; Pre- and post- bicep repair; Subluxation caused by hyperextension




Size	Chest Circum.	Bicep Measurement	Straps	Product No.
X-Small	22"-25"	10"-13"	2" x 16" 2" x 18"	S5 XS
Small	25"-28"	11"-14"	2" x 18" 2" x 20"	S5 S
Medium	28"-31"	11½"-14½"	2" x 20" 2" x 22"	S5 M
Large	32"-36"	12"-14¾"	2" x 21" 2" x 23"	S5 L
X-Large	36"-40"	12½"-15½"	2" x 22" 2" x 24"	S5 XL
2X-Large	39"-43"	13"-16"	2" x 23" 2" x 25"	S5 XXL
3X-Large	42"-48"	13½"-16½"	2" x 25" 2" x 27"	S5 XXXL
4X-Large	47"-52"	14½"-17"	2" x 27" 2" x 29"	S5 XXXXL

Orthomerica - Highland™ Shoulder Sling

- Designed to limit chronic shoulder subluxation/dislocation
- Lightweight, low-profile design is anatomically contoured for maximum patient comfort and function
- Easily trimmable components for customized fitting
- Includes adjustable arm sling





Indications: CVA; Ligamentous rupture; Peripheral nerve injury








Size	Circumference		Product No.	Replacement Pads
	Distal	Proximal		
Small	7"-9"	9"-11"	1220	1228.01
Medium	9"-11"	11"-13"	1221	1228.02
Large	11"-13"	13"-15"	1222	1228.03
Medial Shell Option			1228.05	

Clavicle Straps

- Indications: Clavicular fractures; Posture correction; Acromioclavicular separations

CoreLINE	Clavicle Strap <ul style="list-style-type: none"> • Used to immobilize the clavicle, it is an ideal choice for clavicular fractures • Features contoured urethane foam for better anatomical fit and patient comfort • Fully adjustable, with rear buckles for tension adjustment 				
	Size	Chest Circum.	Product No.		
	X-Small	20"-24"	CL0290		
	Small	24"-30"	CL0292		
	Medium	30"-36"	CL0294		
	Large	36"-42"	CL0296		
	X-Large	42"-48"	CL0298		
2X-Large	48"-52"	CL0299			
Bledsoe	Clavicle Support <ul style="list-style-type: none"> • Urethane foam contoured for fit and comfort • Easy touch closures are fully adjustable • Easily adjust tension with rear buckles 				
	Size	Chest Circum.	Product No.		
	X-Small	20"-24"	SA301001		
	Small	24"-30"	SA301003		
	Medium	30"-36"	SA301005		
	Large	36"-42"	SA301007		
X-Large	42"-48"	SA301009			
DeRoyal	Premium Clavicle Strap <ul style="list-style-type: none"> • Synthetic stockinette over foam padding does not absorb water and can be worn while bathing • Serpentine stitch design to prevent migration • Measurement taken at widest part of chest for men and below the bustline for women 				
	Size	Chest Circum.	Product No.		
			Buckle	Hook-and-Loop	
	Infant	18"-20"	3001-01	3002-01	
	Child	20"-24"	3001-02	3002-02	
	Small	24"-30"	3001-06	3002-06	
	Medium	30"-36"	3001-03	3002-03	
	Large	36"-42"	3001-04	3002-04	
	X-Large	42"-48"	3001-05	3002-05	
	Heavy-Duty Clavicle Strap <ul style="list-style-type: none"> • Foam padding with stockinette covering • Extra long straps for easy adjustment • Measurement taken at widest part of chest for men and below the bustline for women 				
	Size	Chest Circum.	Product No.		
			Buckle	Hook-and-Loop	
	Infant	18"-20"	3004-01	3004-11	
	Child	20"-24"	3004-02	3004-12	
Small	24"-30"	3004-06	3004-16		
Medium	30"-36"	3004-03	3004-13		
Large	36"-42"	3004-04	3004-14		
X-Large	42"-48"	3004-05	3004-15		

Clavicle Straps (continued)

Freeman	<p>Clavicle Splint</p> <ul style="list-style-type: none"> Adjustable crisscross buckle strap construction of 1" wide webbing is covered with natural cotton stockinette foam cushioned shoulder pads Complete buckle assembly is attached to felt back pad for cushioning 			
	Size	Product No.		
	Universal	996-1		
Ossur Americas	<p>Front Closure Clavicle Support</p> <ul style="list-style-type: none"> Unique supportive clavicle splint Front closure design for convenience in application Comfortable padding over clavicle area Measure the chest circumference at the sternum area 			
	Size	Chest Circum.		Product No.
	X-Small	20"-24"		205032
	Small	24"-30"		205033
	Medium	30"-36"		205035
	Large	36"-42"		205037
	X-Large	42"-48"	205038	
	2X-Large	48"-54"	205039	
	<p>Figure 8 Clavicle Splint</p> <ul style="list-style-type: none"> Prong buckle closure ensures correct fit without slippage Measure chest circumference at sternum area 			
	Size	Chest Circum.		Product No.
X-Small	20"-24"	205002		
Small	24"-30"	205003		
Medium	30"-36"	205005		
Large	36"-42"	205007		
X-Large	42"-48"	205008		
ProCare (DJO)	<p>Clavicle Strap</p> <ul style="list-style-type: none"> Foam padded shoulder straps covered with stockinette for added patient comfort 3-way adjustment buckle closure helps prevent strap slippage 			
	Size	Chest Circum.		Product No.
	X-Small	20"-24"		79-85012
	Small	24"-30"		79-85013
	Medium	30"-36"		79-85015
	Large	36"-42"		79-85017
	X-Large	42"-48"	79-85018	
	<p>Heavy Padded Clavicle</p> <ul style="list-style-type: none"> Combination felt and foam padded straps for added comfort Contact closure for easier adjustment Chest circumference measurement taken at bottom of sternum 			
	Size	Chest Circum.		Product No.
	2X-Small	18"-20"		79-85021
	X-Small	20"-24"		79-85022
	Small	24"-30"		79-85023
Medium	30"-36"	79-85025		
Large	36"-42"	79-85027		
X-Large	42"-48"	79-85028		

Clavicle Straps (continued)

Trulife	Contour Clavicle Strap		
	<ul style="list-style-type: none"> • Vertical T-pattern, non-stretch straps provide support and maintain position • Contoured padding protects brachial artery and nerve 		
	Size	Chest Circum.	Product No.
	2X-Small	15"-20"	1640000
	X-Small	20"-27"	1640001
	Small	27"-34"	1640002
	Medium	30"-36"	1640003
Large	36"-42"	1640004	
X-Large	42"-48"	1640005	



AliMed - Elastic Torso Belt Shoulder Immobilizer

- Body band is 4" wide firm knit elastic
- Three wide plastic stays prevent rolling
- Vinyl arm and wrist band fits both left and right
- Felt-lined, front contact closure

Size	Fits Chest Size	Product No.
Small	24"-30"	510642
Medium	30"-36"	510643
Large	36"-40"	510644
X-Large	40"-44"	510645
2X-Large	44"-48"	510646

Females measure just below bust line; Males measure chest at widest point



DeRoyal - Deluxe Sling and Swathe

- Lightweight perforated foam material provides patient comfort
- Adjustable shoulder straps attach to hook receivable sling to ensure accurate fit
- Open or closed elbow designs
- Swathe attaches at any point on sling

Indications: S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Pouch Size	Strap Length	Product No.
Universal	5" x 54"	13" x 19"	9006-00
Universal w/ X-Long Swathe	5" x 58"	13" x 19"	9006-50
Universal Enclosed Elbow	5" x 54"	17" x 19"	9014-00
Universal Pediatric	4 ¹ / ₂ " x 44"	9" x 16 ¹ / ₂ "	9007-00
Swathe Only	5" x 54"	—	9006-10



Universal Foam Sling and Swathe

- Adjustable shoulder straps attach to hook receivable sling to ensure accurate fit
- Swathe attaches at any point on sling
- Perforated foam material for patient comfort
- Extra-long swathe size available

Size	Swathe Size	Pouch Size	Product No.
Universal	5" x 54"	13" x 19"	4442-00
Universal with X-Long Swathe	5" x 58"	13" x 19"	4442-10



Canvas Shoulder Immobilizer with Waist Strap

- Navy canvas sling with web shoulder and waist straps
- Shoulder strap includes O-ring/hook and loop closure
- Thumb loop helps prevent wrist drop
- Universal left or right

Indications: S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Pouch Size	Product No.
X-Small	6" x 12 ¹ / ₂ "	9002-01
Small	7" x 13 ¹ / ₂ "	9002-02
Medium	8" x 16 ¹ / ₂ "	9002-03
Large	9" x 19 ¹ / ₂ "	9002-04
X-Large	9" x 22 ¹ / ₂ "	9002-05



Universal Foam Shoulder Immobilizer

- Contoured perforated foam ensures patient comfort
- Wrist cuff, humeral cuff and shoulder strap are adjustable and can be placed anywhere on belt to achieve desired fit
- Provides shoulder and upper arm immobilization following dislocation or post-surgical application

Indications: S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Belt Size	Wrist Cuff Size	Product No.
Universal	7 ¹ / ₂ " x 48"	5" x 12"	9009-00
Universal*	8" x 54"	3" x 9"	9009-50

*Available in light or dark foam



Velpau Shoulder Immobilizer

- Heavy duty elastic material
- Durable foam padded wrist and humeral cuffs
- Contoured design with hook and loop closure for patient comfort
- Universal left or right

Indications: S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Belt Size	Product No.
X-Small	6" x 29"	9008-00
Small	6" x 32"	9008-01
Medium	6" x 35"	9008-02
Large	6" x 38"	9008-03
X-Large	6" x 41"	9008-04
2X-Large	6" x 44"	9008-05
3X-Large	6" x 47"	9008-06



Shoulder Immobilizer with Foam Straps

- Canvas or lightweight tietex material with foam shoulder and waist straps
- Adjustable straps feature O-ring/hook and loop closure
- Thumb loop helps prevent wrist drop
- Universal left or right

Indications: S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Pouch Size	Product No.	
		Canvas	Tietex
X-Small	6" x 12 ¹ / ₂ "	9023-00	TX9023-00
Small	7" x 13 ¹ / ₂ "	9023-01	TX9023-01
Medium	8" x 16 ¹ / ₂ "	9023-02	TX9023-02
Large	9" x 19 ¹ / ₂ "	9023-03	TX9023-03
X-Large	9" x 22 ¹ / ₂ "	9023-04	TX9023-04



Elastic Shoulder Immobilizer

- Canvas or lightweight elastic material
- Foam wrist and humeral cuffs
- Hook and loop closure; Universal left or right

Indications: S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

Size	Pouch Size	Product No.
X-Small	6" x 29"	1000159
Small	6" x 32"	1010159
Medium	6" x 35"	1020159
Large	6" x 38"	1030159
X-Large	6" x 41"	1040159
2X-Large	6" x 44"	1050159



Hely & Weber - Shoulder Immobilizer

- Padded pressure sensitive straps are totally adjustable and comfortable
- Body strap is easily adjusted or can be removed
- Deep pocket design
- Two adjustment points for each strap allows for better fit

Size	Height x Length	Product No.
X-Small	7 ¹ / ₂ " x 11"	501XS
Small	7 ¹ / ₂ " x 14"	501S
Medium	7 ¹ / ₂ " x 15 ¹ / ₂ "	501M
Large	7 ¹ / ₂ " x 18"	501L
X-Large	7 ¹ / ₂ " x 20"	501XL
*Measure from elbow to metacarpals		



Ossur Americas - Premium Shoulder Immobilizer

- Comfort fit
- Easily adjustable
- Universal fit

Product No. **204100**



Elastic Shoulder Immobilizer

- Padded wrist support
- Contact closure
- Elastic chest band

Size	Waist Circumference		Product No.	
	Male	Female	Male	Female
Small	26"-30"	22"-26"	204033	204043
Medium	30"-34"	26"-32"	204035	204045
Large	36"-40"	32"-36"	204037	204047
X-Large	30"-46"	36"-42"	204038	204048



ProCare (DJO) - Sling and Swathe

- Durable cotton/duck material
- Slide buckle closure on soft cotton woven shoulder strap

Indications: Immobilization of the shoulder after surgery, injury or cast removal

Size	Swathe Size	Envelope Size	Product No.
Small	4" x 44"	9" x 14"	79-84243
Medium	4" x 49"	10" x 16"	79-84245
Large	4" x 54"	11" x 17"	79-84247
X-Large	4" x 59"	12" x 19"	79-84248



Shoulder Immobilizer with Foam Straps

- Soft cotton/poly for patient comfort
- Designed to exert upward pressure and immobilize arm and shoulder
- Waist and shoulder straps are a plush foam laminate for contact closure
- Fits left or right arm

Indications: Immobilization and support of the shoulder and elbow joints

Size	Envelope Size Depth & Length	Product No.
X-Small	6" x 11"	79-84162
Small	7" x 13"	79-84163
Medium	8" x 15"	79-84165
Large	9" x 18"	79-84167
X-Large	10" x 20"	79-84168



Elastic Shoulder Immobilizer

- Elastic waist band construction for compressive support
- Durable yet comfortable foam wrist and nylon humeral cuff with contact closure for ease of application
- Fits left or right arm

Indications: Dislocation or post-surgical immobilization

Size	Rib Measurement		Product No.	
	Male	Female	Male	Female
X-Small	24" - 32"	—	79-84032	—
Small	32" - 36"	24" - 28"	79-84033	79-84043
Medium	36" - 40"	28" - 32"	79-84035	79-84045
Large	40" - 44"	32" - 36"	79-84037	79-84047
X-Large	44" - 48"	36" - 42"	79-84038	79-84048
2X-Large	48" - 52"	42" - 48"	79-84039	79-84049



Acromioclavicular Brace

- Cotton/duck construction and felt shoulder straps provide additional comfort on shoulder area
- Splint helps stabilize shoulder, humerus and elbow following injuries
- Universal sizing fits most patients
- Sizing: 21"-42", chest circumference measurement taken at bottom of sternum

Product No. **79-85010**



Deluxe Shoulder Immobilizer

- Adjustable wrist cuff and shoulder strap may be positioned as needed for proper control and custom fit
- Fits right or left and may be worn without shoulder strap

Indications: Treatment of traumatic or post-surgical immobilization

Size	Waist Measurement	Product No.
Child	up to 38 ¹ / ₂ "	79-84101
Adult/Standard	up to 52"	79-84100
Adult/Standard*	up to 52"	79-96820
Adult/Contoured	up to 52"	79-84100C

*with detachable humeral strap



Professional Products - Semi Universal Shoulder Immobilizer

- 1¹/₂" shoulder strap and ³/₈" foam pad
- Humerus strap is extended to strap to attach posteriorly
- Shoulder strap may be adjusted for use as a hand loop
- Hook and loop, loop-lock closure

Product No. **01460**

Size	Chest Circumference
Small	up to 44"
Large	up to 54"

*Must specify size at time of order



Tapered Arm Sling/Shoulder Immobilizer

- Tapered arm sling for use with a cast
- Durable blue twill pocket has 2" comfort foam shoulder strap for maximum comfort
- Thumb loop maintains hand in neutral position and prevents sling from migrating
- Hook and loop closure

Product No. **02351**

Size	Forearm Measurement
Pediatric	6"-8"
Small	8"-12"
Medium	12"-15"
Large	15"-18"
X-Large	18"-20"
*Must specify size at time of order	



Trulife - FAST-WRAP Sling and Swathe

- Restricts mobility of shoulder and upper arm
- Easy to position swathe
- Criss-cross shoulder straps fasten in front and provide extra arm and wrist support

Indications: Separations; Dislocations; Strains

Size	Sling Dimensions	Product No.
Small	9" x 11 ¹ / ₂ "	1420002
Medium	9 ¹ / ₂ " x 12 ¹ / ₂ "	1420003
Large	10" x 14"	1420004
X-Large	10 ¹ / ₂ " x 15"	1420005



AliMed - Humeral Fracture Orthosis (Over the Shoulder)

- Two-part medial/lateral increases circumferential compression
- Plastic deltoid extension adds protection and helps prevent distal migration
- Full elbow ROM
- Includes two stockinette sleeves

Indications: Diaphyseal humeral fractures; Support after removal of internal or external fixation; As an adjunct to internal fixation



Size	Circumference	Product No.
Small	8"-11"	51312/NA/SM
Medium	11"-14"	51312/NA/MD
Large	12"-15"	51312/NA/LG
X-Large	14"-17"	51312/NA/XL

Measure circumference at largest part of biceps

Corflex - Extended Boxer Splint

- Anatomically designed to provide comfortable support, protection and immobilization of the 4th and 5th metacarpals; extended length provides added wrist support
- Lightweight low density polyethylene shell with soft foam liner and palmer pad and closure straps

Indications: Treatment of acute or sub-acute fractures of the 4th and 5th metacarpals, also known as Boxer's fracture



front view

side view

Size	Wrist Circum.	4th & 5th MCP	Product No.	
			Left	Right
Small	4"-5"	2 ³ / ₄ "-3 ¹ / ₄ "	37-1231	37-1221
Medium	5"-6"	3 ¹ / ₄ "-3 ³ / ₄ "	37-1232	37-1222
Large	6"-7"	3 ³ / ₄ "-4 ¹ / ₄ "	37-1233	37-1223
X-Large	7"-8"	4 ¹ / ₄ "-4 ³ / ₄ "	37-1234	37-1224

Measurement taken from wrist circumference and 4th & 5th MCP joint

Extended and Mid-Length Humeral Splint

- Controls humeral segments through circumferential soft tissue compression of the humerus
- Features deltoid cap extension to prevent migration
- Rigid polyethylene laminated with closed cell foam lining and contact closure
- Includes two stockinette undersleeves and one arm sling




Size	Proximal Circum. (widest part of bicep)	Distal Circum. (above elbow)	Product No.	
			Extended	Mid Length
X-Small	7"-9"	6"-7"	37-2160	37-2180
Small	9"-11"	7"-9"	37-2161	37-2181
Medium	11"-13"	9"-11"	37-2162	37-2182
Large	13"-15"	11"-13"	37-2163	37-2183
X-Large	15"-17"	13"-15"	37-2164	37-2184
2X-Large	17"-19"	15"-17"	37-2165	37-2185

DeRoyal - Colles' Splint

- Pre-molded aluminum design for anatomical fit and proper support
- Perforated design for air circulation
- X-ray lucent

Indications: Colles' injury; Wrist immobilization

Size	Product No.	
	Left	Right
Infant	9105-08	9105-07
Child	9105-04	9105-01
Youth	9105-05	9105-02
Adult	9105-06	9105-03



Plastic Colles' Splint

- Pre-molded plastic design for anatomical fit and proper support
- Perforated design for air circulation
- X-ray lucent

Indications: Colles' injury; Wrist immobilization

Size	Product No.	
	Left	Right
Child	9105-20	9105-24
Youth	9105-21	9105-25
Adult	9105-22	9105-26



Fracture Bracing

- Orthopedic grade thermoplastic construction
- Durable, lightweight, adjustable and easy to remove for proper hygiene
- Allows for movement of muscles surrounding fracture site
- Hook and loop closures help maintain soft tissue compression

Indications: Treatment and management of diaphyseal fractures; Post surgical immobilization

Size	Circumference		Product No.
	Distal	Proximal	
Small	7"-9"	9"-11"	FB2000-10
Medium	9"-11"	11"-13"	FB2000-12
Large	11"-13"	13"-15"	FB2000-14
X-Large	13"-15"	15"-17"	FB2000-16



DonJoy (DJO) - Humeral Stabilizing System

- 3 point stabilization
- 2 lb weight for axial traction
- Full circumference shells and padding
- Breakdown capability converting to upper arm cuff only

Indications: Stable humeral fractures; Stable supracondylar fractures; Biceps/triceps repairs


Size	Product No.	
	Left	Right
Small	11-0113-9-13066	11-0112-9-13066



Orthomerica - Wrist-Hand Orthosis

- Adjustable hook and loop straps allow for swelling and facilitate wound management
- Easily trimmed with cast scissors and reshaped with heat gun
- Large padded thumb component, 2 socks (with thumb)
- Perforated for ventilation


Size	Palm Width	Forearm Circum.	Brace Length	Product No.	
				Left	Right
Small	2 ¹ / ₂ "-3"	7"-9"	7 ¹ / ₄ "	1010	1015
Medium	3"-3 ¹ / ₂ "	9"-10 ¹ / ₂ "	7 ¹ / ₂ "	1011	1016
Large	3 ¹ / ₂ "-4"	10"-12"	8"	1012	1017



Wrist-Hand-Thumb Orthosis


- Tongue is rigid yet comfortable for most traumatic, post-surgical or longer term rehabilitation applications
- Adjustable hook and loop straps allow for swelling and facilitate wound management
- Large thumb component, 2 socks (with thumb)
- Perforated for ventilation

Size	Palm Width	Forearm Circum.	Brace Length	Product No.	
				Left	Right
Small	2 ¹ / ₂ "-3"	7"-9"	7 ¹ / ₂ "	1000	1005
Medium	3"-3 ¹ / ₂ "	9"-10 ¹ / ₂ "	7 ³ / ₄ "	1001	1006
Large	3 ¹ / ₂ "-4"	10"-12"	8 ¹ / ₂ "	1002	1007



Soft Humeral Fracture Brace


- Compression molded foam provides flexibility at the edges with a rigid structural core
- “Shoulder Cap” limits distal migration
- Dual tongue design molded into medial shell for ease of adjustability and compression
- Universal, fits right or left



Size	Proximal Circum.	Distal Circum.	Brace Length	Product No.
X-Small	8"-10"	6"-8"	11"	1229
Small	9"-11"	7"-9"	12"	1230
Medium	11"-13"	9"-11"	13"	1231
Large	13"-15"	11"-13"	13"	1232
X-Large	15"-17"	13"-15"	14"	1233

Humeral Fracture Brace Kit (Shoulder Type)


- Polyethylene shell for comfortable lightweight control
- Soft hook and pile suspension strap passes across chest, back, and under opposite arm
- Shoulder cap limits distal migration yet allows full shoulder movement
- Universal fits left or right



Size	Proximal Circum.	Distal Circum.	Lateral Brace Length	Approx. Medial Length	Product No.
X-Small	8"-10"	6"-8"	11 ¹ / ₂ "	6 ¹ / ₂ "	1089
Small	9"-11"	7"-9"	12"	7"	1090
Medium	11"-13"	9"-11"	13"	7 ¹ / ₂ "	1091
Large	13"-15"	11"-13"	14"	8"	1092
X-Large	15"-17"	13"-15"	14 ¹ / ₄ "	8 ¹ / ₂ "	1093

Humeral Fracture Brace Kit - Bi-Valved Model

- Adjustable hook and loop strap system maintains compression and limits distal migration
- Lined with closed-cell foam
- Perforated for ventilation
- Socks and sling provide comfort and hygiene



Size	Overall Length		Mid Upper Arm Circum.	Product No.	
	Medial	Lateral		Left	Right
Small	5 ¹ / ₂ "	8 ¹ / ₂ "	8"-11"	1350	1355
Medium	5 ¹ / ₂ "	8 ³ / ₄ "	11"-14"	1351	1356
Large	6"	9 ³ / ₄ "	12"-15"	1352	1357
X-Large	6 ³ / ₄ "	10 ³ / ₈ "	12"-15"	1353	1358

Humeral Fracture Brace Kit - Universal Model (Sarmiento Type)


- Allows shoulder range-of-motion
- Adjustable hook and loop strap system maintains compression and limits distal migration
- Lined with closed cell foam; Perforated for ventilation
- Kit includes 2 socks and one premium padded arm sling



Size	Proximal Circum.	Distal Circum.	Brace Length	Product No.	
				Kit	Sleeve Only
X-Small	7"-9"	5"-7"	5 ¹ / ₂ "	1079	1069
Small	9"-11"	7"-9"	7 ¹ / ₄ "	1080	1070
Medium	11"-13"	9"-11"	8"	1081	1071
Large	13"-15"	11"-13"	8"	1082	1072
Large	15"-17"	13"-15"	10 ³ / ₄ "	1083	1073

Ulnar Fracture Brace Kit - Bi-Valved Model

- Maintains reduction of isolated ulnar fracture while allowing adjacent joint movement
- Lined with perforated closed-cell foam for comfort
- Bi-valved anatomical, side specific design (left or right)
- Includes two socks




Size	Proximal Circum.	Distal Circum.	Brace Length	Product No.	
				Left	Right
Small	7"-9"	6"-7"	7"	1030	1035
Medium	9"-11"	7"-8"	8 ¹ / ₄ "	1031	1036
Large	11"-13"	8"-9"	9 ¹ / ₄ "	1032	1037

For brace, only .01 to end of product number

Ulnar Fracture Brace Kit - Universal Model (Sarmiento Type)

- Designed to allow more comfortable elbow motion, wrist flexion and extension; Universal design fits right or left
- Adjustable hook and loop straps maintain compression for fracture alignment and comfort; Lined with closed cell foam
- Kit Includes 2 socks




Size	Proximal Circum.	Distal Circum.	Brace Length	Product No.	
				Kit	Sleeve Only
Small	7"-9"	6"-7"	6"	1050	1040
Medium	9"-11"	7"-8"	6 ³ / ₄ "	1051	1041
Large	11"-13"	8"-9"	8"	1052	1042

Prime™ Elbow System with Wrist Component

- Ulnar sleeve or wrist-hand module to control rotation
- Perforated for ventilation
- Adjustable hook and loop straps facilitate application and removal for wound inspection, dressing and sock changes
- Both kits include two Colles' socks

Indications: Where limited range-of-motion and rotation control are desired



With Ulnar Component						
Size	Upper Arm Circum.		Forearm Circum.	Ulnar Comp. Length	Product No.	
	Distal	Proximal			Left	Right
Small	7"-9"	9"-11"	7"-9"	6"	1120	1125
Medium	9"-11"	11"-13"	9"-11"	7"	1121	1126
Large	11"-13"	13"-15"	11"-13"	8"	1122	1127


With Wrist-Hand Component							
Size	Upper Arm Circum.		Forearm Circum.	Palm Width	Wrist-Hand Comp. Length	Product No.	
	Distal	Proximal				Left	Right
Small	7"-9"	9"-11"	7"-9"	2 ¹ / ₂ "-3"	10 ¹ / ₂ "	1130	1135
Medium	9"-11"	11"-13"	9"-11"	3"-3 ¹ / ₂ "	11 ³ / ₄ "	1131	1136
Large	11"-13"	13"-15"	11"-13"	3 ¹ / ₂ "-4"	12 ¹ / ₄ "	1132	1137
Accessories							
Standard Elbow Joint						1144	

Elbow Orthosis

- Effective in limiting flexion and extension
- Upper and lower arm modular components are inter-changeable
- Perforated and lined with perforated closed cell foam
- Two socks included with each orthosis

Indications: Elbow Arthroplasty; Post-surgical; Post-injury

Size	Mid Upper Arm		Forearm		Product No.
	Circum.	Brace Length	Circum.	Brace Length	
Small	8"-10"	7 ¹ / ₄ "	7"-9"	6"	1110
Medium	10"-12"	7 ¹ / ₂ "	9"-11"	6 ¹ / ₂ "	1111
Large	12"-14"	8"	11"-13"	7 ³ / ₄ "	1112
X-Large	14"-16"	10 ⁷ / ₈ "	13"-15"	8"	1113



ProCare (DJO) - Humeral Fracture Brace/Over the Shoulder

- Lightweight, foam lined brace is designed for the management of mid-shaft and distal ulnar fractures
- Circumferential soft tissue compression provides cast-like rigidity for optimal control of ulnar segments
- Each brace comes with two (2) double thick cotton stockinettes
- Universal fit for left or right

Size	Bicep Circumference	Length	Product No.
Small	8"-10"	10"	79-97953
Medium	10"-13"	10 ¹ / ₄ "	79-97955
Large	12"-15"	11 ¹ / ₄ "	79-97957
X-Large	14"-17"	12 ³ / ₄ "	79-97958



Humeral Cuff

- Prefabricated foam lined orthosis designed for the management of midshaft and distal third humeral fractures
- Provides total contact circumferential soft tissue compression for optimal control of humeral segments while allowing range-of-motion of both shoulder and elbow joints
- Each brace comes with two (2) double thick cotton stockinettes
- Universal fit for left or right

Size	Bicep Circumference	Length	Product No.
Small	8"-11"	6"	79-99213
Medium	10"-13"	6 ¹ / ₄ "	79-99215
Large	12"-15"	7 ¹ / ₄ "	79-99217
X-Large	14"-17"	8"	79-99218



Colles Splint

- Premolded, padded aluminum splint with ventilation holes to provide air circulation for patient comfort

Indications: Wrist support and immobilization prior to surgery or casting

Size	Length x Width	Product No.	
		Left	Right
X-Small	5 ³ / ₄ " x 2 ¹ / ₂ "	79-72112	79-72122
Small	7" x 3 ¹ / ₂ "	79-72113	79-72123
Medium	8 ¹ / ₂ " x 3 ¹ / ₂ "	79-72115	79-72125
Large	9" x 4 ¹ / ₂ "	79-72117	79-72127



Colles Splint

Ulnar Fracture Brace

- Lightweight, foam lined brace designed for the management of mid-shaft and distal ulnar fractures
- Circumferential soft tissue compression provides cast like rigidity for optimal control of ulnar segments
- Full range of wrist and elbow motion is permitted to help prevent atrophy, increase vascularity and encourage bone growth
- Each brace comes complete with two double thick cotton stockinettes

Size	Circumference		Product No.
	Distal	Proximal	
Small	6"-7"	8"-11"	79-97943
Medium	6 ³ / ₄ "-7 ³ / ₄ "	9 ¹ / ₂ "-12 ¹ / ₂ "	79-97945
Large	7 ¹ / ₂ "-8 ¹ / ₂ "	11"-14"	79-97947



Ulnar Fracture Brace

Wrist Hand Orthosis

- Lightweight orthosis designed to comfortably restrict motion of the wrist and hand while holding hand in neutral position
- Fabricated with a foam lined orthopedic grade thermoplastic
- Splint may be trimmed with cast scissors and formed with a heat gun
- Each brace comes complete with two double thick cotton stockinettes

Size	Palm Width	Length	Product No.	
			Left	Right
Small	2 ¹ / ₂ "-3"	7 ¹ / ₂ "	79-97973	79-97963
Medium	3"-3 ¹ / ₂ "	7 ¹ / ₂ "	79-97975	79-97965
Large	3 ¹ / ₂ "-4"	8"	79-97977	79-97967



Wrist Hand Orthosis

Wrist Hand Thumb Orthosis

- Designed to restrict motion of the wrist and thumb while holding the hand in a neutral position
- Fabricated with a foam lined orthopedic grade thermoplastic
- Orthosis provides rigid cast-like immobilization of the hand, wrist and thumb without obstructing web space between the thumb and fore finger
- Each brace comes complete with two double thick cotton stockinettes

Size	Palm Width	Length	Product No.	
			Left	Right
Small	2 ¹ / ₂ "-3"	7 ¹ / ₂ "	79-97713	79-97723
Medium	3"-3 ¹ / ₂ "	7 ¹ / ₂ "	79-97715	79-97725
Large	3 ¹ / ₂ "-4"	8"	79-97717	79-97727



Wrist Hand Thumb Orthosis

RCAL - Over the Shoulder Humeral Fracture Brace

- A universal design for humeral diaphyseal fractures and deep humeral bone bruises
- Allows range-of-motion at shoulder and elbow
- Trimmable polyethylene and lined with closed-cell foam
- Lightweight, durable, and cannot absorb moisture



Size	Mid-Humeral Circumference	Product No.
Toddler	3 ¹ / ₂ "-6 ¹ / ₂ "	54HFB-OS-T
Child	4"-6 ¹ / ₂ "	54HFB-OS-C
Youth	5"-8"	54HFB-OS-Y
X-Small	5"-7 ¹ / ₂ "	54HFB-OS-ES
Small	7 ¹ / ₂ "-9 ¹ / ₂ "	54HFB-OS-S
Medium	9 ¹ / ₂ "-11 ¹ / ₂ "	54HFB-OS-M
Large	11"-15"	54HFB-OS-L
X-Large	13"-18"	54HFB-OS-XL

Trulife - Ulnar Fracture Brace

- Stabilizes ulnar while allowing free range-of-motion of wrist and elbow
- Easy, single-handed donning and doffing

Indications: Ulnar fractures



Size	Circumference		Product No.	
	Distal	Proximal	Left	Right
Small	5 ¹ / ₂ " -6 ¹ / ₂ "	6 ¹ / ₂ " -8 ¹ / ₂ "	JS-1612	JS-1611
Medium	6 ¹ / ₂ " -7 ¹ / ₂ "	8 ¹ / ₂ " -11 ¹ / ₂ "	JS-1614	JS-1613
Large	7 ¹ / ₂ " -8 ¹ / ₂ "	11 ¹ / ₂ " -14 ¹ / ₂ "	JS-1616	JS-1615
X-Large	8 ¹ / ₂ " -9 ¹ / ₂ "	14 ¹ / ₂ " -16 ¹ / ₂ "	JS-1618	JS-1617

Humeral Fracture Orthosis

- Stabilizes humerus
- Single-handed donning and doffing

Indications: Non-displaced mid-shaft to distal third humeral fractures



Size	Circumference			Product No.	
	Distal	Proximal	Length	Left	Right
Small	7"-9"	9"-11"	7 ¹ / ₂ "	JS-1512	JS-1511
Medium	9"-11"	11"-13"	8"	JS-1514	JS-1513
Large	11"-13"	13"-15"	8"	JS-1516	JS-1515
X-Large	13"-15"	15"-17"	11"	JS-1518	JS-1517
Accessories					
Arm Sling	7178D				

Universal Extended Length Humeral Fracture Brace

- Humeral brace with shoulder cap
 - Cross-shoulder strap prevents distal migration
 - Low density, trimmable polyethylene
- Indications:** Non-displaced mid-shaft to distal third humeral fractures



Size	Circumference		Product No.
	Bicep	Proximal	
X-Small	5"-7"	7"-9"	1182090000
Small	7"-9"	9"-11"	1182090001
Medium	9"-11"	11"-13"	1182090002
Large	11"-13"	13"-15"	1182090003
X-Large	13"-15"	15"-17"	1182090004

Over Shoulder Extension - Humeral Fracture Orthosis

- Stabilizes humerus while allowing free range-of-motion of shoulder and elbow joint
 - Cross-shoulder strap prevents distal migration
 - One piece, moldable thermoplastic
- Indications:** Non-displaced mid-shaft to distal third humeral fractures



Size	Side	Circumference		Product No.
		Distal	Proximal	
Small	Right	7"-9"	8"-11"	JS-1521
	Left	7"-9"	8"-11"	JS-1522
Medium	Right	9"-11"	11"-14"	JS-1523
	Left	9"-11"	11"-14"	JS-1524
Large	Right	11"-14"	14"-16"	JS-1525
	Left	11"-14"	14"-16"	JS-1526
X-Large	Right	14"-16"	16"-18"	JS-1527
	Left	14"-16"	16"-18"	JS-1528

Static Wrist Hand Finger Positioning Splint

- Finger channels maintain position
 - Easy, single-handed donning and doffing
- Indications:** Treatment of burns; Hemiplegia; Post-fracture; Post-surgery; Arthritis of wrist, hand and fingers



Size	Metacarpal Width	Wrist Circum.	Forearm Circum.	Product No.	
				Left	Right
X-Small	1 1/2"-2 1/2"	4"-6"	5"-6"	JS-1002	JS-1001
Small	2 1/2"-3"	5"-6"	7 1/2"-9"	JS-1004	JS-1003
Medium	3"-3 1/2"	6"-7"	9"-10 1/2"	JS-1006	JS-1005
Large	3 1/2"-4"	7-8"	10"-12"	JS-1008	JS-1007
X-Large	4"-4 1/2"	8"-9"	12"-14"	JS-1010	JS-1009