

**BioQuest Prosthetics**

- BioStride..... 2
- PerfectStride II-X3..... 2

**DYCOR Prosthetics**

- AFP (Original ADL foot)..... 3

**Endolite®**

- Echelon..... 4
- Elite Blade..... 5
- Elite VT..... 6
- Elite 2, Epirus..... 7
- Javelin..... 8
- Esprit, Multiflex..... 9
- Dynamic Response 2, Navigator..... 10
- Senior, Stellar..... 11

**Fillauer Companies**

- Element..... 12
- Element DS..... 13
- IBEX..... 14
- Wave Sport, Centri MCV Foot Shell..... 15

**Foresee Orthopedic**

- S.A.F.E..... 16

**Freedom Innovations®**

- Renegade..... 17
- Renegade LP..... 18
- Renegade MX..... 19
- Renegade LP-MX..... 20
- Thrive VS<sup>2</sup>..... 21
- Silhouette..... 22
- Silhouette LP..... 23
- Silhouette VS..... 24
- Silhouette LP-VS, Silhouette Adapters..... 25
- Runway..... 26
- Runway HX..... 27
- Promenade..... 28
- Sierra..... 29
- Highlander, Highlander Max..... 30
- Senator..... 31
- Pacifica..... 32
- Pacifica LP..... 33
- LP Symes..... 34
- WalkTek..... 35
- Freestyle Swim..... 36
- Slalom Ski, Nitro Running..... 37

**Kingsley Manufacturing**

- SACH..... 38-39
- SACH Lite..... 39-40
- SACH Lite, Symes..... 40
- STEN, Dynamic Response..... 41

**Motion Control**

- Motion Foot..... 42

**Ossur® Americas**

- PROPRIO FOOT with EVO..... 43
- Re-Flex Rotate with EVO..... 44
- Re-Flex Shock with EVO..... 45
- LP Ceterus..... 46
- Modular III..... 47
- Vari-Flex with EVO..... 48
- LP Vari-Flex with EVO..... 49
- Talux..... 50
- Flex-Foot Axia..... 51
- Flex-Foot Assure..... 52
- Sure-Flex..... 53
- Flex-Foot Balance, Elation..... 54
- Flex-Symes..... 55
- Chopart, Flex-Run..... 56

**Trulife**

- Seattle Catalyst 9..... 57
- Seattle Catalyst..... 58
- Seattle Cadence..... 59
- Seattle Turas..... 60
- Seattle LP..... 61
- Seattle LightFoot 2 with Titanium Pyramid..... 62
- Seattle LightFoot 2..... 63
- Seattle Carbon LightFoot..... 64
- Seattle LightFoot..... 65
- Seattle Lifecast..... 66
- Kinetic, Seattle Ladies..... 67
- Natural, Seattle SACH, Seattle Foamer..... 68
- Seattle Cosmesis..... 69
- Foot Components and Accessories..... 70

**WillowWood®**

- Duralite..... 71
- Fusion..... 72
- Pathfinder II, Pathfinder II LP..... 73
- Trailblazer, Trailblazer MA (Multi-axial)..... 74
- Impulse, Impulse Symes..... 75
- Magnum Heavy Duty, Carbon Copy 2..... 76
- Carbon Copy High Performance..... 77
- Earthwalk 2, SACH..... 78
- Single Axis, Components..... 79

**Shower Feet**

- Endolite Aqualimb..... 80
- KISS Dynamic Flexible Keel “Stubbie”..... 80
- Stomper, Stomper Jr..... 80

**BioQuest Prosthetics - BioStride Series**

- The BioStride Series includes three models designed to achieve multiplane and multi-axial motion for K3 and K4 amputees
- The cornerstone of the design is a multi-layer anterior facing, convexly curved carbon graphite composite calf shank
- Exclusive use of carbon composite and titanium components reduces weight and optimizes kinetic energy as the body's center of mass moves forward from mid stance to push off
- Weight rating up to 300 lbs
- Available sizes: 22-31 cm
- 36 month warranty; 45 day trial period



Description	Model No.
UltraStride (Tall)	<b>MX3</b>
ActiveStride (Tall)	<b>MX2</b>
NaturalStride (Short)	<b>MX1</b>

Description	Weight Rating		Clearance	
	K3	K4	Male	Female
UltraStride	300 lbs	275 lbs	7 <sup>5</sup> / <sub>16</sub> "	7 <sup>11</sup> / <sub>16</sub> "
ActiveStride	300 lbs	275 lbs	7 <sup>5</sup> / <sub>16</sub> "	7 <sup>11</sup> / <sub>16</sub> "
NaturalStride	320 lbs	295 lbs	5 <sup>7</sup> / <sub>8</sub> "	6 <sup>1</sup> / <sub>4</sub> "

**PerfectStride II-X3**

- Optimizes multi plane and multi-axial motion
- Anterior Facing Convexly Curved Titanium Spring Calf Shank
- Asymmetrical Carbon Graphite Foot Keel
- Aluminum Ankle Coupler; Titanium Fasteners
- Carbon Graphite Shank Support Spring
- Adjustable Heel Height Feature
- Male Or Female Pyramid Connectors (with or without transverse plane rotation)
- Available sizes: 22-31 cm
- Weight rating for K3 up to 300 lbs, K4 up to 285 lbs
- 36 month warranty; 45 day trial period



Description	Clearance	Model No.
Short	7.25"	<b>PX3-S</b>
Tall	8.25"	<b>PX3-T</b>

Replacement Items for PerfectStride and Biostride	
Description	Model No.
Light Flesh Tone Foot shell - Right	<b>RFS</b>
Light Flesh Tone Foot shell - Left	<b>LFS</b>
Dark Flesh Tone Foot shell - Right (Negroid)	<b>RFS-N</b>
Dark Flesh Tone Foot shell - Left (Negroid)	<b>LFS-N</b>
Spectra Sock	<b>SS-R</b>
Ankle/Shank Foam Cover	<b>CFC</b>

### DYCOR - AFP Flexible Keel Single Axis

- Originally known as the ADL (Active Daily Living) or Uniaxial
- Designed for individuals exhibiting functional levels 1 and 2
- Remarkable lightness due to design simplicity, coupled with thermoplastic internal components and an expanded polyethylene foot shell
- 12 month warranty

Model No. 1A



### AFP Flexible Multi-axial Keel

- Also known as Triaxial and FMA
- Designed specifically for the geriatric “community ambulator” patients
- Available in two versions, STM and AGM

**Ordering Example**  
 Side ..... L  
 Size ..... 26  
 Keel ..... M  
 Single Axis ..... 1A  
**L6M1A**

### STM (Soft Tissue Management)

- Designed for individuals exhibiting functional levels 1 and 2
- 12 month warranty

Model No. STM

**STM Ordering Example**  
 Side ..... L  
 Size ..... 26  
 Keel ..... M  
 Type ..... STM  
**L6MSTM**

### AGM (Aggressive Gait Management)

- Designed for functional levels 2 and 3, including extra-ambulatory activities such as golf
- 12 month warranty

Model No. AGM

**AGM Ordering Example**  
 Side ..... L  
 Size ..... 26  
 Keel ..... M  
 Type ..... AGM  
**L6MAGM**

Dycor AFP Foot					
cm	Shoe Size	Foot Size	Soft lbs	Medium lbs	Max lbs
22	5	22	105	124	147
23	6	23	105	124	147
24	7	24	115	136	160
25	8	25	125	148	174
26	9	26	135	160	189
27	10	27	145	173	203
28	11	28	156	185	218
29	12	29	156	185	218
Dycor Extra					
29	12	29	168	204	233
30	13	30	168	204	233



### Endolite - Echelon Foot

- Independent hydraulically controlled plantar and dorsiflexion
- Self leveling of heel height for varied terrain and heel heights
- Provides comfort while traversing varied terrain
- Hydraulic ankle allows for greater toe clearance in swing phase
- Controlled plantar flexion for stability and comfort
- K3 activity level
- Heel height: 3/8"
- Clearance: 4 1/2" - 4 7/8"
- Weight rating up to 275 lbs
- Foot shell and spectra sock included
- Available sizes: 22-30 cm
- 36 month warranty

Model No.

EC



#### Ordering Example

Model ..... **EC**  
 Size ..... **22**  
 Side ..... **L**  
 Spring Kit ..... **1**

**EC22L1**

add a D suffix for dark foot shell

**Echelon Foot Spring Selection Chart**

Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-275
<b>Impact Level</b>								
Walking and Occasional Sports (golf or hiking)	1	1	2	3	4	5	6*	7*
Frequent or Daily Sports (jogging)	1	2	3	4	5	6	7*	8*
Spring set recommendations are for transtibial amputees. For trans-femoral, select a spring set one level lower								
*Size 25-30 only								


### Elite Blade

- Maximum energy return
- Low profile option available
- Sliding and rotating proximal male pyramid for ease of alignment
- Simple lightweight design is easy to cosmetically finish
- Promotes dynamic tibial progression
- K3-K4 activity level
- Heel height: 3/8"
- Clearance: 14.37"
- Low Profile Clearance: 7.87"
- Sliding and rotating proximal male pyramid, foot shell, and spectra sock included
- Available sizes: 24-30 cm
- Weight rating up to 365 lbs
- 36 month warranty

Description	Clearance	Model No.
Elite Blade	14.37	<b>EB</b>
Low Profile	7.87	<b>EBLP</b>

**Ordering Example**

Model ..... **EB**  
 Size ..... **24**  
 Side ..... **L**  
 Spring Kit ..... **1**  
**EB24L1**  
 add a D suffix for dark foot shell



Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
<b>Impact Level</b>										
Walking/Occasional Sports	1	1	2	3	4	5	6	7	8	9
Frequent or Daily Sports (jogging)	1	2	3	4	5	6	7	8	9	—
Daily High Impact Sports (running)	2	3	4	5	6	7	8	9	special order	

**Elite VT**

- Dynamic precision engineered titanium spring incorporates vertical compression and axial rotation
- Progressive response to axial loading
- Truly independent heel and toe
- Compact, energy efficient design is easy to cosmetically finish
- K3-K4 activity level
- Heel height: 3/8"
- Clearance: 6 5/8"
- Weight rating up to 365 lbs
- Foot shell and spectra sock included
- Available sizes: 24-30 cm
- 36 month warranty

Model No. **EVT**



**Ordering Example**

Model ..... **EVT**  
 Size ..... **27**  
 Side ..... **L**  
 Spring Kit ..... **2**  
 Axial Sprint ..... **1**

**EVT27L21**  
 add a D suffix for dark foot shell

Axial Spring Kits	
Category	Axial Spring
1-3	1
4-5	2
6-7	3
8-9	4

*Standard prematched combinations will provide most active transtibial amputees a suitable range of vertical compression and axial rotation. You may achieve a softer or stiffer action by selecting one of the alternative springs*

Elite VT Transtibial Spring Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Impact Level										
Walking/Occasional Sports	1	1	2	3	4	5	6	7	8	9
Frequent or Daily Sports (jogging)	1	2	3	4	5	6	7	8	9	—
Daily High Impact Sports (running)	2	3	4	5	6	7	8	9	—	—
Spring set recommendations are for transtibial amputees. For trans-femoral, select a spring set one level lower										

### Elite2 Foot

- Enhanced E-carbon heel provides 20% more vertical shock absorption
- Improved glide sock for greater wear and durability
- 13 mm hex head bolts for easy exchanging of toe spring
- Low build height with light weight simple design that is easy to cosmetically finish
- K3-K4 activity level
- Heel height: 3/8"
- Clearance: 5 1/8"
- Weight rating up to 365 lbs
- Foot shell and spectra sock included
- Available sizes: 24-30 cm
- 36 month warranty



**Ordering Example**  
 Model ..... **EL**  
 Size..... **24**  
 Side ..... **L**  
 Ankle..... **1**  
**EL24L1**  
 add a D suffix for dark foot shell

Model No. **EL**

Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
<b>Impact Level</b>										
Walking/Occasional Sports	1	1	2	3	4	5	6	7	8	9
Frequent or Daily Sports (jogging)	1	2	3	4	5	6	7	8	9	—
Daily High Impact Sports (running)	2	3	4	5	6	7	8	9	special order	

Spring set recommendations are for transtibial amputees. For trans-femoral, select a spring set one level lower

### Epirus Foot

- Multi-axial rotation with tripod stability system
- Low profile energy storing and return design
- Biomechanically engineered to provide more natural ankle motion
- K3 activity level
- Heel height: 3/8"
- Clearance: 3 5/16"
- Weight rating up to 275 lbs
- Foot shell and spectra sock included
- Available sizes: 22-30 cm
- 36 month warranty



**Ordering Example**  
 Model ..... **EP**  
 Size..... **22**  
 Side..... **L**  
 Spring Kit ..... **1**  
**EP22L1**  
 add a D suffix for dark foot shell

Model No. **EP**

Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-275
<b>Impact Level</b>								
Walking and Occasional Sports ( golf or hiking)	1	1	2	3	4	5	6	7
Frequent or Daily Sports (jogging)	1	2	3	4	5	6	7	8

Spring set recommendations are for transtibial amputees. For trans-femoral, select a spring set one level lower

\*Size 25-30 only

**Javelin**

- New foot design has been streamlined to focus on the essential mechanisms which govern walking, flexibility, stability and terrain compliance
- Slimline continuous shin and toe spring flexes for greater energy return and reduced strain on the back, knee, and leg muscles
- e-Carbon lay up; easy cosmetic finishing of minimalist profile
- Fitting jig allows accurate assessment of limb performance and alignment before final cutting and drilling
- Full length or low profile available
- K3 activity level
- Heel height: 3/8"
- Clearance: 14.37
- Low Profile Clearance: 7.87
- Foot shell, spectra sock and cutting jig included
- Available sizes: 22-30 cm
- Weight rating up to 365 lbs
- 36 month warranty

Description	Clearance	Model No.
Javelin	14.37	<b>JAV</b>
Low Profile	7.87	<b>JAVLP</b>



**Ordering Example**  
 Model ..... **JAV**  
 Size ..... **22**  
 Side ..... **L**  
 Spring Kit ..... **1**  
**JAV22L1**  
 add a D suffix for dark foot shell

Javelin Foot Spring Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Impact Level										
Moderate Activity Level (K3)	1	1	2	3	4	5	6	7	8*	9*

\*Size 25-30 only



### Esprit Foot

- Low profile energy storing and return design
- Tripod stability system
- e-Carbon lay up
- Truly independent heel & toe
- Easy to cosmetically finish
- K3 activity level
- Heel height: 3/8"
- Clearance: 2<sup>5</sup>/<sub>8</sub>" - 2<sup>7</sup>/<sub>8</sub>"
- Weight rating up to 275 lbs
- Foot shell and spectra sock included
- Available sizes: 22-30 cm
- 36 month warranty



Model No. **ESP**

†Refer to the Spring Selection Chart on Page 7 for the Esprit and Epirus feet.

**Ordering Example**

Model ..... **ESP**  
 Size ..... **26**  
 Side ..... **R**  
 Spring Kit† ..... **4**  
**ESP26R4**  
 add a D suffix for dark foot shell

### Multiflex Foot

- Multi-axial ankle motion; available in standard or slim widths
- K2 activity level
- Heel height: 1/4" - 1 1/2"; Clearance: 1 3/4"
- Weight rating up to 275 lbs
- Available sizes: 22-30 cm
- 24 month warranty



\*Ankle and Snubber must be ordered with foot (see chart below)

Size	Standard Width Foot		Slim Width Foot	
	Left	Right	Left	Right
	(20 mm heel height)		(25 mm heel height)	
22 cm	<b>519119</b>	<b>519120</b>	<b>519141</b>	<b>519142</b>
23 cm	<b>519121</b>	<b>519122</b>	<b>519143</b>	<b>519144</b>
24 cm	<b>519123</b>	<b>519124</b>	<b>519145</b>	<b>519146</b>
25 cm	<b>519125</b>	<b>519126</b>	<b>519147</b>	<b>519148</b>
26 cm	<b>519127</b>	<b>519128</b>	—	—
27 cm	<b>519129</b>	<b>519130</b>	—	—
28 cm	<b>519131</b>	<b>519132</b>	—	—
29 cm	<b>519133</b>	<b>519134</b>	—	—
30 cm	<b>519135</b>	<b>519136</b>	—	—

**K2 Ankle Selection Chart\***

Weight/ lbs	up to 132	133-176	177-220	221-275
	Ankle/ Snubber	W	W	O
	W	O	O	B

**Standard Snubber Kit for DR2 and Multiflex**

Description	Model No.
White 60	<b>409007</b>
Orange 70	<b>409008</b>
Blue 80	<b>409009</b>

Ankle	30 mm Clamp	Pyramid	30 mm Pylon
White (W) 60	<b>379342</b>	<b>409353</b>	<b>409241</b>
Orange (O) 70	<b>379343</b>	<b>409354</b>	<b>409242</b>
Blue (B) 80	<b>379344</b>	<b>409355</b>	<b>409243</b>

**Dynamic Response 2 Foot**

- Amputee adjustable response
- Multi-axial ankle motion; full length toe lever
- Heel height adjustment based on ankle position on foot
- K3 activity level
- Heel height: 1/4" - 1 3/8"; Clearance: 2 3/8"
- Weight rating up to 275 lbs
- Available sizes: 22-30 cm
- 24 month warranty



\*Ankle and Snubber must be ordered with foot (see chart below)

Size	Model No.	
	Left	Right
22 cm	<b>529120</b>	<b>529121</b>
23 cm	<b>529122</b>	<b>529123</b>
24 cm	<b>529124</b>	<b>529125</b>
25 cm	<b>529126</b>	<b>529127</b>
26 cm	<b>529128</b>	<b>529129</b>
27 cm	<b>529130</b>	<b>529131</b>
28 cm	<b>529132</b>	<b>529133</b>
29 cm	<b>529134</b>	<b>529135</b>
30 cm	<b>529136</b>	<b>529137</b>

K3 Ankle Selection Chart*				
Weight/ lbs	up to 132	133-176	177-220	221-275
	Ankle/ Snubber	W	O	O
W		O	B	B

Ankle	30 mm Clamp	Pyramid	30 mm Pylon
White (W) 60	<b>379342</b>	<b>409353</b>	<b>409241</b>
Orange (O) 70	<b>379343</b>	<b>409354</b>	<b>409242</b>
Blue (B) 80	<b>379344</b>	<b>409355</b>	<b>409243</b>

**Navigator Foot**

- Multi-axial with flexible keel
- Durable replaceable foot shell
- Integrated male pyramid
- Full length toe lever; Sandal toe
- K2 activity level
- Heel height: 3/8"
- Clearance: 3 5/16"
- Weight rating up to 275 lbs
- Available sizes: 22-30 cm
- 24 month warranty



Model No. **NAV**

Ankle Selection Chart						
Weight/lbs	100-130	131-150	151-170	171-195	196-220	221-275
Limited Indoor*	WW	WW	WO	OO	OB	BB
Limited Outdoor**	WW	WO	OO	OB	BB	BB
*The amputee is able to walk on level surfaces at slow cadence						
**The amputee is able to walk on uneven surfaces at low to medium speed and can manage small obstacles						
Note: When making an ankle selection for trans-femoral amputee, consider reducing the ball and snubber stiffness by one activity level to enhance the smoothness and ease of rollover						

**Ordering Example**  
 Model ..... **NAV**  
 Size ..... **22**  
 Side ..... **L**  
 Ankle ..... **WW**  
**NAV22LWW**  
 add a D suffix for dark foot shell

### Senior Foot

- Cushioned heel design combines new materials and molding techniques, creating a foot that is comfortable, lightweight and simple to use
- A slim line keel and a detachable foot cosmesis allows easy replacement of the cosmesis
- An integral titanium and alloy pyramid provides alignment while keeping weight to a minimum
- 19 mm heel height (plus or minus 7° with angulation of pyramid)
- K1-K2 activity level
- Heel height: 3/4"
- Clearance: 2<sup>5</sup>/<sub>8</sub>-3"
- Weight rating up to 275 lbs
- Available sizes: 22-28 cm
- Foot shell included
- 24 month warranty



Size	Foot Weight	Model No.	
		Left	Right
22 cm	10 oz	<a href="#">SNR-22L</a>	<a href="#">SNR-22R</a>
23 cm	11 oz	<a href="#">SNR-23L</a>	<a href="#">SNR-23R</a>
24 cm	11.5 oz	<a href="#">SNR-24L</a>	<a href="#">SNR-24R</a>
25 cm	12 oz	<a href="#">SNR-25L</a>	<a href="#">SNR-25R</a>
26 cm	13 oz	<a href="#">SNR-26L</a>	<a href="#">SNR-26R</a>
27 cm	14 oz	<a href="#">SNR-27L</a>	<a href="#">SNR-27R</a>
28 cm	15 oz	<a href="#">SNR-28L</a>	<a href="#">SNR-28R</a>

\*Add the letter D at the end of part number for Brown foot shell

### Stellar Foot

- Adjustable heel cushion with Integral pyramid
- Sandal toe
- K1-K2 activity level
- Heel height: 3/4"
- Clearance: 2<sup>5</sup>/<sub>8</sub>-3"
- Weight rating up to 275 lbs
- Available sizes: 22-28 cm
- Foot shell included
- 24 month warranty



Size	Model No.	
	Left	Right
22 cm	<a href="#">529101</a>	<a href="#">529102</a>
23 cm	<a href="#">529103</a>	<a href="#">529104</a>
24 cm	<a href="#">529105</a>	<a href="#">529106</a>
25 cm	<a href="#">529107</a>	<a href="#">529108</a>
26 cm	<a href="#">529109</a>	<a href="#">529110</a>
27 cm	<a href="#">529111</a>	<a href="#">529112</a>
28 cm	<a href="#">529113</a>	<a href="#">529114</a>

**Fillauer - Element™ Foot System**

- Unique 3rd spring Element technology allows smooth rollover from heel to toe
- Smooth rollover and superior stability
- Split heel for enhanced multi-axial function
- Very low profile - only 3.75" clearance (9.5 cm)
- Foot weight: 340 g
- K3-K4 activity level
- Heel height: 3/8"
- Weight rating up to 275 lbs
- Available sizes: 22-30 cm
- 3 year warranty/30 day patient trial
- Foot shell sold separately Centri MCV (6 month warranty)



Description	Model No.
Element Foot	<b>100-10</b>
MCV Foot Shell (sold separately)	<b>??13033</b>
*replace 1st ? with side and 2nd ? with length - refer to ordering example	

**Ordering Example**  
 Model ..... **100-10**  
 Size ..... **24**  
 Category ..... **B43**  
**100-10-24-B43**

**Foot Shell Ordering Example**  
 Side ..... **45**  
 Foot Length ..... **24**  
 Foot Shell ..... **13033**  
**452413033**

*45=Left; 46=Right  
 Length: 22-30 cm; 6 month warranty*

Category Selection Guide for Element Foot System									
Size	22	23	24	25	26	27	28	29	30
100-119 lbs	A11	A11	A11	A32	A32	—	—	—	—
120-139 lbs	A11	A11	A11	A32	A32	A43	A43	—	—
140-159 lbs	A11	A11	A11	A32	A32	A43	A43	C64	C64
160-179 lbs	A22	A22	A22	A43	A43	A43	A43	C64	C64
180-199 lbs	A22	A22	A22	B43	B43	B54	B54	C64	C64
200-224 lbs	A32	A32	A32	B53	B53	B54	B54	D75	D75
225-249 lbs	—	—	—	C54	C54	C65	C65	D75	D75
250-275 lbs	—	—	—	C64	C64	C65	C65	D75	D75

### Element™ DS (with DuraShock Kit)

- Smooth rollover and superior stability
- Vertical shock and 60° of axial rotation
- K3-K4 activity level
- Heel height: 3/8"
- Clearance: 6.5"
- Foot weight: 632 g
- Weight rating up to 275 lbs
- Available sizes: 22-30 cm
- 3 year warranty/30 day patient trial
- Durashock kit colors available: Black (K), Red (R), and Blue (B)
- Foot shell sold separately Centri MCV (6 month warranty)



Description	Model No.
Element DS Foot	<b>105-10</b>
MCV Foot Shell (sold separately)	<b>??13033</b>
*replace 1st ? with side and 2nd ? with length - refer to ordering example	

**Ordering Example**

Model ..... **105-10**  
 Size ..... **27**  
 Category ..... **C65**  
 Kit Color ..... **K**  
**105-10-27-C65K**

**Foot Shell Ordering Example**

Side ..... **45**  
 Foot Length ..... **27**  
 Foot Shell ..... **13033**  
**452713033**

*45=Left; 46=Right  
 Length: 22-30 cm; 6 month warranty*

Category Selection Guide for Element Foot System with DuraShock Kit									
Size	22	23	24	25	26	27	28	29	30
110-119 lbs	A11R	A11R	A11R	A32R	A32R	—	—	—	—
120-139 lbs	A11R	A11R	A11R	A32R	A32B	A43B	A43B	—	—
140-159 lbs	A11R	A11R	A11R	A32B	A32B	A43B	A43K	C64K	C64K
160-179 lbs	A22R	A22B	A22B	A43B	A43B	A43K	A43K	C64K	C64K
180-199 lbs	A22B	A22B	A22B	B43B	B43K	B54K	B54K	C64K	C64K
200-224 lbs	A32B	A32K	A32K	B53K	B53K	B54K	B54K	D75K	D75K
225-249 lbs	—	—	—	C54K	C54K	C65K	C65K	D75K	D75K
250-275 lbs	—	—	—	C64K	C64K	C65K	C65K	D75K	D75K

**Ibex™ Foot System**

- Innovative Micro-slice pylon
- Multi-axis heel/toe for terrain conformance
- Energy optimized pylon/heel plate design
- Weight rating up to 330 lbs
- Foot weight: 16.2 oz (459 g); build height 6.75" (17 cm)
- Ibex XD weight rating up to 440 lbs
- Ibex XD foot weight: 21.2 oz (602 g); build height 7.63" (19 cm)
- K3-K4 activity level
- Heel height: 3/8"
- Available sizes: 22-30 cm; Ibex XD sizes: 25-30 cm
- 36 month warranty/30 day patient trial
- Foot shell sold separately Centri MCV (6 month warranty)



Description	Model No.
IBEX Foot	<b>120-10</b>
MCV Foot Shell (sold separately)	<b>??13033</b>
*replace 1st ? with side and 2nd ? with length - refer to ordering example	

**Ordering Example**  
 Model ..... **120-10**  
 Size ..... **22**  
 Category ..... **A2**  
**120-10-22-A2**

**Foot Shell Ordering Example**  
 Side ..... **45**  
 Foot Length ..... **22**  
 Foot Shell ..... **13033**  
**452213033**

*45=Left; 46=Right  
 Length: 22-30 cm; 6 month warranty*

Category Selection Guide for Ibex Foot System									
Size	22	23	24	25	26	27	28	29	30
110-119 lbs	A2	A2	A2	B4	B4	—	—	—	—
120-139 lbs	A2	A2	A2	B4	B4	C5	C5	—	—
140-159 lbs	B3	B3	B3	B5	B5	C5	C5	D6	D6
160-179 lbs	B4	B4	B4	B5	B5	C6	C6	D6	D6
180-199 lbs	B4	B4	B4	C6	C6	C6	C6	D7	D7
200-224 lbs	—	B4	B4	C6	C6	D7	D7	D7	D7
225-249 lbs	—	—	—	D7	D7	D7	D7	E8	E8
250-275 lbs	—	—	—	D7	D7	E8	E8	E8	E8
276-300 lbs	—	—	—	E8	E8	E8	E8	E9	E9
301-330 lbs	—	—	—	E8	E8	E9	E9	E9	E9
330-440 lbs	—	—	—	E8XD	E8XD	E9XD	E9XD	E9XD	E9XD

### Wave Sport Foot System™

- Unique wave-spring technology
- Excellent shock absorption
- Maintenance free design
- Full length heel gives superior stability and heel dynamics
- Carbon “C-shaped” pylon provides greater dynamic response
- Compatible with standard athletic shoes
- Light weight carbon fiber construction (550 g for 27 cm foot)
- Moderate to very high activity levels
- Standard Sport Foot build height: 9.5” (24 cm)
- Low Profile build height: 8” (20 cm)
- Weight rating up to 275 lbs
- Available sizes: 22-30 cm
- 36 month warranty/30 day patient trial



Description	Model No.
Wave Sport Foot	<b>127-10-SH-CSF</b>
Low Profile Wave Sport Foot	<b>127-10-LP-CSF</b>
MCV Foot Shell (sold separately)	<b>??13033</b>
*replace 1st ? with side and 2nd ? with length - refer to ordering example	

### Centri MCV Foot Shell

- Custom colors available as special order with 3 weeks lead time
- Shell selection requires Left or Right and Foot Length

	Foot Length (cm)									
	22	23	24	25	26	27	28	29	30	
Left	45 22	45 23	45 24	45 25	45 26	45 27	45 28	45 29	45 30	
Right	46 22	46 23	46 24	46 25	46 26	46 27	46 28	46 29	46 30	

**Foot Shell Ordering Example**  
 L (45) or R (46). . . . . **45**  
 Foot Length . . . . . **24**  
 Foot Shell . . . . . **13033**  
**452413033**  
*Indicate Right or Left when ordering*

**Foresee Orthopedic - S.A.F.E.  
(Stationary Attachment Flexible Endoskeletal)**

- Absorbs energy for smooth rollover throughout its length during flexing of the foot
- Keel compresses as the foot bends
- Keel releases its energy as the body moves forward
- Weight rating up to 250 lbs for all S.A.F.E. feet

**S.A.F.E. I Without Toes**

- Smoother roll over to toe
- Available sizes: <sup>3</sup>/<sub>4</sub>" heel height (22-30 cm)

**S.A.F.E. II With Toes<sup>®</sup>**

- Higher energy storage for release at toe-off
- Available sizes: <sup>1</sup>/<sub>4</sub>" heel height (14-17 cm; 24-30 cm wide)  
- <sup>3</sup>/<sub>8</sub>" heel height (24-30 cm wide), (18-29 cm narrow)  
- <sup>3</sup>/<sub>4</sub>" heel height (22-30 cm wide)
- Carbon keel available <sup>3</sup>/<sub>8</sub>" heel 22-30 cm - narrow only
- Split toe 22-30 cm - narrow only
- 12 month warranty on all feet



**Ordering Chart for S.A.F.E. I, S.A.F.E. II, Syme and Pyramid**

Style	Type	Keel	Heel Height	Heel Density	Width	Size cm	Left/Right	Color
S1=SAFE I	S=Standard	O=Orange	38= <sup>3</sup> / <sub>8</sub> "	S=Soft	W=Wide	14-30	L=Left	1=Caucasian
S2=SAFE II	Y=Syme	Standard-Moderate Activity	34= <sup>3</sup> / <sub>4</sub> "	M=Medium	N=Narrow		R=Right	2=Light Negroid
	W=Waterproof	R=Red	14= <sup>1</sup> / <sub>4</sub> "	F=Firm				3=Medium Negroid
	H=Heavy Duty	Firm-High Activity						4=Dark Negroid
	V=Waterproof Heavy Duty	Y=Yellow						5=Light Latin
	Z=Syme Waterproof	Flexible Low Activity						6=Medium Latin
	C=Carbon Keel							7=Dark Latin
	VE=Vector							

\*Some styles not available in all sizes

**S.A.F.E. Natural<sup>®</sup>**

- Natural cosmetic design, yet easily fits in a shoe
- Stores and releases energy throughout its full length
- Same qualities as the original S.A.F.E. II feet
- Heel Height: 38 = <sup>3</sup>/<sub>8</sub>"
- Width: Wide only
- Available sizes: 22-30 cm
- Choose Left or Right
- 12 month warranty



**S.A.F.E. Ordering Example**

Style ..... **SN**  
 Type ..... **S**  
 Keel ..... **O**  
 Heel Height ..... **38**  
 Heel Density ..... **M**  
 Width ..... **W**  
 Size ..... **22**  
 Side ..... **R**  
 Color ..... **1**  
**SNSO38MW22R1**

**Ordering Chart for S.A.F.E. Natural and Split Toe**

Style	Type	Keel	Heel Density	Color
SN=Natural	S=Standard	O=Orange	S=Soft	1=Caucasian
ST=Split Toe	Y=Syme	Standard-Moderate Activity	M=Medium	2=Light Negroid
	W=Waterproof	R=Red	F=Firm	3=Medium Negroid
	H=Heavy Duty	Firm-High Activity		4=Dark Negroid
	V=Waterproof Heavy Duty	Y=Yellow		5=Light Latin
	Z=Syme Waterproof	Flexible Low Activity		6=Medium Latin
				7=Dark Latin



### Freedom Innovations - Renegade®

- Patented Z-Shock Technology achieves unprecedented efficiency, allowing active tibial progression
- Ultra-lightweight design is virtually maintenance free
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Heel Height: 3/8"
- 36 month warranty on foot module; 6 months on foot shell



Description	Model No.
Renegade Foot	<b>RS1</b>
Foot Shell (default)	<b>FTC-3M-0</b>
Foot Shell (default sandal toe)	<b>FTC-2M-1</b>
High Profile Foot Shell	<b>FTC-2H-1</b>

**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M), High (H) or Sandal toe (S)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Toe Regular (R) or Sandal (S)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model .....	<b>FTC</b>
Generation .....	<b>3</b>
Profile .....	<b>M</b>
Cap .....	<b>0</b>
Size .....	<b>23</b>
Skin Tone .....	<b>L</b>
	<b>4</b>
Toe .....	<b>R</b>
Side .....	<b>L</b>
<b>FTC-3M-023L4-RL</b>	
*Generation 3 not available in High Profile	
*Cap option only available for Generation 3	
*Sandal toe only available in sizes 22-28	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	6 <sup>1</sup> / <sub>8</sub> /154	22/625
26-28	6 <sup>3</sup> / <sub>4</sub> /172	27/765
29-31	7 <sup>1</sup> / <sub>2</sub> /189	33/945

Renegade Stiffness Category Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	—
High Impact Level	2	3	4	5	6	7	8	9	—	—

**Renegade® LP**

- Delivers the same features as the Renegade, but designed to accommodate clearances under 5"
- Designed for K3 and K4 amputees
- The lowest profile vertical shock foot
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Heel Height: 3/8"
- 36 month warranty on foot module; 6 months on foot shell



Description	Model No.
Renegade LP Foot	<b>RS6</b>
Foot Shell (default)	<b>FTC-3M-0</b>
Foot Shell (default sandal toe)	<b>FTC-2M-1</b>
High Profile Foot Shell	<b>FTC-2H-1</b>

**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M), High (H) or Sandal toe (S)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Toe Regular (R) or Sandal (S)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model . . . . .	<b>FTC</b>
Generation . . . . .	<b>3</b>
Profile . . . . .	<b>M</b>
Cap. . . . .	<b>0</b>
Size . . . . .	<b>23</b>
Skin Tone . . . . .	<b>L</b>
	<b>4</b>
Toe . . . . .	<b>R</b>
Side . . . . .	<b>L</b>
<b>FTC-3M-023L4-RL</b>	
<i>*Generation 3 not available in High Profile</i>	
<i>*Cap option only available for Generation 3</i>	
<i>*Sandal toe only available in sizes 22-28</i>	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	4 <sup>3</sup> / <sub>4</sub> /120	22/625
26-28	5 <sup>1</sup> / <sub>8</sub> /131	27/760
29-31	5 <sup>5</sup> / <sub>8</sub> /143	33/930

Renegade LP Stiffness Category Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	—
High Impact Level	2	3	4	5	6	7	8	9	—	—

## Renegade® MX

- Utilizing “Z Shock Technology”, the Renegade MX provides a vertical shock device that keeps users stable and secure during a broader spectrum of activities
- 8° of transverse rotation; 8° of inversion/eversion; 6° of plantar flexion; 7.5° of dorsiflexion
- Accommodates clearance of 7.2”
- User Weight Rating: 285 lbs (130 kg)
- Stiffness Categories: 1-7
- Connector: Integral male pyramid
- Heel Height: 3/8”
- 36 month warranty on foot module, 24 months on rotator and 6 months on foot shell



Description	Model No.
Renegade MX Foot	<b>RS1-M1</b>
Foot Shell (default)	<b>FTC-3M-0</b>
Foot Shell (default sandal toe)	<b>FTC-2M-1</b>
High Profile Foot Shell	<b>FTC-2H-1</b>

### Please provide the following when ordering a foot shell:

- Model number
- Generation 2 or 3
- Profile: Mid (M), High (H) or Sandal toe (S)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Toe Regular (R) or Sandal (S)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

**Foot Shell Ordering Example**

Model . . . . . **FTC**  
 Generation . . . . . **3**  
 Profile . . . . . **M**  
 Cap. . . . . **0**  
 Size . . . . . **23**  
 Skin Tone . . . . . **L**  
 Toe . . . . . **R**  
 Side . . . . . **L**

**FTC-3M-023L4-RL**

\*Generation 3 not available in High Profile  
 \*Cap option only available for Generation 3  
 \*Sandal toe only available in sizes 22-28

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	6 <sup>5</sup> / <sub>8</sub> /170	36/1020
26-28	7 <sup>3</sup> / <sub>8</sub> /188	42/1195
29-31	8 <sup>1</sup> / <sub>8</sub> /206	46/1300

Renegade MX Stiffness Category Selection Chart								
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285
Low Impact Level	1	1	2	3	4	5	6	7
Moderate Impact Level	1	2	3	4	5	6	7	—
High Impact Level	2	3	4	5	6	7	—	—

**Renegade® LP MX**

- Provides the same robust features as the Renegade MX in low profile version
- User Weight Rating: 285 lbs (130 kg)
- Heel Height: 3/8"
- 36 month warranty on foot module, 24 months on rotator and 6 months on foot shell

Description	Model No.
Renegade LP Foot	<b>RS6-M1</b>
Foot Shell (default)	<b>FTC-3M-0</b>
Foot Shell (default sandal toe)	<b>FTC-2M-1</b>
High Profile Foot Shell	<b>FTC-2H-1</b>



**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M), High (H) or Sandal toe (S)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Toe Regular (R) or Sandal (S)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model . . . . .	<b>FTC</b>
Generation . . . . .	<b>3</b>
Profile . . . . .	<b>M</b>
Cap. . . . .	<b>0</b>
Size . . . . .	<b>23</b>
Skin Tone . . . . .	<b>L</b>
	<b>4</b>
Toe . . . . .	<b>R</b>
Side . . . . .	<b>L</b>
<b>FTC-3M-023L4-RL</b>	
*Generation 3 not available in High Profile	
*Cap option only available for Generation 3	
*Sandal toe only available in sizes 22-28	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	5 <sup>5</sup> / <sub>8</sub> /143	34/970
26-28	6/153	39/1100
29-31	6 <sup>1</sup> / <sub>2</sub> /166	45/1275

Renegade LP-MX Stiffness Category Selection Chart								
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285
Low Impact Level	1	1	2	3	4	5	6	7
Moderate Impact Level	1	2	3	4	5	6	7	—
High Impact Level	2	3	4	5	6	7	—	—

**Thrive VS<sup>2</sup> Foot**

- World’s first load-activated carbon fiber prosthetic foot
- Accommodates added weight when lifting or carrying heavy objects
- Utilizes a synergistic, dual-keel design
- Low to moderate impact K3 level ambulators
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Connector: Male pyramid
- Heel Height: 3/8”
- 36 month warranty on foot module; 6 months on foot shell



Description	Model No.
Thrive Foot	<b>FS5</b>
Foot Shell (default)	<b>FTC-3M-0</b>
High Profile Foot Shell	<b>FTC-2H-1</b>

**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

**Foot Shell Ordering Example**

Model ..... **FTC**  
 Generation ..... **3**  
 Profile ..... **M**  
 Cap ..... **0**  
 Size ..... **25**  
 Skin Tone ..... **L**  
                                   **4-R**  
 Side ..... **R**

**FTC-3M-025L4-RR**

\*Generation 3 not available in High Profile  
 \*Cap option only available for Generation 3

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	5 <sup>3</sup> / <sub>4</sub> /145	22/630
26-28	6 <sup>1</sup> / <sub>5</sub> /157	27/760
29-31	7/175	31/870

Thrive Stiffness Category Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	—
High Impact Level	2	3	4	5	6	7	8	9	—	—

**Silhouette®**

- Slim profile and advanced materials facilitate easy and refined cosmesis
- Ideal for low to moderate impact K3 amputees
- Strong, lightweight carbon-fiber construction allows user to walk farther and faster, using less energy
- Multi-axial rotation delivers extended range of motion and exceptional stability
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Connector: Male or female pyramid or Versa Clamp
- Heel Height:  $\frac{3}{8}$ " and  $\frac{3}{4}$ "
- 36 month warranty on foot module; 6 months on foot shell

Description	Model No.	
	3/8" Heel	3/4" Heel
Female Tube Clamp	<b>RS4-01</b>	<b>RS8-01</b>
Male Pyramid	<b>RS4-02</b>	<b>RS8-02</b>
Versa Clamp	<b>RS4-03</b>	<b>RS8-03</b>
Foot Shell (default)	<b>FTC-3M-0</b>	
High Profile Foot Shell	<b>FTC-2H-1</b>	
$\frac{3}{8}$ " heel height - Average Clearance 7.2"		
$\frac{3}{4}$ " heel height - Average Clearance 7.7"		



**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model . . . . .	<b>FTC</b>
Generation . . . . .	<b>3</b>
Profile . . . . .	<b>M</b>
Cap . . . . .	<b>0</b>
Size . . . . .	<b>25</b>
Skin Tone . . . . .	<b>L</b>
	<b>4-R</b>
Side . . . . .	<b>R</b>
<b>FTC-3M-025L4-RR</b>	
<i>*Generation 3 not available in High Profile</i>	
<i>*Cap option only available for Generation 3</i>	

Build Height and Weight with Female Adapter		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	7/178	19/525
26-28	7 $\frac{1}{2}$ /191	23/645
29-31	7 $\frac{7}{8}$ /201	28/785

Silhouette Stiffness Category Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	—
High Impact Level	2	3	4	5	6	7	8	9	—	—

## Silhouette® LP

- Provides same multi-axial rotation design of Silhouette in lower profile with integral male pyramid for easy bolt-on assembly
- Ideal for low to moderate impact K3 amputees
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Connector: Male or female pyramid or Versa Clamp
- Heel Height: 3/8" and 3/4"
- 36 month warranty on foot module; 6 months on foot shell

Description	Model No.	
	3/8" Heel	3/4" Heel
Female Tube Clamp	<b>RS5-01</b>	<b>RS9-01</b>
Male Pyramid	<b>RS5-02</b>	<b>RS9-02</b>
Versa Clamp	<b>RS5-03</b>	<b>RS9-03</b>
Foot Shell (default)	<b>FTC-3M-0</b>	
High Profile Foot Shell	<b>FTC-2H-1</b>	
3/8" heel height - Average Clearance 7.2"		
3/4" heel height - Average Clearance 7.7"		



### Please provide the following when ordering a foot shell:

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model . . . . .	<b>FTC</b>
Generation . . . . .	<b>3</b>
Profile . . . . .	<b>M</b>
Cap . . . . .	<b>0</b>
Size . . . . .	<b>25</b>
Skin Tone . . . . .	<b>L</b>
	<b>4-R</b>
Side . . . . .	<b>R</b>
<b>FTC-3M-025L4-RR</b>	
*Generation 3 not available in High Profile	
*Cap option only available for Generation 3	

Build Height and Weight with Male Adapter		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	6 <sup>5</sup> / <sub>8</sub> /168	15/430
26-28	7 <sup>1</sup> / <sub>8</sub> /181	19/550
29-31	7 <sup>1</sup> / <sub>2</sub> /191	24/685

Silhouette LP Stiffness Category Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	—
High Impact Level	2	3	4	5	6	7	8	9	—	—

**Silhouette® VS**

- Vertical shock absorption softens impact at heel strike to provide comfort and preserve residual limb health
- Optimized sole plate with tri-durometer urethane strips transfers load from lateral to medial plane with each step for more natural medial roll-over
- Low to moderate impact K3 level ambulators
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Heel Heights: 3/8" or 3/4"
- 36 month warranty on foot module; 6 months on foot shell

Description	Model No.	
	3/8" Heel	3/4" Heel
Female Tube Clamp	<b>R14-01</b>	<b>R18-01</b>
Male Pyramid	<b>R14-02</b>	<b>R18-02</b>
Versa Clamp	<b>R14-03</b>	<b>R18-03</b>
Foot Shell (default)	<b>FTC-3M-0</b>	
High Profile Foot Shell	<b>FTC-2H-1</b>	
3/8" heel height - Average Clearance 7.2"		
3/4" heel height - Average Clearance 7.7"		



**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Build Height and Weight with Female Adapter		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	7 1/8/178	19/525
26-28	7 1/2/191	23/645
29-31	7 7/8/201	28/785

**Foot Shell Ordering Example**

Model ..... **FTC**  
 Generation ..... **3**  
 Profile ..... **M**  
 Cap. .... **0**  
 Size ..... **25**  
 Skin Tone..... **L**  
 Side ..... **4-R**  
 Side ..... **R**

**FTC-3M-025L4-RR**

*\*Generation 3 not available in High Profile*  
*\*Cap option only available for Generation 3*

Silhouette VS and Silhouette LP-VS Stiffness Category Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	—
High Impact Level	2	3	4	5	6	7	8	9	—	—



### Silhouette® LP-VS

- Pre-cut LP design is pre-assembled with either a female or male adapter to simplify the fitting process
- Low to moderate impact K3 level ambulators
- User Weight Rating: 365 lbs (166 kg)
- Heel Heights: 3/8" or 3/4"
- 36 month warranty on foot module; 6 months on foot shell

Description	Model No.	
	3/8" Heel	3/4" Heel
Female Tube Clamp	<b>R15-01</b>	<b>R19-01</b>
Male Pyramid	<b>R15-02</b>	<b>R19-02</b>
Foot Shell (default)	<b>FTC-3M-0</b>	
High Profile Foot Shell	<b>FTC-2H-1</b>	
3/8" heel height - Average Clearance 7.2"		
3/4" heel height - Average Clearance 7.7"		



#### Please provide the following when ordering a foot shell:

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

**Foot Shell Ordering Example**

Model . . . . . **FTC**

Generation . . . . . **2**

Profile . . . . . **M**

Cap. . . . . **0**

Size . . . . . **25**

Skin Tone. . . . . **L**

Side . . . . . **4-R**

Side . . . . . **R**

**FTC-2M-025L4-RR**

\*Generation 3 not available in High Profile

\*Cap option only available for Generation 3

Build Height and Weight with Male Adapter		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	6 <sup>5</sup> / <sub>8</sub> /168	15/430
26-28	7 <sup>1</sup> / <sub>8</sub> /181	19/550
29-31	7 <sup>1</sup> / <sub>2</sub> /191	24/685

### Versa-Clamp Adapter for Silhouette™

- Allows use with pin locking system or provides easy fix when a Silhouette pylon has been cut too short
- Extended 34 mm tube clamp can be paired with pre-cut spacer ring to build back pylon height
- Weight rated up to 365 lbs
- 24 month warranty

Model No. **KIT-00-123V0-00**



### Knee Pyramid Adapter for Silhouette™

- Provides additional fitting options for lower extremity amputees
- Compatible with all prosthetic knees that accept a 30 mm pylon
- Warranty: 24 months for user weights up to 300 lbs


Model No. **ACC-00-19500-00**



**Runway® Adjustable Prosthetic Foot**

- User-adjustable heel height foot designed for active K3 amputees
- Allows the amputee to elevate the heel up to 2" with the push of a button
- Designed for active K3 amputees
- User Weight Rating: 255 lbs (116 kg)
- Stiffness Categories: 1-6
- Connector: Integral male pyramid
- Heel Height: 3/8"
- 36 month warranty on foot module; 18 months on ankle and 6 months on foot shell

Description	Model No.
Runway Foot	<b>RS2</b>
Foot Shell (default sandal toe)	<b>FTC-2M-1</b>
Foot Shell (default regular toe)	<b>FTC-3M-1</b>



Runway heel height adjusts via push button

**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M), High (H) or Sandal toe (S)
- Cap: Yes (0) or No (1)
- Available sizes: 22-28 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Toe Regular (R) or Sandal (S)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model .....	<b>FTC</b>
Generation .....	<b>2</b>
Profile .....	<b>M</b>
Cap .....	<b>1</b>
Size .....	<b>23</b>
Skin Tone .....	<b>L</b>
Toe .....	<b>4</b>
Side .....	<b>S</b>
	<b>L</b>
<b>FTC-2M-123L4-SL</b>	
*Generation 3 not available in High Profile	
*Cap option only available for Generation 3	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	3 <sup>5</sup> / <sub>8</sub> /92	22/620
26-28	3 <sup>7</sup> / <sub>8</sub> /98	26/735

Runway Stiffness Category Selection Chart							
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255
Low Impact Level	1	1	2	3	4	5	6
Moderate Impact Level	1	2	3	4	5	6	—
High Impact Level	2	3	4	5	6	—	—

## Runway® HX

- Heel height adjustability from 0-2", utilizing a 4 mm hex key
- Designed for active K3 amputees
- Available sizes: 22-28 cm
- User Weight Rating: 255 lbs (116 kg)
- Stiffness Categories: 1-6
- Average Clearance: 3"
- Heel Height: 3/8"
- Uses Sandal Toe foot shell
- 36 month warranty on foot module; 18 months on ankle and 6 months on foot shell



Description	Model No.
Runway HX Foot	RS2-K1
Foot Shell (default sandal toe)	FTC-2M-1
Foot Shell (default regular toe)	FTC-3M-1

### Please provide the following when ordering a foot shell:

- Model number
- Generation 2 or 3
- Profile: Mid (M), High (H) or Sandal toe (S)
- Cap: Yes (0) or No (1)
- Available sizes: 22-28 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Toe Regular (R) or Sandal (S)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

**Foot Shell Ordering Example**

Model ..... **FTC**  
 Generation ..... **2**  
 Profile ..... **M**  
 Cap ..... **1**  
 Size ..... **23**  
 Skin Tone ..... **L**  
 Toe ..... **S**  
 Side ..... **L**

**FTC-2M-123L4-SL**

\*Generation 3 not available in High Profile  
 \*Cap option only available for Generation 3

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	3 <sup>5</sup> / <sub>8</sub> /92	22/625
26-28	3 <sup>7</sup> / <sub>8</sub> /98	26/735

Runway HX Stiffness Category Selection Chart							
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255
Low Impact Level	1	1	2	3	4	5	6
Moderate Impact Level	1	2	3	4	5	6	—
High Impact Level	2	3	4	5	6	—	—

**Promenade™**

- Designed specifically for low impact amputees
- Partial split toe offers controlled inversion/eversion
- Full length sole plate facilitates a smooth transition from heel strike all the way through toe off
- User Weight Rating: 325 lbs
- Stiffness Categories: 1-5
- Heel Height: 3/8"
- 36 month warranty on foot module; 6 months on foot shell



Description	Model No.
Promenade Foot	<b>VS2</b>
Foot Shell (default)	<b>FTC-3M-0</b>
High Profile Foot Shell	<b>FTC-2H-1</b>

**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model . . . . .	<b>FTC</b>
Generation . . . . .	<b>3</b>
Profile . . . . .	<b>M</b>
Cap . . . . .	<b>0</b>
Size . . . . .	<b>25</b>
Skin Tone . . . . .	<b>L</b>
	<b>4-R</b>
Side . . . . .	<b>R</b>
<b>FTC-3M-025L4-RR</b>	
<i>*Generation 3 not available in High Profile</i>	
<i>*Cap option only available for Generation 3</i>	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	4 <sup>5</sup> / <sub>8</sub> /116	17/490
26-28	4 <sup>7</sup> / <sub>8</sub> /124	21/595
29-31	5 <sup>3</sup> / <sub>8</sub> /135	25/710

Promenade Stiffness Category Selection Chart					
Patient Weight (lbs)	100-126	127-160	161-203	204-258	259-325
Low Impact Level	1	2	3	4	5

**Sierra®**

- Angle top design increases the length of the carbon fiber spring providing a smoother gait and ultimate flexibility
- Full-length heel enables foot flat early in the gait cycle, providing a stable base of support
- Split toe delivers excellent inversion-eversion
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Heel Height: 3/8"
- 36 month warranty on foot module; 6 months on foot shell

Description	Model No.
Sierra Foot	<b>FS1</b>
Foot Shell (default)	<b>FTC-3M-0</b>
High Profile Foot Shell	<b>FTC-2H-1</b>



**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model .....	<b>FTC</b>
Generation .....	<b>3</b>
Profile .....	<b>M</b>
Cap .....	<b>0</b>
Size .....	<b>23</b>
Skin Tone .....	<b>L</b>
Side .....	<b>4-R</b>
Side .....	<b>L</b>
<b>FTC-3M-023L4-RL</b>	
*Generation 3 not available in High Profile	
*Cap option only available for Generation 3	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	5 <sup>1</sup> / <sub>2</sub> /140	18/155
26-28	6/154	23/640
29-31	6 <sup>3</sup> / <sub>4</sub> /172	26/750

Sierra Stiffness Category Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	—
High Impact Level	2	3	4	5	6	7	8	9	—	—

**Highlander®**

- Designed for use with longer residual limbs or knee frames
- Low maintenance design ensures longer life with high durability
- Split toe delivers excellent inversion-eversion
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Connector: Integral male pyramid
- Heel Height: 3/8"
- 36 month warranty on foot module; 6 months on foot shell



**Highlander® MAX**

- Additional carbon fiber layers provide durability and energy return
- User Weight Rating: 500 lbs (227 kg)
- Build Height (average): 3.7 inches
- Custom made and non-returnable
- Warranty: 36 months on foot module; 6 months on foot shell

Description	Model No.
Highlander Foot	<b>FS3</b>
Highlander MAX Foot	<b>FS3-H5</b>
Foot Shell (default)	<b>FTC-3M-0</b>
High Profile Foot Shell	<b>FTC-2H-1</b>

**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Highlander sizes: 22-31 cm;  
Highlander Max sizes: 26-31 cm - custom order
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

**Foot Shell Ordering Example**

Model ..... **FTC**  
 Generation ..... **3**  
 Profile ..... **M**  
 Cap ..... **0**  
 Size ..... **23**  
 Skin Tone ..... **D**  
 Side ..... **4-R**  
 Side ..... **L**

**FTC-3M-023D4-RL**

\*Generation 3 not available in High Profile  
 \*Cap option only available for Generation 3

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	3 <sup>7</sup> / <sub>8</sub> /98	17/485
26-28	4 <sup>1</sup> / <sub>4</sub> /108	21/590
29-31	4 <sup>3</sup> / <sub>4</sub> /121	26/750

Highlander Stiffness Category Selection Chart													
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365	366-420	421-480	481-500
Low Impact Level	1	1	2	3	4	5	6	7	8	9	10	11	12
Moderate Impact Level	1	2	3	4	5	6	7	8	9	10	—	—	—
High Impact Level	2	3	4	5	6	7	8	9	10	11	—	—	—

## Senator®

- Lightweight, energy-return design facilitates smooth, easy stride and durability
- Foot is not pre-bonded to foot shell, simplifying the use of heel wedges to improve alignment
- Simulated ankle motion provides a stable base of support throughout the stance phase
- Low to moderately active K3 level amputees
- User Weight Rating: 300 lbs (136 kg)
- Stiffness Categories: 1-6
- Heel Height: 3/8"
- 36 month warranty on foot module; 6 months on foot shell

Description	Model No.
Senator Foot	<b>VS1</b>
Foot Shell (default)	<b>FTC-3M-0</b>
High Profile Foot Shell	<b>FTC-2H-1</b>



### Please provide the following when ordering a foot shell:

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Available sizes: 22-30 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

**Foot Shell Ordering Example**

Model ..... **FTC**  
 Generation ..... **3**  
 Profile ..... **M**  
 Cap ..... **0**  
 Size ..... **25**  
 Skin Tone ..... **L**  
 Side ..... **4-R**  
 Side ..... **R**

**FTC-3M-025L4-RR**

\*Generation 3 not available in High Profile  
 \*Cap option only available for Generation 3

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	4 <sup>3</sup> / <sub>4</sub> /120	18/520
26-28	5 <sup>3</sup> / <sub>8</sub> /136	23/655
29-31	5 <sup>7</sup> / <sub>8</sub> /149	28/780

Senator Stiffness Category Selection Chart						
Patient Weight (lbs)	100-120	121-145	146-175	176-210	211-250	251-300
Low Impact Level	1	2	3	4	5	6
Moderate Impact Level	2	3	4	5	6	—

**Pacifica®**

- Designed for use with longer residual limbs or those needing more build space
- Split toe design provides stability enhancing inversion and eversion
- Lightweight, durable, suited to all activity levels
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Connector: Integral male pyramid
- Heel Height:  $\frac{3}{8}$ "
- 36 month warranty on foot module; 6 months on foot shell



Description	Model No.
Pacifica Foot	<b>FS2</b>
Foot Shell (default)	<b>FTC-3M-0</b>
High Profile Foot Shell	<b>FTC-2H-1</b>

**Please provide the following when ordering a foot shell:**

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Available sizes: 24-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model . . . . .	<b>FTC</b>
Generation . . . . .	<b>3</b>
Profile . . . . .	<b>M</b>
Cap . . . . .	<b>0</b>
Size . . . . .	<b>25</b>
Skin Tone . . . . .	<b>M</b>
	<b>4-R</b>
Side . . . . .	<b>L</b>
<b>FTC-3M-025M4-RL</b>	
<i>*Generation 3 not available in High Profile</i>	
<i>*Cap option only available for Generation 3</i>	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
24-25	2 <sup>3</sup> / <sub>8</sub> /59	17/485
26-28	2 <sup>3</sup> / <sub>8</sub> /61	22/610
29-31	2 <sup>1</sup> / <sub>2</sub> /63	25/720

Pacifica Stiffness Category Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	—
High Impact Level	2	3	4	5	6	7	8	9	—	—



## Pacifica® LP

- Requires less than 2" of clearance for most commonly used sizes
- Lightweight and durable
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Connector: Integral male pyramid
- Heel Height: 3/8"
- 36 month warranty on foot module; 6 months on foot shell

Description	Model No.
Pacifica LP Foot	<b>FS4</b>
Foot Shell (default)	<b>FTC-3M-0</b>
High Profile Foot Shell	<b>FTC-2H-1</b>



### Please provide the following when ordering a foot shell:

- Model number
- Generation 2 or 3
- Profile: Mid (M) or High (H)
- Cap: Yes (0) or No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model . . . . .	<b>FTC</b>
Generation . . . . .	<b>3</b>
Profile . . . . .	<b>M</b>
Cap . . . . .	<b>0</b>
Size . . . . .	<b>23</b>
Skin Tone . . . . .	<b>L</b>
Side . . . . .	<b>4-R</b>
Side . . . . .	<b>L</b>
<b>FTC-3M-023L4-RL</b>	
*Generation 3 not available in High Profile	
*Cap option only available for Generation 3	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	1 <sup>7</sup> / <sub>8</sub> /47	16/465
26-28	2/50	20/570
29-31	2 <sup>3</sup> / <sub>8</sub> /60	27/760

Pacifica LP Stiffness Category Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	—
High Impact Level	2	3	4	5	6	7	8	9	—	—

**LP Symes**

- Carbon fiber construction delivers energy return and durability
- Direct lamination cup allows for less than 2.1" of clearance
- Ideal for K3/K4 amputees
- User Weight Rating: 365 lbs (166 kg)
- Stiffness Categories: 1-9
- Heel Height: 3/8"
- Uses Mid Profile foot shell
- 36 month warranty on foot module; 6 months on foot shell



Description	Model No.
LP Symes Foot	<b>LP2</b>
Foot Shell (default)	<b>FTC-3M-1</b>

**Please provide the following when ordering a foot shell:**

- Model number
- Generation 3
- Profile: Mid (M)
- Cap: No (1)
- Available sizes: 22-31 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model .....	<b>FTC</b>
Generation .....	<b>3</b>
Profile .....	<b>M</b>
Cap .....	<b>1</b>
Size .....	<b>23</b>
Skin Tone .....	<b>L</b>
	<b>4-R</b>
Side .....	<b>L</b>
<b>FTC-3M-123L4-RL</b>	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	1 <sup>3</sup> / <sub>4</sub> /43	15/420
26-28	1 <sup>7</sup> / <sub>8</sub> /48	19/545
29-31	2 <sup>1</sup> / <sub>8</sub> /53	24/670

LP Symes Stiffness Category Selection Chart										
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	—
High Impact Level	2	3	4	5	6	7	8	9	—	—

### WalkTek™

- Designed specifically for limited community ambulators, K2 level single-speed walkers
- iMX Technology (independent multi-axial suspension) provides controlled multi-axial movement and independent heel and toe compliance for excellent stability
- Ankle height pyramid placement simplifies the fitting process by making adjustments quick and easy
- User Weight Rating: 300 lbs (136 kg)
- Connector: Male pyramid
- Heel Height: 3/8"
- 12 month warranty on foot module; 6 months on foot shell



Description	Model No.
WalkTek Foot	<b>VS3</b>
Foot Shell (default)	<b>FTC-3M-0</b>
Foot Shell (Sandal Toe)	<b>FTC-2M-1</b>
High Profile Foot Shell	<b>FTC-2H-1</b>

**Please provide the following when ordering a foot shell:**

- Model number
- Generation 3
- Profile: Mid (M) or Sandal toe (S)
- Cap: Yes (0) or No (1)
- Regular toe sizes with cap: 21-30 cm; Sandal toe sizes without cap: 22-28 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Toe Regular (R) or Sandal (S)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model . . . . .	<b>FTC</b>
Generation . . . . .	<b>3</b>
Profile . . . . .	<b>M</b>
Cap . . . . .	<b>0</b>
Size . . . . .	<b>23</b>
Skin Tone . . . . .	<b>L</b>
	<b>4</b>
Toe . . . . .	<b>R</b>
Side . . . . .	<b>L</b>
<b>FTC-3M-023L4-RL</b>	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
21-25	3 <sup>3</sup> / <sub>8</sub> /85	17/490
26-30	3 <sup>3</sup> / <sub>4</sub> /97	22/630

WalkTek Stiffness Category Selection Chart				
Patient Weight (lbs)	100-130	131-180	181-240	241-300
Category	1	2	3	4

**Freestyle Swim Foot**

- Push Button technology allows movement between 70° plantar flexion for swimming and a 3/8" heel height setting for walking
- Protective cap snaps over the ankle
- Utilizes LP Symes as base
- 36 month warranty on carbon foot module, 12 months on ankle and 6 months on shell



Description	Model No.
Freestyle Swim Foot	<b>LP2-W2</b>
Foot Shell (default)	<b>FTC-3M-1</b>

**Please provide the following when ordering a foot shell:**

- Model number
- Generation 3
- Profile: Mid (M)
- Cap: No (1)
- Available sizes: 22-30 cm
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

Foot Shell Ordering Example	
Model .....	<b>FTC</b>
Generation .....	<b>3</b>
Profile .....	<b>M</b>
Cap .....	<b>1</b>
Size .....	<b>25</b>
Skin Tone .....	<b>L</b>
	<b>4-R</b>
Side .....	<b>R</b>
<b>FTC-3M-125L4-RR</b>	

Build Height and Weight		
Size (cm)	Height (in/mm)	Weight (oz/g)
22-25	4/102	24/680
26-28	4 <sup>1</sup> / <sub>4</sub> /107	28/790
29-31	4 <sup>5</sup> / <sub>8</sub> /117	33/945

Freestyle Swim Foot Stiffness Category Selection Chart							
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255
Low Impact Level	1	1	2	3	4	5	6
Moderate Impact Level	1	2	3	4	5	6	—
High Impact Level	2	3	4	5	6	—	—

### Slalom Ski Foot

- Designed to be quickly and easily inserted directly into a standard ski binding
- Eliminates the need for a ski boot
- Provides necessary flex and shock absorption for fluid motion
- Better control and range-of-motion increases endurance and allows more time on the slopes
- Build Height: 4"/102 mm; Build Weight: 24 oz/680 g
- Lead time: 10-14 days
- 12 month warranty non-returnable, custom product



Model No. **SKI**

Slalom Ski Stiffness Category Selection Chart			
Patient Weight (lbs)	100-160	161-200	201-250
Category	11	12	13

### Nitro Running Foot

- Custom foot module engineered for jogging and running
- Individually handcrafted carbon-epoxy lay-up techniques ensure optimal vertical displacement and energy return
- Endoskeletal attachment options facilitate adjustment and fabrication
- 9 stiffness categories tailored to amputee weight and activity preference
- Non-returnable, custom product
- Choice of short profile (FS9) or tall profile (FS6)
- 12 month warranty for user weights up to 365 lbs
















Description	Clearance	Model No.
Tall Profile	9.5"	<b>FS6</b>
Short Profile	7.25"	<b>FS9</b>

Spectra Socks		
Sock Length and Width	Foot Sizes	Model No.
12" x 2 <sup>3</sup> / <sub>4</sub> "	21-25 cm	<b>ACC-00-10112-00</b>
14" x 2 <sup>3</sup> / <sub>4</sub> "	26-28 cm	<b>ACC-00-10114-00</b>
18" x 2 <sup>3</sup> / <sub>4</sub> "	29-31 cm	<b>ACC-00-10118-00</b>

**Kingsley - SACH Feet**




*All Kingsley feet are available in Medium Brown 4A and Dark Brown 7A. Latin is also available.*

Description	Heel Cushions	Heel Height	Sizes	Model No.	
<b>Low Profile*</b>	Soft, Medium, Firm	3/4" (18 mm)	5-16 (22-33 cm)	<b>K01</b>	
<b>High Profile</b>	Soft, Medium, Firm	3/4" (18 mm)	6-13 (23-30 cm)	<b>K02</b>	
<b>Cowboy Hi-Heel</b> • Male SACH foot for high heel boots or shoes	Soft, Medium, Firm	1 3/8" (35 mm)	5-13 (22-30 cm)	<b>K03</b>	
<b>Continental</b> • Molded SACH foot for stacked heels	Soft, Medium, Firm	1" (25 mm)	5-11 (22-28 cm)	<b>K041</b>	
<b>Casual</b> • Ladies SACH for stacked heel	Soft, Medium	1 3/8" (35 mm)	5-10 (22-27 cm)	<b>K042</b>	
<b>Fashion</b> • Ladies SACH for medium heel heights, pumps	Soft, Medium, Firm	1 3/4" (45 mm)	5-10 (22-27 cm)	<b>K043</b>	
<b>Hi-Fashion</b> • Ladies SACH for mid to high heels, wedges, and boots	Medium	2 3/8" (60 mm)	5-10 (22-27 cm)	<b>K044</b>	
<b>Flattie</b> • Molded Sach foot for flats, sandals, sneakers, and ski boots	Soft, Medium, Firm	3/8" (10 mm)	5-11 (22-28 cm)	<b>K051</b>	
<b>Flattie Zero</b> • Molded Sach foot for flats, sandals, sneakers, and ski boots	Medium	Flat	5-11 (22-28 cm)	<b>K051Z</b>	
<b>Hi-Style</b> • Ladies SACH foot for high heels spikes	Medium	3 1/2" (89 mm)	5-10 (22-27 cm)	<b>K054</b>	
<b>Strider</b> • Molded SACH foot for tennis shoes and all mens flat heel shoes	Soft, Medium, Firm	1/4"-3/8" (8-10 mm)	5-16 (22-33 cm)	<b>K08</b>	
<b>Strider Zero</b> • Molded SACH foot for tennis shoes and all mens flat heel shoes	Soft, Medium, Firm	Flat	5-16 (22-33 cm)	<b>K08Z</b>	



SACH Feet - Continued					
<b>Beachcomber</b> <ul style="list-style-type: none"> <li>Exposed Keel Molded SACH foot with rigid foam ankle block for sports, recreation and shower wear</li> </ul>	Medium	Flat	5-13 (22-30 cm)	<b>K16</b>	
<b>Wayfarer</b> <ul style="list-style-type: none"> <li>Molded SACH foot for Post-Op; Can be used with sandals and limited barefoot walking</li> </ul>	Medium	Flat	Small (24 cm) Medium (26 cm) Large (28 cm)	<b>K10T</b>	

SACH Lite Feet					
Description	Heel Cushions	Heel Height	Sizes	Model No.	
<b>Litefoot Low Profile</b> <ul style="list-style-type: none"> <li>For geriatric and light activity amputees</li> </ul>	Soft, Medium, Firm	3/4" (18 mm)	5-16 (22-33 cm)	<b>K09</b>	
<b>Continental Litefoot</b> <ul style="list-style-type: none"> <li>For stacked heels</li> </ul>	Soft, Medium, Firm	1" (25 mm)	5-11 (22-28 cm)	<b>K0941</b>	
<b>Casual Litefoot</b> <ul style="list-style-type: none"> <li>Ladies Sach for stacked heels</li> </ul>	Soft, Medium	1 3/8" (35 mm)	5-10 (22-27 cm)	<b>K0942</b>	
<b>Litefoot High Profile</b> <ul style="list-style-type: none"> <li>High profile molded SACH foot</li> </ul>	Soft, Medium, Firm	3/4" (18 mm)	6-13 (23-30 cm)	<b>K0902</b>	
<b>Litefoot Cowboy Hi-Heel</b> <ul style="list-style-type: none"> <li>Male SACH foot for high heel boots or shoes</li> </ul>	Soft, Medium, Firm	1 3/8" (35 mm)	5-13 (22-30 cm)	<b>K0903</b>	
<b>Litefoot Strider</b> <ul style="list-style-type: none"> <li>Molded SACH foot for tennis shoes and all men's flatter heel shoes</li> </ul>	Soft, Medium, Firm	1/4"-3/8" (8-10 mm)	5-16 (22-33 cm)	<b>K0908</b>	
<b>Litefoot Strider Zero</b> <ul style="list-style-type: none"> <li>SACH foot for men's flat heel shoes</li> </ul>	Soft, Medium, Firm	Flat	5-16 (22-33 cm)	<b>K0908Z</b>	
<b>Litefoot Wayfarer - Flat</b> <ul style="list-style-type: none"> <li>Molded SACH foot for post-op</li> </ul>	Medium	Flat	Small (24 cm) Medium (26 cm) Large (28 cm)	<b>K0910T</b>	
<b>Litefoot Fashion</b> <ul style="list-style-type: none"> <li>Ladies molded SACH foot for medium heel heights</li> </ul>	Soft, Medium, Firm	1 3/4" (45 mm)	5-10 (22-27 cm)	<b>K0943</b>	
<b>Litefoot Hi-Fashion</b> <ul style="list-style-type: none"> <li>Ladies molded SACH foot for mid-high heels, wedges and fashion boots</li> </ul>	Medium	2 3/8" (60 mm)	5-10 (22-27 cm)	<b>K0944</b>	


**SACH Lite Feet - Continued**

Description	Heel Cushions	Heel Height	Sizes	Model No.	
<b>Litefoot Flattie</b> <ul style="list-style-type: none"> <li>Molded SACH foot for flats, sandals and ski boots</li> </ul>	Soft, Medium, Firm	3/8" (10 mm)	5-11 (22-28 cm)	<b>K0951</b>	
<b>Litefoot Flattie Zero</b> <ul style="list-style-type: none"> <li>Molded SACH foot for flats and sandals</li> </ul>	Medium	Flat	5-11 (22-28 cm)	<b>K0951Z</b>	
<b>Litefoot Post-Op Flatfoot</b> <ul style="list-style-type: none"> <li>Molded SACH foot for light duty bare foot walking</li> </ul>	None	Flat	Small (24 cm) Medium (26 cm) Large (28 cm)	<b>K10</b>	

**SACH Feet Rated to 500 lbs**

	Description	Heel Height	Heel Cushions	Sizes	Model No.	
<ul style="list-style-type: none"> <li>TiMed integrated 500 lb rated SACH adapter nonreturnable</li> <li>Extra Firm heel cushion and K20 Sturdy toe standard</li> <li>14 mm integrated foot bolt</li> <li>TiMed BK-500 XHD Kit located in Endoskeletal section of catalog</li> </ul>	Low Profile	3/4" (18 mm)	Firm or Extra Firm	5-16 (22-33 cm)	<b>KT1</b>	
	Strider	1/4"- 3/8" (8-10 mm)	Firm or Extra Firm	5-16 (22-33 cm)	<b>KT8</b>	


**Symes Feet**


	Description	Heel Height	Sizes	Model No.	
<ul style="list-style-type: none"> <li>Molded Symes Foot allowing clearances to 1"</li> <li>Medium heel cushion</li> </ul>	Adult Symes	3/4" (18 mm)	5-14 (22-31 cm)	<b>K07</b>	
	Adult Symes Strider	1/4"- 3/8" (8-10 mm)	5-14 (22-31 cm)	<b>K0708</b>	
	Steplite Symes	3/4" (18 mm)	5-13 (22-31 cm)	<b>KC07</b>	
	Steplite Strider Symes	1/4"- 3/8" (8-10 mm)	5-16 (22-33 cm)	<b>KC08</b>	

Description	Measurements	Model No.
Standard foot bolt only	3/8" -16x3"	<b>K20103</b>
Foot bolt assembly with Large Adapter/HD	3/8" -16x3"	<b>K20103AH</b>
Small T-Nut	3/8" -16x3"	<b>K2012</b>
Juvenile Symes Nut Assembly	3/8" -16x3"	<b>K2017A</b>
Heavy Duty Symes Nut	—	<b>K2016</b>





STEN Foot						
	Description	Heel Height	Heel Cushions	Sizes	Model No.	
<ul style="list-style-type: none"> <li>Flexible keel foot ideally suited for the limited community ambulator</li> <li>Unique keel structure offers enhanced stability to better accommodate uneven surfaces while retaining smooth gait</li> <li>Enhanced inversion characteristics</li> <li>Heel cushions available in Soft, Medium and Firm</li> </ul>	STEN Low Profile	3/4" (18 mm)	Soft, Medium, Firm	5-16 (22-33 cm)	<b>K01ST</b>	
	STEN Strider	1/4" - 3/8" (8-10 mm)	Soft, Medium, Firm	5-16 (22-33 cm)	<b>K08ST</b>	
	STEN Strider Zero	Flat	Soft, Medium, Firm	5-16 (22-33 cm)	<b>K08ZST</b>	
	STEN Flattie	3/8" (10 mm)	Soft, Medium, Firm	5-11 (22-28 cm)	<b>K051ST</b>	
	STEN Casual	1 3/8" (35 mm)	Soft, Medium, Firm	5-10 (22-27 cm)	<b>K042ST</b>	

Dynamic Response Step Lite Foot					
	Description	Heel Height	Sizes	Model No.	
<b>Dynamic Response Litefoot</b> <ul style="list-style-type: none"> <li>Designed for the moderately active amputee</li> <li>Offers smooth rollover characteristics, while the composite spring plate gives a forward movement during the transition from mid-stance through toe-off</li> <li>Indicated for the patient that wants to walk comfortably rather than run, the spring plate is fabricated in the low to mid flexion resistance range for a smooth gait</li> <li>Heel Cushions: Medium</li> </ul>	Steplite Low Profile	3/4" (18 mm)	5-16 (22-33 cm)	<b>KC01</b>	
	Steplite Strider	1/4" - 3/8" (8-10 mm)	5-16 (22-32 cm)	<b>K08</b>	
	Steplite Flattie	3/8" (10 mm)	5-11 (22-28 cm)	<b>KC051</b>	
	Steplite Flattie Zero	Flat	5-10 (22-27 cm)	<b>KC051Z</b>	
	Steplite Strider Zero	Flat	5-16 (22-33 cm)	<b>KC08Z</b>	
	Steplite Flattie Litefoot	3/8" (10 mm)	5-11 (22-28 cm)	<b>KC0951</b>	
	Steplite Wayfarer	Flat	Small (24 cm) Medium (26 cm) Large (28 cm)	<b>KC10T</b>	

Heel Cushion Material	
Description	30° Wedge
Extra Soft	<b>K20HW-XS</b>
Soft	<b>K20HW-S</b>
Medium	<b>K20HW-M</b>
Firm	<b>K20HW-F</b>
Extra Firm	<b>K20HW-XF</b>



Double Heel Cushion Blocks		
Description	4" x 6" x 3 1/2"	4" x 6" x 2 1/2"
Extra Soft	<b>K20HCBS-XS</b>	<b>K20HCB-XS</b>
Soft	<b>K20HCBS-S</b>	<b>K20HCB-S</b>
Medium	<b>K20HCBS-M</b>	<b>K20HCB-M</b>
Firm	<b>K20HCBS-F</b>	<b>K20HCB-F</b>
Extra Firm	<b>K20HCBS-XF</b>	<b>K20HCB-XF</b>



**Motion Control - Motion Foot™**

- Very low ground reaction force provides unmatched comfort
- Unequaled range of motion
- 50° plantar and dorsiflexion with hydraulic ankle enables total foot relaxation while sitting and incomparable ability to climb up/down steep slopes
- Both plantar flexion and dorsiflexion resistance may be adjusted independently for comfort or user activities
- K3 activity level
- Heel height: 3/8"
- Clearance: 9.7"
- Available sizes: 22-30 cm
- Weight rating up to 220 lbs
- Foot shell included\* (6 month warranty)
- When ordering foot shell, specify right or left
- Standard foot shell color: Caucasian, lead time for Brown foot shell
- 24 month warranty/60 day patient trial



Size (cm)	Model No.
22	<a href="#">4050195</a>
23	<a href="#">4050196</a>
24	<a href="#">4050197</a>
25	<a href="#">4050198</a>
26	<a href="#">4050199</a>
27	<a href="#">4050200</a>
28	<a href="#">4050201</a>
29	<a href="#">4050202</a>
30	<a href="#">4050203</a>

MotionFoot Foot Shells		
Size (cm)	Model No.	
	Left	Right
22	<a href="#">1702024</a>	<a href="#">1702025</a>
23	<a href="#">1702026</a>	<a href="#">1702027</a>
24	<a href="#">1702028</a>	<a href="#">1702029</a>
25	<a href="#">1702030</a>	<a href="#">1702031</a>
26	<a href="#">1702032</a>	<a href="#">1702033</a>
27	<a href="#">1702034</a>	<a href="#">1702035</a>
28	<a href="#">1702036</a>	<a href="#">1702037</a>
29	<a href="#">1702038</a>	<a href="#">1702039</a>
30	<a href="#">1702040</a>	<a href="#">1702041</a>

### Ossur - PROPRIO FOOT® with EVO™

- The newly redesigned PROPRIO FOOT with EVO utilizes intelligent ankle flexion to help users walk confidently and naturally
- Includes Energy Vector Optimization (EVO) technology, closely mimicking the movement of the human foot
- The resultant lifelike gait reduces stress on back, hips, and knees, enabling users to walk longer
- Detects toe-off, triggering dorsiflexion during swing phase to safely clear the ground
- Low to moderate impact level
- Heel height: adjustable up to 2" (50 mm)
- Clearance: 7"
- Available sizes: 25-30 cm; categories 1-6
- Specify left or right when ordering foot shell
- 12 month limited warranty

Description	Model No.
PROPRIO Foot	<b>PRPE</b>
Foot Cover	<b>FCPE0</b>



**Ordering Example**

Model . . . . . **PRPE**  
 Category . . . . . **5**  
 Foot Size . . . . . **27**  
 Heel Height . . . . . **0**  
**PRPE5270**  
 Foot Cover . . . . . **FCPE0**  
 Foot Size . . . . . **27**  
 Side . . . . . **L**  
**FCPE027L**

\*add a BR suffix for Brown foot cover

PROPRIO Category Selection Guide							
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-194	195-220	221-256
Low Impact Level	1	1	2	3	4	5	6
Moderate Impact Level	1	2	3	4	5	6	—

**Re-Flex Rotate™ with EVO™**

- Features a sleek profile and maintenance-free design
- Integrated shock absorption diminishes stress on the body
- The EVO feature enables users to experience natural gait progression by utilizing the foot plate to its fullest potential
- Low to high impact level
- Available sizes: 22-30 cm; categories 1-8
- Weight rating up to 324 lbs (147 kg)
- Heel heights available in  $\frac{3}{8}$ " (0) or  $\frac{3}{4}$ " (2)
- Clearance: Standard,  $8\frac{1}{2}$ " -  $10\frac{7}{8}$ "; Tall Profile,  $10\frac{1}{2}$ " -  $12\frac{9}{16}$ "
- Standard profile with shock and without shock only available with  $\frac{3}{8}$ " heel height;  $\frac{3}{4}$ " for these are special order
- Add BR to model number when ordering Brown\* foot cover
- 36 month limited warranty; foot covers carry a 2 month warranty

Description	Model No.
with Male Pyramid Adapter (pictured)	<b>RSPE</b>
with Female Pyramid Adapter	<b>RSFE</b>
Tall Profile with Male Pyramid Adapter	<b>RHPE</b>
Tall Profile with Female Pyramid Adapter	<b>RHFE</b>
Standard Profile without Shock	<b>RSXE</b>
Tall Profile without Shock	<b>RHXE</b>
Foot Cover	<b>FCE0</b>



**Ordering Example**

Model..... **RSPE**  
 Category..... **3**  
 Size..... **23**  
 $\frac{3}{8}$ " Heel Height..... **0**

**RSPE3230**

Foot Cover..... **FCE0**  
 Foot Size..... **23**  
 Side..... **L**

**FCE023L**

\*add a BR suffix for Brown foot cover

Re-Flex Rotate with EVO Category Selection Guide									
Patient Weight (lbs)	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324
Low Impact Level	1	1	2	3	4	5	6	7	8
Moderate Impact Level	1	2	3	4	5	6	7	8	—
High Impact Level	2	3	4	5	6	7	8	—	—

### Re-Flex Shock® with EVO™

- Unparalleled vertical shock absorption and energy return for the active user
- Improved attachment clamp for easy alignment and height adjustments
- EVO feature for improved rollover and stability in stance
- Low to high impact level
- Available sizes: 22-30 cm; categories 1-9
- Weight rating up to 366 lbs (166 kg)
- Heel heights available in 3/8" (0) or 3/4" (2)
- Add BR to model number when ordering Brown\* foot cover

Description	Model No.
with Male Pyramid Adapter (pictured)	<b>RSP0</b>
Standard Profile with Shock	<b>SSPE</b>
Tall Profile with Shock	<b>SHPE</b>
Foot Cover	<b>FCE0</b>



**Ordering Example**

Model ..... **RSP0**  
 Category ..... **2**  
 Size ..... **25**  
 3/8" Heel Height ..... **0**  
**RSP02250**  
 Foot Cover ..... **FCE0**  
 Foot Size ..... **25**  
 Side ..... **L**  
**FCE025L**


\*add a BR suffix for Brown foot cover

Re-Flex Shock Category Selection Guide										
Patient Weight (lbs)	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324	325-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	—	—
High Impact Level	2	3	4	5	6	7	8	—	—	—

**LP Ceterus®**

- An extra low profile option for long residual limbs
- Integrated adjustable shock absorption unit
- Rotation with progressive stiffening and no abrupt stops
- Maintenance free with slim profile for easy finishing
- Available in sizes 22-30 cm, categories 1-8
- Clearance: 6" with pyramid and foot cover
- Brown foot shells available, add BR to part no. when ordering
- Maximum patient weight is 325 lbs
- Heel height -  $\frac{3}{8}$ " = 0 or  $\frac{3}{4}$ " = 2
- Male Pyramid (1) or Female Pyramid (2)
- 36 month warranty

Description	Model No.
LP Vari-Flex	<b>VLX0</b>
LP Ceterus Shock Module	<b>CTL0</b>
Foot Cover	<b>FCA0</b>



<b>Shock Kits</b> replace ? with category (see below)	
Description	Model No.
Male Pyramid; Foot size 22-24	<b>CTL0?11</b>
Male Pyramid; Foot size 25-30	<b>CTL0?21</b>
Female Pyramid; Foot size 22-24	<b>CTL0?12</b>
Female Pyramid; Foot size 25-30	<b>CTL0?22</b>

**Ordering Example**

Model. . . . . **VLX0**  
 Category . . . . . **4**  
 Size . . . . . **24**  
 Heel Height . . . . . **0**  
**VLX04240**  
 Shock Kit . . . . . **CTL0**  
 Category . . . . . **2**  
 Pyramid . . . . . **11**  
**CTL0211**  
  
 Foot Cover. . . . . **FCA0**  
 Foot Size . . . . . **24**  
 Side . . . . . **L**  
**FCA024L**  
 \*add a BR suffix for Brown foot cover

<b>LP Ceterus Category Selection Guide</b>									
Patient Weight (lbs)	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325
Low Impact Level	1	1	2	3	4	5	6	7	8
Mod. Impact Level	1	2	3	4	5	6	7	8	—
High Impact Level	2	3	4	5	6	7	8	—	—

### Modular III™

- High energy storage and release
- Lightweight yet extremely strong for maximum performance
- Low to high impact level
- Available sizes: 22-30 cm; categories 1-9
- Clearance: 7<sup>1</sup>/<sub>4</sub>" (or 15<sup>1</sup>/<sub>16</sub>" with pyramid and foot cover)
- Foot covers available in Brown\* and Beige
- Weight rating up to 365 lbs (166 kg)
- Heel heights available in <sup>3</sup>/<sub>8</sub>" (0) or <sup>3</sup>/<sub>4</sub>" (2)
- 36 month warranty; foot covers carry a 2 month warranty

Description	Model No.
Modular III	<b>MDX0</b>
Foot Cover	<b>FCA0</b>

Bolt-On Male Pyramid Adapters			
Description	Weight	Category	Model No.
Male Pyramid Adapter	180 g (6.4 oz)	1-6	<b>CM100007</b>
Male Pyramid Adapter	220 g (7.8 oz)	7-9	<b>CM100009</b>



**Ordering Example**

Model ..... **MDX0**  
 Category ..... **5**  
 Size ..... **27**  
<sup>3</sup>/<sub>4</sub>" Heel Height..... **2**  
**MDX05272**

Foot Cover ..... **FCA0**  
 Foot Size ..... **27**  
 Side..... **L**  
**FCA027L**

\*add a BR suffix for Brown foot cover

Modular III Category Selection Guide										
Patient Weight (lbs)	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324	325-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	*
High Impact Level	2	3	4	5	6	7	8	9	*	*
<small>*Special Order</small>										

**Vari-Flex® with EVO™**

- Vari-Flex with EVO is lightweight, easy to assemble and has a slender profile which makes it easy to cosmetically cover
- The EVO™ feature enables users to experience natural gait progression by utilizing the foot plate to its fullest potential
- Offers natural gait with less fatigue and less strain on the back
- Low to high impact level
- Weight rating up to 365 lbs (166 kg)
- Available sizes: 22-30 cm; categories 1-9
- Clearance: 6<sup>11</sup>/<sub>16</sub>"
- Heel heights available in 3/8" (0) or 3/4" (2)
- Foot covers available in Brown\* and Beige
- 36 month warranty; foot covers carry a 2 month warranty

Description	Model No.
With Male Pyramid (pictured)	<b>VFPE</b>
Without Male Pyramid	<b>VFXE</b>
Foot Cover	<b>FCE0</b>



**30 mm (13/16") Carbon Fiber Tube Kit**

- Choose according to the category and size of the foot module
- Hole Spacing: 38 mm (1<sup>1</sup>/<sub>2</sub>"

Category	Size	Length	Weight	Model No.
1-6	22-30	9" (229 mm)	7.3 oz (216 g)	<b>VFX00003</b>
7-9	22-30	9" (229 mm)	9.7 oz (275 g)	<b>VFX00004</b>
1-6	22-30	15" (381 mm)	10.3 oz (293 g)	<b>VFX00006</b>
7-9	22-30	15" (381 mm)	13 oz (369 g)	<b>VFX00007</b>

**Ordering Example**

Model ..... **VFPE**  
 Category 1-9. .... **6**  
 Foot Size ..... **26**  
 3/8 Heel Height ..... **0**  
**VFPE6260**  
 Foot Cover ..... **FCE0**  
 Foot Size ..... **26**  
 Side. .... **L**  
**FCE026L**  
 \*add a BR suffix for Brown foot cover

Vari-Flex with EVO Category Selection Guide										
Patient Weight (lbs)	99-115	116-130	131-150	151-170	171-194	195-221	222-256	257-287	288-324	325-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Mod. Impact Level	1	2	3	4	5	6	7	8	9	*
High Impact Level	2	3	4	5	6	7	8	9	*	*
*Special Order										



**LP Vari-Flex® with EVO™**

- Integrated male pyramid or EXO-lamination adapter
- Available sizes: 22-24 cm for categories 1-5 and sizes 25-30 for categories 1-9
- Low to high impact level
- Weight rating up to 365 lbs (166 kg)
- Clearance: 2<sup>11</sup>/<sub>16</sub>" with pyramid and foot cover
- Heel height: <sup>3</sup>/<sub>8</sub>"
- Foot covers available in Brown\* and Beige
- 36 month warranty



Description	Model No.
Low Profile with Pyramid** (pictured)	<b>VLPE</b>
Low Profile without Pyramid**	<b>VLXE</b>
Foot Cover	<b>FCE0</b>

**Ordering Example**

Model . . . . . **VLPE**  
 Category 1-9. . . . . **6**  
 Foot Size . . . . . **26**  
<sup>3</sup>/<sub>8</sub> Heel Height . . . . . **0**  
**VLPE6260**  
 Foot Cover . . . . . **FCE0**  
 Foot Size . . . . . **26**  
 Side. . . . . **L**  
**FCE026L**

\*add a BR suffix for Brown foot cover  
 \*\*heel height does not apply to low profile

LP Vari-Flex with EVO Category Selection Guide										
Patient Weight (lbs)	99-115	116-130	131-150	151-170	171-194	195-221	222-256	257-287	288-324	325-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Mod. Impact Level	1	2	3	4	5	6	7	8	9	*
High Impact Level	2	3	4	5	6	7	8	9	*	*
*Special Order										

**Talux®**

- Provides fluid, natural walking motion for a variety of terrain
- Dual “J” springs and carbon fiber keel ensure energy storage and release
- Tarsal core offers multi-axial function and Achilles Strap enhances forward propulsion
- Add BR to model number when ordering Brown\* foot cover
- Available sizes: 23-30 cm; categories 1-8
- Low to moderate impact level
- Weight rating up to 325 lbs (147 kg)
- Heel height: 3/8”
- Clearance: 7”
- 36 month warranty for foot module;
- 2 month warranty for foot cover



Description	Model No.
With Male Pyramid	<b>TLP0</b>
With Tube Clamp (option sold separately)	<b>TLX0</b>
Foot Cover	<b>FCT0</b>

30 mm Carbon Fiber Tube Kit			
Category	Length	Weight	Model No.
1-6	9" (229 mm)	7.60 oz/216 g	<b>TLX01011</b>
7-8	9" (229 mm)	9.70 oz/275 g	<b>TLX01012</b>
1-6	15" (381 mm)	10.30 oz/293 g	<b>TLX01013</b>
7-8	15" (381 mm)	13 oz/369 g	<b>TLX01014</b>

**Ordering Example**  
 Model ..... **TLP0**  
 Category 1-8..... **3**  
 Foot Size ..... **24**  
 Side..... **L**  
**TLP0324L**  
 Foot Cover ..... **FCT0**  
 Foot Size ..... **24**  
 Side..... **L**  
**FCT024L**  
 \*add a BR suffix for Brown foot cover

Low Impact Category Selection Chart									
Patient Weight lbs	97-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324
<b>Foot Size</b>									
23-24	1	1	2	3	4	5	—	—	—
25-26	—	—	2	3	4	5	6	7	—
27-28	—	—	—	3	4	5	6	7	8
29-30	—	—	—	—	4	5	6	7	8
Moderate Impact Category Selection Chart									
23-24	1	2	3	4	5	—	—	—	—
25-26	—	2	3	4	5	6	7	8	—
27-28	—	—	3	4	5	6	7	8	—
29-30	—	—	—	4	5	6	7	8	—

### Flex-Foot<sup>®</sup> Axia<sup>™</sup>

- Multi-axial foot with proportional energy response and sandal toe design
- Guided Roll-Over enhances stability and gait dynamics
- Rated for low to moderate impact levels
- Available sizes: 23-30 cm; <sup>3</sup>/<sub>8</sub>" heel height
- Add BR to model number when ordering Brown\* foot cover
- Weight rating up to 309 lbs (140 kg)
- Clearance: 2<sup>15</sup>/<sub>16</sub>"
- 36 month warranty; foot covers carry a 2 month warranty

Description	Model No.
Axia <sup>™</sup>	LAP0
Foot Cover	FCT0

**Ordering Example**

Model . . . . . **LAP0**

Category 1-4 . . . . . **1**

Foot Size . . . . . **24**

Side . . . . . **L**

**LAP0124L**

Foot Cover . . . . . **FCT0**

Foot Size . . . . . **24**

Side . . . . . **L**

**FCT024L**

\*add a BR suffix for Brown foot cover



**Low Impact Category Selection Chart**

Foot Size (cm) Weight Range	23	24	25	26	27	28	29	30
99-150 lbs	0	0	0	0	1	1	—	—
151-194 lbs	1	1	1	1	1	1	—	—
195-220 lbs	1	1	2	2	2	2	4	4
221-256 lbs	—	—	2	2	2	2	4	4
257-276 lbs	—	—	3	3	3	3	4	4
277-309 lbs	—	—	—	—	—	—	4	4

**Moderate Impact Category Selection Chart**

Foot Size (cm) Weight Range	23	24	25	26	27	28	29	30
99-132 lbs	0	0	0	0	1	1	—	—
133-170 lbs	1	1	1	1	1	1	—	—
171-198 lbs	1	1	2	2	2	2	4	4
199-220 lbs	—	—	2	2	2	2	4	4
221-256 lbs	—	—	3	3	3	3	4	4
257-276 lbs	—	—	—	—	—	—	4	4

**Flex-Foot® Assure**

- Lightweight foot with full length toe lever
- Low build height accommodates long residual limbs
- Low impact level
- Available sizes 19-30 cm; categories 1-5
- Weight rating up to 300 lbs (136 kg)
- Heel height: 3/8"
- Clearance: 5"
- Add BR to model number when ordering Brown\* foot cover
- 36 month warranty
- 2 month warranty for foot cover

Description	Model No.
Flex Foot with clamp (option sold separately)	<b>FAX0</b>
Flex Foot with pyramid	<b>FAP0</b>
Foot Cover (included with foot)	<b>FCA0</b>



30 mm Fiberglass Tube Kit					
Category	Wall	Wall Thickness	Length	Weight	Model No.
Red	1-3	1/8" (3 mm)	11" (279 mm)	10.30 oz/292 g	<b>SFX00001</b>
Red	1-3	1/8" (3 mm)	15" (381 mm)	11.7 oz/331 g	<b>SFX00002</b>
Green	4-5	3/16" (4 mm)	11" (279 mm)	11.7 oz/331 g	<b>SFX00003</b>
Green	4-5	3/16" (4 mm)	15" (381 mm)	14 oz/401 g	<b>SFX00004</b>

Flex-Foot Assure Category Selection Guide					
Patient Weight (lbs)	99-135	136-176	177-214	215-251	252-300
Impact Level	1	2	3	4	5

**Ordering Example**

Model ..... **FAX0**  
 Category ..... **1**  
 Size ..... **19**  
 Side ..... **L**  
**FAX0119L**

Foot Cover ..... **FCA0**  
 Foot Size ..... **19**  
 Side ..... **L**  
**FCA019L**

\*add a BR suffix for Brown foot cover

### Sure-Flex®

- Lightweight, carbon fiber foot with moderate energy return
- Male pyramid or 30 mm tube option; bonded foot cover
- Add BR to model number when ordering Brown\* foot cover
- Low impact level
- Available sizes: 22-30 cm; 1/2" heel height
- Weight rating up to 300 lbs (136 kg)
- Clearance: 5 1/8"
- 24 month warranty for foot module
- 2 month warranty for foot cover

Description	Model No.
Male Pyramid	<b>SFP0</b>
Tube Clamp (sold separately)	<b>SFX0</b>

30 mm Fiberglass Tube Kit				
Category	Wall Thickness	Length	Weight	Model No.
1-4	1/8" (3 mm)	11" (279 mm)	10.20 oz/292 g	<b>SFX00001</b>
1-4	1/8" (3 mm)	15" (381 mm)	11.70 oz/333 g	<b>SFX00002</b>
5-6	3/16" (4 mm)	11" (279 mm)	11.60 oz/331 g	<b>SFX00003</b>
5-6	3/16" (4 mm)	15" (381 mm)	14.10 oz/401 g	<b>SFX00004</b>

Sure-Flex Category Selection Chart	
Category	Weight lbs/kg
1	100-119/45-54
2	120-143/55-65
3	144-174/66-79
4	175-209/80-95
5	210-249/96-113
6	250-300/114-136



**Ordering Example**  
 Model . . . . . **SFP0**  
 Category 0-3 . . . . . **1**  
 Foot Size . . . . . **22**  
 Side . . . . . **L**  
**SFP0122L**

**Flex-Foot® Balance**

- A comfortable multi-axial foot that combines proven durability and stability with comfort
- Full length toe lever and integrated multi-axial design offer a smooth transition from heel strike to toe off
- Low impact level
- Available sizes: 21-28 cm; categories 0-5
- Weight rating up to 300 lbs (136 kg)
- Heel height: 3/8"
- Clearance: 2 7/16"
- 12 month warranty for foot module; 2 month warranty for foot cover



Description	Model No.
Flex-Foot Balance	<b>FBP0</b>
Foot Cover (included with foot; replacement only)	<b>FCT0</b>

**Flex-Foot Balance Category Selection Guide**

Weight (lbs)	79-97	98-134	135-176	177-214	215-251	252-300
Impact Level	0	1	2	3	4	5

**Ordering Example**

Model ..... **FBP0**  
 Category 0-3 ..... **1**  
 Foot Size ..... **23**  
 Side ..... **L**  
**FBP0123L**  
 Foot Cover ..... **FCT0**  
 Foot Size ..... **23**  
 Side ..... **L**  
**FCT023L**

\*add a BR suffix for Brown foot cover

**Elation™**

- Energy storing with progressive stiffening and sandal toe
- Rated for low to moderate impact levels
- Available sizes: 22-28 cm
- Add BR to model number when ordering Brown\* foot cover
- Weight rating up to 220 lbs (100 kg)
- Heel height adjustable up to 50 cm (2")
- Clearance: 3 5/8"
- 12 month warranty



Description	Model No.
Male Pyramid	<b>ELP00</b>
Tube Clamp (included)	<b>ELT00</b>
Foot Cover (included with foot; replacement only)	<b>ELX00</b>

**Ordering Example**

Model ..... **ELP00**  
 Foot Size ..... **24**  
 Side ..... **L**  
**ELP0024L**  
 Foot Cover ..... **ELX00**  
 Foot Size ..... **24**  
 Side ..... **L**  
**ELX0024L**

\*add a BR suffix for Brown foot cover

### Flex-Symes™

- Specifically designed for Symes users; CarbonX active heel
- Build Height with foot cover:  
22-24 cm - 2<sup>1</sup>/<sub>16</sub>”; 25-27 cm - 2<sup>3</sup>/<sub>16</sub>”; 28-30 cm - 2<sup>1</sup>/<sub>4</sub>”
- Low to high impact level
- Add BR to model number when ordering Brown\* foot cover
- Available sizes: 22-24 cm for categories 1-6  
and sizes 25-30 for categories 1-9
- Heel height: 3/8”
- Weight rating up to 365 lbs (166 kg)
- Clearance: 2<sup>3</sup>/<sub>8</sub>”
- 36 month warranty



Description	Model No.
Circular Alignment Wedge	<b>SLP0</b>
Without Adapter	<b>SLX0</b>
Foot Cover	<b>FCA0</b>

Ordering Example	
Model .....	<b>SLP0</b>
Category 1-9.....	<b>1</b>
Foot Size .....	<b>30</b>
Split Toe .....	<b>0</b>
	<b>SLP01300</b>
Foot Cover .....	<b>FCA0</b>
Foot Size .....	<b>30</b>
Side.....	<b>L</b>
	<b>FCA030L</b>
*add a BR suffix for Brown foot cover	

Flex-Symes Category Selection Guide										
Weight (lbs)	97-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324	325-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Mod. Impact Level	1	2	3	4	5	6	7	8	9	*
High Impact Level	2	3	4	5	6	7	8	9	*	*
*Special Order										

**Chopart™**

- Specifically designed for those with chopart amputations as well as lisfranc, pirogoff, boyd and partial foot prostheses.
- Build Height: 1/4" (size 27 cm)
- Low to high impact level
- Standard order for:  
22-24 cm for categories 1-5  
25-26 cm for categories 2-6  
27-28 cm for categories 3-7  
29-30 cm for categories 4-8
- Heel height: 3/8"
- Weight rating up to 324 lbs (147 kg)
- 24 month warranty



**Ordering Example**  
 Model ..... **CHPOO**  
 Category 1-9..... **2**  
 Foot Size ..... **27**  
**CHPOO27**  
 Foot Cover ..... **FCA0**  
 Foot Size ..... **30**  
 Side..... **L**  
**FCA030L**  
 \*add a BR suffix for Brown foot cover

Description	Model No.
Chopart Kit	<b>CHPOO</b>
Foot Cover	<b>FCA0</b>

Chopart Category Selection Guide									
Weight (lbs)	97-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324
Low Impact Level	1	1	2	3	4	5	6	7	8
Mod. Impact Level	1	2	3	4	5	6	7	8	9
High Impact Level	2	3	4	5	6	7	8	9	*
*Special Order									

**Flex-Run™**

- Increased cushioning, lower frequency dynamic make this foot ideal for transfemoral and transtibial users serious about fitness
- Impact Level: High to Extreme
- Weight rating up to 325 lbs (130 kg)
- Weight of Foot: 11.9 oz (336 g)
- Clearance: 11" (279 mm): Flex-Run Junior - 8 1/4" (209 mm)
- One size per category
- Adapter Options: Male or Female bolt on pyramid
- 12 month warranty



Description	Model No.
Flex-Run Kit	<b>FSX001</b>
Flex-Run Junior Kit	<b>FSX107</b>
<b>Components</b>	
Attachment Plate	<b>FWX40001</b>
Male Pyramid Foot Adapter Set (75 g, M10 x 30 mm)	<b>FWX00030</b>
Female Pyramid Foot Adapter Set (100 g, M10 x 35 mm)	<b>CM150161</b>

Flex-Run Category Selection Guide									
Patient Weight (lbs)	97-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	
High to Extreme Impact Level	2	3	4	5	6	7	8	9	



## Trulife - Seattle Catalyst 9

- Split keel provides inversion/eversion for stability on uneven surfaces
- Shock absorption starting at heel strike
- Smooth rollover; controlled plantar flexion
- Narrow profile, solid toe foot
- 9 activity/weight categories to fit broad range of amputee needs
- Adjustable heel height with heel wedges
- Medium to high activity level (K3-K4)
- Weight rating up to 366 lbs
- Available sizes: 22-30 cm; minimum clearance 7.1"
- Includes: foot shell, spectra sock, shoe horn and 3 heel wedges
- 36 month warranty on all hardware, shank and keel components; 6 month warranty on foot shell

### Keel/Shank Selection Instructions

To optimize the selection and ensure amputee's safety, follow the three steps below to determine the appropriate keel/shank.

1. In the appropriate Keel/Shank Selection Table, locate the column that corresponds with the amputee's activity level.
2. With the selected column, locate the amputee's weight.
3. If the amputee has a long BK, carries heavy loads, or will reach a higher activity level within a year, choose the next higher category.



Description	Replacement Foot Shell	Model No.
Caucasian	<b>SFC290</b>	<b>SHF390</b>
Light Brown	<b>SFC291</b>	<b>SHF391</b>
Dark Brown	<b>SFC293</b>	<b>SHF393</b>

**Ordering Example**

Model . . . . . **SHF390**  
 Foot Size . . . . . **24**  
 Side . . . . . **L**  
 Category . . . . . **3**  
 Size Category\* . . . . . **S**  
 Size\*\* . . . . . **4**  
**SHF390-24-L-3S4**

Foot Cover . . . . . **SFC290**  
 Foot Size . . . . . **24**  
 Side . . . . . **L**  
**SFC290-24-L**

\*22-25=S and 26-30=L  
 \*\*Use last digit of size

Activity/Weight Category	Low	Medium-Low	Medium	Medium-High	High
	Walking with Aid	Limited Walking	Walking on Uneven Surfaces	Light Sports	Running, Basketball
9	—	—	310-366 lb 141-166 kg	271-320 lb 123-145 kg	232-274 lb 106-125 kg
8	—	305-366 lb 139-166 kg	266-309 lb 121-140 kg	232-270 lb 106-122 kg	199-231 lb 91-105 kg
7	288-366 lb 131-166 kg	255-304 lb 116-138 kg	221-265 lb 100-120 kg	194-231 lb 89-105 kg	166-198 lb 75-90 kg
6	259-287 lb 118-130 kg	229-254 lb 104-115 kg	199-220 lb 91-99 kg	175-193 lb 80-88 kg	150-165 lb 68-74 kg
5	230-258 lb 105-117 kg	204-228 lb 93-103 kg	177-198 lb 81-90 kg	155-174 lb 70-79 kg	133-149 lb 61-67 kg
4	199-229 lb 91-104 kg	180-203 lb 82-92 kg	155-176 lb 70-80 kg	136-154 lb 62-69 kg	117-132 lb 53-60 kg
3	177-198 lb 81-90 kg	155-179 lb 70-81 kg	144-154 lb 66-69 kg	126-135 lb 57-61 kg	107-116 lb 49-52 kg
2	133-176 lb 61-80 kg	133-154 lb 61-69 kg	132-143 lb 60-65 kg	117-125 lb 53-56 kg	101-106 lb 46-48 kg
1	<132 lb <60 kg	<132 lb <60 kg	<132 lb <60 kg	<117 lb <53 kg	<101 lb <46 kg

**Seattle Catalyst™**

- The “S” shank of the Seattle Catalyst dynamic foot increases performance and versatility over the popular “C” shaped shanks
- Provides enhanced shock absorption, controlled plantar flexion, smooth rollover and split keel stability
- Solid toe
- Built for medium to high activity amputees K3-K4
- Weight rating up to 366 lbs
- Adjustable heel height with heel wedges
- Available sizes: 22-30 cm
- Includes: foot shell, spectra sock, shoe horn and 3 heel wedges
- 36 month warranty on all hardware, shank and keel components; 6 month warranty on foot shell

**Keel/Shank Selection Instructions**

To optimize the selection and ensure amputee’s safety, follow the three steps below to determine the appropriate keel/shank.

1. In the appropriate Keel/Shank Selection Table, locate the column that corresponds with the amputee’s activity level.
2. With the selected column, locate the amputee’s weight.
3. If the amputee has a long BK, carries heavy loads, or will reach a higher activity level within a year, choose the next higher category.



Description	Replacement Foot Shell	Model No.
Caucasian	<b>SFC290</b>	<b>SHF290</b>
Light Brown	<b>SFC291</b>	<b>SHF291</b>
Dark Brown	<b>SFC293</b>	<b>SHF293</b>

**Ordering Example**

Model . . . . . **SHF290**  
 Foot Size . . . . . **23**  
 Side . . . . . **R**  
 Category . . . . . **3**  
 Size Category\* . . . . . **S**  
 Size\*\* . . . . . **3**  
**SHF290-23-R-3S3**

Foot Cover . . . . . **SFC290**  
 Foot Size . . . . . **23**  
 Side . . . . . **R**  
**SFC290-23-R**

\*22-24=S, 25-27=M, 28-30=L  
 \*\*Use last digit of size

Activity/Weight Category	Low	Medium-Low	Medium	Medium-High	High
	Walking with Aid	Limited Walking	Walking on Uneven Surfaces	Light Sports	Running, Basketball
5	296-366 lb 134-166 kg	264-366 lb 120-166 kg	226-366 lb 103-166 kg	181-300 lb 82-136 kg	131-220 lb 59-99 kg
4	231-295 lb 105-133 kg	207-263 lb 94-119 kg	151-225 lb 69-102 kg	69-180 lb 31-81 kg	<131 lb <59 kg
3	136-230 lb 62-104 kg	118-206 lb 53-93 kg	<151 lb <69 kg	<69 lb <31 kg	—
2	<136 lb <62 kg	<118 lb <53 kg	—	—	—

## Seattle Cadence™

- C-Shaped ankle
- Multi-axial
- Solid toe
- Simulates biomechanics of natural foot
- Full length, split, carbon keel for stability on uneven surfaces
- Medium to high activity, K3-K4
- Heel height adjustable with wedges
- Weight rating up to 366 lbs (166 kg)
- Available sizes: 22-30 cm
- Includes: foot shell, spectra sock, shoe horn and 3 heel wedges
- 24 month warranty keel; 6 month warranty on foot shell and spectra sock

### Category Selection Instructions

To optimize the selection and ensure amputee's safety, follow the two steps below to determine the appropriate category.

1. Locate the column that corresponds with the amputee's activity level
2. Within the selected column, locate the amputee's weight

**Important Note:** If the amputee has a long BK, carries heavy loads or will reach a higher activity level within a year, choose the next category higher.

**Warning:** Choosing a lower strength category than what is suggested based on the above procedure and patient data will void the warranty.



**Ordering Example**

Model . . . . . **SHF190**  
 Foot Size . . . . . **22**  
 Side . . . . . **L**  
 Category . . . . . **2**  
 Size Category\* . . . . . **S**  
 Size\*\* . . . . . **2**  
**SHF190-22-L-2S2**

Foot Cover . . . . . **SFC290**  
 Foot Size . . . . . **22**  
 Side . . . . . **L**  
**SFC290-22-L**

\*22-25=S and 26-30=L  
 \*\*Use last digit of size

Description	Replacement Foot Shell	Model No.
Caucasian	<b>SFC290</b>	<b>SHF190</b>
Light Brown	<b>SFC291</b>	<b>SHF191</b>
Dark Brown	<b>SFC293</b>	<b>SHF193</b>

Activity/Weight Category	Low	Medium-Low	Medium	Medium-High	High
	Walking with Aid	Limited Walking	Walking on Uneven Surfaces	Light Sports	Running, Basketball
5	296-366 lb	264-366 lb	226-366 lb	181-300 lb	131-220 lb
	134-166 kg	120-166 kg	103-166 kg	82-136 kg	59-99 kg
4	231-295 lb	207-263 lb	151-225 lb	69-180 lb	<131 lb
	105-133 kg	94-119 kg	69-102 kg	31-81 kg	<59 kg
3	136-230 lb	118-206 lb	<151 lb	<69 lb	—
	62-104 kg	53-93 kg	<69 kg	<31 kg	—
2	<136 lb	<118 lb	—	—	—
	<62 kg	<53 kg	—	—	—

## High Performance Foot Service Tool

- Kit includes: Tool Base, Pin Driver, Pin Setter, Drift Pin
- Available for Cadence, Catalyst 9 and Catalyst

Model No. **SHF-STK**

### Seattle Turas Foot

- Integrated multi-axial composite ankle
- Dynamic split keel design
- Enhanced plantar flexion at heel strike provides early foot stability
- Removable, split toe foot shell included
- Indications: Medium activity, K3
- Heel height: 3/8"
- Weight rating up to 250 lbs
- Minimum clearance: 3.5"
- Available sizes: 22-29 cm
- 24 month warranty keel; 6 month warranty on foot shell

#### Keel Selection Instructions

To optimize the selection and ensure amputee's safety, follow the two steps below to determine the appropriate keel.

1. In the chart below, select amputee's foot size
2. Within the selected column, locate the patient's weight.

**Important Note:** If the amputee has a long BK, carries heavy loads or will reach a higher activity level within a year, choose the next category higher.

Description	Replacement Foot Shell	Model No.
Caucasian	<b>SFC300</b>	<b>SLP300</b>
Light Brown	<b>SFC301</b>	<b>SLP301</b>
Dark Brown	<b>SFC303</b>	<b>SLP303</b>



**Ordering Example**

Model . . . . . **SLP300**

Size . . . . . **22**

Side . . . . . **L**

Keel . . . . . **J7**

Size\* . . . . . **2**

**SLP300-22-L-J72**

Foot Cover . . . . . **SFC300**

Foot Size . . . . . **22**

Side . . . . . **L**

**SFC300-22-L**

\*Use last digit of size



*Removable, split toe foot shell*

Keel	Low	Medium-Low	Medium	Medium-High	High
	Walking with Aid	Limited Walking	Walking on Uneven Surfaces	Light Sports	Running, Basketball
J8	—	—	201-250 lb 91-113 kg	161-200 lb 73-90 kg	136-165 lb 62-74 kg
J7	—	191-250 lb 87-113 kg	151-200 lb 69-90 kg	121-160 lb 55-72 kg	101-135 lb 46-61 kg
J6	180-250 lb 81-113 kg	130-190 lb 59-86 kg	100-150 lb 45-68 kg	80-120 lb 36-54 kg	65-100 lb 29-45 kg
J5	<180 lb <81 kg	130 lb 59 kg	<100 lb <45 kg	<80 lb <36 kg	<65 lb <29 kg

### Seattle LP Foot

- Multi-axial, low profile solid toe foot and foot shell
- C-shape replicates true, anatomical plantar flexion
- Full length, split composite keel for stability
- Removable foot shell
- Ideal for amputees with long residual limbs
- Weight rating up to 300 lbs
- Minimum clearance: 3.6"
- Medium to high activity level (K3-K4)
- Available sizes: 22-30 cm available in categories 1-3; sizes 26-30 cm only available in category 4
- Heel height: 3/8"
- 36 month warranty; 6 month warranty on foot shell



Description	Replacement Foot Shell	Model No.
Caucasian	<b>SFC350</b>	<b>SLP350</b>
Light Brown	<b>SFC351</b>	<b>SLP351</b>
Dark Brown	<b>SFC353</b>	<b>SLP353</b>

**Ordering Example**

Model . . . . . **SLP350**  
 Foot Size . . . . . **22**  
 Side . . . . . **L**  
 Category . . . . . **1**  
 Size Category\* . . . . . **S**  
 Size\*\* . . . . . **2**  
**SLP350-22-L-1S2**

Foot Cover . . . . . **SFC350**  
 Foot Size . . . . . **22**  
 Side . . . . . **L**  
**SFC350-22-L**

\*22-25=S and 26-30=L  
 \*\*Use last digit of size



*Removable low profile foot shell*

Activity/ Weight Category	Low	Medium-Low	Medium	Medium-High	High
	Walking with Aid	Limited Walking	Walking on Uneven Surfaces	Light Sports	Running, Basketball
4	223-300 lb 101-136 kg	223-300 lb 101-136 kg	223-300 lb 101-136 kg	194-260 lb 88-118 kg	179-240 lb 81-109 kg
3	179-220 lb 81-100 kg	179-220 lb 81-100 kg	179-220 lb 81-100 kg	163-192 lb 74-87 kg	150-176 lb 68-80 kg
2	134-176 lb 61-80 kg	134-176 lb 61-80 kg	134-176 lb 61-80 kg	123-161 lb 56-73 kg	112-148 lb 51-67 kg
1	<132 lb <60 kg	<132 lb <60 kg	<132 lb <60 kg	<121 lb <55 kg	<110 lb <50 kg

**Seattle LightFoot2 with Titanium Pyramid**

- Effortless alignment
- Split toe cosmesis
- Medium activity, K3
- Available sizes: 22-30 cm for 3/4" heel height feet; 22-29 cm for 3/8" heel height feet
- Weight rating up to 300 lbs
- Minimum clearance: 3.5"
- 24 month warranty on all hardware, shank and keel components; 12 month warranty on foot shell

**Keel Selection Instructions**

To optimize the selection and ensure amputee's safety, follow the two steps below to determine the appropriate keel.

1. Locate the column that corresponds with the amputee's activity level.
2. Within the selected column, locate the amputee's weight.

**Important Note:** If the amputee has a long BK, carries heavy loads or will reach a higher activity level within a year, choose the next category higher.

**Warning:** Choosing a lower strength keel than what is suggested based on the above procedure and patient data will void the warranty. If your patient's weight exceeds the limits of the chart, please browse the catalog for other foot options or call SPS Customer Service.



Description	Model No.	
	3/8" Heel	3/4" Heel
Caucasian	SLF195	SLF190
Light Brown	SLF196	SLF191
Dark Brown	SLF198	SLF193

Ordering Example	
Model . . . . .	SLF190
Size . . . . .	22
Side . . . . .	L
Keel . . . . .	G9
<b>SLF190-22-L-G9</b>	

Keel*	Low	Medium-Low	Medium	Medium-High	High
	Walking with Aid	Limited Walking	Walking on Uneven Surfaces	Light Sports	Running, Basketball
G9/H9	—	256-300 lb 116-136 kg	251-300 lb 114-136 kg	201-235 lb 91-106 kg	170-200 lb 77-90 kg
G8/H8	—	191-255 lb 87-115 kg	201-250 lb 91-113 kg	161-200 lb 73-90 kg	136-169 lb 62-76 kg
G7/H7	271-300 lb 123-136 kg	166-190 lb 75-86 kg	151-200 lb 69-90 kg	121-160 lb 55-72 kg	101-135 lb 46-61 kg
G6/H6	181-270 lb 82-122 kg	131-165 lb 59-74 kg	101-150 lb 45-68 kg	81-120 lb 36-54 kg	66-100 lb 29-45 kg
G5/H5	<181 lb <82 kg	<131 lb <59 kg	<101 lb <45 kg	<81 lb <36 kg	<66 lb <29 kg

\*Note: "G" keel is for the 3/4" heel height; "H" keel is for the 3/8" heel height.

## Seattle LightFoot2

- Enhanced rollover
- Durable, natural looking, split toe cosmesis
- Appropriate for all lower extremity amputees, including Symes
- Medium activity, K3
- Available sizes: 22-30 cm for 3/4" heel height feet  
22-29 cm for 3/8" heel height feet
- Weight rating up to 300 lbs
- Minimum clearance: 1.5"
- 24 month warranty on all hardware, shank and keel components;  
12 month warranty on foot shell



### Keel Selection Instructions

To optimize the selection and ensure amputee's safety, follow the two steps below to determine the appropriate keel.

1. Locate the column that corresponds with the amputee's activity level.
2. Within the selected column, locate the amputee's weight.

**Important Note:** If the amputee has a long BK, carries heavy loads or will reach a higher activity level within a year, choose the next category higher.

**Warning:** Choosing a lower strength keel than what is suggested based on the above procedure and patient data will void the warranty. If your patient's weight exceeds the limits of the chart, please browse the catalog for other foot options or call SPS Customer Service.

Description	Model No.	
	3/8" Heel	3/4" Heel
Caucasian	SLF165	SLF160
Light Brown	SLF166	SLF161
Dark Brown	SLF168	SLF163

Ordering Example	
Model .....	SLF165
Size .....	24
Side .....	R
Keel .....	H5
<b>SLF165-24-R-H5</b>	

Keel*	Low	Medium-Low	Medium	Medium-High	High
	Walking with Aid	Limited Walking	Walking on Uneven Surfaces	Light Sports	Running, Basketball
G9/H9	—	256-300 lb 116-136 kg	251-300 lb 114-136 kg	201-235 lb 91-106 kg	170-200 lb 77-90 kg
G8/H8	—	191-255 lb 87-115 kg	201-250 lb 91-113 kg	161-200 lb 73-90 kg	136-169 lb 62-76 kg
G7/H7	271-300 lb 123-136 kg	166-190 lb 75-86 kg	151-200 lb 69-90 kg	121-160 lb 55-72 kg	101-135 lb 46-61 kg
G6/H6	181-270 lb 82-122 kg	131-165 lb 59-74 kg	101-150 lb 45-68 kg	81-120 lb 36-54 kg	66-100 lb 29-45 kg
G5/H5	<181 lb <82 kg	<131 lb <59 kg	<101 lb <45 kg	<81 lb <36 kg	<66 lb <29 kg

\*Note: "G" keel is for the 3/4" heel height; "H" keel is for the 3/8" heel height.

### Seattle Carbon LightFoot

- Provides increased dynamic response and improved ML compliance with its split composite keel
- Lightweight design with dynamic response
- Split, composite keel
- Natural looking, split toe cosmesis
- Medium activity, K3
- Heel height: 3/8"
- Weight rating up to 275 lbs
- Minimum clearance: 1.5"
- 24 month warranty on all hardware, shank and keel components; 12 month warranty on foot shell



#### Keel Selection Instructions

To optimize the selection and ensure amputee's safety, follow the two steps below to determine the appropriate keel.

1. Locate the column that corresponds with the amputee's activity level.
2. Within the selected column, locate the amputee's weight.

**Important Note:** If the amputee has a long BK, carries heavy loads or will reach a higher activity level within a year, choose the next category higher.

**Warning:** Choosing a lower strength keel than what is suggested based on the above procedure and patient data will void the warranty. If your patient's weight exceeds the limits of the chart, please browse the catalog for other foot options or call SPS Customer Service.

Description	Foot Size	Model No
Caucasian	22-29 cm	<b>SCF185</b>
Light Brown	22-29 cm	<b>SCF186</b>
Dark Brown	22-29 cm	<b>SCF188</b>
Ladies Caucasian	22-28 cm	<b>SCF195</b>
Ladies Light Brown	22-28 cm	<b>SCF196</b>
Ladies Dark Brown	22-28 cm	<b>SCF198</b>

Ordering Example	
Model . . . . .	<b>SCF185</b>
Size . . . . .	<b>23</b>
Side . . . . .	<b>R</b>
Keel . . . . .	<b>N1</b>
<b>SCF185-23-R-N1</b>	

Keel	Low	Medium-Low	Medium	Medium-High	High
	Walking with Aid	Limited Walking	Walking on Uneven Surfaces	Light Sports	Running, Basketball
N3	251-275 lb	211-275 lb	176-270 lb	146-235 lb	126-200 lb
	114-125 kg	96-125 kg	80-122 kg	66-106 kg	57-90 kg
N2	150-250 lb	120-210 lb	100-175 lb	85-145 lb	75-125 lb
	68-113 kg	54-95 kg	45-79 kg	38-65 kg	34-56 kg
N1	<150 lb	<120 lb	<100 lb	<85 lb	<75 lb
	<68 kg	<54 kg	45 kg	<38 kg	<34 kg



## Seattle LightFoot

- Durable, natural looking
- Split toe cosmesis
- Appropriate for all lower extremity amputees, including Symes
- Medium activity, K3
- Weight rating up to 300 lbs
- Heel height:  $\frac{3}{4}$ " or  $\frac{3}{8}$ "
- Available sizes: 22-30 cm
- Minimum clearance: 1.5"
- 24 month warranty on all hardware, shank and keel components;  
12 month warranty on foot shell



### Keel Selection Instructions

To optimize the selection and ensure amputee's safety, follow the two steps below to determine the appropriate keel.

1. Locate the column that corresponds with the amputee's activity level.
2. Within the selected column, locate the amputee's weight.

**Important Note:** If the amputee has a long BK, carries heavy loads or will reach a higher activity level within a year, choose the next category higher.

**Warning:** Choosing a lower strength keel than what is suggested based on the above procedure and patient data will void the warranty. If your patient's weight exceeds the limits of the chart, please browse the catalog for other foot options or call SPS Customer Service.

Description	Model No.	
	$\frac{3}{8}$ " Heel	$\frac{3}{4}$ " Heel
Caucasian	SLF135	SLF130
Light Brown	SLF136	SLF139
Dark Brown	SLF138	SLF141

Ordering Example	
Model .....	SLF130
Size .....	28
Side .....	L
Keel .....	C5
<b>SLF130-28-L-C5</b>	

Keel*	Low	Medium-Low	Medium	Medium-High	High
	Walking with Aid	Limited Walking	Walking on Uneven Surfaces	Light Sports	Running, Basketball
F9	—	256-300 lb 116-136 kg	251-300 lb 114-136 kg	201-235 lb 91-106 kg	170-200 lb 77-90 kg
F8	—	191-255 lb 87-115 kg	201-250 lb 91-113 kg	161-200 lb 73-90 kg	136-169 lb 62-76 kg
C7/F7	271-300 lb 123-136 kg	166-190 lb 75-86 kg	151-200 lb 69-90 kg	121-160 lb 55-72 kg	101-135 lb 46-61 kg
C6/F6	181-270 lb 82-122 kg	131-165 lb 59-74 kg	101-150 lb 45-68 kg	81-120 lb 36-54 kg	66-100 lb 29-45 kg
C5/F5	<181 lb <82 kg	<131 lb <59 kg	<101 lb <45 kg	<81 lb <36 kg	<66 lb <29 kg

\*Note: "C" keel is for the  $\frac{3}{4}$ " heel height; "F" keel is for the  $\frac{3}{8}$ " heel height.

### Seattle Lifecast

- Original energy storing, dynamic response foot; offers high levels of function and cosmesis to amputees
- “J” design keel
- Split toe
- Life-like cosmesis developed from molds of real feet
- Medium activity, K3
- Weight rating up to 252 lb
- Heel Rise: 3/4” (19 mm)
- Available sizes: 22-29 cm
- Minimum clearance: 3.1” (22-26 cm); 3.4” (27-29 cm)
- 24 month warranty on all hardware, shank and keel components; 12 month warranty on foot shell



#### Keel Selection Instructions

To optimize the selection and ensure amputee’s safety, follow the two steps below to determine the appropriate keel.

1. In the chart below, select amputee’s foot size
2. Within the selected row, locate the patient’s weight

**Important Note:** Seattle Lifecast Keel Selection chart is for medium activity level; choose the next higher keel if amputee is a higher activity level, has a long BK, carries heavy loads or will reach a higher activity level within a year

**Warning:** Choosing a lower strength keel than what is suggested based on the above procedure and patient data will void the warranty. If your patient’s weight exceeds the limits of the chart, please browse the catalog for other foot options or call SPS Customer Service.

Description	Model No.
Caucasian	SFH110
Light Brown	SFH111
Dark Brown	SFH113

Ordering Example	
Model .....	SFH110
Size .....	23
Side .....	L
Keel .....	A4S
<b>SFH110-23-L-A4S</b>	

Size	Weight	Keel
22 cm	<111 lb	A0S
	111-130 lb	A1S
	131-150 lb	A2S
23-26 cm	<155 lb	A3S
	155-185 lb	A4S
	186-220 lb	A5S
	221-252 lb	A6S
27-29 cm	<151 lb	A3L
	151-180 lb	A4L
	181-205 lb	A5L
	206-245 lb	A6L

### Kinetic Foot

- K2 activity level with integrated multi-axial ankle
- Split toe cosmesis (non-removable)
- Female proximal adapter is integrated for ease of adjustment
- Heel Rise: 3/8" (9.5 mm); available sizes: 22-30 cm
- Minimum clearance: 2.8"
- Weight rating up to 300 lbs for medium activity
- 24 month warranty on keel; 12 month warranty on foot shell

Description	Model No.
Caucasian	SKF200
Light Brown	SKF201
Dark Brown	SKF203



Seattle foot shell is not removable

Bumper Selection Chart										
Foot Length Weight	22	23	24	25	26	27	28	29	30	
221-300 lb/100-136 kg	3	3	4	4	4	4	4	4	4	
177-220 lb/81-99 kg	2	2	3	3	3	3	3	3	3	
132-176 lb/60-80 kg	1	1	2	2	2	2	2	2	2	
<132 lb/<60 kg	1	1	1	1	1	1	1	1	1	

Ordering Example	
Model .....	SKF200
Size .....	22
Side .....	L
Bumper .....	1
<b>SKF200-22-L-1</b>	

### Seattle Ladies Foot

- Slim design and low profile, accommodates women's and other narrow shoe styles
- Split toe
- Weight rating up to 245 lb for K3 activity level
- Heel Rise: 3/8" (9.5 mm)
- Available sizes: 22-28 cm
- Minimum clearance: 3"
- 24 month warranty on all hardware, shank and keel components; 12 month warranty on foot shell

Description	Model No.
Caucasian	SFH120
Light Brown	SFH121
Dark Brown	SFH123

Weight lbs	<115	115-135	136-165	166-195	196-245
Foot Size	Keel Selection				
22-24 cm	B2S	B3S	B4S	B5S	B6S
25-28 cm	B2L	B3L	B4L	B5L	B6L

Ordering Example	
Model .....	SFH120
Size .....	22
Side .....	L
Keel .....	B3S
<b>SFH120-22-L-B3S</b>	



#### Keel Selection Instructions

To optimize selection and ensure amputee's safety, follow the steps below to determine the appropriate keel.

1. In the chart left, select amputee's foot size
2. Within the selected column, locate the patient's weight.

**Important Note:** Seattle Ladies Keel Selection chart is for medium activity level; choose the next higher keel if amputee is a higher activity level, has a long BK, carries heavy loads or will reach a higher activity level within a year

**Warning:** Choosing a lower strength keel than what is suggested based on the above procedure and patient data will void the warranty. If your patient's weight exceeds the limits of the chart, please browse the catalog for other foot options or call SPS Customer Service.

**Natural Foot**

- Performance characteristics of a compliant foot with the traditional LightFoot Cosmesis
- Split toe
- Appropriate low to medium-low activity levels
- Available sizes: 22-30 cm; heel rise: 3/8" (9.5 mm)
- Minimum clearance: 3"
- Weight rating up to 225 lbs (102 kg) for medium activity
- 24 month warranty on keel; 12 month warranty on foot shell



Description	Model No.
Caucasian	<b>SNF150</b>
Light Brown	<b>SNF151</b>
Dark Brown	<b>SNF153</b>

**Ordering Example**  
 Model ..... **SNF150**  
 Size ..... **24**  
 Side ..... **R**  
**SNF150-24R**

**Seattle SACH Foot**

- Lightweight; split toe
- Indications: Low to medium-low activity level, K1
- Available sizes: 22-30 cm
- Minimum clearance: 3"
- Weight rating up to 225 lbs (102 kg) for medium activity
- 24 month warranty on keel; 12 month warranty on foot shell



Description	Model No.	
	3/8" Heel	3/4" Heel
Caucasian	<b>SCH100</b>	<b>SCH105</b>
Light Brown	<b>SCH101</b>	<b>SCH106</b>
Dark Brown	<b>SCH103</b>	<b>SCH108</b>

**Ordering Example**  
 Model ..... **SCH100**  
 Size ..... **25**  
 Side ..... **L**  
**SCH100-25L**

**Seattle Foamer Feet**

- Solid foam foot without a keel, designed to be incorporated into a prosthesis for patients with partial foot amputations
- Available sizes: 22-30 cm
- Minimum clearance: 3"
- 12 month warranty



Foot Type	Color	Model No.	
		3/8" Heel	3/4" Heel
LightFoot	Caucasian	—	<b>SFF130</b>
	Light Brown	—	<b>SFF139</b>
	Dark Brown	—	<b>SFF141</b>
Low Profile LightFoot	Caucasian	<b>SFF135</b>	—
	Light Brown	<b>SFF136</b>	—
	Dark Brown	<b>SFF138</b>	—
Lifecast	Caucasian	—	<b>SFF110</b>
Ladies Foot	Caucasian	<b>SFF120</b>	—

**Ordering Example**  
 Model ..... **SFF130**  
 Size ..... **22**  
 Side ..... **L**  
**SFF130-22-L**

### Seattle Endoskeletal Cosmesis

- The Endoskeletal Cosmesis is appropriate for transtibial amputations
- For use with the following Seattle feet: Turas, Carbon LightFoot, LightFoot2, LightFoot2 with Titanium Pyramid, LightFoot, Natural Foot, SACH or Child's Play Foot
- 24 month warranty

Foot Size	Height	Caucasian	Model No.	
			Light Brown	Dark Brown
22-26 cm	10 <sup>1</sup> / <sub>4</sub> "	SEC201	SEC206	SEC210
26-29 cm	15 <sup>1</sup> / <sub>2</sub> "	SEC202	SEC207	SEC211



**Ordering Example**  
 Model . . . . . SEC201  
 Foot Size . . . . . 23  
 Side . . . . . R  
**SEC201-23R**

### Standard Ankle Block

- Designed for use with the Seattle feet in exoskeletal construction
- High strength rigid polyurethane foam surrounding a nylon core
- Appropriate Feet: LightFoot, LightFoot2, Carbon LightFoot
- Product Weight: 426 g
- 24 month warranty

Model No. **SAB310**



### Low Profile Ankle Block

- Appropriate For Lifecast, Ladies Foot, Natural Foot, SACH Foot
- Product Weight: 369 g
- 24 month warranty

Model No. **SAB320**



### Seattle Adult Laminating Core

- High-strength, low profile alternative to the Seattle ankle block
- Requiring minimal keel to socket clearance, the Laminating Core is laminated directly to the socket,
- Appropriate Feet: Lifecast, Ladies Foot, Natural Foot, SACH Foot, LightFoot2, Carbon LightFoot, Seattle LightFoot
- Color: Black
- Product Weight: 160 g
- Patient Weight limit: 300 lbs
- 24 month warranty



**Model No.** **SLC300**

### Seattle Foot Spacers

- Designed for use with Carbon LightFoot, Lightfoot2, Lightfoot
- May be used temporarily to permit easier access to alignment devices or permanently to allow setup on existing limb systems and to accommodate patient growth
- Weight rating up to 300 lbs (136 kg) medium activity level
- 24 month warranty



Product Height	Model No.
1.5"	<b>SFS210</b>
1.65"	<b>SFS220</b>
1.8"	<b>SFS221</b>
0.25"	<b>SFS230</b>
0.5"	<b>SFS231</b>
1.0"	<b>SFS232</b>

SFS-KIT includes SFS230, SFS231, and SFS232

### Seattle Symes Nut

- Provides a broad, solid anterior support surface for the keel
- Appropriate Feet: LightFoot, LightFoot2, Carbon LightFoot
- Product Weight: 42 g; Patient Weight limit: 300 lbs
- Minimum Clearance: 2.0" (5.1 cm), including foot
- 24 month warranty

**Model No.** **SSY300**



Accessories	
Description	Model No.
Seattle Foot Bolt 10 mm x 55 mm	<b>SFB540</b>
Seattle Foot Bolt 10 mm x 32 mm	<b>SFB545</b>
Seattle Foot Bolt 10 mm x 80 mm	<b>SFB542</b>

## WillowWood DuraLite™

- Lightweight (356 g) with high energy return\*
- Proximal adapters allow for easy field height adjustment of foot clearance without drilling holes; Choose Pyramid Adapter (1), Pyramid Receiver (2), or Torsion Receiver (3)
- Moderate to high activity
- Patient maximum weight limit: 300 lbs\*\*
- Available sizes: 23-31 cm
- Heel height: 3/8" (L) or 3/4" (M)
- Maximum clearance: 13"
- Minimum DuraLite clearance:
  - Pyramid adapter: 6.15" for 3/4" heel; 6.55" for 3/8" heel
  - Pyramid receiver: 6.0" for 3/4" heel; 6.4" for 3/8" heel
  - Torsion receiver: 7.8" for 3/4" heel; 8.2" for 3/8" heel
- Option for foot to be pre-cut to minimum clearance available
- Cosmesis: Buff (B), Medium Brown (M) or Tan (T)
- 36-month warranty for foot module; 24-month warranty for Torsion Receiver, and 9-month warranty for the foot cover, all warranties are from the invoice date

### General Guidelines to Select Proper Composite Toe Spring Stiffness:

1. Calculate adjusted body weight. (Adjusted body weight = patient body weight + loads normally or routinely carried)
2. Using adjusted body weight and patient activity level, select the toe spring stiffness from the chart below
3. If the amputee's weight is near the border between two categories:
  - a. choose the higher of the two categories if the amputee is very active or has a longer residual limb
  - b. choose the lower of the two categories if the amputee is moderately active or has a very short residual limb



\*Based on size 26 cm foot, without proximal adapter, foot cover, or Spectra sock

\*\*Body weight plus any loads normally or routinely carried cannot exceed this weight limit

Description	Model No.
DuraLite	<b>DLF</b>
DuraLite w/ precut option	<b>DLF-PRECU</b>
Replacement Foot Cover	<b>FSHT</b>

DuraLite Foot Toe Spring Stiffness Chart				
Patient Weight	Activity Level	Foot Size		
		23-24 cm	25-26 cm	27-31 cm
100-120 lbs	Medium	1	1	1
	High	2	2	2
121-150 lbs	Low	1	1	1
	Medium	2	2	2
151-190 lbs	High	3	3	3
	Low	2	2	2
191-230 lbs	Medium	3	3	3
	High	4	4	4
231-270 lbs	Low	3	3	3
	Medium	4	4	4
	High	5	5	5
271-300 lbs	Low	—	4	4
	Medium	—	5	5
	High	—	—	6

Ordering Example	
Model .....	<b>DLF</b>
Size .....	<b>23</b>
Side .....	<b>L</b>
Color .....	<b>B</b>
Toe Resistance .....	<b>4</b>
Heel Height .....	<b>L</b>
Proximal Adapter .....	<b>1</b>
<b>DLF-23LB-4L-1</b>	
Foot Cover .....	<b>FSHT</b>
Size .....	<b>23</b>
Side .....	<b>L</b>
Color .....	<b>B</b>
Heel Height .....	<b>L</b>
<b>FSHT-23LB-L</b>	

**Fusion™ Foot**

- Provides multi-axial movement and vertical shock absorption
- Allows the composite shank, heel and foot plate to work as one
- Moderate to high activity
- Smooth heel-to-toe transition with inversion/eversion
- Weight rating up to 230 lbs (foot size 23-24 cm); 270 lbs (foot size 25-26 cm); 300 lbs (foot size 27-31 cm)
- Proximal Adapters/Weights: Pyramid Adapter, 205 g (1); Pyramid Receiver, 190 g (2); Torsion Receiver, 420 g (3)
- Heel height: 3/8" (L) or 3/4" (M)
- Available sizes: 23-31 cm
- Clearance for Fusion Foot:
  - Pyramid Adapter: 6.45" for 3/4" heel; 6.85" for 3/8" heel
  - Pyramid Receiver: 6.3" for 3/4" heel; 6.7" for 3/8" heel
  - Torsion Receiver: 8.1" for 3/4" heel; 8.5" for 3/8" heel
- Option for foot to be pre-cut to minimum clearance available
- Cosmesis: Buff (B), Medium Brown (M) or Tan (T)
- 36-month warranty for foot module; 24-month warranty for Torsion Receiver, and 9-month warranty for the foot cover, all warranties are from the invoice date

**General Guidelines to Select Proper Composite Toe Spring Stiffness:**

1. Calculate adjusted body weight. (Adjusted body weight = patient body weight + loads normally or routinely carried)
2. Using adjusted body weight and patient activity level, select the toe spring stiffness from the chart below
3. If the amputee's weight is near the border between two categories:
  - a. choose the higher of the two categories if the amputee is very active or has a longer residual limb
  - b. choose the lower of the two categories if the amputee is moderately active or has a very short residual limb



Description	Model No.
Fusion Foot	<b>FFT</b>
Fusion w/ precut option	<b>FFT-PRECURT</b>
Replacement Foot Cover	<b>FSHT</b>
High Activity Spectra Sock	<b>700-SS001</b>

Fusion Foot Composite Toe Spring Stiffness Chart				
Patient Weight	Activity Level	23-24 cm	25-26 cm	27-31 cm
		100-120 lbs	Medium	1
121-150 lbs	High	2	2	2
	Low	1	1	1
	Medium	2	2	2
151-190 lbs	High	3	3	3
	Low	2	2	2
	Medium	3	3	3
191-230 lbs	High	4	4	4
	Low	3	3	3
	Medium	4	4	4
231-270 lbs	High	5	5	5
	Low	—	4	4
	Medium	—	5	5
271-300 lbs	High	—	—	6
	Low	—	—	5
	Medium	—	—	6

Ordering Example	
Model .....	<b>FFT</b>
Foot size .....	<b>26</b>
Side .....	<b>L</b>
Color .....	<b>B</b>
Toe Spring Stiffness .....	<b>4</b>
Heel Height .....	<b>L</b>
Proximal Adapter .....	<b>1</b>
<b>FFT-26LB-4L-1</b>	
Foot Cover .....	<b>FSHT</b>
Size .....	<b>26</b>
Side .....	<b>L</b>
Color .....	<b>B</b>
Heel Height .....	<b>L</b>
<b>FSHT-26LB-L</b>	



## Pathfinder II® Foot

- Exceptional inversion & eversion due to revised triangular design
- Moderate to high activity
- Weight limit: 350 lbs (activity levels 1-3), 300 lbs (activity level 4)
- Low Profile accommodates amputees with longer residual limbs
- Available sizes: 23-31 cm
- Cosmesis: Buff (B), Medium Brown (M) or Tan (T)
- Heel Height: 3/8" (L), 3/4" (M)
- Proximal Adapter options: Pyramid Adapter (1), Pyramid Receiver Adapter (2), Threaded 4-hole Adapter (3), Through Hole Adapter with Valve (4), Without Slide Adjustment (5)
- 36-month warranty for foot module; 9-month warranty for the foot cover, all warranties are from the invoice date

### General Guidelines to Select Proper Composite Toe Spring Stiffness:

1. Calculate adjusted body weight. (Adjusted body weight = patient body weight + loads normally or routinely carried)
2. Using adjusted body weight and patient activity level, select the toe spring stiffness from the chart below
3. If the amputee's weight is near the border between two categories:
  - a. choose the higher of the two categories if the amputee is very active or has a longer residual limb
  - b. choose the lower of the two categories if the amputee is moderately active or has a very short residual limb



Description	Model No.
Pathfinder II Foot	PFDR
Low Profile Pathfinder II® (not pictured)	PFLP
Replacement Foot Cover	FSHA

Ordering Example	
Model . . . . .	PFDR
Size . . . . .	23
Side . . . . .	L
Color . . . . .	B
Toe Resistance . . . . .	4
Heel Height . . . . .	L
Proximal Adapter . . . . .	1
<b>PFDR-23LB-4L-1</b>	
Foot Cover . . . . .	FSHA
Size . . . . .	23
Side . . . . .	L
Color . . . . .	B
Heel Height . . . . .	L
<b>FSHA-23LB-L</b>	

Composite Toe Spring Stiffness Chart							
Toe Spring Category	Adjusted Body Weight (lbs)						
	100-135	136-170	171-205	206-235	236-280	281-315	316-350
	3	4	5	6	7	8*	9*

\*Not available in sizes 23 and 24

Clearance Range Chart for Pathfinder II Foot (from bottom of heel to top of Proximal Adapter)				
	Pathfinder II with Slide Adjustment		Pathfinder II without Slide Adjustment	
	18 mm (3/4") Heel	10 mm (3/8") Heel	18 mm (3/4") Heel	10 mm (3/8") Heel
with Threaded 4-hole adapter	8.7" (22.1 cm)	9.1" (23.1 cm)	8.3" (21.1 cm)	8.7" (22.1 cm)
with Through-hole adapter with valve	8.7" (22.1 cm)	9.1" (23.1 cm)	—	—
with Pyramid adapter	9.5" (24.1 cm)	9.8" (24.9 cm)	9.3" (23.6 cm)*	9.7" (24.6 cm)*
with Pyramid receiver adapter	9.2" (23.4 cm)	9.5" (24.1 cm)	8.9" (22.6 cm)*	9.3" (23.6 cm)*
Clearance Range Chart for Low Profile Pathfinder II Foot				
with Threaded 4-hole adapter	7.7" (19.6 cm)	8.1" (20.6 cm)	7.3" (18.5 cm)	7.7" (19.6 cm)
with Through-hole adapter with valve	7.7" (19.6 cm)	8.1" (20.6 cm)	—	—
with Pyramid adapter	8.5" (21.6 cm)	8.8" (22.4 cm)	8.3" (21.1 cm)*	8.7" (22.1 cm)*
with Pyramid receiver adapter	8.2" (20.8 cm)	8.5" (21.6 cm)	7.9" (20.1 cm)*	8.3" (21.1 cm)*

\*These heights are based on bolting a WillowWood pyramid adapter or pyramid receiver adapter on the fixed threaded 4-hole adapter

**Trailblazer™**

- Low profile design accommodates active and athletic lifestyles
- Innovative heel design allows for in office customization of foot
- Moderate to high activity
- Weight rating up to 250 lbs\*
- Heel Height: <sup>3</sup>/<sub>8</sub>" (L), clearance 2.78"; <sup>3</sup>/<sub>4</sub>" (M), clearance 2.41"
- Heel Density: (S) Soft, (M) Medium, or (F) Firm
- Available sizes: 23-31 cm
- Cosmesis: Buff (B), Medium Brown (M) or Tan (T)
- 36-month warranty for foot module; 9-month warranty for the foot cover, all warranties are from the invoice date

\*Body weight plus any loads normally or routinely carried cannot exceed this weight limit



Description	Model No.
Trailblazer Fixed Pyramid	<b>TBFP</b>
Trailblazer Rotatable Pyramid	<b>TBRP</b>
Replacement Foot Cover	<b>FSHA</b>

**Ordering Example**

Model ..... **TBFP**  
 Size ..... **24**  
 Side ..... **R**  
 Color ..... **M**  
 Toe Resistance ..... **3**  
 Heel Height ..... **L**  
 Heel Density ..... **M**

**TBFP-24RM-3L-M**

Foot Cover ..... **FSHA**  
 Size ..... **24**  
 Side ..... **R**  
 Color ..... **M**  
 Heel Height ..... **L**

**FSHA-24RM-L**

**Trailblazer™ MA**

- Low profile, lightweight foot with a high activity multi-axial ankle
- Customize amount of ankle movement with unique flexion ring; larger diameter yields increased motion (3 included)
- Compound spring system provides dramatic energy return with long term durability
- Moderate to high activity
- Weight rating up to 250 lbs\*
- Heel Height: <sup>3</sup>/<sub>8</sub>" (L) clearance 3.89"; <sup>3</sup>/<sub>4</sub>" (M) clearance 3.52"
- Available sizes: 23-31 cm
- Cosmesis: Buff (B), Medium Brown (M) or Tan (T)
- 36 month warranty from the invoice date

\*Body weight plus any loads normally or routinely carried cannot exceed this weight limit



Description	Model No.
Trailblazer MA	<b>TBMA</b>
Replacement Foot Cover	<b>FSHA</b>

**Ordering Example**

Model ..... **TBMA**  
 Size ..... **25**  
 Side ..... **L**  
 Color ..... **T**  
 Toe Resistance ..... **3**  
 Heel Height ..... **L**  
 Heel density ..... **M**

**TBMA-25LT-3L-M**

Foot Cover ..... **FSHA**  
 Size ..... **25**  
 Side ..... **L**  
 Color ..... **T**  
 Heel Height ..... **L**

**FSHA-2LT-L**

**Trailblazer Composite Toe Spring Stiffness Chart**

Adjusted Body Weight	Toe Spring Category
Up to 110 lbs	1
111-150 lbs	2
151-190 lbs	3
191-230 lbs	4
231-250 lbs	5

### Impulse® Foot

- Pathfinder II composite material foot plate attached to a Kevlar® nylon keel
- Outstanding flexibility and toe response
- Recessed keel for Endo, Exo or Symes application
- Keel insert densities (Soft, Medium, Firm) included with foot
- Moderate to high activity
- Weight rating up to 350 lbs\*
- Heel Height: Women: 3/8" (L) or 1 1/8" (H); Men: 3/4" (M) or 3/8" (L)
- Clearance: 1.38"
- Available sizes: Women: 22-27 cm, Men: 22-31 cm
- Cosmesis: Buff (B), Medium Brown (M) or Tan (T)
- 12 month warranty, (three weeks allowed for installation)

\*Body weight plus any loads normally or routinely carried cannot exceed this weight limit



Description	Model No.
Men's Impulse	<b>IMPM</b>
Women's Impulse	<b>IMPW</b>
Spacer Block	<b>2LTC-1210</b>
Symes Adapter Kit	<b>2LSC-1400</b>
Shaped Foam Ankle - Mens	<b>2LAM (size, side, color)</b>
Shaped Foam Ankle - Womens	<b>2LAW (size, side, color)</b>

Ordering Example	
Model . . . . .	<b>IMPM</b>
Size . . . . .	<b>26</b>
Side . . . . .	<b>L</b>
Color. . . . .	<b>B</b>
Toe Resistance . . . . .	<b>4</b>
Heel Height . . . . .	<b>M</b>
<b>IMPM-26LB-4M</b>	

Toe Resistance						
Very Low (1)	Low (2)	Reg. (3)	High (4)	Very High (5)	HD Toe 1 (6)	HD Toe 2 (7)
Less than 110 lbs	111-150 lbs	151-190 lbs	191-230 lbs	231-270 lbs	271-310 lbs	311-350 lbs

### Impulse® Symes Foot

- Moderate to high activity
- Symes socket adapter kit included
- Allows for angular and length adjustments
- Keel insert densities (Soft, Medium, Firm) included with foot
- Weight rating up to 250 lbs\*
- Heel Height: Women: 3/8" (L) or 1 1/8" (H); Men: 3/4" (M) or 3/8" (L)
- Available sizes: Women: 22-27 cm, Men: 22-31 cm
- Clearance: 1.75"
- Cosmesis: Buff (B), Medium Brown (M) or Tan (T)
- Toe Resistance: See chart
- 12 month warranty, (three weeks allowed for installation)

\*Body weight plus any loads normally or routinely carried cannot exceed this weight limit



Description	Model No.
Men's Impulse Symes	<b>IMSM</b>
Women's Impulse Symes	<b>IMSW</b>

Ordering Example	
Model . . . . .	<b>IMSM</b>
Size . . . . .	<b>22</b>
Side . . . . .	<b>L</b>
Color. . . . .	<b>B</b>
Toe Resistance . . . . .	<b>1</b>
Heel Height . . . . .	<b>L</b>
<b>IMSM-22LB-1L</b>	

Toe Resistance				
Very Low (1)	Low (2)	Reg. (3)	High (4)	Very High (5)
Less than 110 lbs	111-145 lbs	146-180 lbs	181-215 lbs	216-250 lbs

**Magnum Heavy Duty Foot**

- Heavy Duty System for amputees weighing up to 350 lbs
- Proven liability through laboratory and field testing
- Moderate to high activity
- Lightweight, sleek design
- Available separately or as part of three different endoskeletal kits
- Weight rating up to\*:
  - 350 lbs (160 kg) for U.S. Activity Levels 1-3
  - 300 lbs (135 kg) for U.S. Activity Level 4
- Heel Height: Women: 3/8" (L), Men: 3/4" (M) or 3/8" (L)
- Available sizes: Women, 22-27 cm; Men, 23-31 cm
- Clearance: Foot, 3.5"; Complete system, 8"
- Cosmesis: Buff (B), Medium Brown (M) or Tan (T)
- Heel Insert Densities: Medium
- 12 month warranty, (three weeks allowed for installation)

\*Body weight plus any loads normally or routinely carried cannot exceed this weight limit



Description	Model No.
Mens	<b>MHDM</b>
Womens	<b>MHDW</b>

Toe Resistance		
Low (6)	Regular (7)	High (8)
251-275 lbs	276-300 lbs	301-350 lbs

Ordering Example	
Model . . . . .	<b>MHDM</b>
Size . . . . .	<b>26</b>
Side . . . . .	<b>L</b>
Color . . . . .	<b>B</b>
Toe Resistance . . . . .	<b>6</b>
Heel Height . . . . .	<b>L</b>
<b>MHDM-26LB-6L</b>	

**Carbon Copy 2<sup>®</sup> Foot**

- Energy storing; incorporates a Kevlar<sup>®</sup>/Nylon keel for increased durability and longevity
- Dual spring design provides a smooth transition from heel to toe
- Moderate to high activity
- Weight rating up to 250 lbs\*
- Heel Height: Men: 3/4" (M) or 3/8" (L), Women: 3/8" (L)
- Available sizes: Men, 23-31 cm; Women, 22-27 cm
- Clearance: Men, 3.4"; Women, 3"
- Color: Buff (B), Medium Brown (M), or Tan (T)
- Heel Densities: Soft (S), Medium (M) and Firm (F)
- 12 month warranty, (three weeks allowed for installation)

\*Body weight plus any loads normally or routinely carried cannot exceed this weight limit



Description	Model No.
Mens	<b>CC2M</b>
Womens	<b>CC2W</b>
Endoskeletal Interface Plate	<b>3F-6CC-002</b>
Exoskeletal Interface Plate	<b>3F-6CC-003</b>
Titanium Foot Pyramid	<b>2LTC-1208</b>

Toe Resistance			
Foot Size	Very Low (2)	Low (3)	Regular (4)
Women's 22-27 cm	110 lbs and below	111-150 lbs	151-200 lbs
Men's 23-25 cm	120 lbs and below	121-160 lbs	161-230 lbs
Men's 26-30 cm	130 lbs and below	131-170 lbs	171-250 lbs

Ordering Example	
Model . . . . .	<b>CC2M</b>
Size . . . . .	<b>27</b>
Side . . . . .	<b>R</b>
Color . . . . .	<b>B</b>
Toe Resistance . . . . .	<b>4</b>
Heel Height . . . . .	<b>M</b>
Heel Density . . . . .	<b>M</b>
<b>CC2M-27RB-4M-M</b>	

### Carbon Copy® High Performance Foot

- Designed for highly active amputees
- Keel/Plate attachment allows energy storage at Mid Stance
- Moderate to high activity
- Weight rating up to 250 lbs\*
- Heel Height: Women: 3/8" (L) or 1 1/8" (H);  
Men: 3/4" (M) or 3/8" (L)
- Available sizes: Women, 22-27 cm; Men, 23-31 cm
- Clearance: 3.5"
- Cosmesis: Buff (B), Medium Brown (M) or Tan (T)
- 12 month warranty, (three weeks allowed for installation)

\*Body weight plus any loads normally or routinely carried cannot exceed this weight limit



Description	Model No.
Mens	<b>CHPM</b>
Womens	<b>CHPW</b>

Ordering Example	
Model . . . . .	<b>CHPW</b>
Size . . . . .	<b>22</b>
Side . . . . .	<b>L</b>
Color . . . . .	<b>B</b>
Toe Resistance . . . . .	<b>4</b>
Heel Height . . . . .	<b>L</b>
<b>CHPW-22LB-4L</b>	

Foot Size	Toe Resistance				
	Very Low (1)	Low (2)	Regular (3)	High (4)	Very High (5)
Women's 22-23 cm	70 lbs and below	71-115 lbs	116-160 lbs	161-195 lbs	195-250 lbs
Women's 24-25 cm	80 lbs and below	81-125 lbs	126-165 lbs	166-200 lbs	201-250 lbs
Women's 26-27 cm	90 lbs and below	91-130 lbs	131-170 lbs	171-205 lbs	206-250 lbs
Men's 23-24 cm	80 lbs and below	81-125 lbs	126-165 lbs	166-200 lbs	201-250 lbs
Men's 25-26 cm	90 lbs and below	91-130 lbs	131-170 lbs	171-205 lbs	206-250 lbs
Men's 27-28 cm	100 lbs and below	101-135 lbs	136-175 lbs	176-210 lbs	211-250 lbs
Men's 29-31 cm	110 lbs and below	111-145 lbs	146-180 lbs	181-215 lbs	216-250 lbs

### Foam Foot (Foot Blank)

- The Carbon Copy cosmesis is available as a solid foam foot with no Kevlar reinforcement
- Sizes available: 22-31
- Colors available: Buff (B) or Medium Brown (M)

Ordering Example	
Mens. . . . .	<b>FFTM</b>
Womens. . . . .	<b>FFTW</b>
Size . . . . .	<b>24</b>
Side . . . . .	<b>R</b>
Color. . . . .	<b>M</b>
<b>FFTM-24RM</b>	

**Earthwalk™ 2**

- Low activity; flexible keel foot/ankle system
- “Late stop” allows for easy roll-over so foot can maintain foot flat longer, this action is important for less active amputees who require more stability and less dynamic response
- Ankle provides smooth motion and more stability
- Weight rating up to 250 lbs\*
- Clearance: Women: 3.9” for 1-1/8” heel; 4.7” for 3/8” heel  
Men: 4.3” for 3/4” heel; 4.7” for 3/8” heel
- Heel Height: Women: 3/8” (L) or 1 1/8” (H); Men: 3/4” (M) or 3/8” (L)
- Available Sizes: Women: 22 cm - 27 cm; Men: 23 cm to 31 cm
- Cosmesis: Buff (B), Medium Brown (M) or Tan (T)
- 12 month warranty, (three weeks allowed for installation)



**\*Weight Limit and Guide for Selecting Toe Resistance** --The patient’s body weight plus any loads carried must not exceed 250 lbs (115 kg). This “modified body weight” is used to determine toe resistance. If the “modified body weight” exceeds 250 lbs (115 kg), the Magnum System is recommended. Use the patient’s foot size and “modified body weight” to determine the recommended toe resistance. If the patient is CMS level 3 or 4 increase the toe resistance to the next higher resistance.

Description	Model No.
Mens	<b>EWF2M</b>
Womens	<b>EWF2W</b>

Amputee Weight	Toe Resistance	Ankle Selection	Model No.
Less than 150 lbs	1	Red	<b>EWA2-R</b>
151-175 lbs	2	Blue	<b>EWA2-B</b>
176-220 lbs	3	Yellow	<b>EWA2-Y</b>
221-250 lbs	4	Green	<b>EWA2-G</b>

**Ordering Example**

Model . . . . . **EWF2M**  
 Size . . . . . **26**  
 Side . . . . . **L**  
 Color . . . . . **B**  
 Toe Resistance . . . . . **3**  
 Heel Height . . . . . **M**  
 Ankle Stiffness . . . . . **Y**

**EWF2M-26LB-3M-Y**

**SACH Foot**

- Low activity; lightweight and water-resistant
- Flexible toe belting with three toe resistances; lightweight composite keel
- 3 interchangeable heel inserts
- SACH with titanium pyramid features a permanently molded in pyramid to reduce weight and fabrication time
- Weight rating up to 250 lbs\*
- Heel Height: 3/4” (M) or 3/8” (L); Clearance: 3.4”
- Available sizes: 22-31 cm
- Cosmesis: Unisex; Buff (B), Medium Brown (M) or Tan (T)
- 12 month warranty, (three weeks allowed for installation)



*SACH Foot with Titanium Pyramid*

Description	Model No.
SACH Foot	<b>SFT</b>
SACH Foot with Titanium Pyramid	<b>SFP</b>

Foot Size	Toe Resistance		
	Low (2)	Regular (3)	High (4)
22-31 cm	130 lbs and below	131-180 lbs	181-250 lbs

**Ordering Example**

Model . . . . . **SFT**  
 Size . . . . . **22**  
 Side . . . . . **R**  
 Color . . . . . **B**  
 Toe Resistance . . . . . **3**  
 Heel Height . . . . . **M**

**SFT-22RB-3M**

### Single Axis Foot

- Composite keel and superior design makes this foot 10% to 20% lighter than conventional Single Axis Feet
- Interchangeable anterior and posterior deflection bumpers provide smooth transition from plantar flexion to dorsiflexion
- Weight rating up to 250 lbs adjusted body weight
- Low activity
- Heel Height: 3/4" (M) or 3/8" (L)
- Available sizes: 22-31 cm
- Clearance: Foot only, 3.4"; Assembly, 3.5"
- Cosmesis: Unisex Buff (B), Medium Brown (M) or Tan (T)
- 12 month warranty, (three weeks allowed for installation)
- **Note: Both Posterior and Anterior Bumpers, as well as an Endoskeletal Ankle, is needed to complete the Single Axis Foot**



Model No. SAF

Toe Resistance			
Foot Size	Low (2)	Regular (3)	High (4)
22-27 cm	130 lbs and below	131-180 lbs	181-250 lbs

Ordering Example	
Model .....	SAF
Size .....	26
Side .....	L
Color .....	B
Toe Resistance .....	3
Heel Height .....	M
<b>SAF-26LB-3M</b>	

### Anterior and Posterior Bumpers

- Bumpers are durable polyurethane
- Easily changed to optimize the patient's gait

Model No.	Anterior Bumper Selection Guidelines			Posterior Bumper Selection Guidelines		
	Soft (blue)	Medium (black)	Firm (white)	Soft (yellow)	Medium (green)	Firm (red)
	SAF-00128	SAF-00132	SAF-00134	SAF-00140	SAF-00142	SAF-00144
<b>Foot Size (cm)</b>						
22-23 cm	120 lbs & below 55 kg & below	120-250 lbs 55-115 kg	—	165 lbs & below 75 kg & below	166-250 lbs 76-115 kg	—
24-25 cm	—	250 lbs & below 115 kg & below	—	175 lbs & below 79 kg & below	176-250 lbs 80-115 kg	—
26-27 cm	—	185 lbs & below 84 kg & below	186-250 lbs 85-115 kg	—	250 lbs & below 115 kg & below	—
28-29 cm	—	195 lbs & below 89 kg & below	196-250 lbs 90-115 kg	—	195 lbs & below 89 kg & below	196-250 lbs 90-115 kg
30-31 cm	—	160 lbs & below 73 kg & below	161-250 lbs 73-115 kg	—	160 lbs & below 73 kg & below	161-250 lbs 73-115 kg

Note: Complete set of Anterior and Posterior deflection bumpers - Part No. **SAF-00130**

### Single Axis Ankle Endoskeletal

- Titanium pyramid molded in the ankle
- Pyramid reduces weight and eliminates the need for a foot bolt
- Ankle pin secures the ankle to the Single Axis Foot

Model No. SAF-00100



**Endolite Aqualimb**

- The Aqualimb is a unique, lightweight shower limb
- Water-resistant leg comes with a full length, hollow cosmesis
- Realistically shaped foot is designed for barefoot use and incorporates an anti-slip tread pattern on the sole to maximize grip on wet surfaces
- Structure made from long fiber injection molded components which will not corrode
- Compatible with all socket designs (provided that materials can withstand water contact)
- K1-K2 activity level
- Heel height: 3/16"
- Maximum clearance: 17"
- 24 month warranty



Foot Size	185 lb Weight Limit		220 lb Weight Limit	
	Left Foot	Right Foot	Left Foot	Right Foot
22 cm	629405L	629406L	629705L	629706L
23 cm	629407L	629408L	629707L	629708L
24 cm	629409L	629410L	629709L	629710L
25 cm	629411L	629412L	629711L	629712L
26 cm	629413L	629414L	629713L	629714L
27 cm	629415L	629416L	629715L	629716L

**KISS Dynamic Flexible Keel "Stubbie" Feet**

- Durable, Flexible keel
- Dynamically stores energy

Model No. **CMP27/A**



**Stomper and Stomper Jr.**

- Designed for use in and around water, Stomper and Stomper Jr. make an excellent choice for a prosthesis for shower, beach or pool
- May also be used for static fitting of test sockets
- Easily attaches to 30 mm and 34 mm pylon systems, or 2" PVC
- Stomper Jr. fits 30 mm and 34 mm pylons
- Stomper fits 2" outside diameter PVC pipe

Color	Stomper (2" PVC Pipe)	Stomper Jr (30 mm Pylon)	Stomper Jr (34 mm Pylon)
Black	STOMPBLK	STOMP30BLK	STOMP34BLK
Blue	STOMPBLU	STOMP30BLU	—
Red	STOMPRED	STOMP30RED	—

