LOWER EXTREMITY COMPONENTS

Feet

•	Foresee Orthopedic
•	Freedom Innovations
•	Kingsley
•	Ossur Americas
•	Trulife
•	WillowWood

Endoskeletal and Knees

664-665
665-667
668
668-669
669
670-672
672-674
675-676

Liners

•	ALPS	677-680
•	ESP	680-681
•	Euro International	682-683
•	Evolution	683
•	Medi USA	684-685
•	Ossur	
•	WillowWood	

Sleeves

•	ALPS
•	Comfort Products
•	ESP
•	Euro International
•	New Options Sports
•	P.O.R. Inc
•	Silipos
•	SPS
•	Syncor

Suspension

•	Knit-Rite	694-695
•	Syncor	696

Valves

•	ALPS
•	ESP
•	Evolution Industries
•	Fillauer Companies
•	Trulife

Locks

•	Coyote Designs	. 699-700
٠	Fillauer Companies	701
٠	KISS Suspension	702
•	Ossur Americas.	703
•	WillowWood	703

Socks and Sheaths

•	Comfort	 704

Shrinkers

•	Knit-Rite.																																.7	70)Ç)
---	------------	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	----	----	----	---

UPPER EXTREMITY COMPONENTS

Hosmer

٠	Hooks
٠	Terminal Device
•	Baby Mitt
•	Hands
٠	Wrist and Elbows715
•	Elbow Hinges and Harnesses716

Prosthetic 658

Foresee Orthopedic - S.A.F.E. II

- Stores and releases energy throughout its full length
- When the foot bends, the keel compresses, much like a spring
- As the body moves forward, the keel releases its energy, thrusting the leg forward
- Size Availability:

 - ¹/4" heel wide, 14-17 cm
 ³/8" heel narrow, 18-21 cm
- 12 month warranty



Ordering Example	
Style	S2
Туре	S
Keel	0
Heel Height	38
Heel Density	Μ
Width	Ν
Size	18
Side	L 👘
Color	1
S2SO38MN18L1	

		Orderin	ng Chart	for S.A.F.E	E. 11			
Style	Туре	Keel	Heel Height	Heel Density	Width	Size cm	Left/Right	Color
S2=SAFE II	S=Standard W=Waterproof (¹ 4" wide only)	O=Orange Standard-Moderate Activity R=Red Firm-High Activity Y=Yellow Flexible Low Activity	38= ³ /8" 14= ¹ /4"	S=Soft M=Medium F=Firm	W=Wide N=Narrow	13-24	L=Left R=Right	1=Caucasian 2=Light Negroid 3=Medium Negroid 4=Dark Negroid 5=Light Latin 6=Medium Latin 7=Dark Latin

*Some styles not available in all sizes



Freedom Innovations - Defender

- Smooth rollover with exceptional ground compliance in a foot that is lightweight, robust, and maintenance free
- Increased toe lever length provides the highest degree of energy storage-return available in pediatric foot
- Integrated pediatric pyramid attachment and patented technology eliminates bolt holes in the forefoot; ensuring optimal strength and performance
- User Weight Rating: 100 lbs
- 18 month warranty; 60 day trial period



Description	Model No.
Defender Foot	FS3-P
Foot Shell Sizes 16-18	FTC-1M-1
Foot Shell Sizes 19-21	FTC-2M-1

Please provide the following when ordering a foot shell:

- Model number
- Generation: 1 or 2
- Profile: Mid (M)
- Cap: 1
- Size: 16-21
- Skin Tone: Light (L), Medium (M) or Dark (D)
- Side: Left (L) or Right (R)
- Please refer to Foot Shell Ordering Example

	Build Height and Weigh	t				
Sizes Height Weight						
16-18 cm	2 ³ /8" (61 cm)	6 oz (175 g)				
19-21 cm	3" (76 cm)	8 oz (225 g)				

Foot Shell					
Ordering Exam	Ordering Example				
Model	FTC				
Generation	2				
Profile	Μ				
Сар	1				
Size	20				
Skin Tone	L				
	4-R				
Side	R				
FTC-2M-120L4-RR					

J	
Ð	
Q	
0)	
¥	
⊐.	
Ω	
0	

Defender Stiffness Category Chart							
Weight (lbs)	20-30	31-40	41-50	51-60	61-70	71-80	81-100
Category	1P	2P	3P	4P	5P	6P	1A



Kingsley - Pediatric Feet

Description	Heel Cushions	Heel Height	Sizes	Part No.	
Juvenile SACH Foot	Medium	10-13 cm Flat 14-17 cm ¹ /4" 18-22 cm ³ /8"	10-22 cm	K06	
Juvenile Smooth SACH Foot • Excellent for active children	Medium	14-17 cm ¹ /4" 18-21 cm ³ /8"	14-21 cm	KJ6	•
Juvenile Symes Aluminum plate keel 1" clearance 	Medium	10-13 cm Flat 14-17 cm ¹ /4" 18-22 cm ³ /8"	10-22 cm	K07J	

Ossur - Flex-Foot[®] Junior

- Low Profile energy storing foot with integrated male pyramid
- Pyramid shape contributes to progressive stiffening of the foot
- ³/8" heel height
- Rated to 100 lbs
- Build Height: 1¹³/16"
- Weight: 10.2 oz
- 8 month limited warranty

Description	Model No.
Male Pyramid	JRP0
Foot Shell	JRX50

Category Selection Chart						
Weight (kg)	Weight (lbs)	Category				
15-24 kg	33-53 lbs	1				
25-34 kg	54-75 lbs	2				
35-45 kg	76-99 lbs	3				



Flex Foot Ordering Example Model JRP0 Category 2 Size 19 Side R JRP0219R
Foot Cover Example ModelJRX50 Size19 SideR JRX5019R

Pediatrics

- Custom designed solution for children
- Vertical compliance and efficient energy return are must-haves for high active children
- Designed to be the most flexible of running feet
- Increased cushioning, lower frequency dynamic make this foot ideal for transfemoral and transtibial users
- Impact Level: High to Extreme
- Weight rating up to 325 lbs (130 kg)
- Weight of Foot: 11.9 oz (336 g)
- Build height: 8¹/4" (209 mm)
- One size available
- Adapter Options: Male or Female bolt on pyramid

Description	Model No.
Flex-Run Kit	FSX107
Components	
Attachment Plate	FWX40001
Male Pyramid Foot Adapter Set (75 g, M10 x 30 mm)	FWX00030
Female Pyramid Foot Adapter Set (100 g, M10 x 35 mm)	CM150161

Flex-Run Junior Category Selection Guide								
Patient Weight (Ibs)	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287
High to Extreme Impact Level	2	3	4	5	6	7	8	9

Trulife - Child's Play Foamer Foot

- Solid foam foot without a keel
- Designed to be incorporated into a prosthesis for patients with partial foot amputations
- Minimum Clearance: 3" (7.6 cm)
- Available in sizes 13-21 cm
 - Solid toe sizes 13-15 cm
 - Split toe sizes: 16-21 cm
- 1 year warranty

Foot Size	Heel Height
13-15	¹ /4" (6 mm)
16-18	³ /8" (9.5 mm)
19-21	¹ /2" (13 mm)

Description	Model No.
Caucasian	SFF132
Light Brown	SFF142
Dark Brown	SFF144



Foamer Foot Ordering Example
Model SFF132 Foot Size 13 Side L SFF132-13-L

Prosthetic 662

hetic **PEDIATRICS**

Child's Play Foot

- · Sculpted cosmesis reflect natural appearance of child's foot
- · Delrin keel offers shock absorption and a smooth roll over
- Appropriate for all lower limb pediatric amputees, including Symes
- Weight limit 65 kg/144 lbs
- Available in sizes 13-21 cm
- Solid Toe sizes 13-15; Split Toe sizes 16-21
- 24 month warranty on Keel and 12 month warranty on shell Keel Selection Instructions:
 - 1. In the charts below, find the appropriate weight,
 - activity level and limb length factors.
 - 2. Sum these three factors.
 - 3. With the factor sum, find the appropriate keel in the
- Keel Selection Table and check for availability by foot size. **Note:** Sizes 13-15 cm do not have split toes cosmesis.

Warning: Choosing a lower strength keel than what is suggested based on the above procedure and patient data will void the warranty.

	Та	ble I		Table II	
Weight	Factor	Weight	Factor	Activity Level	Factor
13-14 kg	0.8	31-34 kg	3.8	Low	-0.6
30-33 lbs		70-75 lbs		Medium-Low	-0.4
15-16 kg 34-37 lbs	1.0	35-37 kg 76-83 lbs	4.0	Medium	0
17-18 kg		38-40 kg		Medium-High	0.3
38-41 lbs	1.2	84-89 lbs	4.2	High	0.5
19-20 kg 42-45 lbs	1.8	41-44 kg 90-97 lbs	4.8	Table III	
21-22 kg	2.0	45-48 kg	E O	Limb Length	Factor
46-49 lbs	2.0	98-105 lbs	5.0	AK/Short BK	-0.4
23-24 kg 50-53 lbs	2.2	49-51 kg 106-113 lbs	5.2	Medium BK	0
				Long BK/Symes	0.3
25-26 kg 54-58 lbs	2.8	52-56 kg 114-123 lbs	5.8		
27-28 kg 59-63 lbs	3.0	57-60 kg 124-133 lbs	6.0		
29-30 kg 64-69 lbs	3.2	61-65 kg 134-144 lbs	6.2		

		Keel Selection Table
FactorSum	Keel	Foot Sizes (cm)
0.5 – 1.4	D1	13, 14, 15
1.5 – 2.4	D2	13, 14, 15, 16, 17, 18
2.5 - 3.4	D3	13, 14, 15, 16, 17, 18, 19, 20, 21
3.5 - 4.4	D4	16, 17, 18, 19, 20, 21
4.5 - 5.4	D5	19, 20, 21
5.5 - 6.4	D6	19, 20, 21

WillowWood - P-POD[™] SACH Foot

- · Made of composite, which makes it extremely durable
- Weight limit of 132 lbs (60 kg)
- Available in Buff (B), Medium Brown (M), or Tan (T)
- Heel Height: 9-17 cm are ¹/4"; 18-22 cm are ³/8"
- Heel Inserts: Soft, medium, and firm heel inserts are included with size 15-22 cm; not used with size 9-14 cm
- 12 month warranty

Model No.

PPS

PPI



SACH Ordering Example
Model
Size10
SideR
ColorB
Toe Stiffness2
PPS-10RB-2

P-POD[™] Impulse Feet

- Uses Pathfinder[™] composite material with a Kevlar[®] heel; material combination provides durability and energy return
- Weight limit of 132 lbs (60 kg)
- · Withstands the wear and tear children put it through
- Available in Buff (B), Medium Brown (M), or Tan (T)
- Heel Height: Sizes 13-17 cm are ¹/4", 18-22 cm are ³/8"
- Heel Inserts: Soft, medium, and firm heel inserts are included with size 15-22 cm; not used with size 13-14 cm
- 12 month warranty

Model No.

Foot Size	Low Toe (2)	Regular Toe (3)	High Toe (4)
13-14	<20 lbs	20-35 lbs	>35 lbs
15-16	<30 lbs	30-45 lbs	>45 lbs
17-18	<35 lbs	35-50 lbs	>50 lbs
19-20	<45 lbs	45-70 lbs	>70 lbs
21-22	<60 lbs	60-90 lbs	>90 lbs



Impulse Ordering Example
Model
Size
SideL
Color
Toe Stiffness
PPI-17LM-3

ESP - Haberman Alig	nment Devices	
 High grade, lightweight Titanium construction 		
Description and Part Weights	Product No.	
 Central Core Foundation for attachment of all other modules Build height: .650 in 97 g 	HAD-CC-P	
Locking Bar WedgeOne included with every HAD-CC-P10 g	HAD-W-P	
 Sliding 4-Hole Plate (threaded) Build height: .200 in 63 g 	HAD-4HPT-P	
Sliding 4-Hole Plate (unthreaded)Build height: .200 in60 g	HAD-4HPU-P	
Sliding PyramidBuild height: .400 in45 g	HAD-SP-P	
Sliding Pyramid Receiver • Build height: .600 in • 44 g	HAD-SPR-P	
Rotating Housing Plate Must be used with HAD-RSA-P Build height: .200 in 31 g 	HAD-RHP-P	
 Sliding Rotating Threaded Adapter Must be used with HAD-TP-P or HAD-TPR-P Build height: .600 in 43 g 	HAD-SRTA-P	
Threaded PyramidBuild height: .75 in39 g	HAD-TP-P	
Threaded Pyramid ReceiverBuild height: 1.062 in41 g	HAD-TPR-P	

Prosthetic 665

ESP - Haberman Alignme High grade, lightweight Titanium construction	ent Devices - Continue	a
Description and Part Weights	Product No.	
Rotating Slide Adapter Must be used with HAD-RHP Build height: .400 in 527 g 	HAD-RSA-P	
 Sliding Tube Clamp Adapter Build height: .125 in 50 g 	HAD-STCA-P	A de la de l
Threaded Female Dovetail Build height: .600 in 42 g 	HAD-TFD-P	

Fillauer - MightyMite [®] Components			
Description	Product No.		
 Sach Foot Adapter with Medium Ti Pyramid Includes foot bolt Rated to 132 lbs 	FC100110		
 ¹/4" Sach Foot Growth Compensation Spacer Includes lock washer Rated to 132 lbs 	FC100104		
 ¹/2" Sach Foot Growth Compensation Spacer Includes lock washer Rated to 132 lbs 	FC100102		
4-Hole Connector with Medium Ti PyramidRated to 132 lbs	FC100114		
Rotatable Medium Pyramid ReceiverRated to 132 lbsIncludes screw set	FC100134		
 Rotatable Small Pyramid Receiver Rated to 99 lbs Includes screw set 	FC100144		
Rotatable PyramidRated to 132 lbs	FC100135		

.

Fillauer - MightyMite® Components - Continued			
Description	Product No.		
 AK Socket Adapter with Rotation With lamination dummy Rated to 132 lbs 	FC100130		
 4-Hole Adapter with Rotation Rated to 132 lbs Includes screw set 	FC100132		
 4-Hole Distal Attachment Block Rated to 132 lbs Includes screws, vinyl caps and lamination plate 	FC100150		
 BK Socket Adapter with Medium Titanium Pyramid Includes lamination cap Rated to 132 lbs 	FC100152		
Alignment Adaptor	FC100160		
Tube Clamp - Aluminum • Rated to 132 lbs Tube Clamp - Titanium • Rated to 132 lbs	FC100122 FC100123		
Integrated 22 mm Pylon - Aluminum Medium Pyramid Receiver Rated to 99 lbs 	FC100118		
Integrated 22 mm Pylon - Alloy Steel Medium Titanium Pyramid Receiver Medium Aluminum Pyramid Receiver Rated to 132 lbs 	FC100119 FC100120		

	PEDIATRICS
al	

Fillauer MightyMite [®] Components - Continued			
Description	Product No.		
Integrated 22 mm Pylon - Alloy Steel Adult Pyramid Receiver Rated to 132 lbs 	FC100121		
Pylon Tube - Alloy Steel • Rated to 132 lbs Pylon Tube - Aluminum • Rated to 99 lbs	FC100142 FC100140		
 4-Bar Knee with pyramid attachment Rated to 132 lbs Standard 22 mm pylon Adjustable Extension Assist 	FC100510		
 4-Bar Knee without proximal attachment Rated to 132 lbs Standard 22 mm pylon Adjustable Extension Assist 	FC100500		
 MiniMac Knee with pyramid attachment Automatic locking mechanism Rated to 132 lbs Lightweight (11.1 oz) Adjustable stance flexion 	FC100520		
Proximal Pyramid AdapterRated to 132 lbs	FC100501		
Proximal Threaded AdapterRated to 132 lbs	FC100504		

à

Grace Prosthetics - Pediatric Grace Plate

- A 4-hole thermoplastic or laminated socket adaptor that fits almost all pediatric prosthetic components
- Comes complete with necessary screws for attachment

Product No.

1210-P




Ossur - Components					
Description	Weight Limit	Material	Weight	Product No.	
4-Prong Socket Adapter	100 lbs	Stainless Steel	79 g/2.8 oz	AJ-114040	
Female Pyramid Insert	100 lbs	Aluminum	36 g/1.3 oz	AJ-142300	
Male Pyramid Insert	100 lbs	Stainless Steel/ Aluminum	36 g/1.3 oz	AJ-133300	
4-Hole Socket Adapter	100 lbs	Aluminum	46 g/1.6 oz	AJ-122100	
4-Hole Male Pyramid Adapter	100 lbs	Stainless Steel/ Aluminum	39 g/1.4 oz	AJ-233100	
4-Prong Female Pyramid Adapter	100 lbs	Aluminum	34 g/1.2 oz	AJ-242100	
Female Pyramid Tube Clamp	100 lbs	Aluminum	48 g/1.7 oz	AJ-342100	

Ossur - Components - Continued					
Description	Weight Limit	Material	Weight	Product No.	
22 mm TubeLength: 250 mm	100 lbs	Aluminum	104 g/ 3.7 oz	AJ-712010	
Sach Foot Adapter	100 lbs	Stainless Steel/ Aluminum	54 g/ 1.9 oz	AJ-833100	8
 Total Knee JR Polycentric knee with geometric locking system Adjustable stance flexion; Adjustable extension promoter 160 degrees of knee flexion available 24 month warranty 	100 lbs		395 g/ 0.87 lbs	1100	
Total Knee JR 4-Hole Adapter	100 lbs	Aluminum		AJ-122300	
Total Knee JR Pyramid Adapter	100 lbs	Aluminum/ Titanium		AJ-831300	
Junior Shin Ferrule Protects knee unit and increases ease of foam cover application 				1145	
Junior Knee Cap (small/black)				1150	
Junior Knee Cap (large/beige)				1151	

Otto Bock - Components			
Description	Weight Limit	Product No.	
 Child's Rotatable AK Socket Adapter With (4) 506G3=M6 adjustment set screws Laminated into the socket Connects with components distally 	Up to 45 kg/99 lbs.	4R110	

PEDIATRICS Prosthetic 669



TiMed - Components				
Description	Weight Limit	Product No.		
Titanium Adapter, 25mm to 30mm • 32 g	300 lbs / 136 kg	A-30-25		
Titanium Adapter, 25mm to 30mm • 49 g	176 lbs / 80 kg	A-25-30		
Aluminum Tube, 22mm • 105 g • Color: Red • Length: 13" (330 mm)	99 lbs / 45 kg	A-400B		
 Carbon Fiber Tube, 22mm 68 g Length: 13" (330 mm) 	176 lbs / 80 kg	A-422C		
 4-Hole Connector with Pyramid 59 g 	176 lbs / 80 kg	P-192	a constant of the second secon	
 4-Hole Connector with Pyramid Receiver 63 g 	176 lbs / 80 kg	P-193		
Titanium Tube Clamp, 22mm Pyramid Receiver • 49 g	150 lbs / 68 kg	P-203A	10	
Aluminum Tube Clamp, 22mm • 49 g • Color: Orange	99 lbs / 45 kg	P-203B	Le Maria	
Small Double Ended Adapter, 32mm • 86 g	176 lbs / 80 kg	P-401	100	
Small Double Ended Adapter, 45mm (not pictured) • 80 g	176 lbs / 80 kg	P-402		
Small Double Ended Adapter, 60 mm • 84 g	176 lbs / 80 kg	P-403		

Prosthetic	
671	

TiMed - 22 mm Components - Continued				
Description	Weight Limit	Product No.		
Small Double Ended Adapter, 75 mm 86 g 	176 lbs / 80 kg	P-404		
Titanium Tube with Pyramid Receiver, 22 mm • 72 g	176 lbs / 80 kg	P-422		
Small 4-Hole Rotating Base • 39 g	176 lbs / 80 kg	P-502	Service	
Small Titanium Rotating Pyramid 26 g 	176 lbs / 80 kg	P-505		
Small Titanium Rotating Pyramid Receiver • 38 g	176 lbs / 80 kg	P-506		
Small Titanium SACH Foot Adapter with M8 Bolt • 50 g	176 lbs / 80 kg	P-509	Ì	
Small 4-Ear Titanium Connector • 34 g	176 lbs / 80 kg	P-510		
22 mm Small Flexible 3-Ear Titanium Connector • 45 g	176 lbs / 80 kg	P-511-3		
 25 mm Titanium Tube Clamp (not pictured) 48 g 	176 lbs / 80 kg	P-503		
Titanium Tube Clamp with standard 30 mm Pyramid Receiver, 22mm • 65 g	150 lbs / 68 kg	P-523		



TiMed - 22 mm Components - Continued			
Description	Weight Limit	Product No.	
 Aluminum Pediatric Single Axis Knee 183 g Color: Orange 	99 lbs / 45 kg	PK-100	
 Aluminum Pediatric 4-Bar Knee 176 g Color: Blue & Red 	120 lbs / 55 kg	PK-400	

Trulife - Child's Play [®] System				
Description	Weight Limit	Product No.		
 Child's Play Foot Pyramid Foot bolt not included (SFB544) Titanium 27 g 	120 lbs	SCA297		
 Child's Play SACH Foot Adapter Includes 8 mm foot bolt Aluminum 67 g 	110 lbs	P13035		
 Child's Play 22 mm Clamp Adapter Titanium 35 g 	120 lbs	SCA292		
 Child's Play 22 mm Pyramid Clamp Aluminum body with SS pyramid insert Color: Blue 40 g 	88 lbs	P13025		
 Child's Play 23 mm Variable Height Adapter Titanium 34 g 	120 lbs	SCA293		
 Child's Play 32 mm Variable Height Adapter Titanium 43 g 	120 lbs	SCA294		

Trulife - C	hild's Play® Sys	tem	
Description	Weight Limit	Product No.	
 Child's Play 45 mm Variable Height Adapter Titanium 56 g 	120 lbs	SCA295	Tr
 Child's Play 60 mm Variable Height Adapter Titanium 70 g 	120 lbs	SCA296	00
 Child's Play 22 mm 4-hole Pyramid Adapter For use with laminating block SLB240 Titanium 48 g 	120 lbs	SCA290	
 Child's Play 22 mm Female Pyramid Receiver 4-hole pattern may be mounted to laminating block SLB240 Aluminum 54 g Color: Blue 	88 lbs	P13003	
 Child's Play Knee 4-Bar polycentric design Friction swing control and extension spring assist Height: 4" (11 cm) 140° flexion For long residual limbs, remove pyramid and use Symes Nut lamination adapter SSY303 Color: Blue & Red 55 g 	121 lbs	SSK610A	
 Child's Play 22 mm 4-Ear Socket Pyramid Aluminum with SS pyramid insert 41 g Color: Orange 	110 lbs	P13043	
 Child's Play 22 mm x 12" Tube Aluminum 114 g Color: Blue 	88 lbs	P13010	
Child's Play Composite Flex 22 mm x 16" Pylon • 120 g • Color: Purple	120 lbs	SSL260-P3	Child's Play
 Child's Play Laminating Block Rigid polyurethane foam surrounding a nylon core Height: 1.3"; Diameter: 3.6" 93 g 	144 lbs	SLB240	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1

Pediatrics



Trulife	Trulife - Child's Play [®] System - Continued				
Description					
Child's Play Ankle Block	Product Wt.	Sizes	Weight Limit	Product No.	
 High strength rigid polyurethane foam surrounding a nylon core Height: 3.5" 	123 g 146 g 179 g	13-15 cm 16-18 cm 19-21 cm	69 lbs 89 lbs 144 lbs	SAB311 SAB312 SAB313	
	Height	Sizes	Weight Limit	Product No.	
 Child's Play Foot Spacer Permits access to alignment devices or use permanently to allow setup on existing limb systems 	6 mm (.25") 25 mm (1") 6 mm (.25") 25 mm (1") 6 mm (.25") 30 mm (1.2")	13-15 cm 13-15 cm 16-18 cm 16-18 cm 19-21 cm 19-21 cm	69 lbs 69 lbs 89 lbs 89 lbs 144 lbs 144 lbs	SFS240 SFS241 SFS242 SFS243 SFS244 SFS245	
Child's Play Laminating Core	Height	Sizes	Weight Limit	Product No.	
 Alternative to the ankle block Laminates directly to the socket providing adaptation to the foot 	63 mm 79 mm 75 mm	13-15 cm 16-18 cm 19-21 cm	144 lbs 144 lbs 144 lbs	SLC301 SLC302 SLC303	
	Height	Sizes	Weight Limit	Product No.	
 Child's Play Symes Nut Abrasive etched finish of the aluminum ensures secure socket adhesion during fabrication Minimum clearance 4 cm (1.6") 	13 mm (0.5")	13-15 cm 16-18 cm 19-21 cm	69 lbs 89 lbs 144 lbs	SSY301 SSY302 SSY303	0
	1		Color	Product No.	
 Child's Play Cosmesis Calf circumference 20 cm (10.25") Available for sizes 13 to 21 cm (sp Specify side 		ering)	Caucasian Light Brown Dark Brown	SEC220 SEC221 SEC223	

Prosthetic
675

Description	POD [™] Foundation Con Product No.	
Titanium Symes Nut • Rated to 132 lbs • 30 g	PPF-FS006	
 Laminated Socket Attachment Block Rated to 132 lbs 33 g 	PPF-700250	
 4-Hole Titanium Pyramid Rated to 132 lbs 30 g 	PPF-227014	
 4-Hole Aluminum Pyramid Receiver Rated to 132 lbs 30 g Color: Blue 	PPF-137004	
 51 mm Aluminum Dual Receiver Rated to 132 lbs 65 g Color: Blue 	PPF-138051	
 27 mm Aluminum Dual Receiver Rated to 132 lbs 44 g Color: Blue 	PPF-138027	
Aluminum Foot Pyramid Receiver • Rated to 132 lbs • 35 g • Color: Blue	PPF-1209	
Titanium Foot Pyramid AdapterRated to 132 lbs15 g	PPF-1208	
Aluminum 22 mm Tube Clamp Adapter • Rated to 132 lbs • 45 g	PPF-130000	

WillowWood P-	-POD [™] Founda	tion Component	ts - Continued
Description		Product No.	
 22 mm Aluminum Pylons 300 mm/11.8" (91 g) 400 mm/15.7" (160 g) 		PPF-130030 PPF-130040	
 22 mm Pylons with Fixed Aluminum 300 mm/11.8" (160 g) 400 mm/15.7" (199 g) 	n Adapter	PPF-131330 PPF-131340	
 1/4" 4-Hole Spacer Plates 10 g 		PPF-994114	~
1/8" 4-Hole Spacer Plates ● 5 g		PPF-994118	
 Growth Plate Kit Enables prosthesis to be lengther to use a new pylon or adapters Kit includes ¹/8", ¹/4", ¹/2" and 1 for 	-	PPF-GKIT	
¹ /8" Growth Plate		PPF-G125	
¹ /4" Growth Plate		PPF-G250	
¹ /2" Growth Plate		PPF-G500	
Exoskeletal Laminating Core Color: Blue 		PPF-1540	
Description	Right Product No.	Left Product No.	
P-POD [™] Exoskeletal Ankle Block Sach 15-18 cm (small)	PPS-1518SR	PPS-1518SL	9
P-POD [™] Exoskeletal Ankle Block Sach 19-22 cm (large)	PPS-1922LR	PPS-1922LL	
P-POD [™] Exoskeletal Ankle Block Impulse 15-18 cm (small)	PPI-1518SR	PPI-1518SL	9
P-POD [™] Exoskeletal Ankle Block Impulse 19-22 cm (large)	PPI-1922LR	PPI-1922LL	



ALPS - General Purpose (GP) Locking and Cushion Liners

- Distal end encapsulation retrofits most competitors' liners
- Slightly tacky feel reduces bunching in the popliteal region
- · Warranty: one liner 12 months, two liners 24 months

	Produ	ict No.
Size	3 mm Locking	6 mm Locking
16-19 cm	GPDT16-3	GPDT16-6
20-23 cm	GPDT20-3	GPDT20-6
24-25 cm	GPDT24-3	GPDT24-6
Please see the adult se	ection for larger sizes.	

	Produ	<u>ct No</u> .
Size	3 mm Cushion	6 mm Cushion
16-19 cm	GPFR16-3	GPFR16-6
20-23 cm	GPFR20-3	GPFR20-6
24-25 cm	GPFR24-3	GPFR24-6
Please see the adult section	for larger sizes.	



GP Locking Liner

GP Cushion Liner

Proximal Seal Interface (PSI) Liner

- Unlimited restriction of stretch on the proximal portion of the liner where there is no fabric
- Terry fabric acts as an air wick
- Works well with elevated vacuum systems currently available on the market
- · Excellent redistribution of forces will soothe the skin
- Measurement is taken 6 cm up from the distal end
- Warranty: One liner, 12 months; two liners, 24 months

	Produc	ct No.
Size	Tapered	6 mm
16-19 cm	PSI16-6/3	PSI16-6
20-23 cm	PSI20-6/3	PSI20-6
24-25 cm	PSI24-6/3	PSI24-6
Please see the adult sec	tion for larger sizes.	



EasyLiner[™] Fabric Reinforced Liner

PEDIATRICS


- Combines the softness, strength and elasticity of the ALPS EasyGel offering comfort and durability
- Provides a lower coefficient of friction relieving the skin of any shear force
- · Measurement is taken 6 cm up from the distal end
- 12 mm distal pad protects bony prominence
- · Contains antioxidants that help with the skin healing process
- · Warranty: one liner 12 months, two liners 24 months

	EasyLiner [™] Locking		
Size	3 mm Uniform	6 mm Uniform	6/3 mm Tapered
16-19 cm	ELDT-16-3HP	ELDT-16-6HP	ELDT-16-6/3HP
20-23 cm	ELDT-20-3HP	ELDT-20-6HP	ELDT-20-6/3HP
24-25 cm	ELDT-24-3HP	ELDT-24-6HP	ELDT-24-6/3HP
Please see th	ne adult section for lar	per sizes.	

Please see the adult section for larger sizes

	EasyLiner [™] Cushion		
Size	3 mm Uniform	6 mm Uniform	6/3 mm Tapered
16-19 cm	ELFR-16-3	ELFR-16-6	ELFR-16-6/3
20-23 cm	ELFR-20-3	ELFR-20-6	ELFR-20-6/3
24-25 cm	ELFR-24-3	ELFR-24-6	ELFR-24-6/3
Please see the	e adult section for larg	ger sizes.	





EasyLiner[™] Super Stretch

- Allows complete freedom of movement and at the same time achieves superior total contact with the skin without restricting blood flow
- Beneficial for scar and difficult-to-fit patients
- · Measurement is taken 6 cm up from distal end
- Warranty: 6 months

Size	6 mm Uniform Product No.
10-15 cm	ELPX-10
16-19 cm	ELPX-16
20-23 cm	ELPX-20
24-25 cm	ELPX-24
Please see the adult section for larger sizes.	



Clearpro[™]

- Fabricated from transparent silicone, allowing constant observation of the residual limb
- Medical grade silicone allows freedom of movement with the added comfort of as little as 1.5 mm of material behind the knee
- Designed with additional silicone along the distal end to reduce pistoning and increase distal stabilization
- · Measurement is taken 6 cm up from the distal end
- Warranty: 12 months


	Produ	<u>ict No.</u>	
Size	1.5 mm 3 mm		
14-15 cm	SSA14-1.5	SSA14-3.0	
16-17 cm	SSA16-1.5	SSA16-3.0	
18-19 cm	SSA18-1.5	SSA18-3.0	
20-21 cm	SSA20-1.5	SSA20-3.0	
22-23 cm	SSA22-1.5	SSA22-3.0	
24-25 cm	SSA24-1.5	SSA24-3.0	
Please see the adult section	n for larger sizes.		

Clear Gen[™]

- Fabricated from transparent silicone, allowing constant observation
 of the residual limb
- Available in two thicknesses and 6 sizes
- · Measurement is taken 6 cm up from the distal end
- Warranty: 12 months

	Produ	ict No.
Size	1.5 mm	3 mm
14-15 cm	GL-14	GL-14-T
16-17 cm	GL-16	GL-16-T
18-19 cm	GL-18	GL-18-T
20-21 cm	GL-20	GL-20-T
22-23 cm	GL-22	GL-22-T
24-25 cm	GL-24	GL-24-T
Please see the adult section for	larger sizes.	







Pediatrics

Thermoliner[™] Custom Liner

• Available in Locking or Cushion

PEDIATRICS

- Easily customize and shape a gel interface to the configuration of the user
- Apply the Thermoliner over the positive and place it in the oven for 30 minutes at 194°F or 90°C and customization is complete
- Can be pulled on like a sock, making donning easier
- · Gel durometer and contour provide additional control and support
- Softer version indicated BK's or bony AK's
- · Measurement is taken 6 cm up from the distal end
- · Warranty: one liner 12 months, two liners 24 months

	Product No.	
Size	3 mm Locking	6 mm Locking
16-19 cm	TFDT16-3	TFDT16-6
20-23 cm	TFDT20-3	TFDT20-6
24-25 cm	TFDT24-3	TFDT24-6
Places soo the adult s	action for larger sizes	



Please see the adult section for larger sizes.

	Produ	Product No.	
Size	3 mm Cushion	6 mm Cushion	
16-19 cm	TFFR16-3	TFFR16-6	
20-23 cm	TFFR20-3	TFFR20-6	
24-25 cm	TFFR24-3	TFFR24-6	
Please see the adult section for larger sizes.			

ESP - Aegis Ultimate Liner

- · Anatomically shaped for maximum comfort
- Pre-flexed design at 43°
- Newly formulated pure silicone gel; does not cold-flow; standard or T.S.B. mold modifications may be used
- Seamless fabric cover eases donning/doffing; reduces pistoning while stretching to accommodate wider limbs
- Uniform 6 mm thickness

Aegis Ultimate - Cushion Liner				
Limb Size	Liner Size	4"	6"	8"
14-20 cm	X-Small	BKUXS4C	BKUXS6C	BKUXS8C
20-26 cm	Small	BKUS4C	BKUS6C	BKUS8C
Please see the	e adult section fo	or larger sizes.		

Aegis Ultimate - Locking Liner				
Limb Size	Liner Size	4"	6"	8"
14-20 cm	X-Small	BKUXS4L	BKUXS6L	BKUXS8L
20-26 cm	Small	BKUS4L	BKUS6L	BKUS8L
Please see the adult section for larger sizes.				



Prosthetic 681

Aegis Streamline Liner - Clear

- This silicone liner features the same anatomical shape as the Aegis Ultimate liner but in a uniform 2 mm wall thickness
- Pre-flexed at 43°
- 10-durometer pigment free silicone to minimize adverse skin reactions; 2 mm thickness
- The locking liner comes with threaded end and adapter nut to accept all current pin options
- Available lengths: 4", 6", 8" or 10" (10" only available in size 20 cm)
- Replace ? with desired length

		Produ	ct No.
Limb Size	Liner Size	Cushion	Locking
17-19 cm	14 cm	BK14?G	BK14?S
19-21 cm	16 cm	BK16?G	BK16?S
21-23 cm	18 cm	BK18?G	BK18?S
23-25 cm	20 cm	BK20?G	BK20?S

Aegis Streamline Liner- Covered

- Anatomically Shaped for maximum comfort
- Pre-flexed design at 43°
- Standard or T.S.B. mold modifications may be used; does not cold-flow
- Compatible with any locking system (free adapter included)
- Seamless fabric cover eases donning/doffing; reduces pistoning while stretching to accommodate wider limbs
- Uniform 2-3 mm thickness

Limb Size/		Produ	ct No.
Liner Size	Length	Cushion Liner	Locking Liner
	4"	BK144G-C	BK144S-C
17-19 cm/	6"	BK146G-C	BK146S-C
14 cm	8"	BK148G-C	BK148S-C
	10"	—	—
	4"	BK164G-C	BK164S-C
19-21 cm/	6"	BK166G-C	BK166S-C
16 cm	8"	BK168G-C	BK168S-C
	10"	—	—
	4"	BK184G-C	BK184S-C
21-23 cm/	6"	BK186G-C	BK186S-C
18 cm	8"	BK188G-C	BK188S-C
	10"	_	—
	4"	BK204G-C	BK204S-C
23-25 cm/	6"	BK206G-C	BK206S-C
20 cm	8"	BK208G-C	BK208S-C
	10"	BK2010G-C	BK2010S-C
Please see the adult section for larger sizes.			



To determine proper sizing:

- 1. Take circumference measurement 2" (5-6 cm) proximal to distal end of residual limb.
- 2. Take length measurement from mid-patella tendon to distal end of residual limb.



To determine proper sizing:

- Take circumference measurement 2" (5-6 cm) proximal to distal end of residual limb.
- 2. Take length measurement from mid-patella tendon to distal end of residual limb.

- Medical grade silicone offers skin compatibility and comfort
- Durable bi-elastic single jersey fabric on the outside provides easy handling for your patient
- Offered in both tapered (6 mm distally to 2 mm proximally) and 6 mm uniform version; locking and cushion
- Alleviates pressure over bony prominences
- To select the correct size, measure the circumference of the residual limb approx. 4 cm from the distal end
- Sizes Available: 16, 18, 20, 22, 23.5
- · Replace ? with size; please specify when ordering



	Product No.		
Description	Cushion	Locking	
1st Class 3 mm Tapered	5-EU-1S7-?	5-EU-1S5-?	
1st Class + Uniform 6 mm	5-EU-1S8-?	5-EU-1S6-?	
Please see the adult section for I	arger sizes.		

Pediatrics

Clear-Liner, Matrix 10

- · Made of top quality medical grade silicone
- · The distal matrix reinforcement is essential for preventing pistoning
- To select the correct size, measure the circumference of the residual limb approximately 4 cm from the distal end
- Sizes Available: 12, 14, 16, 18, 20, 21, 22, 23.5
- · See the adult section for larger sizes
- · Replace ? with size; please specify when ordering

		Product No.		
Description	Thickness	Cushion	Locking	
Two-Color	2 mm	—	5-1S18-?	
Clearline	2.5 mm	5-1S2-?	5-1S1-?	





Prosthetic 683

ContexGel Liner

- Streifeneder polymer Gel SPG
- Uniform-styling
- 3 mm, 6 mm, and 9 mm available
- Easy rolling up through hard-wearing cover
- Thermoplastically deformable at 65°C
- To select the correct size, measure the circumference of the residual limb approx. 4 cm from the distal end

3 cm/ 3 m 3 cm/ 3 m 3 cm/ 6 m 3 cm/ 9 m	nm 5-2S10		S11-3S S11-6S
6 m 6 cm 8 cm/ 9 m)-6S 5-23	S11-6S
9 m	5 0044		
3 cm	1m 5-2S10)-9S 5-23	S11-9S
) cm/ 3 m 1 cm	nm 5-2S10)-3M 5-28	611-3M
) cm/ 6 m 1 cm	וויז 5-2S10)-6M 5-28	611-6M
) cm/ 1 cm 9 m	וויז 5-2S10)-9M 5-28	611-9M
	I cm 6 m	6 mm 5-2S10	6 mm 5-2S10-6M 5-2S





Evolution - Chameleon Liner

- Specifically designed for kids sizing and colors have been chosen for kids; electric blue, hot pink, neon green and fluorescent orange
- · Available in both cushion and locking
- Matrix is standard in Chameleon locking liners and as a special order option in the cushion liners
- · Silky surface allows the liner to slide on itself
- Unique medical grade silicone formula was specifically developed to ensure product quality and performance
- No cover to tear or run
- No coating to crack, peel or wear off

Circumference (inches)	Cushion Liner*	Locking Liner*
5 ¹ /2" - 6 ¹ /4"	ECC-1	ECL-1
6 ¹ /4" - 7"	ECC-2	ECL-2
7" - 6 ⁷ /8"	ECC-3	ECL-3
7 ⁷ /8" - 8 ³ /4"	ECC-4	ECL-4
*Specify color when ordering		

CC-1 ECL-1 CC-2 ECL-2 CC-3 ECL-3 CC-4 ECL-4

Prosthetic 684

hetic **PEDIATRICS**

Medi USA - Sensitive Liner

- Designed for the transtibial amputee with standard or conical shaped limbs
- · Soft silicone gel provides the ultimate in comfort and durability
- Available in 3 mm and 6 mm profiles
- Accommodates bony prominences as well as fleshy limbs
- Twelve sizing options
- Comfortable fit for the majority of transtibial amputees

	Product No.		
Size	3 mm	6 mm	
18 cm	M893118000	M893218000	
20 cm	M893120000	M893220000	
22 cm	M893122000	M893222000	
23.5 cm	M893123000	M893223000	
Please see the adult se	ction for larger sizes.		



Every Medi liner comes with a distal plug and a locking pin adapter

Protect Care Liner

- Protect Care silicone includes vitamin E, reduces adhesion, and allows easy donning and doffing
- Wall 2 mm thick, thicker distal wall, 10 cm matrix
- Fabric covered; can be used with and without the pin adapter
- · Indicated for residual limbs with limitations due to skin problems

Size	Product No.
18 cm	MC99118000
20 cm	MC99120000
22 cm	MC99122000
23.5 cm	MC99123000
Please see the adult section for lar	ger sizes.



Prosthetic 685

RELAX and Plus RELAX Liner

- The Medi Liner RELAX is distinguished by its uniquely formulated silicone gel and the special materials used in its cover
- With its Umbrellan textile cover, long term effectiveness against phantom pain, is now possible without negative side effects
- This cover shields amputees who already have a disturbed nerve structure from electromagnetic influences around the limb

	Product No.		
Size	3 mm Relax	6 mm Plus Relax	
18 cm	M892118000	M892218000	
20 cm	M892120000	M892220000	
22 cm	M892122000	M892222000	
23.5 cm	M892123000	M892223000	
25 cm	M892125000	M892225000	
Please see the adul	t section for larger sizes.		

RELAX Night Care Sock

- Umbrellan[®] for phantom pain, phantom sensations and idiopathic residual limb pain
- Clima-Comfort® material wicks moisture away from the skin
- Double-knitted fabric protects the residual limb from cold, slight impacts, friction, and other external influences
- 5 sizes for a wide range of transtibial and transfemoral limbs*

Size	TT-Lower Circ.	TT-Upper Circ.	Length	TF Groin Circ.	Product No.
Small	18-23 cm	20-26 cm	40 cm	N/A	1880322
Medium	23-30 cm	26-33 cm	40 cm	26-33 cm	1880323
Please see	the adult section	n for larger sizes			

Ossur - Iceross Junior Locking Liner

- Active skin care moisturizes and soothes skin
- · Stabilizing matrix and umbrella minimizes pistoning
- Nylon outer cover offers good radial strength combined with comfortable elasticity (strong adhesion between silicone and cover)
- Indicated for young and users with sensitive skin engaging medium impact daily activities
- Children can decorate their liners and personalize them, there are 4 different colored permanent markers included with the liner
- Matrix: 10 cm
- 6 Month Warranty
- Replace ?? with size 12, 14, 16 or 18

Description	Product No.
3 mm Stabilo Junior	I-8501??
3 mm Dermo Junior	I-8502??








. 2 2 2 2 2 2 2

WillowWood - Alpha[®] P-POD[™] Liners

- Surrounds the residual limb in a layer of skin-friendly, mineral oilbased gel that contains Vitamin E and is heat-moldable
- Alpha P-POD Liners are sized with young amputees in mind bringing comfort to young amputees for the first time
- The new MAX fabric is three times more abrasion-resistant than many other liners on the market and helps to increase the durability around the trim lines and other common wear areas
- This liner is made for those who are highly active and require a heavy-duty liner and have proven to be ideal for children
- 12 month warranty

		3 mm Uniform Product No.		
Size	Color	Cushion	Locking	
X-Small	Buff	ALC-5839-E	ALL-5839-E	
X-Small +	Buff	ALC-5832-E	ALL-5832-E	
X-Small	Med. Brown	ALC-5939-E	ALL-5939-E	
X-Small +	Med. Brown	ALC-5932-E	ALL-5932-E	

Distance from Distal End



Alpha[®] P-POD[™] Titanium Suction Pyramid

- Recommended for laminated sockets
- Features a duckbill valve and works well with an Alpha Suction Seal to create a secure fit
- PPF-SP471 includes 1 suction pyramid, 1 porous plastic insert and 1 foam muffler. When buying PPF-SP470, you receive a laminated socket attachment block in addition to the above
- Patient weight limit is 132 lbs

Description	Product No.
Titanium 4-Hole Suction Pyramid	PPF-SP471
Titanium 4-Hole Suction Pyramid Kit	PPF-SP470
Duckbill Valve Only	PPF-484







ALPS - Extreme Sleeve Fabric Reinforced

- Provides superior comfort due to ALPS GripGel[™], which seals to the user's skin without creating shear forces
- Provides secure suspension by sealing against the skin without restricting circulation
- Available in 6 mm and 3 mm uniform thickness
- 3 month warranty

		Product No.		
Size	Measurement (cm)	6 mm	3 mm	
XS	20-24	SFK16-6	SFK16-3	
S	25-30	SFK20-6	SFK20-3	
Please se	e the adult section for larger si	zes.		



V.I.V.A Valve Integrated Vacuum Activated Sleeve


- · Built in valve to evacuate the air in the socket
- Minimizes compression and irritation proximally
- The 6 mm is recommended for very active users
- Measurements taken 10 cm above patella
- Warranty: 3 months

Measurement		Product No.		
Size	(cm)	3 mm	6 mm	
X-Small	20-24	VIVA16-3	VIVA16-6	
Small	25-30	VIVA20-3	VIVA20-6	
Please see	the adult section for lar	ger sizes.		

EasySleeve[™] Fabric Reinforced

- Gel suction suspension sleeve with HP (high performance) fabric
- Elasticity of the fabric cover allows superior freedom of movement throughout knee flexion and while achieving a secure suspension by sealing against the skin without restricting circulation
- Also available in a sticky version
- · Measurements taken 10 cm above patella
- Warranty: 3 months

Measurement		Product No.		
Size	(cm)	3 mm Uniform	6 mm Uniform	
XXS	16-19	SFR-10-3 HP	SFR-10-6 HP	
XS	20-24	SFR-16-3 HP	SFR-16-6 HP	
Please se	e the adult section for la	raer sizes.		





EasySleeve[™] Super Stretch

- Gel suction suspension sleeve no fabric
- Allows complete freedom of movement and at the same time achieves superior secure suspension by sealing with the skin without restricting circulation
- · Separate fabric cover keeps dirt and clothes away from gel
- Measurements taken 10 cm above patella
- Warranty: 3 months

	Measurement	Product No.	
Size	(cm)	6 mm Uniform	Spandex Cover
XX-Small	10-15	SLE-10	—
X-Small	16-19	SLE-16	CSLE-16
Small	20-23	SLE-20	CSLE-20
Please see t	he adult section for la	arger sizes.	

Reinforced Clearline[™]

- Silicone suspension sleeve with fabric
- Elasticity of the fabric cover allows freedom of movement throughout knee flexion; secure suspension achieved by sealing against the skin without restricting circulation
- 1.5 mm Silicone with Lycra® Fabric
- Measurements taken 10 cm above patella
- Sizes: 28-37 cm; larger sizes available in adult section
- Custom sizes available
- Available in Beige, Brown and Black
- Warranty: 3 months

Size	Measurement (cm)	Beige	<u>Product No.</u> Brown	Black
Small	28-32	SLR-S	SLR-S-N	SLR-S-BLK
Medium	33-34	SLR-M	SLR-M-N	SLR-M-BLK
Large	35-37	SLR-L	SLR-L-N	SLR-L-BLK
Please see	e the adult section fo	or larger sizes	δ.	





Comfort Products - Endurance CRITERION

- The Criterion features a smooth skin lining that gives superior comfort against the skin and reduces the chance for popliteal bunching
- For patients without strenuous activity levels, the Criterion is an excellent choice for comfort and durability
- · Smoother feel; tougher and more durable neoprene material
- Rugged outer nylon cover; 15" length
- Manufactured using an intensive blind stitching and gluing process that greatly reduces seam failure

Circ. 4"		Product No.	
Size	above KC	Black	Beige
00	9-10"	SUC00	SUCB00
0	10 ¹ /2-12"	SUC0	SUCB0



- The Sport has all the features of the Criterion with the added benefit of a nylon lined center section to reduce abrasion
- This design allows the sleeve to give excellent suspension in even the most demanding of conditions
- · Smoother feel; tougher and more durable neoprene material
- Rugged outer nylon cover; 15" length
- Manufactured using an intensive blind stitching and gluing process that greatly reduces seam failure

Circ. 4"		Produ	Product No.	
Size	above KC	Black	Beige	
00	9-10"	SUS00	SUSB00	
0	10 ¹ /2-12"	SUS0	SUSB0	

Please see the adult section for larger sizes.

Endurance ZipSleeve

- Zip Sleeves are made of a higher grade of neoprene than most of the neoprene suspension sleeves available on the market
- Has a smoother feel and is a much tougher and more durable neoprene material than others
- Has a rugged outer nylon cover; 15" length
- Manufactured using an intensive blind stitching and gluing process
 that greatly reduces seam failure

	Circ. 4"	Product No.	
Size	above KC	Black	Beige
00	9-10"	SUZ00	SUSZ00
0	10 ¹ /2-12"	SUZ0	SUZB0
Please see	e adult section for larger	sizes.	







Prosthetic 690

hetic **PEDIATRICS**

ESP - FlexiKids™

- · Suction suspension sleeve
- · Maximizes comfort and freedom of movement
- Preflexed at 43°
- ¹/8" latex-free neoprene
- Measure circumference 3"-4" above the patella
- Colors available: Chocolate Brown, Tan, or Black



		Product No.		
Size	Measurement	Tan	Black	Chocolate Brown
X-Small	6" - 6 ¹ /2"	KS-ESP-XS-T	KS-ESP-XS-B	KS-ESP-XS-C
Small	6 ¹ /2" - 7 ¹ /2"	KS-ESP-S-T	KS-ESP-S-B	KS-ESP-S-C
Medium	7 ¹ /2" - 8 ¹ /2"	KS-ESP-M-T	KS-ESP-M-B	KS-ESP-M-C
Large	8 ¹ /2" - 9 ¹ /2"	KS-ESP-L-T	KS-ESP-L-B	KS-ESP-L-C
X-Large	9 ¹ /2" - 10 ¹ /2"	KS-ESP-XL-T	KS-ESP-XL-B	KS-ESP-XL-C
XX-Large	10 ¹ /2" - 11 ¹ /2"	KS-ESP-XXL-T	KS-ESP-XXL-B	KS-ESP-XXL-C

Euro International - Contex Gel Suspension Sleeve

- High quality gel offers maximum strength and stretch
- Excellent memory
- Contoured proximally to distal to mimic anatomical shape
- No bulky seams at front of knee to limit knee flexion
- Upper circumference: 23-35 cm; Lower circumference: 15-23 cm; 36 cm length


	Product No.	
Size	Beige	Gray
Small	3S10/S	3S10/S-G
Please see the adult section for larger sizes.		



New Options Sports - Hybrid[™] Suspension

- Exclusively designed adjustable suspension for all pediatric below the knee amputees
- Circumferential Velcro strap to help prevent migration
- Protective knee pad for comfort and added longevity
- Adjustable Velcro posterior closure; Open popliteal for comfort
- Shark Skin material inside grips prosthesis for superior hold

Size	Dimensions	Product No.
X-Small	6" - 6.6"	PK77-XS
Small	6.5" - 7.5"	PK77-S
Medium	7.5" - 8.5"	PK77-M
Large	8.5" - 9.5"	PK77-L
X-Large	9.5" - 10.5"	PK77-XL
XX-Large	10.5" - 11.5"	PK77-XXL
Circumference meas	surement 3" - 4" above knee	



P.O.R. Inc. - Latex Rubber

- Dip molded, latex rubber
- Provides effective suspension
- Covers the border of the prosthesis and prosthetic sock
- To choose the proper size, measure leg circumference 5" above the knee to establish thigh size; length 13¹/2"

Size	Thigh Size	Product No.
0	10 ¹ /2" - 12 ¹ /2"	por0
1	11 ¹ /2" - 13 ¹ /2"	1
Please see the a	dult section for larger sizes.	

Heavy Duty Latex Rubber

- Dip molded, latex rubber
- Provides effective suspension
- Covers the border of the prosthesis and prosthetic sock
- To choose the proper size, measure leg circumference 5" above the knee to establish thigh size; length 13¹/2"

Size	Thigh Size	Product No.
1	11 ¹ /2" - 13 ¹ /2"	1HD
1.5	13" - 14 ¹ /2"	1.5HD
Please see the a	adult section for larger sizes.	





Prosthetic 691

Silipos - DuraGel Suspension Sleeve

- Easy to donn and doff
- Mineral oil gel moisturizes and conditions the skin
- Made with DuraGel for added durability
- 3 mm gel thickness

Description	Product No
Suspension Sleeve	18181
Reinforced Sleeve	18191

DuraGel Sizing Guidelines				
	Proximal	Distal		
Relaxed	9" (23 cm)	8" (20 cm)		
Fully Stretched 15" (38 cm) 13" (33 cm)				
Please see the adult section for larger sizes.				

HD-BK Heavy Duty Suspension Sleeve[™]

- Helps to prevent rotation while allowing full range of motion
- Does not bunch-up behind the popliteal area, and keeps skin conditioned
- Circumference: 11"-17"
- Length: 12" (30.5 cm)
- Size: 1

Product No.



Arm Suspension Sleeve[™]

- Helps to maintain suspension
- May be used as BK-Sleeve with children and small adult sizes
- Circumference: 8"-13" (20-30 cm)
- Length: 8³/4" (22 cm)

Product No.

13005

1310



SPS - Neoprene

- Constructed of ¹/8" neoprene
- Nylon covering on one side
- Glued seam in front no air leakage
- Contour shaped to knee
- Lasts longer than latex sleeves
- Beige color

Regular Knee Sleeve				
Size	Measurement	Product No.		
0	10 ¹ /2" - 12"	3530075200		
1	12 ¹ /4" - 14"	3530075300		
1.5	13 ¹ /2" - 15"	3530075350		

Please see the adult section for larger sizes.

	Reinforced Knee Sleeve	
Size	Measurement	Product No.
0	10 ¹ /2" - 12"	3530179200
1	12 ¹ /4" - 14"	3530179300
1.5	13 ¹ /2" - 15"	3530179350
	duit a attac factor and a	

Please see the adult section for larger sizes.

	Reinforced Knee Sleeve with	Strap
Size	Measurement	Product No.
0	10 ¹ /2" - 12"	
1	12 ¹ /4" - 14"	3530207300
1.5	13 ¹ /2" - 15"	_
Please see the a	dult section for larger sizes.	



- · Sleeve is airtight
- No seam over patella
- · No seam over stress prone area
- 15° initial flexion produces a contoured area for patella and increases sitting comfort and suspension during full extension
- Lycra[®] Reinforcement Panel reduces wear at socket to sleeve interface

Size	Measurement 4" Above Knee	Product No.		
00	8" - 10"	DUR-A00		
0	10" - 12"	DUR-000		
1	12" - 14"	DUR-100		
Please see the adult section for larger sizes.				









Durasleeve EZ

- Airtight hypobaric suction
- · Friction control panel minimizes prosthetic wear
- · Seamless pre-flexed knee reduces pressure on the patella

Size	Measurement 4" Above Knee	Product No.
X-Small (0)	11" - 12"	DUR-EZ-0
Small (1)	12" - 14"	DUR-EZ-1
Please see the adult	section for larger sizes.	



Knit-Rite - Power Belt Auxiliary Suspension Belt[™]

- Highly porous 1 mm thin material to keep patients cooler and minimize bulk, yet remain dynamically strong and durable
- The benefits of design and material application can expand the range of trans-femoral users who may benefit from improved body integration and prosthesis control from auxiliary suspension, during daily use or for selected activities, while minimizing objections encountered with other belts
- · Adjustable Waist Left and Right with silicone gripper

			Produ	ict No.
Size	Waist	Fits Leg	Right	Left
X-Small	21 - 28"	12 - 15+"	3L928RXS	3L928LXS
Small	23 - 30"	15 - 20+"	3L928RSM	3L928LSM
Please see	the adult section	on for larger siz	es.	





Power Short Auxiliary Suspension Belt[™]

- Power Short Auxiliary Suspension Belt[™] option, with continuous waist configuration, eliminates the bulk of Velcro closures and is preferred by some patients, especially those who are active
- The wearable stretch/compression range of the Power Short Auxiliary Suspension Belt[™] is up to 8-10" in the waist and hip, depending on size and wearer preference
- Power Short is universal for left or right application

Continuous Waist Plain Leg with Silicone Gripper					
Size	Waist/Hip Unstretched	Waist/Hip Circ.	Fits Leg	Product No.	
X-Small	26/28"	28-36+"	12-15"	3L911CXS	
Small	28/30"	30-38+"	15-20"	3L911CSM	

Please see the adult section for larger sizes.

Continuous Waist Adjustable Leg with Silicone Gripper						
Size	Waist/Hip Unstretched	Waist/Hip Circ.	Fits Leg	Product No.		
X-Small	26/28"	28-36+"	14-21"	3P911AXS		
Small 28/30" 30-38+" 16-22" 3P911ASM						
Please see the adult section for larger sizes.						



Cool-Flex AK Suspension Belt

- Gun-Slinger "Hip-Hugger" design of this suspension belt improves fit and comfort
- The Cool-Flex material is breathable, lets moisture out, and lets air in to keep the wearer cool and comfortable
- Universal design for either right or left application
- The Sure-Grip silicone suspension band effectively secures the prosthesis in place
- · Reverse trimable Velcro attachment to waistband
- Proximal Brim measurement = 16-20" (40-50 cm)
- Size: Small (Please see the adult section for larger sizes)

Product No.

3LCFAKSM



Syncor - TES Belt (Total Elastic Suspension)

- Total Elastic Suspension; Made of neoprene with nylon on the inside; Velcro compatible fabric on the outside
- · Features a low profile, anti-roll, anti-bunch panel
- Cuff area has exposed textured neoprene on inside to securely grasp the socket

Unilateral Belts						
Socket Brim Circ. Waist Measurement					Produ	ict No.
Size	Inches	Metric	Inches	Metric	Left	Right
Infant	9"-12"	23-30 cm	18"-22"	45-56 cm	TES-INF-L	TES-INF-R
Juvenile	11"-14"	28-36 cm	21"-26"	53-66 cm	TES-JUV-L	TES-JUV-R
X-Small	13"-17"	33-43 cm	27"-32"	68-81 cm	TES-000-L	TES-000-R
Small	16"-20"	41-51 cm	33"-41"	84-104 cm	TES-100-L	TES-100-R
DISCONTRACTOR	41	- tions for a lower of a				

Please see the adult section for larger sizes.

Unilateral Belts with No Lateral Cutout						
	Socket Brim Circ. Waist Measurement Product No.					
Size	Inches	Metric	Inches	Metric	Left	Right
Infant	9"-12"	23-30 cm	18"-22"	45-56 cm	TES-INF-LN	TES-INF-RN
Juvenile	11"-14"	28-36 cm	21"-26"	53-66 cm	TES-JUV-LN	TES-JUV-RN
X-Small	13"-17"	33-43 cm	27"-32"	68-81 cm	TES-000-LN	TES-000-RN
Small	16"-20"	41-51 cm	33"-41"	84-104 cm	TES-100-LN	TES-100-RN
Please see	Please see the adult section for larger sizes.					

Bilateral Belts						
Socket Brim Circ. Waist Measurement Product No.						
Size	Inches	Metric	Inches	Metric	Regular	No Lateral Cutout
Infant	9"-12"	23-30 cm	18"-22"	45-56 cm	TES-INF-B	TES-INF-BN
Juvenile	11"-14"	28-36 cm	21"-26"	53-66 cm	TES-JUV-B	TES-JUV-BN
X-Small	13"-17"	33-43 cm	27"-32"	68-81 cm	TES-000-B	TES-000-BN
Small	16"-20"	41-51 cm	33"-41"	84-104 cm	TES-100-B	TES-100-BN
Ploaso soc	the adult so	ction for larger s	izoe			

Please see the adult section for larger sizes.

ALPS - Mini-Valve

- The Mini-Valve essentially is a one way expulsion valve releasing at approximately 1.5 psi, used to release any air trapped within the socket and the limb
- The mini-valve is very easy to use in both laminated and thermoplastic sockets
- Supplied with a rubber seal and a dummy plug for fabrication, and upon completion of fabrication the dummy plug is easily pushed out of the rubber seal and replaced with the valve



Product No.

LMV-009

ESP - Lyn RV Slide

- Suction socket valve with reliable auto expulsion with the added feature of maintained, hands free, open operation
- Slide the operator button in any direction and the valve stays open for easy donning and doffing

Product No.

PA0400





Lyn Valve BK

- 1. Drill hole 2. Tap hole 3. Screw in valve
- New or retro-installation outside of socket
- Constructed of lightweight, durable plastic
- Low profile, ⁷/8" diameter
- Releases air at 1/2 psi without adjustment
- Optional threaded housing used for installation into soft thermoplastics



DescriptionProduct No.Lyn Valve BKPA0200Lyn Valve BK with Manual ReleasePA0202BK Valve threaded housingPA0222Disassembly ToolPA-PWSI-T

Pee Wee Valve

- Use with thermoplastic, blister molded sockets
- Manually operated expulsion
- Reusable installation tool, "The Fork"
- Optional threaded housing used for installation into soft thermoplastics

Description	Product No.
Pee Wee Valve	PA0001
"The Fork" for installation	PA0009
Pee Wee SI	PA-PWSI
Pee Wee SI threaded housing	PA-PWSI-H
Disassembly Tool	PA-PWSI-T





Evolution Industries - Pediatric Valve Starter Kit

- · Is similar in content to the Adult Velocity Valve Starter Kit, except everything is sized specifically for Pediatrics
- If you have not purchased Pediatric Valves or Valve Kits from Evolution Industries previously, this should be included in your first purchase order
- The kit contains the Pediatric Valve Mounting Plate, Pediatric Expulsion Valve Housing, Pediatric Duckbill Valve, a Pediatric Valve Lamination Dummy, Pediatric Plate Lamination Dummy and Wrench

Product No.

PVV01-001

Pediatric Expulsion Valve Kit

· Is similar in content, design and function to the Adult Expulsion Valve Kit

Description	Product No.
Pediatric Expulsion Valve Kit	PVV01-002
Pediatric Valve Mounting Plate	PVMP-007
Pediatric Expulsion Valve Housing	PVV03-001
Pediatric Duckbill Valve	PVV03-003

Fillauer - Mini-Valve for Vacuum Sockets

- · Ideal for use with Silicone TEC or Alpha liners together with any type of air proof sleeve for below knee and knee-disarticulated amputees
- · One way air flow

Product No.


721000060

Trulife - Suction Socket Valve

Description	Product No.					
Small Suction Socket Valve-SFS Style (Push)						
Plastic Adjustable Leak Rate	P19000008					
Plastic Non-Adjustable Leak Rate	P19000009					
Housing Only	P20276					
Small Suction Socket Valve for Laminating (Push)						
Aluminum Adjustable Leak Rate	P19000005					
Aluminum Non-Adjustable Leak Rate	P19000006					
Housing Only	P20328					
Small Suction Socket Valve with 2-part House	sing (Push or Pull)					
Plastic Adjustable Leak Rate	P19000002					
Plastic Non-Adjustable Leak Rate	P19000003					
Housing Only	P20326					











699

Coyote Designs - Grommet Lock

- Ideal for pediatrics upper and lower extremity with easy engagement
- 1" cube in size, can be modified for areas with less space available
- · Water resistant, durable, and lightweight
- 100 lb weight limit

Description	Product No.		
Grommet lock with 10 mm threaded pin	CD104TM		
Grommet lock with 6 mm threaded pin	CD104SM		
Grommet Lock with ¹ /4 x 20" threaded pin	CD104QT		
Additional Pins			
Pin with 10 mm threaded pin	CD104P10		
Pin with 6 mm threaded pin	CD104P6		
Pin with ¹ /4 x 20" threaded pin	CD104PQ		



Small Air-Lock[™] with Pin

- Lock with suction valve built in for superior suspension
- Water resistant, lightweight, and very little wear
- One stage fabrication that separates pin and prosthetic alignment
- · Accommodates small limb shapes and upper extremity
- Includes 1.3" standard pin and lamination dummy

Description	Product No.
Air-Lock Lock with Pin	CD103S
1.4" Lock Pin Only	CD103P8
1.6" Lock Pin Only	CD103P9



Air-Lock Wrench

- · Easy to use
- · Makes lock maintenance fast and easy
- · Equivalent to 13 mm deep well socket

Product No.

CD103WH





Aligned 4-Hole Connector

- · Separates pin and prosthetic alignments
- · Improves lock function and eliminates need for offset connectors
- One stage fabrication in thermoformed or laminated sockets
- Pediatric Size

PEDIATRICS

Product No.

CD103PAF



Fast 4-Hole Connector

- · Creates quick and easy sockets with built-in 4-hole connector
- · Fast four can be used multiple times in thermoformed sockets
- Can be fabricated in thermoformed or laminated sockets
- Pediatric Size

Product No.

CD103PFF

Coyote Quick Adhesive

- · Easy to use
- 50 ml

Pediatrics

- · Sets up extremely fast with no sag, or smell
- · Replaces methyl methacrylate for socket repair
- Used for Air-Lock fabrication and reinforcing socks for use with distal pins
- Dispensing gun and mix tips required (see below)

Product No.

CD4150

50 ml Dispensing Gun

- · Easy to use
- · Drive design provides smooth operation and excellent control
- · Easily operated with one hand to start and stop flow

Product No.

CDDMA51

Adhesive Mix Tips

- Quad design mix tip reduces flow pressure and improves consistency of mixture
- Mix tip precisely combines the two component adhesive eliminating the need for hand mixing
- Allows you to dispense the exact amount of adhesive desired, reducing waste

Description	Product No.
6 Pack	CD401
20 Pack	CD402







Fillauer - MightyMite[®] Shuttle Locks • Small, rugged design for active lightweight amputees

- Tested to 132 lbs •
- Removable locking mechanism for easy service and replacement
- Thermoforming or lamination options
- Compatible with 22 mm endoskeletal systems

Description	Product No.
Cylindrical Housing Modular Shuttle Lock (Does not include Plunger)	FC100204
Titanium Pyramid Modular Shuttle Lock (Does not include Plunger)	FC100202
1" Plunger ¹ /4-20 thread	809826
1 ¹ /2" Plunger ¹ /4-20 thread	809827
1" Plunger metric thread	809826mm
1 ¹ /2" Plunger metric thread	809827mm
Fabrication Kit	125420





Small Modular Shuttle Lock with Delrin Housing

- · For pediatrics, upper extremities and select light duty transtibials
- Weight limit: 220 lbs

Description	Product No.
Small Modular Shuttle Lock w./ Delrin Housing (Includes 1" Plunger)	125410
Small Modular Shuttle Lock w./ Delrin Housing (Includes 1 ¹ /2" Plunger)	125412
1" Plunger metric thread	809826mm
1 ¹ /2" Plunger metric thread	809827mm
Fabrication Kit	125420





KISS - Pediatric Delrin KISS[®] Kit

- Features a delrin distal base which can be laminated into a socket or vacuum formed in a flexible inner socket
- · Includes an aluminum tab to assist in fabrication
- Designed for a 1 socket, 2 liner application
- Total build height: 3 mm

Pediatric 4-Hole Endo KISS[®] Kit

- Features an aluminum 4-hole distal base which also serves as an attachment base for endoskeletal componentry
- · Easily retrofits with PDI and Bulldog shuttle locks
- Designed for a 1 socket, 2 liner application
- Total build height: 3 mm



Pediatric Delrin KISS Kit

Description	Product No.
Delrin KISS Kit	PKS1/A
4-Hole Endo KISS Kit	PKS2/A
Kit Includes:	
Delrin Distal Base for PKS1/A	PCMP1/A
4-Hole Endo Distal Base for PKS2/A	PCMP2/A
Pediatric Adhesive Placard Set (1)	PCMP5/A
Pediatric Velcro Placard Set (1) (includes Ped. Proximal Nuts & Screws)	PCMP22/A
Pediatric Proximal Strap Set (1)	PCMP7/A
Pediatric Distal Strap Set (1)	PCMP6/A
Pediatric Distal Screw Set (1)	PCMP8/A
Pediatric KISS [®] Muffler Set (1)	PCMP9/A
10 mm - 6 mm Distal Threaded Insert (1)	PCMP15/A



Pediatric 4-Hole Endo KISS Kit



Ossur - Icelock[®] 700 Series, Small

- Ratchet lock without distal adapter designed to be used in conjunction with Upper Extremity and Pediatric liners
- Extremely low build height and lightweight
- Very strong and durable
- Removable from inside of socket
- Reusable fabrication tools
- 1 Year Warranty

Description	Weight Rating	Product No.		
Icelock Ratchet 721 Small		L-721000		
Lo	ock Adapters			
Stainless Steel Pyramid	99 lbs	L-773000		
Alloy 4-hole Adapter	99 lbs	L-774000		
A	ccessories			
Ratchet Attachment Pin Sma		L-792100		
Fabrication Tooling Kit		L-780000		



WillowWood - Mini G-Lock

- Designed to meet the needs of upper extremity and pediatric amputees
- Weight limit: 132 lbs
- Standard pin locking mechanism with pin seal; secures locking pin at nearest groove
- Use in Exoskeletal and Endoskeletal prostheses; only use in laminated sockets

Product No.

700-MGL480





Comfort Products - Regal Acrylic Prosthetic

Sock

- The Regal[™] Acrylic prosthetic sock is made of 98% acrylic and 2% Lycra
- It is a very soft stretchable sock which provides the patient with the latest advancement in fiber technology and knitting methods
- With the combination of fibers it gives a degree of stretch and conformance for the residual limb not possible with other fibers while wicking away moisture to ensure greater comfort
- Can be machine washed and dried in the dryer on delicate cycles
- To order, replace ? with desired ply
- · Add 3 to product number for sock with hole in distal end
- Available by special order in 2 ply, 4 ply and 7 ply. Special order sizes from XX-Short to XX-Long and XX-Wide available, as well as a Wide/Mid Length size



Regal Sock (3 Ply, 5 Ply, and 6 Ply Weight)									
		<u>Length</u>		<u>Top To</u>		<u>)6</u>			
Size	Fits Current Size	In.	cm	In.	cm	In.	cm	Product No.	
Narrow X Short	00, 0, A0, A1	8-10	20-25	4-6	10-15	3-3 ¹ /2	7.6-9	PSRNAXS?	
Narrow Short	00, 0, A0, A1	10-12	25-30	4-6	10-15	3-3 ¹ /2	7.6-9	PSRNASH?	
Narrow Regular	00, 0, A0, A1	14-16	35-40	4-6	10-15	3-3 ¹ /2	7.6-9	PSRNARG?	
Narrow Long	00, 0, A0, A1	18-20	45-50	4-6	10-15	3-3 ¹ /2	7.6-9	PSRNALG?	
Narrow X Long	00, 0, A0, A1	22-24	55-60	4-6	10-15	3-3 ¹ /2	7.6-9	PSRNAXL?	

Please see the adult section for larger sizes

Knit-Rite - Sterling Stretch Wash & Wear Prosthetic Sock

- The sterling Stretch sock adds Lycra Spandex to the acrylic to allow the sock to conform to the shape of the limb better
- · Fleecing adds softness to the inside next to the skin
- Available in 3, 5 and 6 ply weight
- · Wash and wear stretch
- · Conforms to shape of limb
- Transports moisture from skin
- · Reinforced toe
- To order, replace ? with desired ply



Sterling Stretch Sock									
		Ler	<u>igth</u>	1	Гор	To	<u>e</u>		
Size	Fits Current Size	In.	cm	In.	cm	In.	cm	Product No.	
Narrow X Short	00, 0, AS, AL	8-10	20-25	4-6	10-15	3-3 ¹ /2	7.6-9	1SS?NAXS	
Narrow Short	00, 0, AS, AL	10-12	25-30	4-6	10-15	3-3 ¹ /2	7.6-9	1SS?NASH	
Narrow Regular	00, 0, AS, AL	14-16	35-40	4-6	10-15	3-3 ¹ /2	7.6-9	1SS?NARG	
Narrow Long	00, 0, AS, AL	18-20	45-50	4-6	10-15	3-3 ¹ /2	7.6-9	1SS?NALG	

Knit-Rite - Soft Sock Cuddly N'Soft

- Full fashion knit sock made of CoolMax[®]/Lycra spandex in easy-to-fit sizes
- The fiber wicks moisture away from skin
- The skin remains drier and more comfortable
- · Soft-Sock feels soft and cuddly worn next to the skin



Soft Sock Regular									
	Fits	Ler	Length Product No		ct No.	Produ	ict No.		
Size	Current Size	In.	cm	1 Ply	3 Ply	5 Ply	6 Ply		
Child X Short	A	6	15	1SP1CHXS					
Child Short	A	8	20	1SP1CHSH	1SP3CHSH	1SP5CHSH	1SP6CHSH		
Child Medium	A	12	30	1SP1CH12	1SP3CH12	1SP5CH12	—		
Narrow Short	B, 0, 1	10-14	25-35	1SP1NASH	1SP3NASH	1SP5NASH	1SP6NASH		
Narrow Medium	B, 0, 1	16-20	40-50	1SP1NAMD	1SP3NAMD	1SP5NAMD	1SP6NAMD		
Narrow Long	B, 0, 1	22-28	55-70	1SP1NALG	1SP3NALG	1SP5NALG	1SP6NALG		
Please see adult se	ction for larger size	S							

Soft Sock Hole-In-Toe						
Fits	Len	<u>igth</u>	Produ	ct No.	Produ	ict No.
Current Size	In.	cm	1 Ply	3 Ply	5 Ply	6 Ply
A	6	15	1SPSCHXS	_	_	_
А	8	20	1SPSCHSH	1SP3SCSH	1SP5SCSH	1SP6SCSH
B, 0, 1	10-14	25-35	1SPSNASH	1SP3SNSH	1SP5SNSH	1SP6SNSH
B, 0, 1	16-20	40-50	1SPSNAMD	1SP3SNMD	1SP5SNMD	1SP6SNMD
-	A A B, 0, 1	Fits Ler Current Size In. A 6 A 8 B, 0, 1 10-14 B, 0, 1 16-20	Fits Current Size Length In. cm A 6 15 A 8 20 B, 0, 1 10-14 25-35 B, 0, 1 16-20 40-50	Fits Length Production Ourrent Size In. cm 1 Ply A 6 15 1SPSCHXS A 8 20 1SPSCHSH B, 0, 1 10-14 25-35 1SPSNASH B, 0, 1 16-20 40-50 1SPSNAMD	Fits Current Size Length Product No. In. cm 1 Ply 3 Ply A 6 15 1SPSCHXS — A 6 20 1SPSCHSH 1SP3SCSH B, 0, 1 10-14 25-35 1SPSNASH 1SP3SNSH B, 0, 1 16-20 40-50 1SPSNAMD 1SP3SNMD	Fits Length Product No. Product No. In. cm 1 Ply 3 Ply 5 Ply A 6 15 1SPSCHXS — — A 6 15 1SPSCHSH 1SP3SCSH 1SP5SCSH B, 0, 1 10-14 25-35 1SPSNASH 1SP3SNSH 1SP5SNSH B, 0, 1 16-20 40-50 1SPSNAMD 1SP3SNMD 1SP5SNMD

Please see adult section for larger sizes

Spacer Socks

- This 1-ply sock has a reciprocal toe that has the seam one inch up from the distal end (like a sheath)
- · Toe is rounded and reinforced for good fit and long wear
- Blended cotton polyester fabric is soft, comfortable and absorbent
- · Prewashed so size remains true wash after wash
- Machine wash, tumble dry
- 6 per package

Weight	Fits Sizes	Length	Product No.
Narrow lightweight	#B, #0	15"	6S-8035-LT
Narrow heavyweight	#B, #0	15"	6S-8035-HV
Please see adult section f	or larger sizes		



22222222

Stretch Spacer Socks

- Made of cotton/Lycra spandex/nylon
- Like the original Spacer Socks, they feature a very soft and smooth cotton finish with a durable sheath type nylon toe
- Elasticity allows over 50% greater range than non-elastic version
- 6 per package

	Fits Sock		Product No.	
Size	Sizes	Length	Regular	Hole-In-One
Narrow	A,B	12"	6S-8S35-MD	6S-SS35-MD
Medium	0, 1	15"	6S-8S40-MD	6S-SS40-MD
Please see ac	lult section for la	rger sizes		

X-Wool[™] Prosthetic Sock

- A blend of wool. lycra, and featuring X-Static[™] the Silver Fiber, a performance fiber system made with pure silver
- Silver is naturally antimicrobial, enabling X-Static to safely inhibit the growth of odor causing bacteria
- Silver also has the best possible thermodynamic properties, so the patient stays cooler in the summer and warmer in the winter
- Absorbent virgin wool is thick and cushioning. Lycra allows garment to contour without stretching out and accommodates flexion
- Fleeced to add softness
- Machine washable, resists shrinkage





	Fits Sock			Product No.		
Size	Sizes	Length	3 Ply	5 Ply		
Narrow X-Short	00, 0, AS, AL	8-10"	1WX3NAXS	1WX5NAXS		
Narrow Short	00, 0, AS, AL	10-12"	1WX3NASH	1WX5NASH		
Narrow Regular	00, 0, AS, AL	14-16"	1WX3NARG	1WX5NARG		
Please see adult section for larger sizes						

Power Tube Prosthetic Pull Sock

- A lower-friction donning sock
- Open on both end

Size	Max. Circumference	Product No.
X-Small	16"	3M-PT08-XS
Small	20"	3M-PT10-SM





Prosthetic 707

Power Pull Donning Sock

- · Made of slicker, more supple fabric for easier donning
- Single-piece integral design provides a more direct and effective pull transfer and eliminates stress points for a more durable product

Size	Тор	Proximal Circum.	Product No.
Size A	4"	less than 8"	3MPP04AS
Size B	6"	8-12"	3MPP06BS
X-Small	8"	12-16"	3MPP08XS
Small	10"	16-20"	3MPP10SM



New Style Sheaths

- 100% Stretch Nylon permits a given sheath size to cover a range of prosthetic sock sizes
- Size range accommodates prosthetic sock sizes A (arm) through AK sizes 4, 5, or wider
- 3 pair per box

		Product No.		
Description	Size	White	Dusk	
Child Short	А	2HPWCHSH	2HPBCHSH	
Child 12"	А	2HPWCH12	2HPBCH12	
Narrow Short	B, 0, 1	2HPWNASH	2HPBNASH	
Narrow Medium	B, 0, 1	2HPWNAMD	2HPBNAMD	
Narrow Long	B, 0, 1	2HPWNALG	2HPBNALG	
Please see adult see	tion for larger	sizes		
	Hole-ii	n-Toe Sheath		
Narrow Short	B, 0	2HPWNHSH	2HPBNHSH	
Regular Short	1, 2	2HPWRHSH	2HPBRHSH	





Prosthetic 708

Silipos - Silosheath Original

• 2 ply, 2 mm gel

PEDIATRICS

- Minimizes skin breakdown
- · Reduces friction, shearing and abrasion on residual limb
- Cushions and protects
- · Softens, soothes and improves the appearance of skin
- Product is washable and reusable

Double Cushion Silosheath[™]

- 2 ply, 3 mm gel, extra gel at distal end
- Minimizes skin breakdown
- · Reduces friction, shearing and abrasion on residual limb
- Cushions and protects
- · Softens, soothes and improves the appearance of skin
- · Product is washable and reusable

Original

Silosheath Extra Life[™]

- 3 ply, 2 mm gel, increased durability
- Minimizes skin breakdown
- · Reduces friction, shearing and abrasion on residual limb
- Cushions and protects
- · Softens, soothes and improves the appearance of skin
- · Product is washable and reusable

	Gel	Approx.
)	Length	Length
	8" (20 cm)	10"-13" (25-33 cr

Please see adult section for larger sizes

STS - AFO Casting Socks

- Designed for use in obtaining lower extremity molds of children
- Available in 4 sizes to cover a broader range of anatomy and ages

Approx.

Distal Circ.

- · Impressions are as accurate as plaster of Paris
- · Easily removed with small scissors; Cast cutter not needed
- · Glass-smooth finish is ideal for CAD/CAM
- Can be shipped with little chance of damage to central fab; ideal for custom pediatric AFOs

Size	Foot Length	Heel to Sock Top	Product No.
Small	3 ¹ /2"	4"	904-S
Medium	4"	5"	904-M
Large	4 ¹ /2"	6 ¹ /2"	904-L
X-Large	5 ¹ /2"	71/2"	904-XL

22222222

Approx.

Proximal Circ.

 $6^{1}/2$ "-9" (17-23 cm) $6^{1}/2$ "-10" (17-25 cm)



Double

Cushion

13415

Original

92815

Double

Extra

Life

13815

Knit-Rite - Compressogrip Elastic Shrinkers

- Tubular elastic compression material for control or prevention of edema following amputation surgery
- Speeds fabrication, fitting and delivery of a prosthesis by uncomplicating application of an elastic bandage
- Approximate single layer pressure at 50% stretch is 20-25 mmHg
- At 100% stretch, compression is doubled. Doubled layers will double compression at both 50% and 100% stretch
- Applications of more than two layers should be monitored as compression may exceed the limits of desired therapeutic pressure
- 92% polyester, 8% spandex; Latex free

			Product No.		
Size	Length	Circ.	White	Beige	
2	18"	5"-7"	2H0S0218	2H0B0218	
2	24"	5"-7"	2H0S0224	2H0B0224	
2 ¹ /2	18"	7"-10"	2H0S0S18	2H0B0S18	
2 ¹ /2	24"	7"-10"	2H0S0S24	2H0B0S24	
Please see adult section for larger sizes					

APPLY THE FIRST LAYER



POSITION THE RING SNUGLY



APPLY THE SECOND LAYER TO A POINT OF DESIRED MAXIMUM COMPRESSION.

Compressogrip Tubular Elastic Bandage

- · Proven effective for: Edema reduction and control
- Interim control of hypertrophic scarring following burn injury
- Dressing and monitoring device retention
- · Smoothing irregularities when casting
- Wide range of sizes; cost effective
- All sizes provided in 11 yard (10 meter) rolls
- 92% polyester, 8% spandex
- Machine wash and dry, warm, no bleach
- Fabric will relax 20% in length after first washing

		Produc	ct No.
Size	Circ.	White	Beige
2"	5"-7"	6S762802	6S763102
2 ¹ /2"	7"-10"	6S7628S2	6S7631S2
2 ³ /4"	8"-11"	6S7628T2	6S7631T2
Please see	adult section for larg	er sizes	



222222222

Hosmer - Child Size Hook- Model 12P

- · Stainless steel, canted fingers, plastisol coated
- 2³/4" (7 cm) long; 3 oz (85 g)
- · Available in brown and Caucasian colors
- Standard ¹/2" 20 thread



	Product No.		
Description	Left	Right	
12P Caucasian	55056	55057	
12PN Brown	54281	54282	

Child Size - Model 10X

- · Aluminum, canted fingers, nitrile rubber lined
- 3³/8" (8.6 cm) long; 3 oz (85 g)
- Standard ¹/2" 20 thread

	30	De paramente
F		The second secon
3	80	
	U	

	Product No.				
Description	Left Right				
10X Hook	55053	55054			

Child Size - Model 10P

- Aluminum, canted fingers, plastisol coated
- 2¹/4" (8.6 cm) long; 3.5 oz (99 g)
- Available in brown and caucasian colors
- Standard ¹/2" 20 thread

	Product No.				
Description	Left Right				
10P Caucasian	55050	55051			
10PN Brown	54283	54284			



Small Adult Size - Model 99X

- Aluminum, canted fingers, nitrile rubber lined
- 3⁷/8" (9.8 cm) long; 3 oz (85 g)
- Standard ¹/2" 20 thread

	Product No.			
Description	Left	Right		
99X Hook	55043	55044		



Small Adult Size - Model 99P

- · Aluminum, canted fingers, plastisol coated
- 3⁷/8" (9.8 cm) long; 3.5 oz (99 g)
- Available in brown and caucasian colors
- Standard ¹/2" 20 thread

	Product No.		
Description	Left	Right	
99P Caucasian	52946	53604	
99PN Brown	56871	56872	

UCLA CAPP Terminal Device

- · Child size, sturdy alternative terminal device
- · Wider palmar face and resilient friction cover for secure prehension
- · Fits right or left with center or outside pull; available in Caucasian and Dark brown
- · Soft, regular or heavy duty prehension springs are easily interchanged with an installation tool
- ¹/2" 20 thread stud; length: 2¹⁵/16" (7.47 cm); 3 oz (85 g)
- Base diameter $1^{1}/2^{"}$ (3.8 mm); thumb opening $2^{5}/8^{"}$ (6.7 cm)

Description	Product No.			
Caucasian				
Regular Spring	71620			
Soft Spring	71630			
Heavy Duty Spring	71633			
Dark Brown				
Regular Spring	59321			
Soft Spring	59322			
Heavy Duty Spring	71634			
Replacement Spring Kit	71632			

UCLA CAPP Delrin Body Friction Wrist

• These lightweight wrists incorporate a friction adjustment, which completely encircles the threaded interface of the terminal device and wrist, eliminating frequent changes in the friction setting

Description	Weight	Product No.
2" Diameter	1.5 oz	55110
1 ³ /4" Diameter	1.0 oz	55109
1 ¹ /2" Diameter	.75 oz	55108











Prosthetic 711



Baby Mitt

- Soft foam passive hands
- Passive hand which allows the child to begin bi-manular activity, achieves body symmetry, helps to protect child during a fall, aids in standing or sitting balance, incorporation of the prosthesis in the body image
- Hands are available in two sizes for infants and toddlers, and two colors
- Does not require the addition of a cosmetic glove
- Standard 1/2"-20 thread stud



		ļ	MCP	Prod	uct No.
Hand Size	Length	Circ.	with Thumb	Left	Right
		Cauca	sian		
X-Small	3" (7.6 cm)	3 ³ /4" (9.5 cm)	5 ¹ /8" (13 cm)	50484	50483
Small	3 ¹ /4" (8.3 cm)	4 ¹ /2" (11.4 cm)	5 ³ /4" (14.6 cm)	50486	50485
		Dark B	rown		
X-Small	3" (7.6 cm)	3 ³ /4" (9.5 cm)	5 ¹ /8" (13 cm)	59121	59120
Small	3 ¹ /4" (8.3 cm)	4 ¹ /2" (11.4 cm)	5 ³ /4" (14.6 cm)	59123	59122

Child CAPP Hand

- Soft, durable urethane cosmesis with life-like appearance
- More functional than most hands
- Adjustable spring tension
- · Improved parent acceptance of cosmesis
- Does not require cosmetic glove
- Smooth, voluntary opening performance
- Standard ¹/2"-20 thread stud

Color	Wrist to Finger Tip	Weight	Cable Exit	<u>Produ</u> Left	<u>ict No.</u> Right
Caucasian	3.3"	5.67 oz/ 161 g	Dorsal	61826	61827
Medium Brown	3.3"	5.67 oz/ 161 g	Dorsal	61828	61829
Dark Brown	3.3"	5.67 oz/ 161 g	Dorsal	61830	61831




Child Dorrance 200 Mechanical Hand

- Dorrance hands provide voluntary opening of the thumb and first two fingers simultaneously
- Adjustable pinch force
- Breakaway release protects wearer if excessive force applied
- MCP circumference: 6¹/2" (16.51 cm)
- Wrist to thumb tip: 3⁷/8" (9.8 cm)
- Weight: 9 oz (255 g)


Description	Side	Glove No. (yyyyy)	Product No.
Dorrance Hand	Right	10204	50538
Dorrance Hand	Left	10203	50539
Replacement Spring, regular			55429
Replacement Spring, light		_	55430

Gloves				
Description	Product No.			
Micro Coated Vinyl	1yyyyyzz3			
PVC	1yyyyyzz1			
Expanded MCV	1yyyyyzz4			
Hard PVC (light prosthesis)	2yyyyyzz1			
Micro Coated Silicone	3yyyyyzz3			
Silicone	3yyyyzz1			

Please reference Hosmer swatch kit (56760) for color options (zz)								
01	02	03	04	05	06	07	08	09
10	11	12	13	14	15	16	17	18



Glove Ordering Example Replace yyyyy with glove no. Replace zz with color choice				
Product No 1yyyyyzz3 Glove No 10204 Color 01 110204013				



A. Wrist circumference

- B. Metacarpal circumference
- C. Middle finger length
- D. Glove end circumference
- E. Length of glove



Prosthetic

714

Child Passive Hand

PEDIATRICS

- Soft foam passive hands available for either ¹/2"-20 stud or M12 stud thread
- Passive hands positionable flexible fingers
- · Reinforced hands have wired fingers attached to the wrist area
- Super reinforced hands have steel bars in the fingers that attach to the wrist area
- Hard PVC passive hands are a lightweight, cosmetically finished shell that slips over the residual limb and does not require a glove

Hands		
Description	Product No.	
Passive	5xxxxx000	
Reinforced Passive	5xxxxx001	
Super Reinforced Passive	5xxxxx002	




Hand Ordering Example Replace xxxxx with hand no. Product No.... 5xxxx001 Hand No. 10901 510901001

Gloves		
Description	Product No.	
Micro Coated Vinyl	1yyyyyzz3	
PVC	1yyyyyzz1	
Expanded MCV	1yyyyyzz4	
Hard PVC (light prosthesis)	2yyyyyzz1	
Micro Coated Silicone	3yyyyyzz3	
Silicone	3yyyyyzz1	

Please reference Hosmer swatch kit (56760) for color options (zz)								
01	02	03	04	05	06	07	08	09
10	11	12	13	14	15	16	17	18

Glove Ordering Example				
Replace yyyyy with glove no.				
Replace zz with color choice				
Product No 1yyyyyzz3				
Glove No 10902				
Color 01				
110902013				



- A. Wrist circumference
- B. Metacarpal circumference
- C. Middle finger length
- D. Glove end circumference
- E. Length of glove

N	leasur	ement	s in mr	n	Size in	Approx. Age		Hand No.	Glove No.
Α	В	С	D	Е	inches	for use	Side	(xxxxx)	(ууууу)
100	100	40	140	152	4	0-18 mos	Left	N/A	10900L
100	100	40	133	151	4	0-18 mos	Right	N/A	10900R
115	115	45	150	173	4 ¹ /2	1-3 yrs	Left	10901	10901
113	117	46	150	173	4 ¹ /2	1-3 yrs	Right	10902	10902
113	123	47	150	186	5	3-6 yrs	Left	10903	10903
115	128	48	153	191	5	3-6 yrs	Right	10904	10904
128	150	58	175	218	5 ³ /4	5-9 yrs	Left	10905	10905
123	150	59	165	216	5 ³ /4	5-9 yrs	Right	10906	10906





715

WE Friction Wrist

· Friction wrist control permits rotation of the terminal device to any position

Description	Diameter	Weight	Product No.
WE-200N Child Size	1 ¹ /2"	1.5 oz	52147
WE-100N Infant Size	1 ¹ /4"	1 oz	52142

APRL Child Flexion Wrist (WF-50)

- · Installs between terminal device and wrist allowing manual positioning of the terminal device
- · Flexes at 0 and 36 degrees
- ¹¹/16" (2.7 cm) long, 1.5 oz (42 g)
- · Recommended for use with the WE-100N
- 1/2" 20 thread and stub

Product No.

52176





E-50F Friction Elbow

- · Friction elbows allow manual positioning of the forearm; friction is adjustable
- 2" diameter
- 7 oz weight

Product No.

50558

E-50 Elbow

- · Child size mechanical elbow with 10 locking positions
- 2" (5.1 cm) diameter; 7.5 oz (210 g)
- 1³/4" (4.4 cm) Elbow Axis to Socket End

Description	Product No.
E-50 Elbow with Outside Cable Exit	50556
E-50A Elbow with Inside Cable Exit	50557







PC-50 Polycentric Below Elbow Hinge

- · Useful in compensating for excessive tissue in the cubital fold area
- Hinges are fabricated from special tool steel which is heat treated at the gear tooth ends for strength and then plated to prevent corrosion
- Small, 4.5 oz (126 g) pair

Product No.

51635

MA-50 Step Up Below Elbow Hinge

- Intended for short, below elbow amputees who have a limited range of elbow flexion
- Multiplies the residual limb action of the forearm by a factor of 2
- Hinge is made of tool steel which is heat treated and plated for long life
- Small, 5 oz (140 g) pair

Product No.

51478

Northwestern Ring Harnesses

- Fabricated with high grade dacron web tape, stainless steel anchor ring to prevent binding or bunching
- 4-bar buckles, axilla loop, Y strap and an anterior suspension
- · Furnished ready to fit

Below Elbow Prosthesis Harness

- Child, ¹/2" (1.3 cm)
- Med-Small, ³/4" (1.9 cm)

Description	Product No.
NRH-75 (child)	51586
NRH-85 (small-medium)	51588


Above Elbow Prosthesis Harness

- Child, ¹/2" (1.3 cm)
- Med-Small, ³/4" (1.9 cm)

Description	Product No.
NRH-175 (child)	51602
NRH-185 (small-medium)	51603







222222222