

POST-OP DEVICES

Becker Orthopedic

- Coupler, Plate and Conical Block168

Flo-Tech®

- Transtibial Post-Op Components168-170
- Transfemoral Post-Op Components171-172
- SOFTIE, Waist Belts172

Hosmer

- Post-Op Adjustable BK Leg172

Orthomerica®

- PREP Adjustable AK Socket172
- PAL-Guard Post-Amputation Limb Guard173

Trulife

- Adjustable BK Pylon173

POST-OP SUPPLIES

Distal End Pads

- ALPS Padsil Distal End Pads173
- ESP Silicone Distal End Pads174
- Hermell Softeze Distal End Pads174
- Ossur Iceross Pads174
- Silipos Silopad Distal End Pads175

Socks

- Knit-Rite Post-Operative Socks175

Dressing

- Ossur Rigid Dressing176

Miscellaneous Supplies

- BK Felt Pad Kit176
- Lambs Wool176
- Post-Op Suspension Belt176
- Thermoformed Socket Backup Hardware Kits176

Becker - Coupler™ Post-Op

- Provides a more secure attachment to the socket
- Lightweight aluminum with Euro four hole design
- Consists of the basic plate plus contoured stainless steel attachment straps
- Straps are held in place using perpendicular grooves and the same screws that attach the coupler to the endoskeletal system



Product No. **P403010**

Coupler Plate

- Decreases fabrication time
- Deep undercuts and knurls eliminating the need to vacuum form the plastic under the plate
- Conical pocket significantly improves the interface with the residual limb model
- Can be removed and reused



Product No. **P403000**

Coupler Conical Block

- Incorporates all the features of the Coupler Plate, except it comes with molded urethane on the conical surface
- This model should only be used in wet laminations
- In many cases, the need for filler is eliminated
- Sealing and grinding are rarely necessary



Product No. **P403005**

Flo-Tech - Transtibial FLO-TECH-TOR™

- Protection from injury during the early stages of rehabilitation
- Easy access to bandages/dressings for attentive wound care
- Control of swelling and edema
- Maintenance of knee flexion and extension alignment
- Volume control; prepares the residuum or definitive prosthesis
- Includes thigh band, socket band, MPT strap, fork suspension, two foam distal end pads
- Belt (TTWBL/TTWBR) not included



Size	MPT	Distal	Range	7" Length Product No.		9" Length Product No.	
				Right	Left	Right	Left
X-Small	14-15"	13-14"	11-15"	TTFT7RXS	TTFT7LXS	TTFT9RXS	TTFT9LXS
Small	15-16"	14-15"	12-16"	TTFT7RS	TTFT7LS	TTFT9RS	TTFT9LS
Medium	16-17"	15-16"	13-17"	TTFT7RM	TTFT7LM	TTFT9RM	TTFT9LM
Large	17-18"	16-17"	14-18"	TTFT7RL	TTFT7LL	TTFT9RL	TTFT9LL
X-Large	18-19"	17-18"	15-20"	TTFT7RXL	TTFT7LXL	TTFT9RXL	TTFT9LXL
2X-Large	19-20"	18-19"	16-22"	TTFT7RXXL	TTFT7LXXL	TTFT9RXXL	TTFT9LXXL
3X-Large	21-22"	20-21"	19-23"	TTFT7RXXXL	TTFT7LXXXL	TTFT9RXXXL	TTFT9LXXXL

Transtibial Variable Control Supra Patellar Preparatory Socket (VCSPS™)

- Polyethylene supra condylar socket restores knee flexion
- Flexible distal reduction system reduces AP, ML, and circumference dimensions in the distal 3rd
- Clamshell design allows use as a protective device
- Can be used with the UFOS™ for controlled weight bearing
- Includes two foam distal end pads, one MPT 1" strap, distal strap and socket band



Size	MPT	Distal	Range	7" Length Product No.		9" Length Product No.	
				Right	Left	Right	Left
X-Small	14-15"	13-14"	11-15"	TTVC7RXS	TTVC7LXS	TTVC9RXS	TTVC9LXS
Small	15-16"	14-15"	12-16"	TTVC7RS	TTVC7LS	TTVC9RS	TTVC9LS
Medium	16-17"	15-16"	13-17"	TTVC7RM	TTVC7LM	TTVC9RM	TTVC9LM
Large	17-18"	16-17"	14-18"	TTVC7RL	TTVC7LL	TTVC9RL	TTVC9LL
X-Large	18-19"	17-18"	15-20"	TTVC7RXL	TTVC7LXL	TTVC9RXL	TTVC9LXL
2X-Large	19-20"	18-19"	16-22"	TTVC7RXXL	TTVC7LXXL	TTVC9RXXL	TTVC9LXXL

Post-Operative

Transtibial Universal Frame Outer Socket (UFOS™)

- Semi-rigid polypropylene frame promotes venous return and prepares the residual limb for a permanent prosthesis
- Permits controlled weight bearing
- Designed to resist knee flexion contractures while maintaining knee alignment
- Fits over the Flo-Tech-Tor or the VCSPS
- Includes: 2 ratcheting socket straps on 7", 3 ratcheting socket straps on 9", 4-hole Grace plate with two sets of screws, fork suspension strap and Velcro socket strap
- Belt (TTWBL/TTWBR) is not included



Size	MPT	Distal	Range	7" Length Product No.		9" Length Product No.	
				Right	Left	Right	Left
X-Small	14-15"	13-14"	11-15"	TTUF7RXS	TTUF7LXS	TTUF9RXS	TTUF9LXS
Small	15-16"	14-15"	12-16"	TTUF7RS	TTUF7LS	TTUF9RS	TTUF9LS
Medium	16-17"	15-16"	13-17"	TTUF7RM	TTUF7LM	TTUF9RM	TTUF9LM
Large	17-18"	16-17"	14-18"	TTUF7RL	TTUF7LL	TTUF9RL	TTUF9LL
X-Large	18-19"	17-18"	15-20"	TTUF7RXL	TTUF7LXL	TTUF9RXL	TTUF9LXL
2X-Large	19-20"	18-19"	16-22"	TTUF7RXXL	TTUF7LXXL	TTUF9RXXL	TTUF9LXXL
3X-Large	21-22"	20-21"	19-23"	TTUF7RXXXL	TTUF7LXXXL	TTUF9RXXXL	TTUF9LXXXL

Transtibial Rehab Package

- This kit includes
 - (1) Flo-Tech-Tor and (1) UFOS

Size	7" Length Product No.		9" Length Product No.	
	Right	Left	Right	Left
X-Small	RP7RXS	RP7LXS	RP9RXS	RP9LXS
Small	RP7RS	RP7LS	RP9RS	RP9LS
Medium	RP7RM	RP7LM	RP9RM	RP9LM
Large	RP7RL	RP7LL	RP9RL	RP9LL
X-Large	RP7RXL	RP7LXL	RP9RXL	RP9LXL
2X-Large	RP7RXXL	RP7LXXL	RP9RXXL	RP9LXXL
3X-Large	RP7RXXXL	RP7LXXXL	RP9RXXXL	RP9LXXXL



Transtibial Prep Package

- This kit includes
 - (1) VCSPS and (1) UFOS


Size	7" Length Product No.		9" Length Product No.	
	Right	Left	Right	Left
X-Small	PP7RXS	PP7LXS	PP9RXS	PP9LXS
Small	PP7RS	PP7LS	PP9RS	PP9LS
Medium	PP7RM	PP7LM	PP9RM	PP9LM
Large	PP7RL	PP7LL	PP9RL	PP9LL
X-Large	PP7RXL	PP7LXL	PP9RXL	PP9LXL
2X-Large	PP7RXXL	PP7LXXL	PP9RXXL	PP9LXXL



Transtibial Complete Package

- This kit includes
 - (1) Flo-Tech-Tor, (1) VCSPS, and (1) UFOS
- Belt (TTWBL/TTWBR) not included

Size	7" Length Product No.		9" Length Product No.	
	Right	Left	Right	Left
X-Small	TTCP7RXS	TTCP7LXS	TTCP9RXS	TTCP9LXS
Small	TTCP7RS	TTCP7LS	TTCP9RS	TTCP9LS
Medium	TTCP7RM	TTCP7LM	TTCP9RM	TTCP9LM
Large	TTCP7RL	TTCP7LL	TTCP9RL	TTCP9LL
X-Large	TTCP7RXL	TTCP7LXL	TTCP9RXL	TTCP9LXL
2X-Large	TTCP7RXXL	TTCP7LXXL	TTCP9RXXL	TTCP9LXXL



Transtibial Waist Belts

Side	Product No.
Right	TTWBR
Left	TTWBL

Transfemoral FLO-TECH-TOR-TF™

- Protective device used during the early period of rehabilitation
- Helps control swelling and prevent injury
- Provides non-ambulatory protection for TF amputations (use in conjunction with the UFOS-TF™ for ambulatory training)
- Includes two piece Neoprene® waist band, auxiliary suspension strap, Neoprene® socket band, two foam distal end pads
- Multi-layered belt system; belt not included



Size	Ischial Range	Distal Range	Product No.	
			Right	Left
Small	16-22"	13-17"	TFFTRS	TFFTLS
Medium	19-25"	16-20"	TFFTRM	TFFTLM
Large	22-28"	19-23"	TFFTRL	TFFTLL
X-Large	26-32"	24-28"	TFFTRXL	TFFTLXL

Transfemoral Universal Frame Outer Socket (UFOS-TF™)

- Semi-rigid polypropylene frame promotes venous return and prepares the residual limb for a permanent prosthesis
- Permits controlled weight bearing
- Designed to resist knee flexion contractures while maintaining knee alignment
- Fits over the Flo-Tech-Tor or the VCSPS
- Includes 2 ratcheting socket straps, Grace plate with two sets of screws, and 2" Velcro socket strap



Size	Ischial Range	Distal Range	Product No.	
			Right	Left
Small	16-22"	13-17"	TFUFRS	TFUFLS
Medium	19-25"	16-20"	TFUFRM	TFUFLM
Large	22-28"	19-23"	TFUFRL	TFUFL

Transfemoral Complete Package

- This kit includes (1) Flo-Tech-Tor-TF and (1) UFOS-TF
- Belt not included

Size	Ischial Range	Distal Range	Product No.	
			Right	Left
Small	16-22"	13-17"	TFCPRS	TFCPLS
Medium	19-25"	16-20"	TFCPRM	TFCPLM
Large	22-28"	19-23"	TFCPRL	TFCPLL



Post-Operative

Transfemoral Waist Belts

Size	<u>Product No.</u>	
	Right	Left
Small	SPS-TFCBS-RS	SPS-TFCBS-LS
Medium	SPS-TFCBS-RM	SPS-TFCBS-LM
Large	SPS-TFCBS-RL	SPS-TFCBS-LL



SOFTIE™ and SOFTIE-TF™

- Soft socket; protects incisions
- Non constricting; designed to resist knee flexion contractures while maintaining knee alignment

Size	Accommodates Lengths	Distal Range	<u>Product No.</u>	
			Right	Left
Transtibial	4 ¹ / ₂ -10"	11-23"	Softie-TT R	Softie-TT L
Transtibial	6 ³ / ₄ -10"	20-28"	Softie-TT RXL	Softie-TT LXL


Size	Accommodates Lengths	Distal Range	Product No.
Transfemoral	7 ¹ / ₂ -13"	13-23"	Softie-TF
Transfemoral	6 ³ / ₄ -10"	23-28"	Softie-TF XL



Hosmer - Post Operative Adjustable BK Leg (PPL-100)

- A post-operative leg which is aligned by adjusting a ball-socket component in the proximal portion and a similar component in the distal section
- A simple disconnect feature permits easy removal of the pylon and foot
- 1¹/₂" inner diameter


Product No. **51672**



Orthomerica - PREP™ Adjustable A/K Socket

- Lightweight, adjustable polypropylene and polyethylene AK socket for early fitting of the new amputee
- Narrow M-L (Ischial containment, NSNA) design for optimum fit and ambulation
- Circumference adjustable via Velcro straps
- Length can be trimmed to facilitate residual limbs
- Distal section contains adjustable, proximal component
- Hardware included to attach all conventional pylon systems

Size	Proximal Circ.	Mid Circ.	Distal Circ.	<u>Product No.</u>	
				Left	Right
Small	15-18"	14-16"	12-14"	4530	4535
Medium	17-21"	16-18"	13-15"	4531	4536
Large	20-24"	20-22"	16-18"	4532	4537



PAL-Guard™ - Post-Amputation Limb Guard

- The adjustability of the lightweight system, ease of donning and doffing, and the ability to remove and clean the Coolfoam™ liner makes this product especially effective for persons with recent amputations
- Coolfoam liner for increased comfort and decreased perspiration
- Lateral uprights help to align the limb and prevent knee flexion contracture
- Lightweight, durable copolymer shell
- Waist Size: One size fits all

Size	Inside Diameter of Distal Shell	Thigh Circumference	Product No.
Small	4"	13"	8401
Medium	5"	16"	8402
Large	6"	19"	8403



Trulife - Adjustable (Wedge) BK Pylon

- Includes side straps
- Weight as shown: 708 g
- Weight Limit: 220 lbs (100 kg)

Product No. **P017002000**



ALPS - Padsil™ (Silicone Distal End Pads)

- Soft medical grade silicone helps maintain total contact
- Measurements taken 6 cm from distal

Description	Diameter	Product No.
Small	14-19 cm	PDS-S
Medium	20-23 cm	PDS-M
Large	24-27 cm	PDS-L
X-Large	28-31 cm	PDS-XL
2X-Large	32-35 cm	PDS-XXL



Post-Operative

ESP - Silicone Distal End Pads

- Designed to provide added distal cushioning and support
- Made of pure, soft silicone gel
- 3/8" thick center, gradually tapering to outer edges

Description	Diameter	Product No.
X-Small	1 3/4"	ESP-DPXS
Small	2"	ESP-DPS
Medium	2 1/8"	ESP-DPM
Large	2 5/16"	ESP-DPL
X-Large	2 3/4"	ESP-DPXL
Plus Size 1	3 1/2"	ESP-DP1
Plus Size 2	4 1/2"	ESP-DP2
Plus Size 3	5 1/2"	ESP-DP3
Plus Size 4	6"	ESP-DP4



Hermell - Softeze Distal End Pads

- A functional distal pad that insures total contact
- Made of non-toxic urethane foam

Description	Diameter	Product No.
Small	2 3/4"	TC-2500
Medium	3 1/4"	TC-2501
Large	4"	TC-2502
X-Large	4 1/2"	TC-2503



Iceross® Pads

- Made from SenSil with 10mm thickness
- Ideal for pressure distribution in highly loaded areas, padding over bony prominences

Size	Fits Iceross Size	Product No.
38 mm (1 1/2")	12-20	I-610338
48 mm (1 7/8")	21-23.5	I-610348
58 mm (2 5/16")	25-28	I-610358
68 mm (2 11/16")	30-32	I-610368
78 mm (3 1/16")	34-36	I-610378
88 mm (3 1/2")	38-42	I-610388
98 mm (3 7/8")	45	I-610398



Silipos - Silopad™ Double Layer Distal End Pad™

- Gel cushioning pad dissipates a medical grade mineral oil (USP) to moisturize and soften residual limb skin
- Approximately 1/8" (3 mm) thick by 4" (10 cm) in diameter
- Dermatologist tested, hypoallergenic, washable and reusable

Product No. **12005**



Polyurethane Distal End Pads (not pictured)

- Reticulated, polyester urethane foam
- Intended for use as an interface in the rigid dressing
- May be sterilized; Color: Charcoal

Description	Product No.
3"	3MW29903
4"	3MW29804
5"	3MW29805
6"	3MW29806

Knit-Rite - Post Operative Dressing Retention Sock

- Soft, compressive fit with a belt suspension to cover and retain postoperative dressings
- Soft Lycra and CoolMax body for very low compression
- Universal right or left belt system for secure suspension

Post-Op Dressing Retention Sock			
Length	Size	Fits	Product No.
28"	4"	10-12"	1RP10428
32"	6"	14-18"	1RP10632
32"	8"	18-24"	1RP10832



Post-Operative

Post Operative Socks

- 3-ply weight, 97% Orlon Acrylic, 3% Lycra spandex
- Modified box toe; fleeced inside
- 4 socks to a box (sold individually)

Un-Sterilized Post-Op Socks			
Length	Top Width	Toe Width	Product No.
14"	5"	3 1/2"	1PL3BS14
18"	5"	4"	1PL3BS18
18"	6"	5"	1PL30518
18"	8"	6"	1PL32618
18"	9"	7"	1PL33718
Pre-Sterilized Post-Op Socks			
18"	5"	4"	1PL3S418
18"	6"	5"	1PL3S518
18"	8"	6"	1PL3S618
18"	9"	7"	1PL3S718



Ossur - Rigid Dressing (ORD)

- Immobilizes, protects and controls post-operative oedema in the transtibial residual limb immediately following transtibial amputation surgery
- Immobilization and oedema control are achieved by providing a firm dressing formed to match the volume and shape of the residual limb
- The ORD is made out of soft and pliable material and can thus be formed to the residual limb
- The dressing is made rigid by creating vacuum in the internal chambers of the product



Description	Size	Circ. 4 cm from Distal End of Residual Limb	Product No.
Ossur Rigid Dressing	Small	up to 42 cm	K-239100
Ossur Rigid Dressing	Large	43-65 cm	K-239101
Vacuum Pump for Rigid Dressing	—	—	K-239105

Note: All circumference measurements are taken over dressings and padding material

SPS - BK Felt Pad Kit

- Pre-patella, medial tibial crest, and lateral tibial crest pads
- 1/4" orthopedic felt

Product No. **209L00PAD**



Lambs Wool

- Recommended for use as a compressive dressing in immediate post-surgical prosthetic procedure

Description	Product No.
1 oz	582110
8 oz	582120



Post-Op BK Suspension Belt

- Prosthetic webbing with elastic suspension
- Cast attachment strap with buckle
- One size, specify left or right, up to 44"

Description	Product No.
Left	919050D
Right	919051D



Thermoformed Socket Backup Hardware Kits

- Vacuum forming attachment kit
- Should be used with double-walled socket constructions
- Universal 4-hole pattern

Description	Product No.
Round Kit	9427D
Square Kit	9432D

